
**IPPB Circuit**

LOT : 060101

GMP &amp; ISO-13485 Certified facility

**WARNING & CAUTION**

1. Discard if patient circuit or its components are damaged or if damage is suspected.
2. Obstruction or leakage of the device can result in serious injury or death.
3. Single patient use.
4. Patient should be constantly monitored whenever this device is in use.
5. Use test lung to check circuit and ventilator before apply to patient.
6. This product contains traces of phthalate. Although conclusive evidence on the risk of phthalate on human health has not been established, procedures leading to increased exposure of phthalate(s) in the treatment of children, pregnant or nursing women should be minimized as a precautionary measure. These types of procedures may include: hemodialysis in peripuberal males, male fetus and male infant of pregnant women, lactating women; exchange transfusion and total parental nutrition in neonates; multiple procedures in sick neonates; and massive infusion into trauma patients. It is strongly advised that the product be used according to trained medical personnel's guidance and the instructions for use before the marked expiration date.

**SPECIFICATIONS**

All the IPPB circuit contains the following parts:

1. Mouthpiece, 22mmOD\* 1pc
2. Short corrugated tube, 6 inch 22mmID\* 1PC
3. CV-2 IPPB control valve \*1pc
4. Corrugated tube, 39 inch 22mmID\*1PC
5. nebulizer\* 1PC
6. Short control line, 5.5mmOD,3.5mmID,20cm \*1PC and IPPB valve T- connector
7. Bird connector\*1PC
8. Long control line, 5.5mmOD\*1PC
9. Tube clamp
10. Oxygen tubing

**ASSEMBLY PROCEDURE**

1. Before using the IPPB, please check whether the tubes are properly assembled without any sign of occlusion or leakage. If present, return to the manufacturer.
2. Follow the IPPB instructional manual to check for tubing dysfunction.
3. For the construction with one pressure line, take out the T-connector and connect it on to the nebulizer, connect one end of short control line to IPPB control valve and the other to the T-connector.  
  
Connect the long line to the T-connector and the IPPB machine
4. For the construction with two lines, connect the oxygen supply tube to the nebulizer and to the IPPB machine (big adapter) and connect the long control line to the expiration valve and the IPPB machine (small adapter)
5. Follow your physician's instruction for nebulizer's medications

**REGULATION**


- ISO-5367 certified.
- All of the materials (EVA,PP,PVC,PE, and Silicone rubber) are design for its own purposes, any re-sterilization may damage the device and reduces safety for patient use.
- The nebulizer's aerosol size is approximately 0.5-5  $\mu\text{m}$
- ISO-5356 certified with 15mmID and 22mmOD adapter.

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#### GaleMed Xiamen Co., Ltd.