LABELLING LAYOUT SHEET

Title - IFU - Blue Line Ultra® Tracheostomy Inner Cannula Single Patient Use - (004/005/192) (Form Fill Seal - Soft Pack Only)

INFORMATION SHOWN:- Print layout on IFU. (To scale)

English - en	French - fr	German - de	Italian - it	Spanish - es	Portuguese - pt	Dutch - nl	Swedish - sv	Danish - da
Y	Y	Y	Y	Y	Y	Y	Y	Y
Norwegian - no	Finnish - fi	Greek - el	Japanese - ja	Czech - cs	Slovak - sk	Hungarian - hu	Polish - pl	Slovenian - sl
Y	Y	Y	N	Y	N	Y	Y	N
Estonian - et	Latvian - lv	Lithuanian - lt	Turkish - tr	Bulgarian - bg	Romanian - ro	Russian - ru	Croatian - hr	Chinese - zh
Y	Y	N	Y	Y	N	N	N	N

CCN. 004/005/192 Print Revision level - REV. 000

Size: Internal pages: 210mm x 148mm (A5), Folded to 105mm x 148mm (A6).

External pages: 315mm x 148mm, Folded to 105mm x 148mm (A6). (All dimensions have a tolerance of ± 2 mm)

Print Details: PMS Pantone Process Black U

Design Authority: Hythe

Material: 60gsm White Bond Paper


Description: Instruction Booklet, 16 leafs (66 pages), stitched twice.

NOTES:

- The layouts comply with the requirements of Corporate Policy CP30 Issue 3
- This artwork replaces 004/000/007 Rev 002 (L2089 Issue 2)
- Copy approval ref: L5306 Issue 000.
- Approval of this document is for final Format and Layout only, following the addition of the Certified Translations. It is the document originators responsibility to include and check that the correct Certified Translation have been used.

Last Change; PWR 4882. MIN H8376

- Addition of Estonian, Latvian, Turkish and Bulgarian.

SMI Labelling  Name: BILLIE BURVILL Date: 19-OCT-2012 <i>Signature held on master</i>	Design Authority RA Name: Date: <i>Signature held on master</i>	Marketing  Name: S. K. M. Date: 22/10/2012 <i>Signature held on master</i>
Compiled By - Robin Hibbitt  <i>Signature held on master</i>	Design Authority R & D / Operations Technology Name: Date: <i>Signature held on master</i>	Date Effective 23 OCT 2012
Compilation Date - 18 Oct 12		Sheet 1 of 33

smiths medical

bringing technology to life

Ref. L5306

ISSUE - 001

LABELLING LAYOUT SHEET

Title - IFU - Blue Line Ultra® Tracheostomy Inner Cannula Single Patient Use - (004/005/192) (Form Fill Seal - Soft Pack Only)

INFORMATION SHOWN:- Print layout on IFU. (To scale)

English - en	French - fr	German - de	Italian - it	Spanish - es	Portuguese - pt	Dutch - nl	Swedish - sv	Danish - da
Y	Y	Y	Y	Y	Y	Y	Y	Y
Norwegian - no	Finnish - fi	Greek - el	Japanese - ja	Czech - cs	Slovak - sk	Hungarian - hu	Polish - pl	Slovenian - sl
Y	Y	Y	N	Y	N	Y	Y	N
Estonian - et	Latvian - lv	Lithuanian - lt	Turkish - tr	Bulgarian - bg	Romanian - ro	Russian - ru	Croatian - hr	Chinese - zh
Y	Y	N	Y	Y	N	N	N	N

CCN. 004/005/192 Print Revision level - REV. 000

Size: Internal pages: 210mm x 148mm (A5), Folded to 105mm x 148mm (A6).

External pages: 315mm x 148mm, Folded to 105mm x 148mm (A6). (All dimensions have a tolerance of ± 2 mm)

Print Details: PMS Pantone Process Black U

Design Authority: Hythe

Material: 60gsm White Bond Paper


Description: Instruction Booklet, 16 leaves (66 pages), stitched twice.

NOTES:




- The layouts comply with the requirements of Corporate Policy CP30 Issue 3
- This artwork replaces 004/000/007 Rev 002 (L2089 Issue 2)
- Copy approval ref: L5306 Issue 000.
- Approval of this document is for final Format and Layout only, following the addition of the Certified Translations. It is the document originators responsibility to include and check that the correct Certified Translation have been used.

Last Change; PWR 4882. MIN H8376

- Addition of Estonian, Latvian, Turkish and Bulgarian.

SMI Labelling  Name: Bill Burvill Date: 19 OCT 2012 Signature held on master	Design Authority RA  Name: S J BURGESS Date: 23 OCT 2012 Signature held on master	Marketing Name: Date: Signature held on master
Compiled By - Robin Hibbitt  Signature held on master	Design Authority R & D / Operations Technology  Name: A HOWARD-GRIGG Date: 19 Oct 2012 Signature held on master	Date Effective 23 OCT 2012
Compilation Date - 18 Oct 12		Sheet 1 of 33

Folding Diagram

Finished Booklet A6

External pages: 315mm x 148mm, Folded to 105mm x 148mm (A6).


66

CCN.004/005/192 REV. 000 10/12

C.C.N. 004/005/192

Print Revision level - REV.000

Size: Internal pages: 210mm x 148mm (A5), Folded to 105mm x 148mm (A6).
External pages: 315mm x 148mm, Folded to 105mm x 148mm (A6).

English

en

Portex® Blue Line Ultra® Tracheostomy Inner Cannula, Single Patient Use

INSTRUCTIONS FOR USE

These instructions are intended for use with the following Portex® products:
100/850/060-100 Replacement smooth inner cannula x 2
100/851/060-100 Replacement fenestrated smooth inner cannula x 2

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

Please note that the code number of the Portex® Single Patient Use Inner Cannula indicates the size of Portex® Tracheostomy tube with which it should be used e.g.

Inner Cannula Code No.	For use with Tracheostomy Tube	Inner Cannula I.D.	Jackson Size (approx.)	Tracheostomy Tube Length 'C'
100/850/060 100/851/060	6.0 mm	5.0 mm	5	64.5 mm
100/850/070 100/851/070	7.0 mm	5.5 mm	–	70.0 mm
100/850/075 100/851/075	7.5 mm	6.0 mm	7	73.0 mm
100/850/080 100/851/080	8.0 mm	6.5 mm	8	75.5 mm
100/850/085 100/851/085	8.5 mm	7.0 mm	–	78.0 mm
100/850/090 100/851/090	9.0 mm	7.5 mm	–	81.0 mm
100/850/100 100/851/100	10.0 mm	8.5 mm	10	87.5 mm

NOTE:
1 Use only Blue Line Ultra® Inner Cannula with Blue Line Ultra® Tracheostomy Tubes.
2 The use of an Inner Cannula reduces the airway diameter by 1.5mm (excluding size 6.0mm tracheostomy tube which reduces the airway diameter by 1.0mm).

Description

A range of sterile single patient use replacement inner cannula intended to assist in the provision of a patent airway in the tracheostomised patient, available in both fenestrated (coloured red) and unfenestrated (clear) versions for use with the following Portex® Tracheostomy tubes.

100/800/060-100	100/811/060-100	100/831/070-090	100/860/060-100
100/802/060-100	100/812/060-100	100/832/070-090	100/870/060-100
100/810/060-100	100/813/060-100	100/833/070-090	

Instructions for use

- Following insertion of the tracheostomy tube and removal of the obturator, insert the inner cannula fully until a click is felt indicating correct positioning.

- If the inner cannula needs to be removed for cleaning disconnect from breathing circuit if attached.
- Remove the inner cannula by means of the ring pull on the inner cannula.
- Replace the soiled inner cannula with a clean inner cannula if required ensuring the inner cannula is correctly positioned by means of the 'click'.
- Re-connect the tracheostomy tube to breathing circuit (if used) and check for security.
- Clean the soiled inner cannula by means of the recommended cleaning instructions.
- Check the inner cannula prior to insertion/re-insertion and discard if kinked or damaged. Replacement inner cannulae are available (100/850 or 100/851).

Cleaning instructions

Daily cleaning of the inner cannula is recommended or whenever the inner cannula becomes contaminated by mucous or secretions. This period of time will vary according to individual patient needs.

Always remove the inner cannula from the tracheostomy tube for cleaning.

- Soak the inner cannula in sterile saline or mild detergent solution.
- Agitate the inner cannula in the solution to loosen any secretions. Using the cleaning brush provided with all Portex® arc shaped tracheostomy tubes (except 100/800 and 100/802) to help remove any dried or sticky secretions.
- After soaking and cleaning, rinse the inner cannula and cleaning brush with fresh sterile saline solution.
- Allow to dry naturally. Once clean and dry the inner cannula should be stored clean and dry and free from particulate matter.

Cleaning instructions – Tracheostomy Tube

The tracheostomy tube should only be cleaned with a sterile saline solution.

Precautions

- Unit packs must remain in the shelf carton until use to prevent damage to the inner cannula.
- Check the inner cannula prior to insertion and discard if kinked or damaged.
- The security of all breathing system connectors should be checked when established and frequently thereafter. Disconnection may be facilitated with the use of a disconnection wedge (100/555/100).
- Patients should be adequately humidified (e.g. using the Portex® Thermovent® T2 HME 100/570/022 or Thermovent® T 100/570/015) to minimise encrustation of the tracheostomy tube and/or inner cannula.
- The patency of the tracheostomy tube lumen must be assured by regular suctioning and regular cleaning of the inner cannula. Check routinely and replace as required to maintain a patent airway. Maximum recommended period of use 30 days.
- Care should be taken to ensure that the inner cannulae do not become kinked or damaged during cleaning, and that no kinked or damaged inner cannulae are re-inserted into the tracheostomy tube.
- The position and patency of the fenestration in the Portex® Fenestrated Tracheostomy tube and inner cannulae should be checked immediately after insertion to determine optimal position of fenestra and frequently thereafter to ensure that tissue granulation has not caused an obstruction.
- If Portex® Tracheostomy tubes and inner cannulae are used outside the hospital, the patient must be instructed by a healthcare professional in the safe use and handling of the product.
- Do not use solutions other than those indicated in the cleaning instructions to clean any part of the tracheostomy tube.
- Do not use any abrasive cleaning implements to clean the inner cannula. Use only the cleaning brush provided. Do not use this brush to clean the tracheostomy tube.
- Blue Line Ultra® Inner Cannula is designed and intended for single patient use. It should not be re-sterilised by any method.

4

5

Sheet 3 of 33

LA003/1-1

C.C.N. 004/005/192

Print Revision level - REV.000

Size: Internal pages: 210mm x 148mm (A5), Folded to 105mm x 148mm (A6).

External pages: 315mm x 148mm, Folded to 105mm x 148mm (A6).

Warnings

- 1 Do not lubricate the inner cannula as this may occlude the inner cannula and prevent the inner cannula being retained within the tracheostomy tube.
- 2 Contact with electrosurgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air, or ignite in an enriched oxygen environment (e.g. anaesthesia).
- 3 Never us a Portex® Fenestrated Tracheostomy tube for manual or mechanical ventilation unless a Portex® non fenestrated inner cannula is in place to prevent gas leakage through the fenestrations in the tracheostomy tube.
- 4 If a Portex® Fenestrated Tracheostomy tube is used with a Decannulation cap, either a fenestrated inner cannula must be used or if a non fenestrated inner cannula is present, it should be removed, and in both cases the tracheostomy tube cuff must be deflated. Failure to do so may lead to partial or complete airway obstruction.
- 5 Caution must be exercised when using a Decannulation Cap in conjunction with a cuffed fenestrated tracheostomy tube. The user must ensure that the cuff is completely deflated and that a fenestrated inner cannula is used, or that the inner cannula has been removed. Failure to do so may result in total or partial obstruction of the airway. Observe the patient for signs of respiratory distress and immediately remove the Decannulation Cap if required.
- 6 Do not lubricate the inner cannula as this may occlude the inner cannula-and prevent the inner cannula being retained within the tracheostomy tube.
- 7 If it is not possible to remove the inner cannula from the tracheostomy tube, do not try to forcibly remove. Both the inner cannula and tracheostomy tube should be removed together and replaced with a new tracheostomy tube and inner cannula.
- 8 Ensure that the inner cannula is both the appropriate diameter and length for the tracheostomy tube in use. This can be checked by comparing size markings on the tracheostomy tube flange and patient labels with those on the inner cannula unit pack label. Use this inner cannula only with tracheostomy tubes from the Portex® Blue Line Ultra® range.
- 9 Use of an inner cannula of the wrong diameter may cause difficulties in introduction or unnecessarily restrict gas flow. An inner cannula which is too long may protrude excessively from the outer tube, leading to tracheal damage or occlusion. Use of an inner cannula which is too short may lead to build-up of secretions which could cause infection or-blockage.
- 10 The use of an Inner Cannula reduces the airway diameter by 1.5mm (excluding size 6.0mm tracheostomy tube which reduces the airway diameter by 1.0mm).
- 11 Do not fit an inner cannula when using a tracheostomy tube in conjunction with a right-angled connector that incorporates a co-axial fresh gas tube which protrudes into the connector lumen (e.g. Norman elbow). The protruding fresh gas tube could occlude the inner cannula, restricting expiratory flow of gas, resulting in barotrauma/hypoxia.
- 12 During and after attachment of the breathing system to the tracheostomy tube connector, avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion.

Portex, Blue Line Ultra, Thermovent, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

©2012 Smiths Medical. All rights reserved.

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Portex, Blue Line Ultra, Thermovent designmerke, Portex designmerke, samt Smiths Medical designmerke er varemerker som tilhører Smiths Medical. Symbolet ® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.