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Instructions for Use

Asthma Monitor[®] AM1+

781090

Version 08.00

für/for Firmware \geq 7.1 - AM1+

für/for Firmware \geq 7.3 - AM1+ (BT Option)

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This manual has been composed with the assistance of Adobe® Creative Suite® programs Adobe InDesign®, Adobe Photoshop® and Adobe Acrobat®. The complete pages have been produced with Adobe InDesign® CS4 for PC and are also available as PDF files. Acrobat Reader is required in order to open these PDF files.

Trademarks

Adobe, Creative Suite, InDesign, Photoshop, and Acrobat are registered trademarks of Adobe Systems Incorporated in the USA and/or other countries.

JAEGER is a registered trademark. The original company JAEGER does no longer exist. eResearchTechnology GmbH is entitled to use the trademark JAEGER.

Copyright

All rights, also those of translations, are reserved. Written permission of **eResearchTechnology GmbH (ERT)** shall be obtained for the reproduction of the manual or an excerpt thereof by any means which includes printing, photostatic duplication, microfilms or any other process.

eResearchTechnology GmbH
Sieboldstrasse 3
97230 Estenfeld, Germany
Tel. +49 9305 720-9891
Fax +49 9305 720-7891
CustomerCare@ert.com
www.ert.com

Note

This document contains copyright information. All rights are reserved. It is not allowed to copy, duplicate or translate this manual in any other language without having the written approval of **ERT** before. **ERT** reserves the right for itself to alter the information included in this document without notice.

Names of individuals mentioned in the context of this manual are fictitious - any resemblance to living or deceased persons is purely incidental and not intended.

This instruction manual describes the standard device "**AM1+**" and the Bluetooth option "**AM1+BT**".

Whenever you see this symbol, the respective description will only be valid for the "**AM1+BT**".

Subject to technical modifications.

Copyright © 2014 by eResearchTechnology GmbH, Estenfeld4

Table of Contents

Declaration of Conformity.....	73
Notes on Safety in this Instruction Manual	74
Asthma Monitor® AM1+	75
What is the Main Difference?.....	76
How to Adjust your Asthma Monitor?.....	77
When to Measure?	77
Preparing for the Measurement	78
Switching on.....	79
Main Menu.....	80
Measurement	81
Results.....	83
Trend.....	84
Additional Parameter	84
Medication	85
Events (not available for AM1+BT)	87
Symptoms.....	88
Display List of Existing Measurements	90
Evaluate Measurements	91
Transfer of Measuring Data.....	92
Data Transfer via Bluetooth (not available for AM1+)	95
Set Display Language and Phone Numbers.....	100
Display Date and Time	102
Display Settings and Information.....	103
Cleaning.....	105
Disposal of Sensor.....	106
Batteries	107
Set Display Contrast.....	109
Error Checklist	110
General Safety Precautions.....	111
Literature.....	120
Notes on EMC according to EN60601-1-2.....	121
Return of Goods in Medical Institutions.....	125
Accessories	131
Technical Data	132
Warranty	134

Indications for Use

The Asthma Monitor AM1+ is an electronic measurement device to monitor the lung function (determination of the respiratory flows and volume) with high reproducibility wherever and whenever is a need of. The AM1+ measures the flow during expiration serving for the calculation of further parameters as FEV1.

The AM1+ is used to monitor the respiratory status of human beings in the areas asthma, chronic obstructive pulmonary disorder and in areas like occupational medicine, clinical trials and disease management.

The patient is informed of the results by numeric values for selected parameters (e.g. PEF, FEV1). Furthermore a visual control unit, displayed as a kind of traffic lights, allows an immediate indication of the measurement based on criteria defined by the patient's physician.

The device saves the results of a measurement (always with date and time) automatically in an internal database. In addition, a questionnaire functionality can be called up by the use of a software package (AMOS) to record e.g. the "Quality of Life" status. When enabled, the AM1+ can be programmed with a couple of questions, where the patient can select then from a couple of different answers. This information is also stored in the internal database and can be transmitted for evaluation to a standard PC using the software package AMOS.

The AM1+ is designed to replace ordinary peak flow meter, diary and pencil by a single system. Easy handling, sturdy and handy design allow the Asthma Monitor AM1+ being used almost everywhere: at work, at home, in school, for experts opinion, research or clinical trial purposes and in occupational medicine.


The AM1+ operates with serial data transmission whereas the AM1+ with BT option can be operated with serial and Bluetooth data transmission.

Caution:


Federal law restricts this device to sale by or on the order of a physician.

Dear Patients,


Asthma and COPD are obstructive diseases of the respiratory tract caused by airway inflammation and bronchial hyperreactivity. The diseases may be characterized by paroxysmal dyspnoea. In order to make an effective treatment plan, your doctor recommended to control your lung function parameters regularly. Regular measurements allow you to assess your condition and to take the correct measures in case of an emergency.


Bronchioles in patients suffering from asthma and COPD

1. Spasm of the smooth muscles

- plus
2. swelling of mucous membrane

- plus
3. production of viscous mucus


Declaration of Conformity

The original document of the Declaration of Conformity is deposited at Accompanying Documents.

Notes on Safety in this Instruction Manual

Following the ANSI recommendations (American National Standards Institute) for safety notes, specific passages of this instruction manual are clearly marked as safety notes.

Degree of Danger	Injury to persons	Damages to property	Use in Case of:
	X		DANGER indicates an immediate hazardous situation, which, if not avoided, may result in serious injury or death. Limited to extremely dangerous situations.
	X		WARNING indicates a potential hazardous situation, which, if not avoided, may result in serious injury or death.
	X	(X)	CAUTION indicates a potential hazardous situation, which, if not avoided, may result in minor or slight injury. Also used to indicate precarious procedures.


Additional icons used in this manual:

			Important information on operation and other useful information. Does not warn of dangerous or harmful situations.
			Tips, general information and information on operation.

Asthma Monitor® AM1+

The Asthma Monitor AM1+ is an electronic peak flow meter and allows for the measurement of your peak flow wherever and whenever you are in need of. Use your individual companion to measure your peak flow at home, at work or during leisure time. AM1+ measures peak flow during expiration (PEF) in litres per minute. The greater the lumen is, the quicker air passes through and the higher the measured value is. If the airways are constricted, the value is lower.

The traffic light system of the AM1+, which has to be set up by your doctor, allows you to quickly assess the measurement.

Mouthpiece with sensor*
(Type Rotary Flow Sensor)

Control keys
The symbols of these multi-functional keys appear in the display.

Display
Symbols, measurement results or
messages are displayed.

* in the following called sensor

What is the Main Difference?

With AM1+ data will be transmitted via cable to any available analog modem.

With AM1+BT data can be transmitted **wireless** via a **bluetooth** modem. In addition, the **bluetooth** version provides the possibility of transmitting data via mobile phone as **SMS** to the ERT Data Center.

How to Adjust your Asthma Monitor

Your doctor will adjust your AM to your **personal best value**. Your personal best value is defined as the value achieved under optimal (symptom-free) conditions.

 Everything OK, no danger!

 Attention!

 Danger! Emergency!

Depending on this value your doctor will enter values for **"Attention!"** and **"Danger! Emergency!"** and provide you with the respective instructions.

When to Measure?

Observe the instructions of your doctor.

Regular measurement is recommended:

- directly after getting up in the morning, at noon, and in the evening shortly before going to sleep
- always when you suffer from the following symptoms: cough, wheeze, chest tightness, shortness of breath
- after inhalation of your sympathomimetic agent, in order to control its effect
- repeatedly, if the values differ greatly
- repeatedly, if medication is changed
- repeatedly, in case of a bronchial infection

Preparing for the Measurement

Prior to starting the measurement, put the sensor into the AM as shown below. Following that remove the sealing cap from the sensor.

If the sensor (taken from the delivered "Set of three Replacement Rotary Flow Sensors" item no: 141355) is inserted in the medical product by a medical assistant, the medical assistant has to disinfect his/her hands before (the general hygiene standards valid for private practices and hospitals have to be observed).

If the sensor is inserted in the medical product by the patient, the notes on cleaning as described under "Cleaning" must be observed.


**The sensor is intended for single patient use, only.
A disinfection of the sensor makes it useless for operation.**

Switching on

There are two different ways to switch on the **AM**:

1. If a measurement shall be performed **immediately**, press  to switch on the device. The following appears on the display:


Now perform the measurement.

See "**Measurement**". If the measurement is not performed within the next 10 seconds, the Main menu will be displayed.

2. You can also switch on the device with pressing . A start screen appears and the results of the **last** measurement are displayed.

After five seconds the Main menu will be displayed.

Press  and you can "freeze" the display of results. The Main menu will only appear again if  or  is pressed. With  additional information on the last and saved measurements can be displayed. See "Show list of available measurements".

Press  to switch off the device.

Main Menu


The programs

"Measurement", **"Medication"**, **"Events"** and **"Symptoms"** can be selected in the Main menu of the **AM1+**.


The programs

"Measurement", **"Medication"**, **"Symptoms"** and **"Bluetooth"** can be selected in the Main menu of the **AM1+BT**.

Main menu **AM1+**:

Main menu **AM1+BT**:

Press **◀** or **▶** to select one of these programs.

Press **OK** to confirm your selection. The selected program will be displayed.

Press **OFF** to switch off the device.

Measurement

It is recommended NOT to perform measurements in direct sunlight, as the sensor could be damaged.

1. The following is displayed:

Take a very deep breath
and exhale as hard and
fast as possible

- Please note: The measurement has to be performed within the next 10 seconds.
2. Take the AM1+ in your right hand so that you can press the keys with your thumb.
 3. Inhale deeply. Press your lips **firmly** around the mouthpiece and exhale **as quickly and powerfully as possible**. **Attention:** Keep hands off the opening of the mouthpiece!
 4. Wait until you hear a beep. The value is displayed and the arrow points to one of the color-coded markers.

It is recommended to repeat the measurement with pressing **OK**, as it greatly depends on your cooperation.


The best value measured within a period of 10 minutes will be displayed.

The best value will be saved automatically.

The device must not be switched off.

An invalid trial is indicated as follows:

Press **OK** and repeat the measurement.

Bluetooth Version in SMS Mode:

If your Asthma Monitor is connected to a mobile phone which is reachable via bluetooth when the measurement is being performed, data of the current measurement will be automatically transmitted immediately after the measurement.

Please note: The Asthma Monitor only provides SMS. It is not checked whether the SMS was delivered successfully.

Results

The peak flow value "PEF*" is displayed in liters per minute. The "traffic light function" indicates the difference between the measured value and your personal best value and possibly necessary measures. An arrow points to the respective color.

-  **Arrow points to red: Danger! Emergency!**
 -  Repeat the measurement and follow the instructions of your physician or case manager.
 - 

-  **Arrow points to yellow: Attention!**
 -  Repeat the measurement and follow the instructions of your physician or case manager.
 - 

-  **Arrow points to green: Everything OK - no danger!**
 -  Follow the instructions of your physician or case manager.
 -  Press **OK** to perform another measurement.


*Depends on the measuring mode set by your physician.

Further available parameter:

FEV1 = Forced Expiratory Volume after one second

Trend

In addition to the displayed peak flow value the morning values of the past seven days are displayed.

Additional Parameter

Press  and  and an additional parameter will be displayed:


FEV1 Forced Expiratory Volume after one second.

Medication

Take the drugs prescribed by your doctor and observe your therapy plan!

1. Select "**Medication**" from the Main menu. The following is displayed.

2. Press the keys  and  to document the administration of up to four different drugs (indicated by "0" to "3"). Press **ESC** to abort the program.

0 =

1 =

2 =


3 =

Select one of the above-mentioned figures and confirm your entry with **OK**. The following is displayed:

Following that the Main menu will be displayed.

3. In order to control the efficacy of the applied drug, you can perform a post-measurement after about 10 minutes. Perform this measurement as described under "Measurement". The difference between pre- and post-measurement is displayed as follows:

Negative sign "-": Post-measurement shows worsened results.

Positive sign "+": Post-measurement shows improved results.

If the results of the post-measurement are not clearly improved, wait a few minutes and perform another measurement. Follow the instructions of your physician or case manager.

Events (not available for AM1+BT)

Almost every patient suffering from bronchial asthma or COPD knows specific situations where his/her condition normally gets worse, for example, under great effort (sports) or if he/she gets in contact with allergens such as pollen or wheats.

1. Select **"Events"** from the Main menu. The following is displayed:

2. Press the keys  and  to document up to four different events. (For example: "0" = none, "1" = sports, "2" = at work, "3" = pollen). Press **ESC** to abort the program.

0 =

1 =

2 =

3 =

Select one of the above-mentioned figures and confirm your entry with **OK**. The following is displayed:

Following that the Main menu will be displayed.

Symptoms


During the previously described events various symptoms such as cough, dyspnoea or sputum may appear. They can also be documented by indicating a degree of severity.

1. Select **"Symptoms"** from the Main menu. The following is displayed:


2. Press the keys  and  to enter one of four degrees of severity: "0" = none, "1" = slight, "2" = moderate, "3" = severe.
Press **ESC** to abort the program.
Select a figure and confirm with **OK**.
3. The following will be displayed:

Again, one of the four degrees of severity can be entered.
Select a figure and confirm with **OK**.

4. The following will be displayed:

Again, one of the four degrees of severity can be entered.

Select a figure and confirm with **OK**.

As soon as you have confirmed the degree of severity for "Sputum", the following will be displayed:


If any further questions should arise regarding operation of your AM, please refer to your physician or case manager!

Display List of Existing Measurements

As soon as your AM is switched on and  is pressed, PEF is displayed as result of the **last** measurement.

Press  and additional information on the last measurement will be shown.


Press  again and another previous measurement will be displayed.

 shows the previous measurement.  shows the next measurement. **ESC** and the results of the **last** measurement will be displayed again.

Evaluate Measurements


Visit your doctor at regular intervals (for example every three month). He will transmit the measuring data from the AM to his personal computer and will evaluate your measurements.

Memory almost full!
Please contact your
doctor

The indicator consists of a vertical bar on the left with three colored squares: green at the top, yellow in the middle, and red at the bottom. To the right of the text is a vertical bar with seven horizontal segments: five green, two yellow, and one red at the bottom.

The AM memory is almost full (80%). Visit your physician or case manager as soon as possible to have your measurements evaluated.

Memory full!
Please contact your
doctor immediately

The indicator consists of a vertical bar on the left with three colored squares: green at the top, yellow in the middle, and red at the bottom. To the right of the text is a vertical bar with seven horizontal segments: five green, two yellow, and one red at the bottom.

The AM memory is full. It is impossible to save further measurements. Visit your physician or case manager as soon as possible to have your measurements evaluated.

Transfer of Measuring Data

Measuring data can be transferred to a PC in your physician's office via modem.

Please contact your physician or case manager.

For safety reasons it is very important that while the device is connected to the modem by a cable, the subject must not approach the mouthpiece nor perform any measurements. Only IEC 60950-1 compliant computers are allowed to be connected to the device.

Data Transfer via Modem

Connect the modem to the telephone network as described in the Instruction Manual of the modem.

Connect the modem to the AM1+ which is switched off. The modem menu appears on the display of the AM1+. Press **◀** and **▶** to select either **"Data Transfer via modem"** or **"Modem Test/Initialization"**.

Data Transfer via Modem

"Modem Test/Initialization"

Confirm your selection with **OK**. The selected program will be executed.

In the program **"Data Transfer via Modem"** the following steps will be performed and the corresponding messages displayed:

1. "Initialize modem"
2. "Dial telephone no."
3. "Connection to host"
"Transferring data"
"Please wait"
4. "Transfer OK"

If an error occurred during data transfer, the following message will be displayed:

5. "An error occurred!"

In such a case perform an initialization by starting the **"Modem Test/Initialization"** program.

In the "**Modem Test/Initialization**" program a maximum of 10 preselected telephone numbers can be dialed, until a telephone connection is made. See "Set Display Language and Phone Numbers".

The following steps will be performed and the corresponding messages displayed:

1. "Initialize modem"
2. "Testing no. 01.....10"
3. "Connection to host"
"Transferring data"
"Please wait"
4. "Transfer OK"

If a connection is not possible, the following will be displayed:

5. "An error occurred!"

In this case contact your physician or case manager.

Data Transfer via Bluetooth (not available for AM1+)

The AM1+BT provides the possibility of transferring data wireless

- a) as **SMS** to a bluetooth (BT) mobile phone or
- b) to a BT modem.

Data Communication via SMS


Before data can be transferred via **SMS**, the Asthma Monitor and the mobile phone must be **paired**.

Step 1

Activate the "**Bluetooth**" option on your mobile phone to ensure that it is identified by other bluetooth devices like the AM1+BT.

Step 2

Simultaneously press the keys  and  of the Asthma Monitor for about 2 seconds.


Step 3

Press the **ESC** key **twice**.


The Bluetooth menu appears on the Asthma Monitor.

Step 4

Activate the Bluetooth mode by pressing **I/O**.

Step 5

Press **OK** to start the search for bluetooth devices in your vicinity. The message "Search for BT Device" appears on the screen.


The maximum distance between Bluetooth receiver and Asthma Monitor is 5 m.

Step 6

The following screen displays the devices found and their name. Select the name of your mobile phone and confirm with **OK**.

Example:


Step 7

The Asthma Monitor (including serial number) appears on the display of your mobile phone. Confirm the connection on your mobile phone and enter there the following PIN number for authentication: **1609**

Step 8

In order to start the SMS mode, select the option **"SMS"** on your Asthma Monitor and confirm with **OK**.

Pairing was successful, if this symbol is displayed in the Main menu.

Data Communication via Bluetooth Modem

If data are to be transferred via bluetooth modem, Asthma Monitor and bluetooth modem must be paired - analogous to data transfer via SMS.

Follow Step 1 to 7 as described before. However, select "**Modem**" as data transfer mode (Step 8).


Please note that the destination number of the data center must be saved on the Asthma Monitor for data transfer via **SMS** and **Modem**. This number is saved via AMOS Software on the personal computer.

Disconnect an Existing Connection

If you want to disconnect the currently existing connection, press **I/O** .

Example:


Main Menu of "AM1+BT", if the AM is Paired with a Mobile Phone (SMS Mode)


Select this symbol to start the transfer of data which have not yet been sent via SMS (1 SMS per measurement).

Main Menu of "AM1+BT", if the AM is Paired with a Modem (Modem Mode)


Select this symbol to start data transfer via modem (analog data transfer protocol) to the data center.

Set Display Language and Phone Numbers

The languages German, English, French or Spanish can be set for the display.


When the AM1+ is switched off, press both keys simultaneously.

The currently set language appears on the display.

Press **▼** or **▲** to select the desired language and confirm with **OK**.
Press **ESC** to abort the language setting.

Up to 10 telephone numbers which shall be dialed via a modem can be displayed and changed. See "Transfer of Measuring Data".

Identification digit 01 to 10


The telephone number can consist of 25 digits at the maximum. The digits 0123456789!*# as well as blanks can be used.

Press **▼** or **▲** to select an identification digit. Confirm your selection with **OK**. Now the telephone number can be changed by pressing **▼**, **▲**, **OK** and **ESC**. A "?" indicates the last position which can be edited. If the "?" is confirm with **OK**, the dial mode for the just displayed/changed telephone number can be defined.

Selected dial mode


Press **▼** and **▲** to select either the dial mode "Tone" or "Pulse". Confirm your selection with **OK**. The following may be displayed:

Press **ESC** and the BT menu will appear.


Display Date and Time

Date and time can be displayed.

When the AM1+ is switched off, press both keys simultaneously.


Date and time are displayed.

Press **ESC** to abort the date/time display.

Date and time can be set when your AM1+ is connected to the PC of your physician or case manager.

Display Settings and Information


While date/time is displayed, press both keys simultaneously.


Battery voltage
Number of recorded measurements
Number of recorded medications, events and symptoms


Press ▼ and subsequent displays can be shown. Press **ESC** to abort the program.


Patient ID


Threshold yellow/green (L/min)
 Threshold red/yellow (L/min)
 Information for service engineer, only
 Daylight-saving date

Alarm times to remind the patient to perform a measurement.

Version and date (software)
 Serial number (device)
 Indication of measurement mode set for this device.
 PEF = Peak Flow
 Your physician can also set a different measurement mode:
 FEV1 = Forced Expiratory Volume after one second

Current number, date and time of the last measurement
 PEF of the last measurement
 FEV1 of the last measurement

[Down Arrow] to display the previous measurement.
 [Up Arrow] to display the next measurement.

Cleaning

It is recommended to clean the sensor of the AM at regular intervals.

Sensor Cleaning:

Release and remove sensor, rinse it with distilled water in which a mild cleansing agent is soluted.

Get rid of the remaining water by shaking the sensor. Air-dry and reinsert the sensor.

The sensor is intended for single patient use, only. If the AM is passed on to another patient, the AM surfaces have to be cleaned and disinfected with a disinfectant*¹. See next page for information on disposal of sensor.

*¹ We recommend combined cleansing and disinfecting agents with the basic substance quaternary ammonium compound, for example "Cleanisept Wipes" by Dr. Schumacher GmbH.

Disposal of Sensor

The sensor is intended for single patient use, only. Dispose of used sensors if the AM is passed on to another patient.

Disposal:

It is absolutely to avoid that patient or medical assistant or the sensor becomes contaminated with sputum during disassembly of the disposable sensor. Therefore, release and remove the sensor by pulling the disposable sensor downwards (see picture).

Batteries

To operate the AM, insert the three 1.5 V batteries (Micro AAA type) into the battery box at the rear side of the device.

Pay attention to polarity (+/-)!

Due to functional and safety reasons, batteries must not be removed if the AM is connected to the PC.

Warning
Change batteries!

If this message is displayed, change the batteries of your AM within the next 8 days. The data remain saved. If the batteries are not changed within the next 8 days, the AM will then automatically be switched off as soon as the batteries are exhausted.

Remove the batteries, if the AM is not used for a longer period of time. All of the data saved will be kept even if the batteries are removed. The internal clock will stop after about 5 days. As stated in the Error Check-list, please visit your doctor if the date and/or time on the device is incorrect.

After having the batteries exchanged, close the battery compartment with the lid and keep it closed.

Set Display Contrast

The display contrast can be set when the battery compartment is open. It is not necessary to remove the batteries.

Screw to adjust display contrast

Error Checklist

Error Description		Reason	Action
No response during power ON	A.	AM batteries are empty	Replace batteries
	B.	Batteries are not inserted correctly	Insert batteries correctly
Clock is not working correctly			Consult your physician or case manager
Result of measurement is questionable	A.	Flow Sensor is not inserted correctly	Insert Flow Sensor correctly
	B.	Flow Sensor is dirty	Clean Flow Sensor according to cleaning instructions
	C.	Flow Sensor is faulty	Replace Flow Sensor

If the proposed actions do not lead to perfect recovery of the AM's normal functionality, please contact your doctor or case manager or the ERT Hotline.

General Safety Precautions

This instruction manual is regarded to be a part of the instrument, and should always be kept on hand.

The instruction manual describes the present state of the device/system including software and accessories with regard to the fundamental requirements of the MDD 93/42/EEC. Exact adherence to the instructions issued is a prerequisite for perfect and intended functioning of **ERT** instruments.

Deviation from Intended Use

Any non-observance of the procedures (such as preparing for the measurement and methods, disinfecting procedures, use of accessories and replacement parts etc.) described in the Instructions for Use results in a deviation from intended use.

In case of a deviation from intended use the operator/user has to supply proof of meeting all corresponding fundamental requirements. This is possible by performing a corresponding conformity assessment procedure within in-house manufacture (see § 12, paragraph 1 last sentence of MPG (= Medizinproduktegesetz/ Medical Devices Act)).

The operator/user is, however, not only responsible for performing the conformity assessment correctly but is also completely liable for defective products - i.e. the operator/user is not only liable for his/her modification of the medical product.

ERT only guarantees for the safety, reliability and functioning of the device if:

- installation, extension, modifications, and repairs are exclusively carried out by personnel authorized for these tasks by ERT.
- the ambient conditions at the place of installation are suitable for the device.
- the device is used according to the instruction manual.
- the device is not damaged. After unpacking, the device should be checked for damages. If the device is damaged, do not use it, but return it to get a replacement device.

Patient Safety

The AM is a device with an internal power supply (battery operation - no direct connection to mains power supply).

The physician must not touch the patient and the battery contacts of the AM or the computer at the same time.

Additional equipment connected to medical electrical equipment must comply with the respective EN or ISO standards (e.g. EN 60950 for data processing equipment). Furthermore, all configurations shall comply with the requirements for medical electrical systems.

Anybody connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems.

Attention is drawn to the fact that local laws take priority over the above mentioned requirements. If in doubt, consult your local representative or the technical service department.

Radiated Interference

The device meets the regulations according to EN 60601-1-2 regarding the interference radiated and received. The device should not be installed in the vicinity of high-frequency devices, X-ray equipment and motors or transformers with high installed power rating since electric or magnetic interference fields may falsify the result of measurements or make taking measurements impossible. Due to this, the vicinity of power lines is to be avoided as well. Existing environmental interferences may cause deviations of the measuring values without impairing the device's function. Therefore, it is recommended to keep a distance of about 2 meters from possible error sources when using the device.

More detailed information can be found in the EMC tables of the instruction manual of your device/system.

Ambient Conditions

The AM must not be operated in rooms or in the presence of flammable anaesthetic mixture with air or flammable anaesthetic mixture with oxygen or nitrous oxide.

The AM has to be effectively protected against moisture. Therefore, it is required that the AM is always stored in the black bag. The device corresponds to IP20 degree of protection. Measurements in the rain or in the shower are not allowed.

CAUTION

Measuring Mode

As the combination with an IEC 60950-1 proofed PC or modem can lead to a summation of the leakage current, the AM must not be connected to a PC or modem during the measurement. When the measuring mode is activated, the battery cover has to be closed.

Should the measuring values of the AM be changed after a longer period of use, a new sensor should be used.

During operation at an ambient temperature of 40°C, the display can warm up to 42.8°C.

CAUTION

Interfaces

The AM must only be connected to a PC or to a modem that correspond to the respective standards for data processing devices (e.g. IEC 60950).

If the connection cable is defective, it has to be replaced by a new one. The physician must not touch the patient and live parts at the same time.

CAUTION

Medical Supervision

A qualified physician has to reassess all AM measurements. An interpretation by the AM is only significant if it is considered together with other clinical findings.

CAUTION

Contraindications and Possible Adverse Effects

According to "ATS/ERS TASK FORCE: STANDARDISATION OF LUNG FUNCTION TESTING" (ERS Journals Ltd 2005) performing lung function tests can be physically demanding for minority of patients. It is recommended that patients should not be tested within 1 month of a myocardial infarction. In rare cases spirometry testing can lead to syncope due to extensive exhalation.

Putting the Device into Operation

Temperature changes may give rise to condensation in the device. Consequently, the device has to adapt to the ambient temperature before putting it into operation.

Cleaning and Hygiene

Prior to measuring a patient, his/her anamnesis is to be checked in order to avoid a contamination of the device and a resulting cross contamination of the next patient.

The AM must not be connected to a PC or modem when it is being cleaned.

Referring to humidity and water which may get into the device, AM corresponds to the safety degree IP 20. This means that the AM can be cleaned with a damp (not wet) cloth which does not produce fluff.

More detailed information can be found under "Cleaning". Chemicals required for operation or care of the unit must always be stored, prepared, and made available in specially marked vessels to prevent misuse.

Maintenance

The device doesn't require to perform preventive inspection, maintenance and calibration.

No part of the AM should be replaced by the subject/ Doctor. Damaged and frayed plugs, receptacles, housing and the display glass (if available) should be replaced by an authorized specialist or engineer of **ERT Customer Care**.

If the device is opened by an unauthorized person, the warranty will be void.

Before turning ON the device, always check whether the device is free from defects.

Immediate maintenance is necessary, if:

- the display glass bursts or breaks: **Caution: risk of injury**
- the device has been mechanically stressed (e.g. impact, damage to the housing)
- the connection cable is defective. The connection cable has to be replaced with a new one.
- liquid got into the device
- coverings have fallen off.

Children should not get in contact with disposables, accessories and packing material and cleaning and disinfection substances.


Recycling

Keep the national regulations of your country when disposing the medical product and its accessories. Improper disposal of the device and/or its accessories may result in serious environmental hazard.

Batteries

The batteries must not be exchanged if the AM is connected to a PC or modem.

Graphical Symbols

The typeplate can be found at the rear side of the device. To see the typeplate, open the lid of the battery compartment.

Note Instructions for use

Caution!

General warning sign

Stand-by.
Device ON/OFF switch

Date of production

Manufacturer

Applied Part of Type BF

Disposal in compliance with WEEE

Barometric pressure limits

IP20

Protection against intrusion of solid objects with a diameter $\geq 12,5\text{mm}$; no protection against liquid intrusion

SN

Serial Number

Battery type AAA

CE sign with code number of the Notified Body. The certified quality assurance system of **eResearchTechnology GmbH** corresponds to the international standard of ISO 13485.

Rx only

CAUTION: FEDERAL U.S. LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

Possible source of interference

The radiation intensity of the bluetooth module is below the SAR limits which are demanded by the EC Directive 1999/519/EEC.

Attention: Only devices recommended by the manufacturer must be used for sending SMS via a bluetooth modem.

Attention: Sending SMS via a bluetooth modem is only possible if appropriately configured devices are used.

Attention: For safe data transfer, the distance between AM and receiver must not exceed 5 meters.

“Approved in accordance to R&TTE directive transmitter module marked by CE, manufactured by MITSUMI incorporated OEM product.”

“This device complies with Part 15 of the FCC Rules. Operation is subjected to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including that may cause undesired operation.”

This device contains FCC-ID POOWML-C46.

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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EN 60601-1 Medical Electrical Equipment Part 1: General Requirements for Safety

EN60601-1-2 Medical electrical equipment

Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests

DIN VDE 0100 Part 710 (IEC 60364-7-710) Electrical installations of buildings - Part 7- 710: Requirements for special installations or locations - Medical locations

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The safety precautions and operational procedures indicated in this chapter refer to Germany. Different regulations and standards may apply in other countries.

Notes on EMC according to EN60601-1-2

The use of accessories not recommended by ERT may result in an increased electromagnetic radiation or a reduced interference immunity of the AM.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The AM is intended for use in the electromagnetic environment specified below. The customer or the user of the AM should ensure that it is used in such an environment.		
Emission tests	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1 Class B	The AM uses RF energy only for its internal function. Therefore, its RF emissions are not likely to cause any interference to nearby electronic equipment
RF emissions CISPR 11	Group 1 Class B	The AM is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The AM is intended for use in the electromagnetic environment specified below. The customer or the user of the AM should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV 8 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	Not applicable Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	Not applicable Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_T (> 95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 s	Not applicable Not applicable Not applicable Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AM requires continued operation during power mains interruptions, it is recommended that the AM be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.
NOTE	U_T is the A.C. mains voltage prior to application of the test level.		

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The AM3 is intended for use in the electromagnetic environment specified below. The customer or the user of the AM3 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the AM3, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter: Recommended protection distance:
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	$d = 1.17 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1.17 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2,33 \sqrt{P}$ for 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) acc. to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflections from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AM3 is used exceeds the applicable RF compliance level above, the AM3 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the AM3.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distance between portable and mobile RF Communications Equipment and the AM

The AM is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the AM can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the AM as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17 \sqrt{1/V} * \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{m/V} * \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.33 \sqrt{m/V} * \sqrt{P}$
0.01	Not applicable	0.12	0.23
0.1	Not applicable	0.37	0.74
1	Not applicable	1.17	2.33
10	Not applicable	3.7	7.37
100	Not applicable	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not be apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Return of Goods in Medical Institutions

Recommendations for action -

for all staff members having contact with potentially contaminated returns.

Returns or returned goods are all products returned to the producer or the supplier, irrespectively of whether or not they have been used, e.g. due to complaints, for repair or maintenance. Those products might have had contact with biological substances or highly active pharmaceuticals (e.g. cytostatics, radioactive medicines) and could be contaminated by them. If in doubt, the goods to be returned should be treated as contaminated products.

Due to infectious agents, pathogens or pharmaceuticals contaminated goods pose a potential hygienic risk for all persons having contact with the returns.

This leaflet shall minimize the potential hygienic risk when handling returned goods. Among other things, this information is based on the legal standards of the Biostoff-Verordnung (= Biological Substances Regulation) and the Employment Protection Act (both valid in Germany). A more detailed reference list of applicable rules and regulations can be obtained from the BVMed (info@bvmed.de).

In order to protect your and our employees who handle contaminated parts and to optimally examine such parts, you should consider the following:

1. Assessment of returns before reshipment

Irrespectively of whether or not a contamination risk is known, products which had direct or indirect contact with biological working substances (e.g. blood, secretions or other body fluids) or with highly active pharmaceuticals (e.g. gloves of the clinical or surgical staff) or hazardous substances (provocative substances) should be regarded as potentially dangerous for health.

According to the regulation on biological substances, substances of risk group 3 can cause severe illness in humans and pose a serious risk for staff members (e.g. tuberculosis or hepatitis). Substances of risk group 4 cause severe illness in humans and pose a serious risk for staff members (e.g. Ebola or smallpox). Regarding national and international regulations for the transport of potentially infectious substances (ADR, IATA-DGR) the **risk groups 2** (e.g. staphylococcus aureus) **and 3** defined by the Biological Substances Regulation are classified as **Transport Category B**.

Products which are potentially contaminated by biological working substances of risk group 4 according to the regulation on biological substances as well as products which are potentially contaminated by pharmaceuticals posing a serious risk to health (e.g. X-ray contrast agents and cytostatics) must not be returned to the producer.

The following applies for products which are potentially contaminated by biological working substances of risk group 3 + 2 according to the regulation on biological substances and which are thus classified as Transport Category B: Please contact ERT **before** returning the goods and observe the regulations for the transport of dangerous goods.

2. Cleaning

If the products had contact with biological working substances (for example blood or other body fluids), they have to be cleaned and disinfected in a combined cleaning and disinfection procedure, unless a differing agreement has been made with ERT.

As a rule, the products also have to be cleaned in order to minimize adhesions and bacterial contamination. Unless body fluids, body tissue or contrast agents etc. have caused the product defect, deposits and adhesions should be removed carefully without damaging or altering the product, if possible. For this, refer to the notes on the preparation of products.

Products which have been contaminated by highly active pharmaceuticals have to be cleaned appropriately with tap water.

If in doubt, contact ERT for agreeing on the further actions!

3. Disinfection/Sterilization

After cleaning, the products have to be disinfected and/or sterilized (only if permitted for this medical product) in order to exclude endangering of your and our employees.

If in doubt or in case of suspected material incompatibility please consult ERT.

4. Packaging

To avoid any contamination, the cleaned and disinfected product has to be packed as follows:

a) Put it into a sealable primary packing.

Parts with sharp edges need to be packed particularly safe.

b) Put the primary packing in a waterproof secondary packing (if possible use hard packing material).

c) Pack the secondary packing with a neutral packing material.

For "Packaging and labeling of non-contaminated products": see point 6.

5. Labeling

If a concrete risk of infection (e.g. HIV, Hepatitis B or C) is known to be present, this risk has to be noted on the packing of the returned goods and/or in the accompanying documents.

For "Packaging and labeling of non-contaminated products": see point 6.

6. Packaging and labeling of non-contaminated products

If the procedures described under point 2 and 3 are not applied, the contaminated product has to be returned in a combined packing complying with the packaging instruction P 650 ADR after contacting ERT, if necessary. Proceed as follows:

- a) Put the product into a liquid-tight, sealable packing (e.g. tearproof plastic bag) (primary packing)
- b) Put the primary packing into a (if possible, liquid-tight) protective packing (secondary packing); for liquid materials insert an adequate amount of absorbing material between the primary and the secondary packing.
- c) Pack the secondary packing with an additional outer packaging (padded envelope or cardboard box).
- d) Label the outer packaging with the corresponding UN no. 3373 for diagnostic or clinical samples and add the note:

"Biologischer Stoff, Kategorie B/Biological Substance, Category B"

7. Dispatch

Please note that non-decontaminated returns with suspected pathogens of risk group 3 are excluded from mailing. Diagnostic or clinical samples of UN no. 3373 which have been packed according to packing instruction P 650 are not subject to any further regulations on the transportation of dangerous goods and may be transported by a forwarding agent or a parcel service.

For this purpose use the accompanying shipping documents of the forwarding agent/carrier containing the corresponding valid transportation regulations, e.g. GGVSE (Road and Railway Dangerous Goods Regulation).

The product is then dispatched to the address indicated by the manufacturer.

Address:

eResearchTechnology GmbH

Sieboldstrasse 3

97230 Estenfeld, Germany

Tel. +49 9305 720-9891

Fax +49 9305 720-7891

CustomerCare@ert.com

Certificate of Hygiene

This certificate must be attached to ANY product complaint, ANY return of medical products and accessories, ANY repair order and ANY return of studies.

Name of product: _____

REF (ERT item no.): _____

LOT (batch no.): _____

It is herewith confirmed by signature that (please mark appropriate box):

- the enclosed medical product had **no** contact with blood or other body fluids so that it is **hygienically safe**.
- the enclosed medical product had contact with blood or other body fluids during its use. The product has been cleaned and decontaminated as follows.
 - Disinfection by wiping all accessible surfaces with
Disinfectant: _____
Concentration: _____
Reaction time: _____
- Other procedure (please indicate):
- Steam sterilization (3 minutes at 134 °C or 15 minutes at 121 °C)
- the enclosed medical product could **not** be decontaminated.
Reason:

Send returns to the following address: Sender's signature and address:

eResearchTechnology GmbH
Abteilung Wareneingang
Sieboldstrasse 3
97230 Estenfeld
Germany

Accessories

Rotary flow sensor kit

no. 141 355

Three rotary flow sensors for Asthma Monitor AM1+.

Sealing cap for rotary flow sensor

no. 835 560

Bag

no. 141 360

Handy bag made of durable nylon with shoulder belt and belt loop for the AM1+.

Color: black

Address:

eResearchTechnology GmbH
Sieboldstrasse 3
97230 Estenfeld, Germany
Tel. +49 9305 720-9891
Fax +49 9305 720-7891
CustomerCare@ert.com

Technical Data

Principle:

Determination of flow and volume via exchangeable infrared rotary flow sensor.

Range:

PEF 60 - 840 liters/minute

FEV1 0.5 - 8 liters

Accuracy:

PEF $\pm 5 \%$ or ± 20 L/min

FEV1 $\pm 3 \%$ or ± 0.05 liters

Storage capacity: 1200 measurements and 600 diary entries

Power supply: Three 1.5 V batteries, Type AAA
service life: about 150 days

Dimensions

Length x width x height 112 x 82 x 34 mm

Weight 145 g (batteries included)

Ambient conditions

Temperature $+10^{\circ}$ C to $+40^{\circ}$ C

Relative humidity 15 to 95%, non-condensing

Barometric pressure 700 to 1060 hPa

Transport and storage conditions

Temperature -20° C to $+50^{\circ}$ C

Relative humidity 15 to 95%, non-condensing

Barometric pressure 600 to 1200 hPa

Moisture protection IP20

Medical classification Active Medical Device Class IIa

Applied part: Type BF (whole device)

Protection class Battery-powered device

Mode of operation Continuous operation

Max. resistance

70 Pa/L/s at 14 L/s

Interface

RS-232 (9600 Baud, 8 bit, 1 Stop bit, no parity)

Bluetooth (optional)

The expected operational lifetime of the AM1+ is 5 years.

Asthma Monitor® AM1+ corresponds to the recommendations of ATS/ERS.

Warranty

The device has been conceived carefully and has been tested thoroughly before shipment.

Warranty is given for a period of 12 months from the date of purchase proof of which has to be furnished by the bill and/or warranty registration card.

Within this period ERT is obliged to replace defective parts free-of-charge. The period of warranty is thereby not extended. ERT automatically acquires possession of replaced parts. This warranty excludes extended claims especially those for redhibition, diminution or damages.

This warranty shall not apply if defects are due to misuse, neglect of instructions comprised in the manuals or repairs and/or alterations of any party other than ERT. Batteries and parts subject to wear and tear are not included in this warranty. If a claim is illegitimate, ERT is entitled to charge the costs for testing and shipment.

Remember to ship devices together with the warranty form and/or purchase bill to our authorized distributor.

Address:

eResearchTechnology GmbH
Sieboldstrasse 3
97230 Estenfeld, Germany
Tel. +49 9305 720-9891
Fax +49 9305 720-7891
CustomerCare@ert.com

Warranty Registration Card

Asthma Monitor[®]

AM1+

Date of purchase:

Distributor/manufacturer(stamp/sign):

Stamp and sign of the distributor/manufacturer confirm the validity of this warranty card. If date of purchase, stamp and sign of the distributor are missing, the card is to be regarded as invalid.

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The peak flow value* (PEF) is displayed in liters/minute. The "traffic light function" indicates the difference between the measured value and your personal best value and possibly necessary measures. An arrow points to the respective color.

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Arrow points to red: Danger! Emergency!

Repeat the measurement and follow the instructions of your physician or case manager.

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Arrow points to green: Everything OK

Repeat the measurement and follow the instructions of your physician or case manager.

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*Depending on the measurement mode set by the physician. The following setting is also possible:

FEV1 = Forced expiratory volume after one second

© 2014 eResearchTechnology GmbH or one of its affiliates.
All rights reserved

eResearchTechnology GmbH
Sieboldstrasse 3
97230 Estenfeld, Germany
+49 9305 720-9891 tel
+49 9305 720-7891 fax

www.ert.com

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Copyright © 2014 bei eResearchTechnology GmbH, Estenfeld, Germany