**English**  **en**

## Bivona® Fome-Cuf® Adult Tracheostomy Tube

### Description

The Bivona® Fome-Cuf® Adult Tracheostomy Tube is a radiopaque silicone tube with a foam-filled silicone cuff. The tube has a red wing pilot cuff and an integral 15mm connector that can be rotated on the flange to allow angular adjustment of the breathing circuit.

The B551 series comes with a SidePort™ AutoControl™ airway connector.

The B551 series has a talk attachment and SidePort™ AutoControl™ airway connector.

Each tube is individually packaged with an obturator, twill tape, disconnect wedge and Cuff Maintenance Device (CMD). This product is manufactured with Bivona® SuperSlick® formulation.

The contents of this package are sterile unless the package is opened or damaged.

### Intended Use

This tube is intended to provide direct airway access for a tracheostomized patient for up to 29 days. It may be reinserted up to 10 times for single patient reuse. The talk attachment is available to allow verbal communication for patients with functional vocal cords. This attachment allows an auxiliary source of gas to flow past the vocal cords when the control port is occluded.

### Contraindications

Do not use with lasers or electrocautery devices.

### Warnings

- This device must only be used by clinicians trained in its use.**
- Do not occlude the tracheostomy tube with any device, including one-way speaking valves, or total airway occlusion may occur.**
- Ensure that the lubricant gel used does not occlude the lumen of the tube preventing ventilation of the patient.**
- Verify the tracheostomy tube position by bronchoscopic visualization and palpation of the cuff. Clamping the red wing pilot port can result in accidental neck dissection. Clamping the red wing pilot port can also result in neck dissection. Clamping the red wing pilot port can also result in neck dissection. Do not clamp the red wing pilot port.**
- Remove the CMD™ and plug the red wing pilot port. The foam should remain compressed during cannulation.**
- Insert the obturator into the tube and lubricate the patient end with a water soluble lubricant.
- Insert the tube and obturator into the patient's stoma, directed downwards, and immediately remove the obturator in a curved direction. See Illustration 8.
- Open the red wing pilot port and allow the cuff to passively expand. If the desired seal is not achieved, one of the AutoControl™ set up options should be used to inflate the cuff.
- Confirm proper tube placement and secure with twill tape around the patient's neck.**
- Periodic evaluation of the cuff is recommended to measure residual cuff volume, evaluate any condensate, and verify cuff integrity. Prior to evaluation, remove any secretions that may have pooled above the cuff. Cuff volume is determined by applying a vacuum to remove all air
- Do not use a mechanical device to clamp the red wing pilot port as damage may occur.**

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and then turning the stopcock off to the CMD™. Release the CMD™ and record the volume. Changes in cuff volume alert the practitioner to tracheal changes.

Prior to decannulation, remove any secretions that may have pooled above the inflated cuff. Evacuate the cuff from the cuff using the CMD™ until the red wing pilot port collapses, clamp the pilot port with your fingers and plug it. The port should remain collapsed, indicating complete evaution of the cuff.

### 8551 Series (Talk Tube):

Attach the thumb connector to a dry air line and adjust for flow to approximately 4 lpm. Occlude the thumb port when speech is desired. Increase the flow as required to obtain adequate voice. If the line becomes occluded, flush with saline and suction above the cuff.

**CAUTION:**
Disconnect gas flow when not in use or the patient's upper airway may become dry. Do not use the tube with any obturator.

Additional Fome-Cuf® user's manuals are available on request.

### Reprocessing Instructions

- Make sure the cuff is fully expanded and then plug the red wing pilot port to prevent fluids from wetting the foam.
- Gently wash the tube (inside and out) and the obturator in a mild, detergent solution using a soft bristled brush. Rinse thoroughly. Do not use alcohol-based cleansers or disinfectants in an emergency. Always inspect the product prior to reuse. Do not use a tube that is cut or damaged.
- Use the following instructions to clean the product:
  - CAUTION:** The thumb control adapter and the Cuff Maintenance Device cannot be placed in an autoclave or steam sterilizer.
- Sterilization**
Insert the obturator in the tube with the red wing pilot port open and wrap along with the disconnect wedge in protective clean cloth or place them in a sterilization pouch. Sterilize in a gravity displacement steam autoclave at 121 °C (250°F) for 40 minutes. Do not flush sterilize. Do not use deep vacuum "flash" cycles or pulse vacuum cycles.

### Directions

- Attach the CMD™ to the red wing pilot port and evacuate the cuff until the red wing pilot port collapses. Clamp the red wing pilot port. Remove the CMD™ and plug the red wing pilot port. The foam should remain compressed during cannulation.
- Insert the obturator into the tube and lubricate the patient end with a water soluble lubricant.
- Insert the tube and obturator into the patient's stoma, directed downwards, and immediately remove the obturator in a curved direction. See Illustration 8.
- Open the red wing pilot port and allow the cuff to passively expand. If the desired seal is not achieved, one of the AutoControl™ set up options should be used to inflate the cuff.
- Confirm proper tube placement and secure with twill tape around the patient's neck.**
- Periodic evaluation of the cuff is recommended to measure residual cuff volume, evaluate any condensate, and verify cuff integrity. Prior to evaluation, remove any secretions that may have pooled above the cuff. Cuff volume is determined by applying a vacuum to remove all air

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