

Bivona® Cuffless Adjustable Neck Flange Hyperflex™ Pediatric Tracheostomy Tube

I.D. mm	A mm
2.5	0-55
3.0	0-60
3.5	0-65
4.0	0-70
4.5	0-75
5.0	0-80
5.5	0-85

Portex, Inc.

5700 West 23rd Avenue, Gary, Indiana 46406 USA

European Representative: **Portex Ltd.**, Hythe Kent CT21 6JL, UK
 PKG-DFU-AHF-2 1997-11 Rev. 11

Bivona®, Portex®, Smiths®, Hyperflex™, SuperSlick®, and XL-PCN™ are trademarks of Portex, Inc.

smiths
 Smiths Medical
 A part of Smiths Group plc

ENGLISH

Description

The Bivona Cuffless Adjustable Neck Flange Hyperflex Pediatric Tracheostomy Tube is a radiopaque wire reinforced silicone tube with an integral 15 mm ISO swivel, and a sliding lockable neck flange. Each tube is individually packaged with twill tape and a disconnect wedge. This product is manufactured with Bivona SuperSlick® XL-PCN™ formulation.

The contents of this package are sterile unless the package is opened or damaged.

Intended Use

Not intended for long term use.

This tube provides **temporary airway access** for a tracheostomized patient when determining the optimal tube length for a patient. This tube must be replaced with a fixed neck flange tracheostomy tube when the optimal length is determined.

Contraindications

- Do not use in home care settings.
- Do not use with lasers, electrosurgical devices or MRI.

Warnings

- Failure to use a disconnect wedge when removing connections to the 15 mm ISO swivel can result in significant damage to the tube. Do not pull on the connection. See *Illustration A*.
- Do not use a tube that is cut or damaged. Use of a damaged tube can result in airway compromise.

Cautions

- **For temporary use only.** Contact Portex, Inc. to obtain a fixed neck flange tracheostomy tube after optimal length is determined.
- For single patient use only. Do not reprocess.
- Do not leave neck flange unlocked when in use.
- This device must be thoroughly inspected for signs of damage or wear prior to each use.
- Do not allow unsupported weight to pull on this tube.
- Use only water soluble lubricants with this tube. Use of petroleum based lubricants can damage this tube.
- This package insert presents general information on how to use this device. However, it is necessarily the responsibility of the healthcare practitioner to assess these instructions and to take appropriate measures as needed to assure the product and method for use are compatible with individual patient needs. It is also the responsibility of the healthcare practitioner to assure that the instructions for use and maintenance are understood by and provided to the care giver.

Directions

1. Insert an introducer into the tube without protruding beyond the patient end and lubricate with a water soluble lubricant.
2. Insert the tube and introducer into the patient's stoma, and immediately remove the introducer.
3. Adjust the tube to the desired length. Align the white stripe on the shaft with the diamond mark on the adjustable neck flange locking ring.
4. Secure the neck flange by pulling the locking tab over the locking pin.
5. Confirm proper tube placement and secure with twill tape around then patient's neck.
6. The depth of the tracheostomy tube shaft can be adjusted by releasing the locking ring. Sterile water makes adjustment easier. Depth markings indicate the length of the shaft as measured from the patient side of the neck flange.

CAUTION: Never leave the tube in place with the locking mechanism released.

Attention, see instructions for use

Sterile, using ethylene oxide

Lot number

Use by


Catalogue Number

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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