Bivona® Uncuffed Neonatal/Pediatric Tracheostomy Tube

These Instructions for use apply to the following Bivona® Cuffless Tracheostomy Tubes with straight or "Y" shaped neck flange (Illustrations show the Straight Neck Flange version):

I.D. mm	Straight Neck Flange		"Y" Neck Flange	
	Neonatal	Pediatric	Neonatal	Pediatric
2.5mm	60SN025	60SP025	60N025	60P025
3.0mm	60SN030	60SP030	60N030	60P030
3.5mm	60SN035	60SP035	60N035	60P035
4.0mm	60SN040	60SP040	60N040	60P040
4.5mm	-	60SP045	-	60P045
5.0mm	-	60SP050	-	60P050
5.5mm	-	60SP055	-	60P055

Description

The Bivona® Uncuffed Neonatal/Pediatric Tracheostomy Tube is a wire reinforced radiopaque silicone tube with a 15mm connection that can be rotated on the flange to allow angular adjustment of the breathing circuit. Each tube is individually packaged with an obturator, twill tape, and a disconnect wedge.

The contents of this package are sterile unless the package is opened or damaged.

Intended Use

This tube is intended to provide direct airway access for a tracheotomized patient for up to 29 days. It may be reprocessed up to 5 times for single patient reuse.

Contraindications

- Use in conjunction with lasers or electrosurgical devices where their output may contact and damage the tube.
- Use with MRI or electrosurgical devices.

Warnings

This device must only be used by clinicians trained in its use.

- Ensure that the lubricant gel used does not occlude the lumen of the tube, preventing ventilation of the patient.
- Verify the tracheostomy tube position by bronchoscopic visualization and/or chest X-ray to ensure correct placement. Incorrect placement could result in trauma to the trachea or respiratory obstruction.
- During and after attachment of the breathing system to the tracheostomy tube connector, avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion.
- Tracheostomy tubes must be cleaned/replaced regularly to suit individual patients needs.
- Failure to use a disconnect wedge when removing connections to the 15mm connector can result in significant damage to the tube. Do not attempt to remove or pull on the 15mm connector. See Illustration A.
- Do not use a tube that is cut or damaged. Use of a damaged tube can result in airway compromise.

Cautions

- This device must be thoroughly inspected for signs of damage or wear prior to each use.
- The security of all connectors should be checked when the circuit is established and frequently thereafter.
- The disconnect wedge supplied should be kept with the product at all times in order to allow easy disconnection in an emergency.
- Use only water soluble lubricants with this tube. Use of petroleum based lubricants can damage this tube.
- Patients breathing gas should be adequately humidified to minimize encrustation of the tracheostomy tube and prevent mucosal damage.
- The patency of the tracheostomy tube lumen must be assured by regular suctioning.
- This package insert presents general information on how to use this device. However, it is necessarily the responsibility of the healthcare practitioner to assess these instructions and to take appropriate measures as needed to assure the product and method for use are compatible with individual patient needs. It is also the responsibility of the healthcare practitioner to assure that the instructions for use and maintenance are understood by and provided to the caregiver.
- Dispose of this product in a safe manner according to local guidelines for disposal of contaminated waste.
- Follow universal precautions as specified by the Center for Disease Control and Prevention, USA.

Directions

- 1 Insert the obturator into the tube and lubricate the patient end with a water soluble lubricant.
- 2 Insert the tube and the obturator into the patient's stoma, curved end down, and immediately remove the obturator in a downward direction. See Illustration B.
- 3 Confirm proper tube placement and secure with twill tape around the patient's neck.

Reprocessing Instructions

- 1 Remove the 15mm connector by inserting disconnect wedge between the plastic connector and base of the silicone neck flange. See Illustration C.
- 2 Gently wash the tube and connector (inside and out) and the obturator in a mild detergent solution with a soft bristled brush. Rinse thoroughly. Excessive scrubbing, scraping, or stretching will damage this device. Always inspect the product prior to reuse. Do not use a tube that is cut or damaged.
- 3 Allow the tube to air dry, then reattach the connector. (Position the flange on the connector toward the patient.)

4A Sterilization

Insert the obturator into the tube and wrap along with the disconnection wedge in protective lint-free cloth or place in a sterilization pouch. Sterilize in a gravity displacement steam autoclave at 121°C (250°F) for 40 minutes. Do not use deep vacuum "flash" cycles or pulse vacuum cycles.

4B Sanitization

Allow the tube to air dry completely and store with the obturator in place in a clean plastic bag. Prior to sanitization, remove the tube from the plastic bag. Place the tube and the obturator side by side in a pan of rapidly boiling distilled water. Cover the pan and REMOVE FROM HEAT. Allow to cool prior to use and handle only by the obturator handle and the neck flange.

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I.D. mm	Straight Neck Flange		"Y" Neck Flange	
	Neugeborene	Kinder	Neugeborene	Kinder
2,5mm	60SN025	60SP025	60N025	60P025
3,0mm	60SN030	60SP030	60N030	60P030
3,5mm	60SN035	60SP035	60N035	60P035
4,0mm	60SN040	60SP040	60N040	60P040
4,5mm	-	60SP045	-	60P045
5,0mm	-	60SP050	-	60P050
5,5mm	-	60SP055	-	60P055

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© 2005, Smiths Medical family of companies
Bivona® is a trademark of the Smiths Medical family of companies.
The symbol ® indicates the trademark is registered in the US Patent and Trademark Office and certain other countries.
Smiths Medical ASD, Inc.
Smiths Medical International Ltd