PORTEX

Bivona® Cuffless Adult Tracheostomy Tube

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Bivona® hettelest Trakeostomör för voksne

English
These instructions contain important information for safe use of the product. Read the entire contents of these instructions for Use, including Warnings and Precautions, before using this product. Follow warnings, precautions and instructions could result in death or serious injury to the patient. The responsibility of the healthcare practitioner is to ensure that the instructions for use and maintenance are understood by and provided to the patient and caregiver. © 2019 Smiths Medical. All rights reserved.

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1. DESCRIPTION:
Bivona® Cuffless Adult Tracheostomy Tube is a radiopaque silicone tube with an integral 15 mm connector. Each tube is individually packaged with an external decannulation inflator 15 cm long. Bivona® Cuffless Adult Tracheostomy Tube is designed with a hollow obturator. This device is intended for single patient use. The contents of this package are sterile unless the package is opened or damaged.

Ref.	LD (mm)
60A150	5.0
60A160	6.0
60A170	7.0
60A180	8.0
60A190	9.0

2. INDICATIONS FOR USE:
This tube is intended to provide direct airway access for a tracheostomized patient for up to 29 days. It may be reinserted into the airway as a permanent tracheostomy. Do not use for patients with tracheomalacia.

3. CONTRAINDICATIONS:
Do not use in conjunction with tubes or electrocaval devices where their output may contact and damage the tube.

4. WARNINGS:
4.1 Use only water soluble lubricants with this tube. Use of petroleum based lubricants can damage this tube. Care must be taken to prevent silicone lubricants from extruding from the obturator. Failure to do so may result in reduced breathing capability.
4.2 Verify use of tracheostomy tube position by bronchoscopic visualization and chest X-ray to ensure correct placement. Incorrect placement could result in trauma to the trachea or respiratory distress.
4.3 The security of all connectors should be checked when the circuit is established and frequently thereafter. Failure to do so may result in reduced ventilation. During and after attachment of the breathing circuit to the tracheostomy tube connector, avoid application of excessive inward linear force to avoid the risk to prevent accidental disconnection or occlusion.
4.4 Care must be taken to avoid excessive force when connecting a circuit to the connector. Excessive force can lead to airway occlusion.
4.5 Do not use a tube that is cut or damaged. Use of a damaged tube can result in airway patency and reduced breathing capabilities.
4.6 Care must be exercised when using the decannulation inflator cap to ensure that the patient does not aspirate. Observe the patient for signs of respiratory distress and immediately remove the decannulation inflator cap if necessary.
4.7 The patient's oxygen saturations must be monitored regularly to avoid restricted breathing capabilities.
4.8 This tube and its connector device must always be fully inserted into the stoma and the neck spacer secured by tightening to ensure that it remains in place. Failure to do this could result in leaking of the tube and/or a disconnected obturator.
4.9 Patient's breathing should be adequately humidified to minimize encrustation of the tracheostomy tube which can result in restricted breathing capability.
4.10 The presence of encrustations and other biological material is vital to removing bacterial contamination from the product. This must be carried out before reinsertion.
4.11 During repositioning, aggressive scrubbing of the tube or use of hard bristled brushes or sharp objects can damage the tube surface. If necessary, use a soft brush.
4.12 Sanitization should be undertaken within the two hour period before tube re-insertion. Once a tube has been sanitized, it should be stored in the original packaging until ready for use. Do not use unless you have disinfected the tube and have re-sterilized it.
4.13 If a patient is to have a tracheostomy tube replacement, ensure that the guidelines concerning tube changing and potential respiratory distress are followed.
4.14 Patients using the Bivona® Cuffless Adult Tracheostomy Tube should be monitored for any signs of respiratory distress.

5. CAUTIONS:
5.1 Federal USA Law restricts this device to use by or on the order of a physician.
5.2 Bivona® Cuffless Adult Tracheostomy Tube is not intended for use with the instructions for re-processing solutions that the product may be re-used on the same patient.
5.3 Care must be taken to ensure that the device is not re-used on the same patient.
5.4 These products have been validated up to 121° C (250° F) for up to 30 days. Do not expose the product to temperatures in excess of 121° C (250° F). Do not attempt to re-sterilize the product by autoclave. Do not use these vacuum 'flash' or vacuum clean cycles.
5.5 Following cleaning, handle the product only by the neck flange or connector to avoid disconnection, contamination or damage to the device.
5.6 Do not reprocess decannulation cap to prevent compromise of cap.
5.7 Do not use the product after the expiration date of the package.
5.8 Do not use the product after the expiration date of the package.
5.9 Do not use the product after the expiration date of the package.
5.10 Do not use the product after the expiration date of the package.
5.11 Do not use the product after the expiration date of the package.

Smiths Medical ASD, Inc.
10000 W. 16th Ave.
Denver, CO 80202 USA
Phone: +1 (303) 733 3500
Fax: +1 (303) 733 5648
www.smiths-medical.com

Smiths Medical International, Ltd.
Hythe, Dorset DT2 6DL
UK
Phone: +44 (0) 1329 444636
Fax: +44 (0) 1329 444642
www.smiths-medical.com

Australian Representative:
Smiths Medical Australia Pty. Ltd.
5/145 Street
Briar Hill Technology Park
Eight Mile Plains
Queensland 4113
Australia
Phone: +61 (0) 7 3271 3000
Fax: +61 (0) 7 3271 3011
www.smiths-medical.com

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Smiths Medical and Portex design trademarks: SuperSilic and Bivona are trademarks of the Smiths Medical family of companies. The symbol in the trademark is registered in the U.S. Patent and Trademark Office and other countries around the world.
Patending.
© 2019 Smiths Medical family of companies. All rights reserved.

Date of Manufacture:
Herstellungsdatum
Feststillingsdato
Fecha de fabricación
Date de Fabrication
Ημερομηνία κατασκευής
Data di fabbricazione
Produktdatum

Smiths Medical and Portex design trademarks: SuperSilic and Bivona are trademarks of the Smiths Medical family of companies. The symbol in the trademark is registered in the U.S. Patent and Trademark Office and other countries around the world.
Patending.
© 2019 Smiths Medical family of companies. All rights reserved.

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Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
Attenzione : la legge federale (U.S.A.) limita la vendita di questo dispositivo ai medici o alla prescrizione medica.