LABELLING LAYOUT SHEET

Title - IFU - Blue Line Ultra Suctionaid Tracheostomy Tube - 004/005/164

INFORMATION SHOWN:- Print layout on IFU. (To scale)

English - en	French - fr	German - de	Italian - it	Spanish - es	Portuguese - pt	Dutch - nl	Swedish - sv	Danish - da
Y	Y	Y	Y	Y	Y	Y	Y	Y
Norwegian - no	Finnish - fi	Greek - el	Japanese - ja	Czech - cs	Slovak - sk	Hungarian - hu	Polish - pl	Slovenian - sl
Y	Y	Y	N	Y	N	Y	Y	N
Estonian - et	Latvian - lv	Lithuanian - lt	Turkish - tr	Bulgarian - bg	Romanian - ro	Russian - ru	Croatian - hr	Chinese - zh
N	N	N	N	N	N	N	N	N

CCN. 004/005/164

Print Revision level - REV .001

Size: 210mm x 148mm (A5) (Folded Size: 105mm x 148mm (A6))

Print Details: PMS Pantone Process Black

Design Authority: Hythe


Material: 60gsm White Bond Paper

Description: Instruction booklet, printed one colour black throughout. Booklet stitched twice. 50 leafs (100 pages). Pages 99 and 100 are blank.

NOTES:

a) The layouts comply with the requirements of Corporate Policy CP30 Issue 3

b) This IFU replaces 004/000/002 (L1987)

Last Change; PWR 4745, CCF H8079; Origination of new IFU to CP30 Issue 3. Updates include new branding, addition of Distribution Statement + Do not reuse text to meet requirements of MDD, translations reorganized for consistency with corporate language order, removal of text relating to "Tracheostomy Change Kit" as advised by R&D, references to "Velcro" updated inline with "Patent & Trademark" requirements. Addition of new "Caution" and "Latex Free" symbols to glossary; existing symbols retained for use with existing packaging. Addition of Phthalate symbol + website.

SMI Labelling / Date	Design Authority RA DA Date	Marketing / Date
<i>Signature held on master</i>	<i>Signature held on master</i>	<i>Signature held on master</i>
Compiled By - Alie Richardson	Product Engineering Date	Date Effective
<i>Signature held on master</i>		26 May 2011
Compilation Date - 08 Mar 11	<i>Signature held on master</i>	

EXAMPLE LAYOUTS (To Scale)

PORTEX[®]

**Blue Line Ultra[®] Suctionaid[®]
Tracheostomy Tube**

en INSTRUCTIONS FOR USE	da BRUGSANVISNING
fr MODE D'EMPLOI	no BRUKSANVISNING
de GEBRAUCHSANLEITUNG	fi KÄYTTÖOHJEET
it ISTRUZIONI PER L'USO	el ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
es INSTRUCCIONES DE USO	cs NÁVOD K POUŽITÍ
pt INSTRUÇÕES DE UTILIZAÇÃO	hu HASZNÁLATI UTASÍTÁS
nl GEBRUIKSAANWIJZING	pl INSTRUKCJA STOSOWANIA
sv BRUKSANVISNING	

Smiths Medical International Ltd.
Hythe, Kent, CT21 6JL, UK
Tel: +44 (0)1303 260551
Distributed by:
Smiths Medical ASD, Inc.
Keene, NH 03431, USA
Tel: +1 603 352 3812
Australian Representative:
Smiths Medical Australasia Pty. Ltd.
Brisbane, QLD 4113, Australia.
Tel: +61 (0)7 3340 1300
www.smiths-medical.com

CE
0473


CCN.004/005/164 REV.001 02/11

smiths medical

EXAMPLE LAYOUTS (To Scale)

English

en

Portex® Blue Line Ultra® Suctionaid® Tracheostomy Tube

INSTRUCTIONS FOR USE

These instructions are intended for use with the following Portex® products:

100/860/060-100 Tracheostomy Tube, Suctionaid® Soft-Seal® Profile® Cuff sizes 6.0-10.0mm, suction control valve, neck tapes, identification labels

100/870/060-100 Tracheostomy Tube, Suctionaid® Smooth Inner Cannula, Soft-Seal® Profile® Cuff, size 6.0-10.0mm, suction control valve, tube holder, cleaning brush, identification labels

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

Tracheostomy Tube O.D. (mm)	Tracheostomy Tube I.D. (mm)	Tracheostomy Tube Length C (mm)	Backdon Size (approx.)	Inner Cannula I.D. (mm)	Cuff Resting Diameter (mm)
9.2	6.0	64.5	5	5.0	20
10.5	7.0	70.0	6	5.5	24
11.3	7.5	73.0	7	6.0	30
11.9	8.0	75.5	8	6.5	30
12.6	8.5	78.0	105°	7.0	30
13.3	9.0	81.0	105°	7.5	30
14.0	10.0	87.5	105°	8.5	30

WARNING:

● READ THE COMPLETE INSTRUCTIONS FOR USE BEFORE USE

● THIS DEVICE MUST ONLY BE USED BY CLINICIANS TRAINED IN ITS USE

NOTE:

1 Use only Blue Line Ultra® Inner Cannula with Blue Line Ultra® Tracheostomy Tubes.

2 The use of an Inner Cannula reduces the airway diameter by 1.5mm (excluding size 6.0mm tracheostomy tube which reduces the airway diameter by 1.0mm).

3 Do not use Portex® Contingated Inner Cannula (100/528; 100/522) with Blue Line Ultra® Tracheostomy Tubes.

Description

A sterile range of single use Tracheostomy tubes manufactured from clear biocompatible Polyvinyl Chloride (PVC) incorporating the following features:

● Smooth arc shaped design with tapered tip to assist placement and removal of both tracheostomy tube and smooth inner cannula with minimal disruption to patient's anatomy.

● Thermosensitive material with sufficient initial rigidity for insertion which then conforms to the individual patient's upper respiratory tract at body temperature ensuring minimum trauma.

● 15mm connector for connection to anaesthetic or ventilator breathing circuit is located on the tube not the Inner Cannula.

● Anatomically shaped flexible flange for improved patient comfort and improved security of the Tracheostomy tube.

● New Profile® Soft-Seal® cuff design with reduced bulk and smooth transition to tracheostomy tube to improve ease of insertion.

● A hollow obturator with atraumatic tip for use with or without guide wire for improved ease of insertion. The obturator retaining clips prevent obturator movement during insertion.

● A blue cuff inflation line and pilot balloon with integral Luer check valve for inflating and monitoring the cuff.

● A clear suction line terminating above the cuff which allows any pooled secretions to be aspirated. The suction line terminates

with a connector for suction devices and a Luer port to connect to either a syringe or the vacuum control valve which is included in the pack.

● Self-adhesive Patient Notes Label providing tracheostomy tube size and manufacturing batch details.

The range is available in a number of pack configurations as listed below:

Tracheostomy Tube:

Available cuffed in sizes 6, 7, 7.5, 8, 8.5, 9 and 10mm I.D.

Complete with hollow obturator and tracheostomy tapes.

Tracheostomy Tube with Smooth Inner Cannula:

Available cuffed in sizes 6, 7, 7.5, 8, 8.5, 9 and 10mm I.D. with two appropriately sized single patient use smooth inner cannulae, cleaning brush, tracheostomy tube holder and hollow obturator.

Inner Cannula:

Two single patient use smooth inner cannulae for sizes 6, 7, 7.5, 8, 8.5, 9 and 10mm Tracheostomy tube.

For full details on kit contents and ordering information consult the relevant product literature

Indications for Use

Portex® Blue Line Ultra® Suctionaid® Tracheostomy Tube is indicated for airway maintenance of tracheostomised patients. Suctionaid® allows aspiration of contaminated mucous and subglottic secretions that collect and build up between the tracheostomy tube cuff and the glottis.

Contraindications

There are no known contraindications for use of this device.

Relative Contraindications

Morbid obesity and/or neck oedema (skin to trachea distance may render tracheostomy tube too short)

Adverse Reactions

Reported adverse reactions associated with tracheostomy tubes are many and diverse. Consult standard textbooks and medical literature for specific adverse reaction information.

Warnings

1 Prior to removal of cuffed Tracheostomy tubes, all air must be completely removed from the cuff to prevent damage to the trachea and stoma.

2 If the Tracheostomy tube is lubricated prior to insertion, ensure the lubricant does not occlude the lumen of the tube preventing ventilation of the patient

3 Do not inflate the cuff with a measured volume of air, or by feel of pressure from the syringe, since little resistance should be felt during inflation.

4 During anaesthesia, nitrous oxide may diffuse into the cuff causing an increase or decrease in cuff pressure.

5 Cuff pressures should be monitored and adjusted routinely. Over-inflation of the cuff may result in permanent damage to the trachea.

6 Contact with electro surgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air or ignite in an enriched oxygen environment (e.g. anaesthesia).

7 Do not use this Tracheostomy tube in patients with abnormal upper airway anatomy or pathology as this may result in partial total airway obstruction.

8 Do not lubricate the inner cannula as this may occlude the inner cannula and prevent the inner cannula being retained within the tracheostomy tube.



9 If it is not possible to remove the inner cannula from the tracheostomy tube, do not try to forcibly remove. Both the inner cannula and tracheostomy tubes should be removed together and replaced with a new tracheostomy tube and inner cannula.

10 Where an inner cannula is used, ensure that it is both the correct diameter and length for the tracheostomy tube in use.

Sheet 3 of 51

LA003/1-2

EXAMPLE LAYOUTS (To Scale)

Size markings are provided on the tube flange, pack labels and additional patient labels to help identify the correct inner cannula to use. Use only an inner cannula that is supplied specifically for use with the Portex® Blue Line Ultra® range. 11 Use of an inner cannula of the wrong diameter may cause difficulties in introduction or unnecessarily restrict gas flow. An inner cannula which is too long may protrude excessively from the outer tube, leading to tracheal damage or occlusion. Use of an inner cannula which is too short may lead to build-up of secretions which could cause infection or blockage. 12 Tracheostomy tubes must be changed regularly to suit individual patient's needs. 13 The use of an Inner Cannula reduces the airway diameter by 1.5mm (excluding size 6.0mm tracheostomy tube which reduces the airway diameter by 1.0mm). 14 During and after attachment of the breathing system to the tracheostomy tube connector, avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion. 15 Never insert the guidewire into a tracheostomy tube with the red 'J' tip straightener in place to avoid aspiration of the straightener. 16 Ensure that secretions are suctioned from above the cuff immediately before deflating the cuff of the tracheostomy tube in order to minimise secretions entering the lungs on cuff deflation. 17 Do not suction/lavage through suction lumen with cuff deflated to avoid solutions entering the lungs. 18 Do not use a decumulation cap with this product. The decumulation cap is indicated for use only with fenestrated tracheostomy tubes. 19 If a serial percutaneous tracheostomy dilator is used as a loading dilator/obturator to insert this tracheostomy tube, ensure the dilator can be inserted and removed from the tracheostomy tube without undue force before use. 20 If the suction lumen is used for vocalisation, avoid using excessive pressure and/or high oxygen flow rates to avoid trauma to the trachea. 21 The suction lumen must not be used for vocalisation in a freshly formed stoma as air intended for vocalisation may leak through the stoma preventing speech or cause surgical emphysema. 22 The suction lumen must never be used for vocalisation in patients with upper airway obstruction as the air supply may result in a build up of pressure in the subglottic space. 23 Suctionaid® does not eliminate the need for routine endobronchial suctioning and oral hygiene procedures. 24 Do not fit an inner cannula when using a tracheostomy tube in conjunction with a right-angled connector that incorporates a co-axial fresh gas tube which protrudes into the connector lumen (e.g. Norman elbow). The protruding fresh gas tube could occlude the inner cannula, restricting expiratory flow of gas, resulting in barotrauma/hypoxia. Precautions 1 The security of all breathing system connectors should be checked when the circuit is established and frequently thereafter. Disconnection may be facilitated with the use of a disconnection wedge (100/555/000). 2 Patients should be adequately humidified (e.g. using the Portex® Thermovent® T100/570/015) to minimise encrustation of the Tracheostomy tube and/or inner cannula lumen and prevent mucosal damage. 3 The patency of the Tracheostomy tube lumen must be assured by regular suctioning and regular cleaning of the inner cannula (if one is used). Check routinely and replace as required to maintain a patent airway. Maximum recommended period of use 30 days. 4 Cuff pressure and volume should be monitored, adjusted and recorded routinely so that damaging overinflation or aspiration from underinflation does not occur. 5 Devices used in or during inflation of the cuff must be clean and free from all foreign matter. The inflation device should be removed from the inflation valve immediately after use, and the dust cap fitted. Test inflate the cuff before insertion. 6 Guard against cuff damage by avoiding contact with sharp edges. 7 The inflation line valve may interfere with Magnetic Resonance Imaging (MRI) picture clarity. Ensure the valve is positioned away from the area being scanned.	8 Care should be taken to ensure that the inner cannulae do not become kinked or damaged during cleaning, and that no kinked or damaged inner cannulae are re-inserted into the Tracheostomy tube. 9 Following insertion, verify the tube position and the location of the fenestra, for example by means of a chest X-ray. 10 If Portex® Tracheostomy tubes are used outside the hospital, the patient must be instructed by a healthcare professional in the safe use and handling of the product. A patient care booklet is also available from Smiths Medical Customer Services. 11 Only sterile saline should be used to clean accessible parts of the tracheostomy tube whilst in-situ. If the tracheostomy tube is removed from the patient, it must be discarded. 12 Do not use any abrasive cleaning implements to clean the inner cannula. Use only the cleaning brush provided. Do not use this brush to clean the tracheostomy tube. 13 Blue Line Ultra® Tracheostomy Tubes and Inner Cannulae are designed and intended for single patient use. They should not be re-sterilised by any method. 14 Repositioning of the in-situ tracheostomy tube while the cuff is inflated should be avoided. 15 When suctioning via the suction line has been completed the suction device including the control valve should be removed, and the suction line capped to minimise the risk of infection. 16 Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician. Instructions for use 1 The integrity of the cuff and inflation system should be checked prior to insertion. 2 Prior to insertion check that the obturator can be removed from the tube having first released the retaining clips by means of a slight twist. Check that the inner cannula (where fitted) can be inserted and removed from the Tracheostomy tube. Re-insert the obturator ensuring that it 'clicks' into position.   3 A small amount of water soluble lubricant can be applied to the outer tube and obturator tip to ease insertion. 4 Insert the Tracheostomy tube in accordance with currently accepted medical techniques and remove the obturator. 5 Inflate the cuff (where fitted) with the minimum amount of air to provide an effective seal. 6 Secure the Tracheostomy tube with the cotton tapes or Portex® Tracheostomy tube holder (100/503/200) if provided. 7 Insert the inner cannula (where fitted) into the tracheostomy tube. Be sure that the inner cannula is securely inserted. A slight snap should be felt indicating proper positioning. NOTE: The inner cannula should be routinely checked, cleaned or replaced at regular intervals. 8 Remove rectangular Patient Notes Label from the backing sheet and add to the patient's notes in order to record details of the tracheostomy tube used. 9 Remove the 'Tracheostomy Flange Tag' from the backing sheet (removing the circular 'cut-out'), fold in half and place over the tracheostomy tube 15mm connector in order to aid identification of the size of tube used and, if appropriate, alert the clinical staff to the presence of fenestrations. 10 Two inner cannulae are provided so that one inner cannula can continue to be used whilst the other inner cannula is cleaned. 11 Check the inner cannula prior to insertion/re-insertion and discard if kinked or damaged. Replacement inner cannulae are available (100/850).
---	--

EXAMPLE LAYOUTS (To Scale)

English

**Portex® Single Use
Tracheostomy Tube Holder**
INSTRUCTIONS FOR USE
REF: 100/503/200

With the tracheostomy tube connector pointing towards you:


- 
- 1 Take the holder and ensure that the Portex™ logo is facing away from you.
 - 2 Open out one of the narrow VELCRO® brand straps (A), and thread it all the way through the back of the tracheostomy flange (B) through the left-hand hole.
 - 3 Fold the narrow strap (A) back and secure to the wide band material.
 - 4 Pass the holder behind the patient's neck and repeat **2** and **3** for the right-hand side.
 - 5 Adjust the holder using the wide VELCRO® brand strap (D) to fit the patient, ensuring the holder is securely joined.

One size fits all adults

Portex, Profile, Soft-Seal, Blue Line Ultra, Thermoport, Portex design mark and Smiths Medical design mark are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

©2011 Smiths Medical family of companies. All rights reserved.

VELCRO® is a registered trademark of Velcro Industries B.V.

VELCRO® is a registered trademark of Velcro Industries B.V.

Using the Suctionaid®

- 1 Check the suction port located above the cuff. If it is in the trachea and not in pre-tracheal tissue following tube insertion to ensure correct function and prevent flushed saline from exiting the stoma site.
- 2 With the cuff inflated attach the suction device, either syringe or low level suction to the suction line using the suction control valve, and slowly aspirate the secretions up to a maximum suction level of 300mm Hg.
- 3 If resistance is met on aspiration, flush the suction line with either air or sterile saline and aspirate again. If resistance to instillation of air or saline is still met and the instilled fluid cannot be withdrawn, the suction line may be blocked and the tracheostomy tube should be changed or an alternative means of suction instituted.
- 4 Once suctioning is complete remove the suction device including the control valve and cap the Luer port at the proximal end of the suction tube. The suction control valve is intended for single patient use and can be cleaned using sterile saline solution. Replacement suction control valves are also available (100/515/203).
- 5 If using continuous low level suction, monitor the suction pressure regularly. If the suction pressure increases disconnect the suction and follow the above procedure for resistance to suction.

Replacing the Inner Cannula

- 1 Disconnect the tracheostomy tube from the breathing system
- 2 Remove the old inner cannula
- 3 Insert a new or spare inner cannula
- 4 Reconnect the breathing system and check for security
- 5 Clean or discard the old inner cannula

Inner Cannula Cleaning Instructions

Daily cleaning of the inner cannula is recommended or whenever the inner cannula becomes contaminated by mucus or secretions. This period of time will vary according to the individual patient's needs.

Always remove the inner cannula from the Tracheostomy tube for cleaning.

- 1 Soak the inner cannula in sterile saline or mild detergent solution for up to 15 minutes.
- 2 Agitate the inner cannula in the solution to loosen any secretions. Using the cleaning brush provided to help remove any dried or sticky secretions.
- 3 After soaking and cleaning, rinse the inner cannula and cleaning brush with fresh sterile saline solution.
- 4 Allow to dry naturally. Once the inner cannula is clean and dry it should be stored clean and dry and free from particulate-matter.

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									
www.smiths-medical.com/pththalates									
Quantity									
Cuff Resting Diameter									
Sterilised using ethylene oxide									
Do not use if package is damaged.									
Keep dry									
Keep away from sunlight									
Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.									

www.smiths-medical.com/pththalates									
Quantity									
Cuff Resting Diameter									
Sterilised using ethylene oxide									
Do not use if package is damaged.									
Keep dry									
Keep away from sunlight									
Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.									

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