DE

Technische Daten (DIN EN ISO 9360 - 1)

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GB

Technical data (DIN EN ISO 9360 - 1)

Indication: For patients' protection and humidifying of respiratory air.
To be placed between ETT / tracheostomy tube and Y-piece of breathing circuit.

Tidal volume (VT), recommended range, (when applied close to patient): 300 – 1500 ml

Pressure drop: **Inspiration:** 0,5 l / s = 1,4 h Pa
1,0 l / s = 3,3 h Pa
1,5 l / s = 5,5 h Pa
Expiration: 0,5 l / s = 1,3 h Pa
1,0 l / s = 3,2 h Pa
1,5 l / s = 5,4 h Pa

Humidity loss:
VT 1000 ml = 13,9 mg/l
VT 750 ml = 12,0 mg/l
VT 500 ml = 11,1 mg/l

Compressible volume: 45 ml

Bacterial / viral filtration: 99,999%

Filtration: Polyurethan foam with CaCl₂ impregnation

Weight: 29 g

CO₂ connector: yes

Connection: 22 mm f / 15 mm m - 22 mm m / 15 mm f

Instruction for use

For patient's protection – Filter should be placed between ETT / tracheostomy tube and Y-piece of breathing circuit. Unit can also be connected directly to ventilator in- or expiratory. Prior to use check proper fixation of connections, including Luer cap, and check for leakage. If necessary, compensate ventilation due to added dead space of filter. After connection of the filter to artificial respiration the air flow is to be supervised up to a reliable adjustment of parameters.

Contraindications

- Filter should be placed near to patient. Do not use outside of recommended tidal volume range, as the added dead space may cause CO₂ retention at too small tidal volumes or increased respiratory resistance and unsatisfactory humidification at too high tidal volumes.
- Do not use on dehydrated patients, patients with heavy secretions from the lungs and airways or on patients with very large air leaks where most of the expiratory gases do not pass through the Filter/HME.
- Do not use in connection with heated humidifier or nebulizer (disconnect filter/HME for intermittent aerosol treatment), and do not add moisture tount, as this might result in increased resistance.

Caution!

- Do not use if package is damaged.
- Filter should be replaced every 24 hours or more often, if required, to avoid resistance increase due to accumulation of secretions.
- Single use only.
- Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.
- Re-sterilisation (for reuse) will affect the function of unit.
- Used filters should be disposed of immediately.

Disposal

Disposal may only be carried out in accordance with the applicable national regulations for waste materials.

Report

All serious incidents associated with the product should be reported to the manufacturer and the competent authority.

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Primed Halberstadt Medizintechnik GmbH