Nellcor™ Adult-Neonatal SpO₂ Sensor with Wraps

- en** Instructions For Use
- fr** Mode d'emploi
- de** Gebrauchsanleitung
- nl** Gebruiksaanwijzing
- es** Instrucciones de uso
- it** Istruzioni per l'uso
- sv** Bruksanvisning
- da** Brugervejledning
- pt** Instruções de uso
- ru** Инструкция по применению
- zh** 使用说明

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Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

IP22 Protected against access to hazardous parts with a finger and protected against vertically falling water drops when sensor tilted up to 15°.

This device is not made with natural rubber latex or DEHP.

Indications/Contraindications

The Nellcor™ Adult-Neonatal SpO₂ Sensor with Wraps, model OXI-A/N, is indicated for use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing less than 3 kg or more than 40 kg.

The OXI-A/N is applied to the sensor site with Nellcor disposable adhesive wraps, model ADH-A/N or model FOAM-A/N. The ADH-A/N and FOAM-A/N wraps are intended for single patient use only.

Use this sensor only with Nellcor instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). This sensor integrates Nellcor OxiMAX technology into its design. When connected to an OxiMAX-enabled instrument, this sensor uses OxiMAX technology to provide additional advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each manufacturer of Nellcor-compatible instruments is responsible for determining whether and under what conditions its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument.

The OXI-A/N is contraindicated for use on patients who exhibit allergic reactions to the pressure-sensitive adhesive on the wraps.

Instructions for Use

Reusable sensors may be used on the same site for a maximum of 4 hours, provided the site is inspected routinely to ensure skin integrity, correct positioning, and adhesion of the sensor wrap. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients.

To apply the OXI-A/N sensor:

- 1) Place the wrap on a flat surface, adhesive (sticky) side up, with the tab on the left **A**.
- 2) Remove the small piece of paper backing and enough of the large piece of backing to expose both alignment holes **B**.
- 3) Position the sensor on the wrap so the alignment bumps protrude through the holes. Press down firmly in the center of the sensor to ensure that the wrap adheres tightly **C**.
- 4) Remove the remainder of the paper backing from the wrap. The sensor and wrap are now ready to be applied to the patient. Select an appropriate site as described below.

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

Adult (> 40 kg): The preferred site is around an index finger with the cable positioned along the top of the finger. Alternatively, use a thumb or another finger, with the cable positioned along the palm; or around a great toe, with the cable positioned along the sole of the foot.

Neonate (< 3 kg): The preferred site is the foot, below the toes, with the cable positioned along the sole of the foot. Alternatively, place the sensor around the palm of a hand, below the fingers, with the cable positioned along the palm.

- 5) Position the OXI-A/N so that the notches are centered on the side of the foot or finger. The cable end must be positioned along the appropriate surface, as described in #4 above.
- 6) Wrap the OXI-A/N around the foot or finger so the optical components (and alignment bumps) oppose each other across the site. Press the OXI-A/N to ensure that it adheres snugly.
- 7) Wrap the remaining portion of the sensor wrap around the site, loosely enough to ensure good circulation.
- 8) Plug the OXI-A/N into the instrument as described in the instrument operator's manual, and verify proper operation.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned — or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

Removal and Cleaning

Remove the disposable wrap from the OXI-A/N by peeling it away from the sensor. The sensor may be surface-cleaned with a solution such as 70% isopropyl alcohol. If low-level disinfection is required, use a 1:10 bleach solution. Do not use undiluted bleach (5%–5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor could occur.

Caution: Do not expose connector pins to cleaning solution as this may damage sensor.

To clean or disinfect the sensor using the recommended wipe method:

- 1) Saturate a clean, dry gauze pad with the cleaning solution. Wipe all surfaces of the sensor and the cable with this gauze pad.
- 2) Saturate another clean, dry gauze pad with sterile or distilled water. Wipe all surfaces of the sensor and cable with this gauze pad.
- 3) Dry the sensor and cable by wiping all surfaces with a clean, dry gauze pad.

To clean or disinfect the sensor using the recommended soak method:

- 1) Place the sensor in the cleaning solution, such that the sensor head(s) and desired length of cable are completely immersed.
- Caution**: Do not immerse the connector end of the cable as this may damage the sensor.
- 2) Dislodge air bubbles by gently shaking the sensor and cable.
- 3) Soak the sensor and the cable for 10 minutes.
- 4) Remove from cleaning solution.
- 5) Place the sensor and the cable in room temperature sterile or distilled water for 10 minutes.
- 6) Remove from the water.
- 7) Dry the sensor and cable by wiping all surfaces with a clean, dry gauze pad.

Caution: Using excessive force when removing the disposable wrap may damage the sensor.

Caution: Do not sterilize by irradiation, steam, or ethylene oxide.

Warnings

- 1) Failure to apply the OXI-A/N properly may cause incorrect measurements.
- 2) High oxygen levels may predispose a premature infant to develop retinopathy. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accord with accepted clinical standards and considering the accuracy range of the oximeter being used.
- 3) Using the OXI-A/N in the presence of bright lights may result in inaccurate measurements. In such cases, cover the sensor site with an opaque material.
- 4) Reusable sensors must be moved to a new site at least every 4 hours. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients. If skin integrity changes, move the sensor to another site.
- 5) Use only Nellcor wraps designed for use with this sensor. Do not use tape. Use of additional tape or various other types of wraps can cause skin damage.
- 6) Intravascular dyes or externally applied coloring such as nail polish, dye, or pigmented cream may lead to inaccurate measurements.
- 7) Excessive motion may compromise performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
- 8) When the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements.
- 9) As with all medical equipment, carefully route cables to reduce the possibility of patient entanglement and strangulation.
- 10) Do not use the OXI-A/N or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the OXI-A/N may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.
- 11) Do not alter or modify the OXI-A/N. Alterations or modifications may affect performance or accuracy.
- 12) For additional warnings, cautions or contraindications when using this sensor with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the manufacturer of the instrument.
- 13) If the sensor is misapplied with excessive pressure for prolonged periods, a pressure injury can occur.

If you have questions regarding any of this information, contact Covidien's Technical Services Department, or your local Covidien representative.

Warranty

To obtain information about a warranty, if any, for this product, contact Covidien's Technical Services or your local Covidien representative.


Accuracy Specifications


For the accuracy specification range when used with Nellcor monitors, refer to information provided with the monitor, or (in the U.S.), contact Covidien's Technical Services Department. Outside the U.S., contact your local Covidien representative.

For the accuracy specification range of this sensor when used with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the instrument manufacturer.

Additional Copies of Instructions

Additional copies of these Instructions are available at no charge by calling Covidien or its authorized distributors. Also, permission is hereby granted under Covidien copyrights to Purchasers of products obtained from Covidien or its authorized distributors to make additional copies of these instructions for use by such purchasers.

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COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.