

Nellcor™ Pediatric-Infant SpO₂ Sensor with Wraps

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
es	Instrucciones de uso
it	Istruzioni per l'uso
sv	Bruksanvisning
da	Brugsanvisning
pt	Instruções de uso
el	Οδηγίες χρήσης
ru	Инструкция по применению
zh	使用说明

IP22 Protected against access to hazardous parts with a finger and protected against vertically falling water drops when sensor tilted up to 15°.

This device is not made with natural rubber latex or DEHP.

Indications/Contraindications

The Nellcor™ Pediatric-Infant SpO₂ Sensor with Wraps, model OXI-P/I, is indicated for use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing between 3 kg and 40 kg.

The OXI-P/I is applied to the sensor site with Nellcor disposable adhesive wraps, model ADH-P/I or model FOAM-P/I. The ADH-P/I and FOAM-P/I wraps are intended for single patient use only.

Use this sensor only with Nellcor instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). This sensor integrates Nellcor OxiMax technology into its design. When connected to an OxiMax-enabled instrument, this sensor uses OxiMax technology to provide additional advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each manufacturer of Nellcor-compatible instruments is responsible for determining whether and under what conditions its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument.

The OXI-P/I is contraindicated for use on patients who exhibit allergic reactions to the pressure-sensitive adhesive on the wraps.

Instructions for Use

Reusable sensors may be used on the same site for a maximum of 4 hours, provided the site is inspected routinely to ensure skin integrity, correct positioning, and adhesion of the sensor wrap. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients.

To apply the OXI-P/I sensor:

- Place the wrap on a flat surface, adhesive (sticky) side up, with the tab on the left **A**.
- Remove the small piece of paper backing and enough of the large piece of backing to expose both alignment holes **B**.
- Position the sensor on the wrap so the alignment bumps protrude through the holes. Press down firmly in the center of the sensor to ensure that the wrap adheres tightly **C**.
- Remove the remainder of the paper backing from the wrap. The sensor and wrap are now ready to be applied to the patient. Select an appropriate site as described below.

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

Pediatric (15–40 kg): The preferred site is around an index finger with the cable positioned along the top of the finger. Alternative sites are around the thumb or another finger with the cable positioned along the top of the finger; or a great toe with the cable positioned along the sole of the foot.

Infant (3–15 kg): The preferred site is around a great toe, with the cable positioned along the sole of the foot.

- Position the OXI-P/I so that the notches are centered on the side of the toe or finger. The cable end must be positioned along the appropriate surface, as described in #4 above.
- Wrap the OXI-P/I around the toe or finger so the optical components (and alignment bumps) oppose each other across the site. Press the OXI-P/I to ensure that the sensor adheres snugly.
- Wrap the remaining portion of the sensor wrap around the site, loosely enough to ensure good circulation.
- Plug the OXI-P/I into the instrument as described in the instrument operator's manual, and verify proper operation.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.



Removal and Cleaning

Remove the disposable wrap from the OXI-P/I by peeling it away from the sensor. The sensor may be surface-cleaned with a solution such as 70% isopropyl alcohol. If low-level disinfection is required, use a 1:10 bleach solution. Do not use undiluted bleach (5%–5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor could occur.

Caution: Do not expose connector pins to cleaning solution as this may damage sensor.

To clean or disinfect the sensor using the recommended wipe method:

- Saturate a clean, dry gauze pad with the cleaning solution. Wipe all surfaces of the sensor and the cable with this gauze pad.

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