


OXIMAX[®] MAX-I

Infant O₂ Sensor

English	3
Français	6
Deutsch	9
Nederlands	13
Italiano	16
Español	19
Svenska	22
Dansk	25
Norsk	28
Suomi	31
Português	34
РУССКИЙ	37
中文	41
Polski	44
Cesky	48
Slovenian	51
Magyar	54
Ελληνικά	57
TÜRKÇE	61

OXIMAX[®]
MAX-IInfant O₂ Sensor

② Single patient use only  Latex free  **Rx ONLY**

Directions for Use

Indications/Contraindications

The Nellcor[®] OxiMAX[®] infant oxygen sensor, model MAX-I, is indicated for single patient use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing between 3 and 20 kg.

The MAX-I is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

Instructions for Use

1. Remove plastic backing from the MAX-I and locate transparent windows (b) on the adhesive side. Windows cover optical components. Note corresponding alignment marks (a) on the nonadhesive side and the dashed line (c) midway between the marks. **(1)**

A great toe is the preferred MAX-I location. Alternatively, apply the sensor to another digit of similar size, for example, a thumb.

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

2. Orient the MAX-I so the window next to the cable is aligned on the bottom of the great toe as shown. The cable should extend towards the heel. **(2)**
3. Wrap the MAX-I firmly, but not too tightly around the toe. Windows must oppose each other. **(3)**
4. Wrap any excess tape loosely around toe. Use additional tape provided to secure cable across bottom of foot, loosely enough to maintain good circulation. **(4)**
5. Plug the MAX-I into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

Reapplication

The MAX-I can be reused on the same patient as long as the adhesive tape attaches without slippage.

WARNINGS

1. **Do not use the MAX-I or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the MAX-I may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.**

CAUTIONS

1. In the event of damage to the sterile packaging, do NOT re-sterilize. Follow local governing ordinances and recycling instructions regarding disposal or recycling of sensors.
2. Failure to apply the MAX-I properly may cause incorrect measurements.
3. While the MAX-I is designed to reduce the effects of ambient light, excessive light may cause inaccurate measurements. In such cases, cover the sensor with opaque material.
4. Circulation distal to the sensor site should be checked routinely. The site must be inspected every 8 hours to ensure adhesion, skin integrity, and correct optical alignment. If skin integrity changes, move the sensor to another site.
5. Intravascular dyes or externally applied coloring such as nail polish, dye, or pigmented cream may lead to inaccurate measurements.
6. Excessive motion may compromise performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
7. Do not immerse in water or cleaning solutions. Do not resterilize. Immersion or resterilization could damage the sensor.
8. If the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements.
9. Do not alter or modify the MAX-I. Alterations or modifications may affect performance or accuracy.
10. For additional warnings, cautions or contraindications when using this sensor with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the manufacturer of the instrument.

Note: High oxygen levels may predispose a premature infant to develop retinopathy. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accordance with accepted clinical standards and considering the accuracy range of the oximeter being used.

Accuracy Specifications

For the accuracy specification range when used with Nellcor monitors, refer to information provided with the monitor, or (in the U.S.), contact Nellcor's Technical Services Department. Outside the U.S., contact your local Nellcor representative.

For the accuracy specification range of this sensor when used with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the instrument manufacturer.

Use this sensor only with Nellcor OxiMAX instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor OxiMAX sensors (Nellcor-compatible instruments). This sensor integrates Nellcor OxiMAX technology into its design. When connected to an OxiMAX-enabled instrument, this sensor uses OxiMAX technology to provide additional advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each Nellcor OxiMAX-compatible instrument manufacturer is responsible for determining optimal compatibility conditions and settings for its instruments to provide safe and effective use of each Nellcor OxiMAX sensor. This includes specifications and/or warnings, cautions, or contraindications. Refer to each instrument operator's manual or consult manufacturer for complete instructions, warnings, cautions, or contraindications on the use of this sensor with their Nellcor OxiMAX compatible instrument.

Additional Copies of Instructions

Additional copies of these Instructions are available at no charge by calling Nellcor or its authorized distributors. Also, permission is hereby granted under Nellcor Puritan Bennett copyrights to purchasers of products obtained from Nellcor or its authorized distributors to make additional copies of these instructions for use by such purchasers.

This sensor is sold under the following U.S. patents and foreign equivalents, for single patient use only: 6,591,123; 6,708,049; and 4,830,014. Any other use of this sensor is not authorized by Nellcor under any patents.

© 2005 Nellcor Puritan Bennett Inc. All rights reserved.

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Extra exemplaren van deze Instructies zijn kosteloos verkrijgbaar bij Nellcor en haar officiële distributeurs. Bovendien wordt hierbij aan kopers van producten bij Nellcor of zijn erkende vertegenwoordigers, onder het Nellcor Puritan Bennett copyright, toestemming gegeven extra kopieën te maken van deze instructies voor eigen gebruik.
© 2005 Nellcor Puritan Bennett Inc. Tüm hakları saklıdır.
Tyco Healthcare UK LTD