

OXIMAX[®] MAX-N

Neonatal/Adult O₂ Sensor


English	3
Français	6
Deutsch	10
Nederlands	14
Italiano	18
Español	22
Svenska	26
Dansk	29
Norsk	32
Suomi	35
Português	38
РУССКИЙ	41
中文	45
Polski	48
Cesky	52
Slovenian	55
Magyar	58
Ελληνικά	62
TÜRKÇE	66


1


2

3

4

OXIMAX[®]
MAX-NNeonatal/Adult O₂ Sensor

② Single patient use only  Latex free  **Rx ONLY**

Directions for Use

Indications/Contraindications

The Nellcor[®] OxiMAX neonatal/adult oxygen sensor, model MAX-N, is indicated for single patient use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for neonates weighing less than 3 kg or adults weighing more than 40 kg.

The MAX-N is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

Instructions for Use

1. Remove plastic backing from the MAX-N and locate transparent windows on the adhesive side. Windows cover optical components. Note corresponding alignment marks (a) on nonadhesive side and dashed line (b) midway between the marks. **(1)**
2. Orient the MAX-N so the dashed line is on the lateral edge of the site.
Neonates: The preferred site is a foot. Alternatively, use a hand. The window next to the cable goes on the sole of the foot as shown. **(2)**
Adults: The preferred site is an index finger. Alternatively, other fingers may be used. The window next to the cable goes on the nail side, distal to the first joint. Do not place on a joint. Note that the cable must be positioned on the top of the hand. **(3)**

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

3. Wrap the MAX-N firmly, but not too tightly around the foot or finger. Windows must oppose each other.
4. Plug the MAX-N into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

Reapplication

The MAX-N can be reused on the same patient as long as the adhesive tape attaches without slippage.

Enclosed adhesive “dots” are provided for reapplication. Place a transparent dot over each window as shown, then remove the protective paper that covers each dot. The sensor is now ready to be reapplied to the same patient. **(4)**

WARNINGS

- 1. Do not use the MAX-N or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the MAX-N may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.**

CAUTIONS

1. In the event of damage to the sterile packaging, do NOT re-sterilize. Follow local governing ordinances and recycling instructions regarding disposal or recycling of sensors.
2. Failure to apply the MAX-N properly may cause incorrect measurements.
3. While the MAX-N is designed to reduce the effects of ambient light, excessive light may cause inaccurate measurements. In such cases, cover the sensor with opaque material.
4. Circulation distal to the sensor site should be checked routinely. The site must be inspected every 8 hours to ensure adhesion, skin integrity, and correct optical alignment. If skin integrity changes, move the sensor to another site.
5. Intravascular dyes or externally applied coloring such as nail polish, dye, or pigmented cream may lead to inaccurate measurements.
6. Excessive motion may compromise performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
7. Do not immerse in water or cleaning solutions. Do not resterilize. Immersion or resterilization could damage the sensor.
8. If the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements.
9. Do not alter or modify the MAX-N. Alterations or modifications may affect performance or accuracy.
10. For additional warnings, cautions or contraindications when using this sensor with Nellcor-compatible instruments, refer to the instrument operator’s manual or contact the manufacturer of the instrument.

Note: High oxygen levels may predispose a premature infant to develop retinopathy. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accordance with accepted clinical standards and considering the accuracy range of the oximeter being used.

Accuracy Specifications

For the accuracy specification range when used with Nellcor monitors, refer to information provided with the monitor, or (in the U.S.), contact Nellcor's Technical Services Department. Outside the U.S., contact your local Nellcor representative.

For the accuracy specification range of this sensor when used with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the instrument manufacturer.

Use this sensor only with Nellcor OxiMAX instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor OxiMAX sensors (Nellcor-compatible instruments). This sensor integrates Nellcor OxiMAX technology into its design. When connected to an OxiMAX-enabled instrument, this sensor uses OxiMAX technology to provide additional advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each Nellcor OxiMAX-compatible instrument manufacturer is responsible for determining optimal compatibility conditions and settings for its instruments to provide safe and effective use of each Nellcor OxiMAX sensor. This includes specifications and/or warnings, cautions, or contraindications. Refer to each instrument operator's manual or consult manufacturer for complete instructions, warnings, cautions, or contraindications on the use of this sensor with their Nellcor OxiMAX compatible instrument.

Additional Copies of Instructions

Additional copies of these Instructions are available at no charge by calling Nellcor or its authorized distributors. Also, permission is hereby granted under Nellcor Puritan Bennett copyrights to purchasers of products obtained from Nellcor or its authorized distributors to make additional copies of these instructions for use by such purchasers.

This sensor is sold under the following U.S. patents and foreign equivalents, for single patient use only: 6,591,123; 6,708,049; and 4,830,014. Any other use of this sensor is not authorized by Nellcor under any patents.

© 2005 Nellcor Puritan Bennett Inc. All rights reserved.

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Extra exemplaren van deze Instructies zijn kosteloos verkrijgbaar bij Nellcor en haar officiële distributeurs. Bovendien wordt hierbij aan kopers van producten bij Nellcor of zijn erkende vertegenwoordigers, onder het Nellcor Puritan Bennett copyright, toestemming gegeven extra kopieën te maken van deze instructies voor eigen gebruik.
© 2005 Nellcor Puritan Bennett Inc. Tüm hakları saklıdır.
Tyco Healthcare UK LTD