Portex® Adjustable Flange Tracheostomy Tube

en INSTRUCTIONS FOR USE
de GEBRAUCHSANLEITUNG
da BRUGSANVISNING
es INSTRUCCIONES DE USO
fr MODE D'EMPLOI
el ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
it ISTRUZIONI PER L'USO
no BRUKSANVISNING
nl GEBRUIKSAANWIJZING

pt INSTRUÇÕES DE UTILIZAÇÃO
sv BRUKSANVISNING
fi KÄYTTÖOHJEET
cs NÁVOD K POUZITÍ
pl INSTRUKCJA STOSOWANIA
hu HASZNÁLATI UTASÍTÁS
tr KULLANIM TALİMATI
et KASUTUSJUHEND
lv LIETOŠANAS INSTRUKCIJA







Smiths Medical International Ltd.
Hythe, Kent, CT21 6JL, UK
Tel: +44 (0)1303 260551
www.smiths-medical.com

Distributed by:
Smiths Medical ASD, Inc
Keene, NH 03431, USA

Australian Representative:
Smiths Medical Australasia Pty. Ltd.
61 Brandl Street, Eight Mile Plains,
QLD 4113, Australia.
Tel: +61 (0) 7 3340 1300

CCN 004/001/151-1
© 2005 Smiths Medical family of companies
All rights reserved.

smiths

Portex® Adjustable Flange Tracheostomy Tube

INSTRUCTIONS FOR USE

- Read the instructions completely before use.
- This device is intended for use by trained clinicians.

These instructions apply to the use of the following Portex® products:

- 100/523/060-100 *Tracheostomy tube, Siliconised PVC, Adjustable Flange, Soft-Seal® cuff, (Low Pressure/High Volume), sizes 6.0mm - 10.0mm I.D.*
- 100/526/060-100 *Tracheostomy tube, Siliconised PVC, Adjustable Flange, Uncuffed, sizes 6.0mm - 10.0mm I.D.*

Description

A range of sterile single-use tracheostomy tubes with an adjustable flange for use on patients with abnormal skin to trachea depth and/or to facilitate stoma care. Manufactured from Siliconised Polyvinyl Chloride (PVC) incorporating the following features:

- Thermosensitive material with sufficient initial rigidity for intubation which then conforms to the individual patient's upper respiratory tract at body temperature ensuring minimum trauma.
- Anatomically shaped design with radio-opaque Blue Line® to confirm correct tube placement by X-Ray.

Instructions For Use

- 1 The integrity of the cuff and inflation system should be checked prior to insertion.
- 2 Prior to use, loosen the locking screw on the adjustable flange by turning counter clockwise [Fig. 1] and rotate the flange 180° into position as shown in [Fig. 2]. With the locking screw still loosened, check that the obturator can be removed from the tube. Re-insert the obturator.
- 3 Insert tracheostomy tube in accordance with currently accepted medical techniques and slide the flange into the correct position [Fig. 3].
- 4 Remove the obturator, slide the flange forward against the skin surface and lock the flange in position by turning the screw clockwise [Fig. 4].

- Inflate the cuff (where fitted) with the minimum amount of air to provide an effective seal.
- Secure the tracheostomy tube with the cotton tapes provided or a Portex® Tracheostomy Tube Holder.

Stoma Site Care

- Access may be gained to the stoma site to facilitate cleaning by loosening the cotton tapes or Portex® Tracheostomy Tube Holder and fully releasing the screw on the flange by turning counter clockwise [Fig. 1].
- Slide the flange forward while securely holding the tracheostomy tube in position, and clean the stoma site and/or the flange with sterile saline [Fig. 5].
- Securely holding the tracheostomy tube in position, slide the flange back and lock into position by turning the screw clockwise [Fig. 4].

Cuff Resting Diameter (mm)

Tracheostomy Tube I.D. (mm)	6.0	7.0	8.0	9.0	10.0
Soft-Seal® Cuff	20	24	30	30	30

Cautions

- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Follow universal precautions as specified by the Centre for Disease Control and Prevention, USA, or relevant advisory bodies in the country where this device is being used.
- Patients should be adequately humidified to minimise encrustation of the tracheostomy tube lumen and prevent mucosal damage.
- The patency of the Tracheostomy Tube lumen must be assured by regular suctioning. Check routinely and replace as required to maintain a patent airway. Maximum recommended period of use 30 days.
- Cuff pressure should be monitored and adjusted routinely.
- Devices used in or during inflation of the cuff must be clean and free from all foreign matter. The inflation device should be removed from the inflation valve immediately after use, and the dust cap fitted.
- Guard against cuff damage by avoiding contact with sharp edges.
- The inflation line valve may interfere with Magnetic Resonance Imaging (MRI) picture clarity. Ensure the valve is positioned away from area being scanned.
- If Portex® tracheostomy tubes are used outside the hospital, the patient must be instructed by a health care professional in the safe use and handling of the product.
- Do not use solutions other than sterile saline to clean any part of this tracheostomy tube.
- Dispose of the tracheostomy tube in a safe manner according to local guidelines for disposal of contaminated medical waste.

Warnings

- 1 Prior to removal of cuffed Tracheostomy Tubes, all air must be completely removed from the cuff to prevent damage to trachea and stoma.
- 2 If the Tracheostomy Tube is lubricated prior to insertion, ensure the lubricant does not occlude the lumen of the tube preventing ventilation of the patient.
- 3 Do not inflate with a measured volume of air, or by feel of pressure from the syringe, since little resistance should be felt during inflation. Over inflation may result in tracheal damage.
- 4 During anaesthesia, nitrous oxide may diffuse into the cuff causing an increase or decrease in cuff pressure.
- 5 Contact with electrosurgery electrodes or laser surgery beams must be avoided because P.V.C. will produce toxic fumes in air or ignite in an enriched oxygen environment (e.g. anaesthesia).
- 6 If the flange adjustment is already at the minimum setting it may be necessary to consider replacing the Adjustable Flange Tracheostomy tube with Portex® tracheostomy tube with a standard geometry to avoid the potential trauma from cuff movement caused by increased leverage [Fig. 6].
- 7 Do not position the flange beyond the maximum setting as this may result in obstruction of the inflation line where fitted.
- 8 Do not use this tracheostomy tube in patients with abnormal airway anatomy or pathology as this may result in partial or total airway obstruction.
- 9 The security of all breathing system connectors should be checked when the circuit is established and frequently thereafter. Failure to check may result in the patient not being adequately ventilated. Disconnection may be facilitated with the use of a disconnection wedge (product code 100/555/100 if not included).

Portex, Soft-Seal, Blue Line and the Smiths design mark are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2005 Smiths Medical family of companies. All rights reserved.

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Portex, Soft-Seal, Blue Line a Smiths design mark jsou ochranné známky společnosti Smiths Medical. Symbol® označuje, že daná ochranná známka je registrována u U.S. Patent and Trademark Office (Patentový a známkový úřad USA) a v některých dalších zemích. Všechny ostatní názvy, ochranné známky a servisní značky zde uvedené jsou vlastnictvím svých příslušných majitelů.