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Copyright © 2017 ASSKEA GmbH, Gebesee.

The safety of the **ASSKEA secretion aspirators M- and S-series** complies with the acknowledged rules of technology and the requirements of the **German Medical Devices Act**.

The **ASSKEA secretion aspirators M- and S-series** bear the **CE marking CE0843** in accordance with EU Council Directive 93/42/EEC concerning medical devices and meet the essential requirements of Annex I of this directive.

The **ASSKEA secretion aspirators M- and S-series** have been tested in accordance with IEC 62353.

The **quality management system** applied by ASSKEA GmbH is certified in compliance with the relevant international standards.

The **ASSKEA secretion aspirators M- and S-series** are medical aspirators classified as class IIa in accordance with EU Council Directive 93/42/EEC, Annex IX.

Errors and omissions excepted.

Contents

1	User Information	4	3.1.3 Connecting the ASSKEA secretion aspirators S-series	18
1.1	Using this Instruction Manual	4	3.1.4 Connecting the ASSKEA secretion aspirators M-series	18
1.2	Icons	5	3.2 Startup	18
1.2.1	Device, Packaging and Accessories	5	3.2.1 Connecting the reusable secretion canister system	19
1.2.2	Operating Elements	6	3.2.2 Connecting the ASSKEA disposable secretion canister system	20
1.3	Symbol Convention	6	3.2.3 Suction tube	22
1.4	Glossary	6	3.3 Operation of the ASSKEA secretion aspirators S-series	22
1.5	Intended Use	8	3.3.1 Setting the vacuum	22
1.5.1	ASSKEA secretion aspirators S-series	8	3.3.2 Aspiration	22
1.5.2	ASSKEA secretion aspirators M-series	8	3.4 Operation of the ASSKEA secretion aspirators M-series	23
1.5.3	Indications	8	3.4.1 Display and status LEDs	23
1.5.4	Contraindications	8	3.4.2 Setting the vacuum	25
1.5.5	Restrictions on Use	8	3.4.3 Aspiration	26
1.6	Basic Safety Instructions – ATTENTION!	8	3.5. Decommissioning	27
1.7	User Requirements	10	4 Maintenance	28
1.8	Information on Product Liability	10	4.1 Cleaning and Care	28
1.9	Material Compatibility	10	4.1.1 General Information	28
2	Product Description	11	4.1.2 Cleaning and disinfection of the surface of the device	28
2.1	Whole view	11	4.1.3 Cleaning and disinfection of the reusable secretion canister system	29
2.1.1	ASSKEA secretion aspirators S-series	11	4.1.4 Cleaning and disinfection of the external canister 'Bag'	29
2.1.2	ASSKEA secretion aspirators M-series	11	4.1.5 Cleaning of the connecting tube	29
2.2	Product Contents	12	4.1.6 Cleaning/disposal of the suction tube	30
2.3	Product Properties	12	4.1.7 Disposable products	30
2.3.1	ASSKEA secretion aspirators S-series	12	4.2 Maintenance and Service	30
2.3.2	ASSKEA secretion aspirators M-series	13	4.3 Testing of the ASSKEA secretion aspirators M- and S-series	31
2.3.3	Information on the double filter system	14	5 Troubleshooting	31
2.4	Warranty	15	5.1 ASSKEA secretion aspirators S-series	31
3	Operation	16	5.2 ASSKEA secretion aspirators M-series	32
3.1	Set-up and Startup	16	6 Transport, Storage and Disposal	34
3.1.1	Operating elements of the ASSKEA secretion aspirators S-series	17	6.1 Decontamination prior to Shipment	34
3.1.2	Operating elements of the ASSKEA secretion aspirators M-series	17	6.2 Storage	34

6.3	Disposal	34
7	Technical Data	35
7.1	Technical Data of the ASSKEA secretion aspirators S-series	35
7.2	Technical Data of the ASSKEA secretion aspirators M-series	36
7.3	EMC Information	37
8	Order Information Accessories	42
9	Publishing Information	43

1 User Information

1.1 Using this Instruction Manual

Please read this entire instruction for use before operating the **ASSKEA secretion aspirators M- and S-series** for the first time. Alongside medical professionals and trained relatives of patients, it is also intended that patients use the devices themselves, providing that they have been trained accordingly. In this case, all functions can be used safely.

Please read the safety instructions (chapter 1.6) to avoid hazards.

This instruction for use is a component of the **ASSKEA secretion aspirators M- and S-series**. Therefore please keep this instruction for use in an easily accessible place.

Include this instruction for use when passing the **ASSKEA secretion aspirators M- and S-series** on to third parties.

1.2 Icons

1.2.1 Device, Packaging and Accessories

Symbol	Meaning
	Attention: possible bodily injury, health risks or possible property damage
	Follow the instruction for use.
	Power supply unit
	Protection class II
	Protection class: Type BF (Body Floating)
	Temperature limitation
	Air pressure limitation
	Humidity limitation
	Protect from moisture
	Do not use, if packaging is damaged!

Symbol	Meaning
	Note with useful information and tips
	This device must not be disposed of in domestic waste.
	Order number
	Serial number
	Lot number
	Date of manufacture
	Manufacturer
	Operating time (power-on time and power-off time)
	Do not reuse.
	Fuse

1.2.2 Operating Elements

Symbol	Meaning
	Adjustment direction for vacuum setting using the vacuum regulator
	LED status indicator for external power supply (only M-series)
	LED status indicator for double filter system (only M20D and M30D)
	LED status indicator for the rechargeable battery (only M-series)
	LED status indicator ON/OFF (only M-series)

1.3 Symbol Convention

Symbol	Meaning
•	Enumeration
1. 2.	Perform the process in the specified order.

1.4 Glossary

A

Aspirate	Aspirate is the generic term for secretions, bodily fluids and liquids used for flushing that are typically accumulated when aspirating the upper respiratory system. It can be easily aspirated using the devices described here.
----------	--

C

Contamination	Contamination means that bacteria and viruses from the aspirate have come into contact with the device.
---------------	---

D

DFS®	Double filter system An external bacterial filter and a bacterial filter integrated into the aspirator make up the double filter system. The double filter system effectively protects the interior of the device from contamination and overflow. It enables safe processing and rapid reuse of the product.
------	---

Degree of protection The degree of protection indicates the protection of applied parts against electric shock. Applied parts of type BF have to be installed separately to the earth and are not suitable for direct application on the heart.

I

incl. Abbreviation for 'inclusive'

IP22 International Protection/Protection Class
The Protection Class defines the degree of protection of the device against contact and ingress of liquids.

The **ASSKEA secretion aspirators M- and S-series** are protected against finger access and falling water drops at an inclination of up to 15°.

M

MRI Abbreviation for magnetic resonance imaging
Here, a very strong magnetic field makes it possible to create sectional images of the human body for the purpose of assessing the organs.

O

Overflow Overflow means that the secretion is sucked into the interior of the device.

P

Processing The processing procedure is required for each new patient. The term processing denotes the process in which parts coming or potentially coming into contact with aspirate are cleaned, disinfected and replaced, if necessary. The processing procedure must only be performed by ASSKEA GmbH or an authorized service partner of ASSKEA GmbH.

R

resp. Abbreviation for 'respectively'

1.5 Intended Use

1.5.1 ASSKEA secretion aspirators S-series

The **ASSKEA secretion aspirators S-series** are mains-operated medical aspirators and serve for temporary aspiration of aspirate from the trachea. The **ASSKEA secretion aspirators S-series** are ideally suited for applications in the stationary area, but can also be used for care and rehabilitation in the private area (homecare) for the respective patient.

1.5.2 ASSKEA secretion aspirators M-series

The **ASSKEA secretion aspirators M-series** are network-independent medical aspirators and serve for temporary and primarily spontaneous aspiration of aspirate from the trachea. Typical areas of application are:

- in outpatient and inpatient care,
- in the homecare area, especially for the aspiration of aspirate from tracheostomy patients.

1.5.3 Indications

- Tracheostomy patients
- Aspiration on impaired respiratory function
- Aspiration of blood, secretion and food residues from the oral cavity, the pharyngeal zone and the bronchial system
- Evacuation of vacuum mattresses and vacuum splints

1.5.4 Contraindications

The **ASSKEA secretion aspirators M- and S-series** are contraindicated for the following applications:

- Liposuction
- Gynecology applications
- Thoracic drainage
- Continuous drainage
- Applications in wound areas

1.5.5 Restrictions on Use

- In medical rooms where potential equalization is necessary (e.g. heart surgery)
- In hazardous areas

1.6 Basic Safety Instructions – ATTENTION!

- Risk of damage due to improper power supply

Improper operation causes overvoltage in the device, which may be transmitted to the operator.

- Ensure prior to startup, that the mains supply for connecting the **ASSKEA secretion aspirators S-series** is suitable for 230V alternating current with a mains frequency of 50Hz and, for connecting the **ASSKEA secretion aspirators M-series**, that it is suitable for 100V to 240V alternating current with a mains frequency of 50–60Hz.
- Only operate the **ASSKEA secretion aspirators M-series** with the provided power supply unit (type GTM91099-6015-3.0-T2 by GlobTek Inc.)

- Risk of application under conditions that are not approved

- The devices are not intended to be used outdoors/during transportation.
- The devices are not suitable for applications within an MRI environment.

- Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspi-

rate cause health risks.

- Never use the device without an external bacterial filter.
(exception: ASSKEA disposable secretion canister system)
- Check that the bacterial filter is dry and clean prior to every application to ensure faultless function!
(exception: ASSKEA disposable secretion canister system)
- Always perform the aspiration with a suitable sterile aspiration set. The suction tube must never come into direct contact with the aspiration site.
- Follow the hygiene, cleaning and decontamination instructions.

- Hazard of persons due to improper handling

- Only use the device for its intended purpose.
- Never use the device for aspiration in the low vacuum range (e.g. for thorax drainage).
- If aspiration procedures are carried out too often, slight bleeding may result.
- When using the devices for the treatment of children, the vacuum must be set within the range of –100 and –200 mbar!
- When using the power supply unit, make sure the power supply unit is connected to the mains supply (100 V - 240 V AC) only after the power cord plug of the power

supply unit has first been connected to the aspirator.

The separation of the power supply unit from the mains supply must occur in exactly the opposite sequence (first separate the power supply unit from the mains supply (100 V - 240 V AC) and then the power cord plug from the aspirator).

- Risk of damage due to electromagnetic phenomena

Medical electrical equipment is subject to special precautionary measures regarding electromagnetic compatibility and must be installed and operated in accordance with the EMC information provided in the accompanying documentation! (see chapter 7.3)

- Hazard of persons due to strangulation

People may strangle themselves on the tubing or the power cord.

- Ensure that no unauthorized/uninvolved personnel is near the device during aspiration.
- Store the device incl. accessories in the shipping carton.

- Safety defects due to improper accessories and spare parts

The use of accessories and spare parts other than those recommended by ASSKEA GmbH may compromise the safety and function of

the device. Damage caused by using non-recommended accessories and spare parts or by improper use is not covered by warranty in any case.

- Only use original accessories and spare parts.

- Damage to the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Do not use the device in case of apparent housing damage.

- Damage to the device by heat

- Do not cover the power supply unit.
- Keep the device as well as the power cord and power supply unit away from other heat sources.

- Damage to the device by ingress of liquids

The **ASSKEA secretion aspirators M- and S-series** have the IP classification IP22 against the ingress of liquids. Nevertheless protect the device from moisture (see symbol in chapter 1.2.1).

- Do not use the device near splashing water.
- Do not use the device in damp rooms or while bathing or showering.
- Do not allow the power supply unit, power switch and the socket for the power supply

- unit, as well as the display film, to get wet.
- Never submerge the device in water or other liquids (also not while not in operation).

- Known or identifiable conditions for medical care within a domestic environment

- Children and pets must be kept away from the device to ensure that the device is not knocked over or dropped.
- Prior to connecting the power supply unit, ensure that the voltage of the device corresponds to the domestic power supply.
- Do not use the device in damp rooms, bathrooms or showers. Do not allow the power supply unit, plug and switch unit to get wet. Never submerge the device in water or other liquids (also not while not in operation).

- Inspection of the internal power supply

Since the internal battery is not maintained automatically in a fully operational condition, the state of charge must be checked periodically and, if necessary, the battery must be replaced by service staff. The battery may only be replaced by trained service staff, since the replacement by inadequately trained personnel could result in a hazard (such as excessive temperatures, fire or explosion)!

1.7 User Requirements

The **ASSKEA secretion aspirators M- and S-series** must only be operated and used by instructed and appropriately trained persons. Familiarize yourself with the mode of operation of the **ASSKEA secretion aspirators M- and S-series** prior to startup.

1.8 Information on Product Liability

The liability for the operation of the device is channeled to the operator in the following cases:

- the device is used outside its intended use,
- the device is not used in accordance with the instruction for use,
- the device is opened by unauthorized personnel,
- installation, settings, enhancements, routine maintenance or repairs are performed by unauthorized personnel,
- original accessories and spare parts have not been used,
- the device is used beyond its lifetime of 5 years.

1.9 Material Compatibility


Aggressive substances may damage the device and the accessories.

- Please follow the cleaning and care instructions (chapter 4.1).
-

2 Product Description

2.1 Whole view


2.1.1 ASSKEA secretion aspirators S-series (S20K, S20, S30)

- A** Secretion canister (reusable) with integrated overflow protection in the lid (see fig. 29)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Vacuum meter
- E** Vacuum regulator
- F** Power switch
- G** Mains connection 230V AC
- H** Mains cable
- I** Suction tube

Fig. 1


2.1.2 ASSKEA secretion aspirators M-series a) Analog (M20)

- A** Secretion canister (reusable) with integrated overflow protection in the lid (see fig. 29)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Vacuum meter
- E** Vacuum regulator
- F** Power switch
- G** Socket for power supply unit (12V DC)
- H** Power supply unit (GlobTek GTM91099-6015-3.0-T2)
- I** Mains cable
- J** Suction tube

Fig. 2

b) Digital (M20D, M30D)

- A** Secretion canister (reusable) with integrated overflow protection in the lid (see fig. 29)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Buttons for vacuum adjustment
- E** ON/OFF button
- F** Max. button
- G** Status indicator
- H** Socket for power supply unit (12V DC)
- I** Power supply unit (GlobTek GTM91099-6015-3.0-T2)
- J** Mains cable
- K** Suction tube

Fig. 3

2.2 Product Contents

- the **ASSKEA secretion aspirator S-series** or the **ASSKEA secretion aspirator M-series**
- the instruction for use of the **ASSKEA secretion aspirators M- and S-series**
- a mains cable (for **ASSKEA secretion aspirators S-series**)
- the power supply unit GlobTek GTM91099-6015-3.0-T2 incl. mains cable (for **ASSKEA secretion aspirators M-series**)
- reusable secretion canister system
- multilingual charging instructions (only for **ASSKEA secretion aspirators M-series**)
- instructions for safe handling of battery packs (only for **ASSKEA secretion aspirators M-series**)
- 'Used Medical Device' label and decontamination certificate
- test report according to EN 62353
- optional accessories (depending on the order)

2.3 Product Properties

Hazards of persons due to improper handling

- Use the device for its intended use only.
- Never use the device for aspiration in the low vacuum range (e.g. for thorax drainage).

Damage to the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Do not use the device in case of apparent housing damage.

2.3.1 ASSKEA secretion aspirators S-series

The **ASSKEA secretion aspirators S-series** are especially handy aspiration devices for stationary use.

The devices are operated via a power supply cable to the 230V mains power supply. For **S20** and **S30**, the desired vacuum is generated by an electric motor driven maintenance-free piston pump and for **S20K**, the desired vacuum is generated by a vacuum pump with shadedpole motor (short-time operation). An overtemperature protection

prevents the aspiration device from overheating and leads to automatic switch-off.

After it is switched on, the vacuum pump creates a vacuum in the tubing system and secretion canister, which is used to aspirate secretion (with a suitable sterile aspiration set). The aspirate is directed away from the patient and collected in the secretion canister. If the secretion canister is full, the aspiration is stopped by means of a mechanical overflow protection integrated into the lid of the reusable secretion canister or a self-sealing bacterial filter in the disposable liner 'OneWay'.

The flow rate is adjusted by using the vacuum regulator.

The reusable secretion canister or resp. the holder for the external canister 'Bag' is attached by directly pushing it onto the guiding rails at the housing of the aspiration device (figs. 17, 18). The vacuum connection between the aspiration device and the secretion canister is made via the connecting tube (fig. 1, 2, 3 (C)).

2.3.2 ASSKEA secretion aspirators M-series

The **ASSKEA secretion aspirators M-series** are especially handy, battery powered aspiration devices for mobile use.

All mobile devices are operated via the internal battery or via the supplied power supply unit, that also can be used to recharge the battery. For **M20** and **M20D**, the desired vacuum is generated by an electric motor driven maintenance-free piston pump and for **M30D**, the desired vacuum is generated by a membrane pump. An overtemperature protection prevents the battery from overheating during operation and while charging under unfavorable environmental conditions.

After it is switched on, the vacuum pump creates a vacuum in the tubing system and secretion canister, which is used to aspirate secretion (with a suitable sterile aspiration set). The aspirate is directed away from the patient and collected in the secretion canister. If the secretion canister is full, the aspiration is stopped by means of a mechanical overflow protection integrated into the lid of the reusable secretion canister or a self-sealing bacterial filter in the disposable liner 'OneWay'.

The flow rate is adjusted by using the vacuum regulator. The reusable secretion canister or the holder for the external canister 'Bag' is attached by directly pushing it onto the

guiding rails at the housing of the aspiration device (figs. 17, 18). The vacuum connection between the aspiration device and the secretion canister is made via the connecting tube (fig. 1, 2, 3 (C)).

The **ASSKEA secretion aspirators M-series** are fundamentally divided according to the configuration versions into the following groups:

a) Analog (M20) (fig. 4)

Just like for the stationary aspiration devices, the vacuum is adjusted by using the vacuum regulator (additional air regulation) in the analog, mobile version. All **ASSKEA secretion aspirators M-series** are operated via the internal battery or via the supplied power supply unit, that also can be used to recharge the battery.

b) Digital (M20D, M30D) (fig. 5)

The digital, mobile versions are equipped with a keyboard foil that can be used to adjust the desired vacuum value and to regulate it using an electronic sensor. All **ASSKEA secretion aspirators M-series** are operated via the internal battery or via the supplied power supply unit, that also can be used to recharge the battery.


Fig. 4


Fig. 5

2.3.3 Information on the double filter system

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never use the device without an external bacterial filter.
(exception: ASSKEA disposable secretion canister system)
- Check that the bacterial filter is dry and clean prior to every application to ensure faultless function!
(exception: ASSKEA disposable secretion canister system)

The ASSKEA double filter system consists of the external bacterial filter and the bacterial filter integrated into the aspirator.

For the reusable secretion canister system, a separate external hydrophobic (liquid-repellent) bacterial filter is used.

For the disposable secretion canister system, an external self-sealing bacterial filter is integrated into the disposable liner 'OneWay'. The external bacterial filters (reusable as well as disposable secretion canister system) and the integrated bacterial filter form the ASSKEA double filter system and in combination they are 99.999% effective against bacteria and viruses.

The ASSKEA double filter system effectively protects against overflow and contamination of the interior of the device. It permits fast, simple and cost-effective processing.

Service life and reuse

The **external bacterial filter** (for reusable secretion canister system) is intended for **single use**. For this reason, it must be replaced for each new patient. When used exclusively for one patient, replace the external bacterial filter every 2 weeks at the latest.

When replacing the external bacterial filter, follow the supplied instruction for use!

The **integrated bacterial filter** as well as the external bacterial filter are **not suitable for reuse**. To ensure consistent performance, the integrated bacterial filter system must be replaced, if it has come into **contact with the aspirate** (blocked) or during **maintenance and repair**. This also applies after the filter service life (DFS LED, see below) has expired.

For the digital versions **M20D** and **M30D** of the **ASSKEA secretion aspirators**, an **LED** additionally helps to indicate the filter status of the double filter system (DFS).

If the **DFS LED** illuminates green, the function standby of the filter system is guaranteed. If it illuminates red, the integrated or the external bacterial filter may be blocked or the aspiration is hindered by a blockage in the tubing system. Before replacing the **DFS**, first check the aspiration system for possible blockages.

If the integrated bacterial filter is blocked, it must be replaced by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH.

2.4 Warranty

ATTENTION: It is not permitted to modify the device!

The **ASSKEA secretion aspirators M- and S-series** are covered by warranty for 2 years. It is neither extended nor renewed by warranty work.

The battery is covered by warranty for 6 months.

Wearing parts are excluded from the warranty.

ASSKEA GmbH is responsible for impacts on safety, reliability and specified performance, only if:

- original ASSKEA accessories and spare parts are used,
- maintenance and repair are performed by professionals authorized by ASSKEA GmbH or by ASSKEA itself,
- the affected product is used and operated in accordance with the instruction for use and within its intended use.

All warranty claims are voided, if:

- the device is opened by unauthorized persons,
- the safety seal is removed/broken,
- repairs are carried out by unauthorized persons,
- structural modifications are made to the device,

since the basic safety of the device can no longer be guaranteed and functional limitations may occur in these cases.

3 Operation

Hazard of persons due to improper handling

- Please read through sections 3.1 and 3.2!
- Only perform aspiration in the respiratory tract, if you have been instructed by specialized personnel!
- Only use suitable sterile aspiration catheters or aspiration sets!
- If aspiration procedures are carried out too often, slight bleeding may result.
- When using the devices for the treatment of children, the vacuum must be set within the range of –100 and –200 mbar!

Malfunction due to aspirated secretion

- Make sure that you empty the reusable secretion canister when it is half-full. If the reusable secretion canister is full, the mechanical overflow protection is triggered. This disrupts the aspiration process.
- Ensure that the disposable liner 'OneWay' is replaced timely. The integrated overflow protection is triggered, if the disposable liner 'OneWay' is full. This disrupts the aspiration process.

- Turn off the device when emptying the reusable secretion canister or when replacing the disposable liner 'OneWay'.
- If secretion has been aspirated into the device, it must be properly processed by ASSKEA GmbH or a service partner authorized by ASSKEA GmbH!

3.1 Set-up and Startup


Known or identifiable conditions for medical care within a domestic environment

- Children and pets must be kept away from the device to ensure that the device is not knocked over or dropped.
- Prior to connecting the mains cable, ensure that the mains voltage stated on the device corresponds to the voltage of the mains power supply.
- Do not use the device in damp rooms, bathrooms or showers. Do not allow the mains cable, plug and switch unit to get wet. Never submerge the device in water or other liquids (also not while not in operation).

The operating elements, connections and the startup of the stationary mains-operated aspiration devices (**ASSKEA secretion aspirators S-series**) and of the mobile network-independent aspiration devices (**ASSKEA secretion aspirators M-series**) are explained to you in the following sections:

3.1.1 Operating elements of the ASSKEA secretion aspirators S-series


- A** Secretion canister (reusable) with integrated overflow protection in the lid (see fig. 29)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Vacuum meter
- E** Vacuum regulator
- F** Power switch
- G** Mains connection 230V AC
- H** Mains cable
- I** Suction tube

Fig. 6

3.1.2 Operating elements of the ASSKEA secretion aspirators M-series


a) Analog (M20)

- A** Secretion canister (reusable) with integrated overflow protection in the lid (see fig. 29)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Vacuum meter
- E** Vacuum regulator
- F** Power switch
- G** Socket for power supply unit (12V DC)
- H** Power supply unit (GlobTek GTM91099-6015-3.0-T2)
- I** Mains cable
- J** Suction tube

Fig. 7

b) Digital (M20D, M30D)

- A** Secretion canister (reusable) with integrated overflow protection in the lid (see fig. 29)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Buttons for vacuum adjustment
- E** ON/OFF button
- F** Max. button
- G** Status indicator
- H** Socket for power supply unit (12V DC)
- I** Power supply unit (GlobTek GTM91099-6015-3.0-T2)
- J** Mains cable
- K** Suction tube

Fig. 8

3.1.3 Connecting the ASSKEA secretion aspirators S-series

Check whether the stated mains cable (see chapter 7 'Technical Data') is present and whether it shows potential damages. Replace the mains cable immediately, if there is any damage.

For the mains-operated aspiration devices, use the mains connection (fig. 6 (G)) to connect the device to the mains power supply with the aid of the mains cable. Make sure that the device is positioned in such a way that it is possible to easily disconnect it later. Information on the permissible environmental conditions during operation can be found in chapter 7 'Technical Data'.

3.1.4 Connecting the ASSKEA secretion aspirators M-series

Check whether the stated power supply unit incl. mains cable (see chapter 7 'Technical Data') is present and whether it shows potential damages. Replace the power supply unit incl. mains cable immediately, if there is any damage.

For the network-independent devices, use the socket for the 12 V power supply unit (fig. 7 (G), fig. 8 (H)) to connect the device to the mains power supply with the aid of the supplied power supply unit to charge or operate it as needed. Make sure that the device is positioned in such a way that it is possible to easily disconnect it later. Information on the permissible environmental conditions during operation can be found in chapter 7 'Technical Data'.

3.2 Startup

The following explanations for startup apply to all versions of the **ASSKEA secretion aspirators M- and S-series**.

- Remove the device and the accessories from the packaging.
- Always place the device on a sturdy and flat surface.
- Prior to the initial startup it is essential that you follow the safety instructions in chapter 1.6.
- Always have one backup external bacterial filter ready, since it is absolutely necessary for safe operation (does not apply for disposable secretion canister system).

3.2.1 Connecting the reusable secretion canister system

1. Screw the lid of the reusable secretion canister sufficiently tight and attach the canister to the device.


Fig. 9

2. Insert the connecting tube with the blue elbow fitting toward the side of the device. Make sure that the side of the external bacterial filter with the lettering 'IN' is facing the canister (directed toward patient)!


Fig. 10

3. Place the grey elbow adapter of the connection tube onto the connection in the middle of the reusable secretion canister lid (above the overflow protection).


Fig. 11

4. Place the suction tube onto the second connection (patient connection). Ensure that the tube is fixed sufficiently tight.

The reusable secretion canister or resp. the holder for the external canister 'Bag' is attached by directly pushing it onto the guiding rails at the housing of the aspiration device (figs. 12, 13).

Please notice that up to four mounting positions are usable for the secretion collection canister (figs. 12, 13) on all devices.


Fig. 12


Fig. 13

3.2.2 Connecting the ASSKEA disposable secretion canister system

Malfunction due to collapsing disposable liner 'OneWay'

A leak in the lid or resp. damage of the canister or the T-piece may cause air to flow into the external canister 'Bag'. This may lead to the collapse of the disposable liner 'OneWay'.

- Inspect the disposable secretion canister system to ensure that the lid of the disposable liner 'OneWay' is firmly connected to the external canister 'Bag'.
- Ensure that all connections are firmly attached and properly connected.
- Ensure that the external canister 'Bag' is undamaged and the T-piece is firmly attached.
- Follow the instruction for use supplied by the manufacturer!

The original ASSKEA disposable secretion canister system consists of the external canister 'Bag', the holder for the external canister 'Bag', the disposable liner 'OneWay', the connecting and the suction tube.


Fig. 14

Connection designation

A Vacuum connection

B Patient connection

Please also follow the instruction for use supplied with the disposable secretion canister system!


Fig. 15

1. Remove the disposable liner 'OneWay' from the packaging and fully extend it.


Fig. 16

2. Place the disposable liner 'OneWay' in the reusable external canister 'Bag'. Press the lid's edges firmly down to ensure proper sealing.


Fig. 17

3. Attach the prefitted connecting tube of the disposable liner 'OneWay' to the bottom end of the T-piece located at the external canister 'Bag'.


Fig. 18

4. Connect the vacuum connection of the device with the corresponding vacuum connection of the external canister 'Bag' (top end of the T-piece). Use the supplied connecting tube to do so.


Fig. 19

5. Connect the patient connection of the disposable liner 'OneWay' (fig. 14 (B)) to the suction tube. Ensure that the tube is firmly attached.

Alternative disposable secretion canister systems and holders are available on request.

3.2.3 Suction tube

For all of the **ASSKEA secretion aspirators M- and S-series** you can connect the suction tube with the aspiration catheter using a fingertip.

3.3 Operation of the ASSKEA secretion aspirators S-series

3.3.1 Setting the vacuum


Fig. 20

1. Set the power switch (fig. 1 (F)) to 'I'.
2. Keep the opening of the suction tube (fig. 1 (I)) and the fingertip (fig. 21) shut.
3. Turn the vacuum regulator (fig. 1 (E)) until the vacuum meter (fig. 1 (D)) displays the desired vacuum. Please note that higher vacuum values are set by


Fig. 21

- turning the vacuum regulator clockwise and that the indicator of the manometer moves counterclockwise at the same time.
4. You can also perform fine adjustments to the set vacuum via the additional air opening on the fingertip.

3.3.2 Aspiration

1. Perform the aspiration procedure as you were instructed by specialized personnel.
2. Control the aspiration procedure and the desired vacuum via the additional air opening on the fingertip (fig. 21).
3. Decommission the device after the aspiration procedure according to chapter 3.5.

Please note that the start-up of the pump is not possible from an existing vacuum of -0.5 bar when switching on again.

Please note that the **S20K** is only designed for **short-term** operation (20 minutes)! Please also note the Technical Data for this.

In any case, please take note of further important information in section 2.3.3 regarding the handling of the ASSKEA double filter system!

3.4 Operation of the ASSKEA secretion aspirators M-series

3.4.1 Display and status LEDs

a) Analog secretion aspirator M20 (see fig. 22)


Fig. 22

Symbol	Description	Meaning of the display
	LED illuminates permanently green	Power supply unit is connected; voltage available
	LED illuminates permanently yellow	Power supply unit is connected, power supply permissible but does not correspond to the intended input voltage of 12V DC. → Battery is being charged!
	LED illuminates permanently red	Input voltage is outside the permissible range, automatic switch-off to protect the device.
	LED illuminates permanently green	Battery 100% to 40%
	LED illuminates permanently yellow	Battery capacity < 40% to 10%
	LED illuminates permanently red	Battery capacity < 10% to 0%, immediately plug in the external power supply. After completion of the aspiration procedure, turn off the device for charging!
	LED illuminates permanently red in spite of external power supply	Battery is faulty, operation only possible via plugged in power supply unit
	LED blinks green	Battery is being charged!
	LED blinks red	Overtemperature protection, interruption of the charging/discharging process until the temperature of the battery is in the permissible range again.
	LED illuminates permanently green	Pump is switched on.

b) Digital secretion aspirators M20D and M30D (see fig. 23)


Fig. 23

Symbol	Description	Meaning of the display
	LED illuminates permanently green	Power supply unit is connected; voltage available
	LED illuminates permanently yellow	Power supply unit is connected, power supply permissible but does not correspond to the intended input voltage of 12V DC. → Battery is being charged!
	LED illuminates permanently red	Input voltage is outside the permissible range, automatic switch-off to protect the device.
	LED illuminates permanently green	Battery 100% to 40%
	LED illuminates permanently yellow	Battery capacity < 40% to 10%
	LED illuminates permanently red	Battery capacity < 10% to 0%, immediately plug in the external power supply. After completion of the aspiration procedure, turn off the device for charging!
	LED illuminates permanently red in spite of external power supply	Battery is faulty, operation only possible via plugged in power supply unit
	LED blinks green	Battery is being charged!
	LED blinks red	Overttemperature protection, interruption of the charging/discharging process until the temperature of the battery is in the permissible range again.
	LED illuminates permanently green	Double filter system operational
	LED illuminates permanently yellow	Filter monitoring time has been reached. The integrated bacterial filter should be replaced and the monitoring time should be reset promptly. For this purpose, return the device to ASSKEA GmbH or to one of the service partners authorized by ASSKEA GmbH.
	LED illuminates permanently red	DFS® blocked, replacement of the double filter system must be performed, especially the integrated bacterial filter has to be replaced. Adherence detection. The tubing system has to be checked for blockages.

Symbol	Description	Meaning of the display
	LED illuminates permanently green	Pump is switched on.
LED -0.1 bar to -0.8 bar	LED illuminates red for 2 seconds	Display of the set or preset vacuum (target value display). Factory-provided default: -0.4 bar.
	LED illuminates permanently green	Display of current vacuum (actual value display).

3.4.2 Setting the vacuum

a) Analog secretion aspirator M20


Fig. 24


Fig. 25

1. Set the power switch (fig. 2 (F)) to 'I'.
2. Keep the opening of the suction tube (fig. 2 (J)) and the fingertip (fig. 25) shut.
3. Turn the vacuum regulator (fig. 2 (E)) until the vacuum meter (fig. 2 (D)) displays the desired vacuum. Please note that higher vacuum values are set by turning the vacuum regula-

tor clockwise and that the indicator of the manometer moves counterclockwise at the same time.

4. You can also perform fine adjustments to the set vacuum via the additional air opening on the fingertip.

b) Digital secretion aspirators M20D and M30D

The digital network-independent versions operate vacuum-controlled. The vacuum is set via the keyboard foil on the front side of the devices, whereat the tolerance for the accuracy of the vacuum display is 10%. This technology enables efficient and economical operation with the energy resources during battery operation.


Fig. 26

1. Switch the device on by pressing the ON/OFF button (fig. 26 (1)). The ON LED signalizes that the pump is switched on by illuminating green.

Immediately after switching on, one of the 8 vacuum display LEDs (fig. 3 (G)) illuminates red for 2 seconds to show the last set or during initial operation preset vacuum value (target value display).

2. Set your desired vacuum via the arrow keys (fig. 27 (2)). When the device is switched off later, the last set value is saved.


Fig. 27

3. The display via the red LED signalizes the selected target value of the vacuum to you. The automatic switch over to the current vacuum via a green LED occurs after 2 seconds.

4. You can also perform fine adjustments to the set vacuum via the additional air opening on the fingertip.


Fig. 28

For immediate maximum vacuum:

- If you require the maximum vacuum quickly, press the 'max' button (fig. 29 (3)). You can also set the maximum vacuum by pressing the 'up' button several times.
- The preselection of the maximum vacuum is signalized via a red LED for 2 seconds (fig. 29 (4)).


Fig. 29

3.4.3 Aspiration

1. Perform the aspiration procedure as you were instructed by specialized personnel.
2. Control the aspiration procedure and the desired vacuum via the additional air opening on the fingertip (example in fig. 28).
3. Decommission the device after the aspiration procedure according to chapter 3.5.

Please note that, when using the analog aspiration device, the start-up of the pump is not possible from an existing vacuum of -0.5 bar when switching on again.

In any case, please take note of further important information in section 2.3.3 regarding the handling of the ASSKEA double filter system!

During the aspiration procedure, ensure that the battery has a sufficient capacity. The display of the battery condition must signalize green or yellow. If the display turns red, the external power supply unit must be plugged in immediately. In this case, the aspiration can be continued without any interruption.

3.5 Decommissioning

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never use the device without an external bacterial filter. (exception: ASSKEA disposable secretion canister system)
- Replace the external bacterial filter according to the used canister system every 2 weeks at the latest when used exclusively for one patient!
- Wear appropriate disposable gloves in doing so.
- The reuse of bacterial filters for multiple patients is prohibited for safety reasons!
- For each new patient, processing by ASSKEA GmbH or an authorized service partner of ASSKEA GmbH is strictly required!
- After each aspiration, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- The disposal of aspirate and contaminated components must be performed in a proper manner.

1. Switch off the device after the aspiration by setting the power switch (fig. 1 (F), fig. 2 (F)) to 'O' on analog secretion aspirators or resp. by pressing the ON/OFF button (fig. 3 (E)) for 3 seconds (confirmed by an acoustic signal) on digital secretion aspirators.
2. Disconnect the power supply unit or the mains cable from the electrical connection (100V to 240V AC) and afterwards disconnect the device plug from the secretion aspirator.
3. Remove the suction tube (fig. 1 (I), 2 (J), 3 (K)) from the aspiration catheter and the secretion canister.
4. Clean (see chapter 4.1) the connecting tube and dispose of (see chapter 6.3) or clean (see chapter 4.1) the suction tube in a proper manner.
5. Empty and clean the reusable secretion canister according to chapter 4.1.3. If you use the disposable canister system, clean the external canister 'Bag' according to the instructions in chapter 4.1.4 and dispose of the disposable liner 'OneWay' in a proper manner according to chapter 6.3.
6. Clean the surface of the device according to chapter 4.1.2.
7. Store the device in the shipping carton until next operation.

4 Maintenance

4.1 Cleaning and Care

4.1.1 General Information

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never use the device without an external bacterial filter (exception: ASSKEA disposable secretion canister system).
- Replace the external bacterial filter according to the used canister system every 2 weeks at the latest when used exclusively for one patient!
- Wear appropriate disposable gloves in doing so.
- The reuse of bacterial filters for multiple patients is prohibited for safety reasons!
- For each new patient, processing by ASSKEA GmbH or an authorized service partner of ASSKEA GmbH is strictly required!
- After each aspiration, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.

- The disposal of aspirate and contaminated components must be performed in a proper manner.

Health risks due to the handling of disinfectants

- The use of appropriate protective clothing during disinfection is recommended.
- Follow the manufacturer's disinfectant instructions.

Risk of damage to the device due to improper cleaning agents

- Do not use disinfectants that contain acetone. These may damage or disfigure the housing components as well as the secretion canister lids.
- Follow the instruction for use provided by the manufacturers of the utilized disinfectants particularly with respect to material and surface compatibility as well as the concentration information.
- ASSKEA GmbH recommends Sekusept® aktiv for immersion disinfection of the accessories and Incidin® Plus resp. Incidin® Liquid for wipe disinfection of the device.

Disinfection is not mandatory, if the device is used for one patient only (especially in in-home care).

Disinfection is mandatory, if used in an inpatient setting!

4.1.2 Cleaning and disinfection of the surface of the device

Clean on a regular basis and disinfect at least once a week the surfaces of the device.


- All of the devices can be wiped with a damp, lint-free cloth.
- Follow the previous chapter 4.1.1 for wipe disinfection.

After prolonged use, minor discolorations may occur on the plastic housing components. These do not affect the function of the device, however.

If liquids come into direct contact with the interior of the device, the device must be inspected by ASSKEA GmbH or an authorized service partner of ASSKEA GmbH.

4.1.3 Cleaning and disinfection of the reusable secretion canister system

1. Empty the reusable secretion canister and dispose of the aspirate in a proper manner.
2. Remove the ball (overflow protection) before performing the following cleaning instructions (fig. 30 (C)).

3. Rinse all components of the reusable secretion canister system under running water.
4. Immerse all components of the reusable secretion canister system in disinfectant solution, considering the specific concentration quantity. The reusable secretion canister system can also be autoclaved at 121°C for 20 minutes.
5. Afterwards, rinse the components thoroughly and allow them to dry. When using the cleaner 'neodisher AN' (manufacturer DR. WEIGERT), the cleaning in a special dishwasher is also possible, considering the manufacturer's instructions.
6. Assemble the dry components according to fig. 30.

If the canister system is used and disinfected frequently, ASSKEA GmbH recommends replacing it **every 4 weeks**!

4.1.4 Cleaning and disinfection of the external canister 'Bag'

1. Rinse the external canister 'Bag' under running water.
2. Immerse the external canister 'Bag' in disinfectant solution, considering the specific concentration quantity. You can also autoclave the external canister 'Bag' at 121°C

for 20 minutes.

3. Afterwards rinse the external canister 'Bag' thoroughly and allow it to dry.

ASSKEA GmbH recommends to replace the external canister 'Bag' **every 4 weeks** at the latest and for each new patient.

4.1.5 Cleaning of the connecting tube

1. Applies only for ASSKEA reusable secretion canister system: Remove the external bacterial filter prior to cleaning the connecting tube.
2. Rinse the silicone tube with clear water and place it in the disinfectant solution recommended by ASSKEA GmbH for immersion disinfection.
3. Afterwards, thoroughly rinse the silicone tube with clear water and allow it to dry!

ASSKEA GmbH recommends to replace the connecting tube **every 4 weeks** at the latest and for each new patient!

4.1.6 Cleaning/disposal of the suction tube

1. Dispose of the suction tube in a proper manner for each new patient!
2. If used for one patient only, rinse the suction tube with clear water after each aspiration and place it once a day in the disinfectant solution recommended by ASSKEA GmbH for immersion disinfection.
3. Afterwards, thoroughly rinse the suction tube with clear water and allow it to dry!

ASSKEA GmbH recommends to replace the suction tube **every 4 weeks** at the latest, even when used for one patient only, since the material can become discolored and brittle due to frequent applications!

4.1.7 Disposable products

Dispose of all articles intended for single use (such as the disposable liner 'OneWay' as well as the external bacterial filter) in a proper manner (see chapter 6.3).

4.2 Maintenance and Service

Risk due to maintenance and service during operation

It is not permitted to perform maintenance and service while the device is in use!

Risk due to the performance of maintenance and service by non-authorized persons

Maintenance and service must only be performed by specialized personnel authorized by ASSKEA GmbH or by ASSKEA GmbH itself.

ATTENTION: It is not permitted to modify the device!

- The **ASSKEA secretion aspirators M- and S-series** are maintenance-free, if used according to the instruction for use, exempt from components with limited durability.
- Perform a visual and functional check prior to each application. Also include the accessories of the correspondent device in doing so.
- Opening of a device and repairs must only be performed by ASSKEA GmbH or by specialized personnel authorized by ASSKEA GmbH in compliance with the service documentation specified by the manufacturer as well as with technical and hygienic precautionary measures.

- The device may be sent back for repair to ASSKEA GmbH directly or via the specialty dealer you purchased the device from.
- Clean and disinfect all accessories prior to returning the device. The device itself must be treated with a surface disinfectant.
- Please notify ASSKEA GmbH prior to returning the device to prevent time lags and mark the shipping carton with the 'Used Medical Device' label.

The corresponding product return form with the 'Used Medical Device' label can be found on the web page of ASSKEA GmbH at www.asskea.de in the service area under product return.

ASSKEA GmbH neither ensures proper functioning of the **ASSKEA secretion aspirators M- and S-series** nor is ASSKEA GmbH liable for property damage or personal injury, if

- original ASSKEA accessories or spare parts are not used,
- the user instructions in this instruction for use are not followed,
- installation, settings, modifications, enhancements and repairs are not performed by ASSKEA GmbH or personnel authorized by ASSKEA GmbH.

4.3 Testing of the ASSKEA secretion aspirators M- and S-series

ASSKEA GmbH offers its partners and customers fast and proper processing as well as required testing services.

5 Troubleshooting

5.1 ASSKEA secretion aspirators S-series

Malfunction	Probable causes	Remedy
Device does not start.	Mains plug not plugged in	Check electrical plug connection.
	No mains voltage	Please contact service! (fuse defective)
	Device is too warm, temperature protection of the motor has been triggered.	Allow the device to cool down, temperature protection automatically resets after cooling.
Insufficient flow rate	External bacterial filter polluted, moistened	Replace the external bacterial filter.
	Integrated bacterial filter polluted, moistened	Please contact service!
	Leak in the tubing system	Check that secretion canister lid and tubing connections are tight and no leakages are apparent.
	Vacuum regulator (additional air valve) not closed	Close vacuum regulator (additional air valve) (fig. 1 (E)).
No flow rate	Manual overflow protection of reusable secretion canister system closed	Check fluid level in the secretion canister, empty the canister.
	Hydrophilic overflow protection of the disposable canister system blocked	Replace the disposable liner 'OneWay'.
	External bacterial filter blocked	Replace the external bacterial filter.
	Integrated bacterial filter blocked	Please contact service!
	Suction tube blocked	Clean/replace suction tube.

5.2 ASSKEA secretion aspirators M-series

Malfunction	LED display	Probable causes	Remedy
Device does not start.	 DFS  ON LEDs do not illuminate	Battery is discharged.	Connect the power supply unit for charging.
	 LED illuminates red	Impermissible input voltage	Connect the correct power supply unit, check the power supply.
Device turns itself off (applies to M20).	 LED blinks red, beeper on after 10s	Temperature of the battery too high	Switch off the device and allow it to cool down. Connect the power supply unit and switch on the device again.
	 ON LED switches off		
Motor turns itself off (applies to M20D, M30D).	 LED blinks red, beeper on after 10s, then the device turns itself off after a further 20s	Temperature of the battery too high	Switch off the device and allow it to cool down. Connect the power supply unit and switch on the device again.
	 ON LED switches off		
Device only operates with maximum vacuum (applies to M20D, M30D).	 LED illuminates red	Battery is discharged.	After the aspiration process, switch off the device for charging.
	 LED illuminates green when power supply unit is connected		
Device does not charge.	 LED blinks red	Temperature of the battery too high for charging process	Allow the device to cool down. Automatic resumption of charging when battery has cooled down.
	 LED illuminates green		
	 LED illuminates red	Battery is defective.	Please contact service!
	 LED illuminates green		

Malfunction	LED display	Probable causes	Remedy
Insufficient flow rate	 LED illuminates red	Battery is weak.	Connect the power supply unit for charging.
	 LED does not illuminate		
	 LED illuminates green/yellow	External bacterial filter blocked	Replace the external bacterial filter.
		Integrated bacterial filter blocked	Please contact service!
No flow rate	 LED illuminates green	Leak in the tubing system or in the canister system	Check that secretion canister lid and tubing connections are tight and no leakages are apparent.
	 LED illuminates green/yellow	Float ball of the overflow protection blocks the aspiration lid inlet.	Check fluid level in the secretion canister, empty the canister.
		Hydrophilic overflow protection of the disposable canister system blocked	Replace the disposable liner 'OneWay'.
	 LED illuminates green	External bacterial filter blocked	Replace the external bacterial filter.
	 DFS	Integrated bacterial filter blocked	Please contact service!
		Suction tube blocked	Clean/replace suction tube.

6 Transport, Storage and Disposal

6.1 Decontamination prior to Shipment

Prior to passing on the **ASSKEA secretion aspirators M- and S-series**, the devices must be properly processed by ASSKEA GmbH or by specialized personnel authorized by ASSKEA GmbH to protect subsequent users. Processing is mandatory in accordance with the German Medical Device Operator Ordinance (MPBetreibV), German Medical Devices Act (MPG) and the manufacturer's instructions.

ASSKEA GmbH offers its partners and customers fast and proper processing as well as required testing services.

The **ASSKEA secretion aspirators M- and S-series** must be cleaned and disinfected prior to shipment to ASSKEA GmbH. Please follow the instructions in chapter 4.1! Please affix the supplied 'Used Medical Device' label to the shipping carton! Please give ASSKEA GmbH advance notice of your product return. The product return form can be found on our web page www.asskea.de in the service area under product return.

6.2 Storage

Store the **ASSKEA secretion aspirators M- and S-series** as indicated in the Technical Data (chapter 7)!

The battery of the **ASSKEA secretion aspirators M-series** must be charged prior to storage of the device. This ensures that the device is operational at all times.

Fully recharge the battery, if the **ASSKEA secretion aspirator M-series** is not used for a longer period of time (approx. 10 months).

6.3 Disposal

- The components of the device must be disposed of in a proper manner at the end of the product's service life.
- Ensure that the disposed components are clean and carefully sorted by material.
- The housing material has a material symbol mark and is fully recyclable.
- Decontaminate the device and the accessories prior to disposal.
- According to EU Directives 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) and 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II), the device must not be disposed of in domestic waste.
- The device and accessories may be disposed of via ASSKEA GmbH or the service partner.
- Outside of the EU: Follow the disposal requirements of your country!

7 Technical Data

7.1 Technical Data of the ASSKEA secretion aspirators S-series

Flow rate* (measuring point at suction tube nozzle)	ASSKEA S20: 18 l/min $\pm$ 2 l/min (low flow) ASSKEA S20K: 18 l/min $\pm$ 2 l/min (low flow) ASSKEA S30: 28 l/min $\pm$ 2 l/min (high flow)
Vacuum	max. -90kPa (~ -900mbar; -675 mmHg) (high vacuum) Conversion factor: 1kPa ~ 7.5 mmHg
Canister	Reusable secretion canister system (1000 ml); Disposable secretion canister system (1000 ml)
Suction tube	PVC, inner diameter 6 mm, wall thickness 2 mm, length 150 cm, vacuum stable up to -0.9 bar
Mains voltage	230 V AC
Mains frequency	50 Hz
Power consumption	ASSKEA S20K: 156 W ASSKEA S20: 80 W ASSKEA S30: 80 W
Current consumption at 230 V	ASSKEA S20K: 0.72 A ASSKEA S20: 0.35 A ASSKEA S30: 0.35 A
Dimensions basic device (H x W x D)	290 mm x 259 mm x 130 mm
Weight (basic device)	ASSKEA S20K: 3.5 kg ASSKEA S20: 3.4 kg ASSKEA S30: 4.1 kg
Operating time	ASSKEA S20K: 20 minutes ON, 40 minutes OFF ASSKEA S20: continuous ASSKEA S30: 45 minutes ON, 10 minutes OFF

Degree of protection pursuant to IEC 60601-1	Type BF
Risk classification pursuant to 93/42/EEC, IX	Ila
Protection class pursuant to IEC 60601-1	II
IP degree of protection	IP22
CE marking	CE0843
Sound emission	without vacuum: ASSKEA S20K: approx. 66 dB(A) ASSKEA S20: approx. 55 dB(A) ASSKEA S30: approx. 57 dB(A)
Ambient conditions	Transport/storage: -20°C to +60°C Operation: 0°C to +35°C Relative humidity: 5 to 80 % non-condensing Air pressure: 800 hPa...1060 hPa (at maximum 2000 m altitude above sea level)
Mains cable	H05VV-F 3G1.0 mm ² , length: 3 m, Plastro Mayer
Item number (REF)	ASSKEA S20K: 100011-3 ASSKEA S20: 100009-3 ASSKEA S30: 100010-3

* The information provided may differ depending on the altitude above sea level, the prevailing air pressure and the air temperature.

7.2 Technical Data of the ASSKEA secretion aspirators M-series

Flow rate* (measuring point at suction tube nozzle)	ASSKEA M20: 18 l/min $\pm$ 2 l/min (low flow) ASSKEA M20D: 18 l/min $\pm$ 2 l/min (low flow) ASSKEA M30D: 26 l/min $\pm$ 2 l/min (high flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg) (high vacuum) Conversion factor: 1 kPa ~ 7.5 mmHg
Canister	Reusable secretion canister system (1000 ml); Disposable secretion canister system (1000 ml)
Suction tube	PVC, inner diameter 6 mm, wall thickness 2 mm, length 150 cm, vacuum stable up to -0.9 bar
Nominal voltage of the power supply unit	In: AC 100 – 240 V ~ 50-60 Hz / 1.5 A Out: DC 12 V / 5A
Maximum load current	5.0 A
Permissible input voltage	12 V
Power consumption at 12 V	ASSKEA M20: 47 W ASSKEA M20D: 57 W ASSKEA M30D: 52 W
Power input at 12 V	ASSKEA M20: 3.5 A ASSKEA M20D: 3.5 A ASSKEA M30D: 3.8 A
Degree of protection pursuant to IEC 60601-1	Type BF
Risk classification pur- suant to 93/42/EEC, IX	Ila
Protection class pursu- ant to IEC 60601-1	II
IP degree of protection	IP22
CE marking	CE0843
Sound emission	without vacuum: approx. 62 dB(A)

Ambient conditions	Transport/storage: -20°C to +60°C Operation: 0°C to +35°C Relative humidity: 5 to 80 % non-condensing Air pressure: 800 hPa...1060 hPa (at maximum 2000 m altitude above sea level)
Battery	10.8 V; 5800 mAh (lithium ion battery)
Charging time, if battery is empty	2.5 h – 3 h
Charging time, if battery is approx. 50% full	1.5 h – 2 h
Power supply unit	GTM91099-6015-3.0-T2, GlobTek Inc.
Mains cable of the power supply unit	H03VVH2-F x0,75 2.5 A 250V~; length: 1.8 m
Dimensions basic device (H x W x D)	290 mm x 259 mm x 130 mm
Weight (basic device)	ASSKEA M20: 3.1 kg ASSKEA M20D: 3.1 kg ASSKEA M30D: 3.7 kg
Operating time	ASSKEA M20: approx. 45 minutes, 10 minutes OFF ASSKEA M20D: approx. 45 minutes, 10 minutes OFF ASSKEA M30D: approx. 40 minutes, 10 minutes OFF
Runtime in battery operation	ASSKEA M20: approx. 45 minutes ASSKEA M20D: approx. 45 minutes ASSKEA M30D: approx. 40 minutes
Item number (REF)	ASSKEA M20: 100003-3 ASSKEA M20D: 100005-3 ASSKEA M30D: 100006-3

* The information provided may differ depending on the altitude above sea level, the prevailing air pressure and the air temperature.

7.3 EMC Information

Medical electrical equipment is subject to special precautionary measures regarding electromagnetic compatibility and must be installed and operated in accordance with the EMC information provided in the accompanying documentation!

Portable and mobile RF communication equipment can affect medical electrical equipment!

The use of accessories and spare parts (incl. transformers and cables) not recommended by ASSKEA GmbH may increase the emission of electromagnetic interference or reduce the electromagnetic immunity of the devices. Damages caused by using non-recommended accessories and spare parts or by improper use are not covered by warranty in any case.

Only use original **ASSKEA** accessories and spare parts!

The use of the recommended accessories and spare parts (incl. transformers and cables) in other devices than the **ASSKEA secretion aspirators M- and S-series** may increase the emission of electromagnetic interference or reduce the electromagnetic immunity. Damages caused by using recommended accessories and spare parts in other devices or by improper use are not covered by warranty in any case.

Only use the accessories and spare parts with the **ASSKEA secretion aspirators M- and S-series!**

The **ASSKEA secretion aspirators M- and S-series** must not be used directly adjacent to or stacked with other devices. If operation adjacent to or stacked with other devices is necessary, monitor the **ASSKEA secretion aspirators M- and S-series** in this configuration to verify proper operation!

The **ASSKEA secretion aspirators M- and S-series** meet the requirements of IEC 60601-1-2/EN 60601-1-2 'Electromagnetic Compatibility – Medical Electrical Equipment'. Electromagnetic interference is therefore reduced to a minimum.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The **ASSKEA secretion aspirators M- and S-series** are intended for operation in the electromagnetic environment specified below. The customer or the user of an **ASSKEA secretion aspirator M- and S-series** must ensure that it is operated in such an environment.

Emissions measurement	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Class B	The ASSKEA secretion aspirators M- and S-series are appropriate for use in all facilities incl. private residences and those directly connected to the public power supply network that also supplies buildings used for residential purposes.
Harmonic emissions pursuant to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions pursuant to IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The **ASSKEA secretion aspirators M- and S-series** are intended for operation in the electromagnetic environment specified below. The customer or user of an **ASSKEA secretion aspirator M- and S-series** must ensure that it is used in such an environment.

Immunity testing	IEC 60601 test level	Compliance level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) pursuant to IEC 61000-4-2	±6kV contact discharge	±6kV contact discharge	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
	±8kV air discharge	±8kV air discharge	
Fast transient electrical disturbances/bursts pursuant to IEC 61000-4-4	±2kV for power supply lines	±2kV for power supply lines	Main supply voltage should be of the same quality as in a typical commercial, living or hospital environment.
	±1kV for input and output lines	±1kV for input and output lines	
Interference voltages/surges pursuant to IEC 61000-4-5	±1kV voltage differential mode	±1kV voltage differential mode	Main supply voltage should be of the same quality as in a typical commercial, living or hospital environment.
	±2kV common mode	±2kV common mode	

Immunity testing	IEC 60601 test level	Compliance level	Electromagnetic Environment – Guidance
Voltage dips, short interruptions and voltage variations on power supply input lines pursuant to IEC 61000-4-11	<5% U_T (>95% dip of U_T) for ½ period 40% U_T (60% dip of U_T) for 5 periods 70% U_T (30% dip of U_T) for 25 periods <5% U_T (>95% dip of U_T) for 5 s	<5% U_T (>95% dip of U_T) for ½ period 40% U_T (60% dip of U_T) for 5 periods 70% U_T (30% dip of U_T) for 25 periods <5% U_T (>95% dip of U_T) for 5 s	Main supply voltage should be of the same quality as in a typical commercial, living or hospital environment. If the user of an ASSKEA secretion aspirator M- and S-series requires continued operation also during mains power outages, it is recommended to power the ASSKEA secretion aspirator M- and S-series from an uninterruptible power supply or a battery.
Magnetic field at power frequency (50/60 Hz) pursuant to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic for a typical location in a typical commercial, living or hospital environment.
NOTE: U_T is the AC mains voltage prior to the application of the test levels.			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The **ASSKEA secretion aspirators M- and S-series** are intended for operation in the electromagnetic environment specified below. The customer or user of an **ASSKEA secretion aspirator M- and S-series** must ensure that it is used in such an environment.

Immunity testing	IEC 60601 test level	Compliance level	Electromagnetic Environment – Guidance
Conducted RF disturbances pursuant to IEC 61000-4-6	3 V _{effective value} 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to the ASSKEA secretion aspirators M- and S-series incl. cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ for 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the data provided by the transmitter manufacturer and d is the recommended separation distance in meters (m). The field strength from fixed RF transmitters should be lower than the compliance level^b in each frequency range as determined by an electromagnetic site survey^a. Interference is possible in the vicinity of devices bearing the symbol shown below.
Radiated RF disturbances pursuant to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment with respect to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength exceeds the compliance level above in the location in which the **ASSKEA secretion aspirator M- and S-series** is used, the **ASSKEA secretion aspirator M- and S-series** should be monitored to verify proper operation. If abnormal performance is observed, additional measures may be necessary, such as e.g. re-orienting or relocating the **ASSKEA secretion aspirator M- and S-series**.
- b For the frequency range of 150 kHz to 80 MHz, the field strength should be lower than 3 V/m.

Recommended protection ratio between portable and mobile RF communications equipment and the ASSKEA secretion aspirators M- and S-series

The **ASSKEA secretion aspirators M- and S-series** are intended for operation in an electromagnetic environment in which the radiated RF disturbances are controlled. The customer or the user of an **ASSKEA secretion aspirator M- and S-series** can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **ASSKEA secretion aspirator M- and S-series** – depending on the output of the communications equipment as indicated below.

Rated power of the transmitter W	Protection ratio based on the frequency of the transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended protection ratio d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W), according to indication of the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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