
DEVILBISS®

DeVilbiss® Humidifier Model 9100D

Instruction Guide

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USE INSTRUCTIONS

CAUTION– Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

This humidifier is intended for single-patient use only. The humidifier may be cleaned and reused by the same person. It must not be used by more than one person.

USE INSTRUCTIONS (FOR SINGLE-PATIENT USE)

This product is intended to be used with continuous positive airway pressure (CPAP), some auto-adjusting CPAP devices, bilevel systems and non-invasive ventilation devices for the treatment of Obstructive Sleep Apnea (OSA).

⚠ WARNING

This product is not recommended for use with auto-adjusting CPAP systems that use snoring for pressure adjustment. The humidifier is known to attenuate snoring, which prevents the auto-adjusting CPAP from detecting most snoring.


⚠ WARNING

This product should not be used with any external heat source not recommended by Sunrise Medical.

IMPORTANT PARTS OF YOUR HUMIDIFIER (Figure 1)

The DeVilbiss Humidifier is supplied fully assembled. After removing the humidifier and tubing from the packaging, inspect the chamber, o-ring seal and tubing for damage before using. Prior to first use, the humidifier should be washed, using the *Care and Cleaning Instructions* on page 6.

The DeVilbiss Model 9100D Humidifier will produce the following humidification levels at the flow rate shown under ambient conditions of 72°F (22°C) and approximately 25% relative humidity.

At the extreme condition of minimum water level and dry ambient air (104°F/40°C, 5% R.H.), the DeVilbiss Humidifier will provide approximately 25% R.H. at the nasal mask.


Figure 1

The DeVilbiss Humidifier consists of the following parts:

1. Chamber Lid
2. Chamber Base
3. O-Ring Seal (Shown assembled in chamber base)
4. Slide Latch (6 latches)
5. 15" (38 cm) 22 mm diameter connector tubing

FILLING INSTRUCTIONS

CAUTION– Use only fresh distilled or sterilized water that is at room temperature. Do not add any medication or other additives to the water.

Ensure that all six slide latches are positioned in the locked (🔒) position and that the o-ring seal is properly positioned. Hold the humidifier vertical with the inlet and discharge ports facing upward. Pour enough water in both ports (**Figure 2**) until the water level in both halves of the chamber is between the minimum and maximum fill lines on the top surface of the humidifier lid. The water level must be between the fill lines across the entire length of the lines (**Figure 3**). The humidifier can also be filled by adding a pre-measured amount of water, between 425ml and 525ml, into both ports.

Before each use, ensure that the humidifier is clean (reference *Care and Cleaning Instructions*) and filled to the appropriate level with fresh distilled or sterilized water that is at room temperature.

⚠️ WARNING

Do not fill the humidifier while it is attached to the CPAP or ventilator. To prevent electrical hazards, do not allow water to contact the CPAP or ventilator.

CAUTION– Inspect the o-ring before each use for damage or deterioration. Replace the o-ring if it is damaged or if the o-ring will not provide a proper seal. Use replacement o-ring part number 9100D-600.

CONNECTION AND SET-UP INSTRUCTIONS

⚠️ WARNING

For safe operation, do not place the humidifier on top of the CPAP or ventilator.

⚠️ WARNING

The humidifier must always be at a lower level than the mask and the CPAP or ventilator.

⚠️ WARNING

The humidifier must be level for proper operation.

If the CPAP or ventilator is small enough to securely fit on the top surface of the humidifier, it can be placed on top of the humidifier (**Figure 4**). Otherwise, the CPAP or ventilator must be placed beside the humidifier (**Figure 5**).

1. Connect one end of the 15" (38 cm) connector tubing to either humidifier port. Connect the free end of the patient circuitry (patient tubing and mask system assembly) to the remaining humidifier port.
2. With the CPAP or ventilator turned off, connect the free end of the 15" (38 cm) connector tubing to the outlet of the CPAP or ventilator.


Figure 2


Figure 3


Figure 4


Figure 5

CARE AND CLEANING INSTRUCTIONS

⚠ WARNING

The humidifier must be cleaned before each use.

1. Disconnect all tubing from the humidifier ports.
2. Disassemble the humidifier chamber by moving each of the slide latches (**Figure 6**) to the unlocked (🔓) position and separate the chamber lid from the chamber base.
3. If the o-ring seal needs to be removed from the chamber base, use a blunt instrument to assist in removal. Use care to prevent cutting or stretching the o-ring during removal.

NOTE– The o-ring seal must be removed if using the dishwasher for cleaning.

4. Wash the tubing, o-ring seal, chamber lid and chamber base in a solution of warm water and mild dishwashing detergent. Rinse the parts thoroughly and allow to air dry before reusing.

Once the humidifier parts are clean and dry, reassemble the humidifier as follows:

1. Replace the o-ring into the o-ring groove along the edge of the chamber base. Ensure that o-ring is positioned in the groove along the entire length.
2. Assemble the chamber lid to the chamber base by first ensuring that each of the slide latches are in the unlocked (🔓) position (**Figure 6**).
3. Place the chamber lid onto the chamber base and slide the six latches into the locked (🔒) position (**Figure 7**). Inspect the o-ring seal and verify that the o-ring is properly seated.
4. Firmly reconnect the 15" connector tubing and the patient circuitry to the humidifier ports.

NOTE– The humidifier chamber lid and base are dishwasher safe if placed in the top rack of the dishwasher. The tubing and o-ring seal are not dishwasher safe. It is recommended that the humidifier be washed separately and that the water temperature be at least 130°F (54°C). Thoroughly rinse the chamber lid and base after removing from the dishwasher and allow to air dry.

CAUTION– Do not use any solutions that contain bleach, chlorine, alcohol, or aromatic-based ingredients, including scented oils, to clean the humidifier assembly or tubing. These solutions may reduce the life of the product. Do not place the tubing or humidifier in direct sunlight, as they may harden and crack.


Figure 6


Figure 7

WARRANTY

The DeVilbiss Humidifier Model 9100D is warranted to be free from defective workmanship and material for a period of one year from date of purchase. Any defective part(s) will be repaired or replaced at Sunrise Medical's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. This warranty does not cover normal wear and tear on the o-ring seal. If repair is necessary, contact your Sunrise Medical provider or Sunrise Medical Service Department for instructions:

U.S.A 800-333-4000 (814-443-4881)

Canada 905-660-2459

Europe 44-138-444-6688

NOTE- This warranty does not cover providing a loaner unit, compensating for costs incurred in rental while said unit is under repair, or costs for labor incurred in repairing or replacing defective part(s).

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

DECLARATION OF CONFORMITY

Manufacturer:

Sunrise Medical

Address:

Respiratory Products Division
100 DeVilbiss Drive
Somerset, Pennsylvania 15501-2125
USA

Product Designation:

CPAP Humidifier

Type, Model:

DeVilbiss® Humidifier, model 9100D

We herewith declare that the above-mentioned product complies with the requirements of EC Directive 93/42/EEC and the following:

Quality System Standards Applied:

ISO9001/EN46001

Notified Body RWTÜV

MDD Annex II Applied

Place:

Somerset, PA

Signature:

Name:

Michael Chellson

Title:

Director of Quality and Regulatory
Affairs

European Contact:

Sunrise Medical Ltd.
Sunrise Business Park
High Street
Wollaston, West Midlands DY8 4PS
ENGLAND
44-138-444-6688

CE 0044

®Registered U.S. Patent and Trademark
Office and other countries.

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®Eingetragen bei dem U.S. Patent and Trademark Office sowie in anderen Ländern.
®Registrerad hos U.S. Patent and Trademark Office och i andra länder.
®Registrert U.S. Patent and Trademark
Sunrise Medical Canada, Inc.
Sunrise Medical Pty. Limited