**EN**  
page 2

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## User Instructions

### Icons in this Manual

In this manual, icons are used to highlight specific information. The meaning of each icon is explained in the table below:

ICON	EXPLANATION
	<b>Warning!</b> Risk of personal injury.

### Intended Use

The patient circuit is intended to be connected to Breas Vivo 50.

The intended use of the patient circuit is to deliver air, with or without added oxygen, invasively or non-invasively, from the Vivo 50 and, in addition to the pressure/volume control regulation provided by the Vivo 50, control the inspiration and expiration phases by opening and closing the exhalation valve.

The patient circuit is intended for home and hospital use.

The patient circuit is intended to be connected to patients and maintained by persons who are trained to use ventilators.

### Safety Information

- Read this manual thoroughly so that you completely understand how the patient circuit is operated before taking it into use, to ensure correct usage and maximum performance.
- Breas Medical AB reserves the right to make changes to this product without any prior notification.
- Do not use a damaged patient circuit.

## Symbols Shown on the Units

The meaning of the symbols shown on the patient circuit label is explained in the table below.

SYMBOL	EXPLANATION
0123	CE marking applies according to the directive for Medical Devices MDD 93/42/EEC.
	Product number
	Batch code
	Read user instructions
	Caution
	Keep away from sunlight
	Does not contain any natural latex rubber.
	Manufacturer

## Material Specifications


PARTS	MATERIALS
Patient tube	Polyester with silicone rubber cuffs
Control pressure tube	Silicone rubber
Exhalation valve body	Polycarbonate
Membrane assembly	Polycarbonate, polypropylene and silicone rubber

## Technical Specifications

MISCELLANEOUS	SPECIFICATIONS
Rated flow for the patient circuit	150 l/min., ≤ 0.2 kPa
Length of the breathing tube, in resting condition (without being held under tension)	1.8 m

## Connect the Patient Circuit to the Vivo 50

- 1 Connect the tube to the patient air outlet on the side panel. If a bacterial filter is used, connect it between the circuit and the patient air outlet.
- 2 Connect the other end of the patient circuit to the exhalation valve.
- 3 Connect the control pressure tube to the membrane nipple of the exhalation valve.
- 4 Connect the control pressure tube to the exhalation valve control pressure outlet on the side panel.

## Cleaning


- Always be careful when cleaning to ensure that you do not damage any equipment.
- Take care not to allow water to enter into the control pressure tube. The inside of this tube does not need to be cleaned.  
Any water in the control pressure tube or the exhalation valve may lead to malfunctions.
- Always clean the parts or use a new set when used by a new patient.  
All parts that come into contact with the respiration gas must be cleaned.

**1** Disconnect the patient circuit from the Vivo 50 and the tube from the exhalation valve.

**2** Unscrew the exhalation valve cover completely together with the membrane assembly. The membrane can be cleaned properly while still fitted to its cover.

Let the control pressure tube remain connected to the nipple to prevent water from entering into the membrane.

**3** Place the dismantled parts in hot water containing mild detergent.

**4** Remove fouling with a brush.

**5** Rinse parts thoroughly under running hot water.

**6** Shake water out of all parts.

**7** Dry the parts completely.

## Disinfection

The patient circuit can be disinfected by the following solutions, provided that the disinfectants are used according to the manufacturer's instructions. Do not disinfect the control pressure tube.

- Gigasept® FF (10% solution), for 15 minutes up to 20 cycles
- Steranios 2%, for 10 minutes up to 20 cycles

Before disinfection, always clean the parts as described above.

After disinfection, thoroughly rinse the parts in running water for 2 minutes, in order to remove any residue of the disinfectant. Make sure to rinse all surfaces. Dry the parts before use.

The effectiveness of this disinfection method has not been validated.  
It is recommended with regard to material durability only.

## **Autoclaving**

The parts can withstand to be autoclaved at 126 °C (258 °F) for 15 minutes. This treatment can be repeated up to 20 times. Do not autoclave the control pressure tube.

The effectiveness of this sterilization method has not been validated.  
It is recommended with regard to material durability only.

## **After cleaning**

Check that there is no visible damage.

## **Recommended Usage Period**

3 to 6 months.

## **Storage**

Store in dust-free location.

The parts shall not be stored where they are exposed directly to sunlight.

## **Disposal**

The patient circuit must be disposed of and recycled in accordance with the local environmental regulations regarding the disposal of used equipment and waste.

## **Ordering Numbers**

### **ORDERING NUMBERS**

Breas Vivo 50 Patient Circuit, Single Limb with Exhalation Valve, Reusable	005055
---	--------

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