





**DURAVENT®**

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EN

## I. Foreword

These instructions are applicable to the Fahl® tracheostomy tubes; they are intended to ensure appropriate use of Fahl® tracheostomy tubes by doctors/physicians, nurses and patients or users.

### **Please read the instructions for use carefully before using the product for the first time!**

Keep them in an easily accessible place for future reference.

Please save the packaging as long as this tracheostomy tube is being used. It contains important product information!

The Fahl® tracheostomy tubes are medical devices as defined by the MPG (German medical devices act) and conform to, among others, the "Basic Requirements" as per Annex I of Directive 93/42 EEC.

## II. Intended Use

Fahl® tracheostomy tubes are intended to stabilise the tracheostoma following laryngectomy or tracheostomy.

The tracheostomy tube is designed to keep the tracheostoma open.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

### **Users must be trained by medical professionals on how to use the Fahl® tracheostomy tubes safely and appropriately.**

The extra long tracheostomy tubes can be necessary specifically in cases of deep tracheal stenoses.

Fahl® LINGO tube variants are intended for tracheostomised patients who retain their larynx or for laryngectomees who use a voice prosthesis (users of shunt valve).

## III. Warnings

### **Patients must be trained by medical professionals on how to use the tube safely.**

Fahl® tracheostomy tubes must never be blocked, e.g. by mucus or encrustation. Blockage can cause asphyxiation!

Mucus in the trachea can be aspirated (removed by suction) through the tracheostomy tube with a tracheal suction catheter.

Damaged tracheostomy tubes must not be used and must be discarded immediately. Use of a damaged tube may result in airway compromise.

When inserting and removing a tracheostomy tube, irritations, coughing or light bleeding can occur. If you experience persisting bleeding consult your treating doctor/physician immediately.

Tracheostomy tubes are single-patient products and intended for use by one patient only.

Do not re-use with another patient. Do not reprocess for re-use by another patient.

Do not use tracheostomy tubes during laser treatment (laser therapy) or with electrosurgical devices! Contact with a laser beam may damage the tube.

In the presence of oxygen-rich mixtures the tube may ignite, which may cause injury from heat and exposure to decomposing and toxic combustion products.

### **CAUTION**

Tracheostomy tubes with speaking function are recommended only for tracheostomised patients with normal secretion and inconspicuous mucus tissue.

Tracheostomy tubes containing metal parts must under no circumstances be used during x-ray examinations or radiotherapy to avoid risk of severe skin damage for example. If it is necessary to wear a tracheostomy tube during such examinations, use a medical grade polymer tracheostomy tube containing no metal parts. In the case of some medical grade polymer tracheostomy tubes with silver speaking valve, for example, the valve and its safety chain can be removed completely by detaching the inner cannula with speaking valve from the outer cannula.

### **CAUTION**

Patients who have heavy secretion outflow and/or tend to develop granulation tissue under radiation therapy or who develop encrustations should only use a sieved cannula variant after having previously consulted their doctor/physician. Regular medical supervision and a short replacement interval (one week) for the tracheostomy tube is highly recommended, because the sieve in the outer cannula can lead to an intensified formation of granulation tissue.

## IV. Complications

The following complications could emerge when using this product:

Impurities (contamination) around the stoma may make it necessary to remove the tube; contamination may cause infections, requiring the use of antibiotics.

Unintentional inhalation of a tube that was not sized correctly may need to be removed by a doctor/physician. If mucus blocks the tube it should be removed and cleaned.

## V. Contraindications

**Do not use if the patient is allergic to the material.**

### CAUTION

**During mechanical ventilation do not use tracheostomy tubes without a cuff**

### CAUTION

**Only use sleeved/fenestrated tracheostomy tubes during mechanical ventilation after consultation with the treating medical doctor/physician.**

### CAUTION

**Laryngectomised (patients without larynx) must not use tracheostomy tubes with a speaking valve. This can lead to severe complications and can cause suffocation!**

## VI. Precautions

The correct tube size should be selected by the doctor/physician or the trained medical professionals.

Use only auxiliary equipment with a 22 mm connection with Fahl® KOMBI adapters to prevent accidental loss of accessories or damage to the tube.

We highly recommend always keeping at least two replacement tubes on hand to ensure continuous supply.

Do not perform any type of adjustment repair or alteration on this tracheostomy tube or compatible accessories. If a product is damaged, it must be discarded immediately.

### CAUTION

**United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).**

## VII. Product description

Duravent® tracheostomy tubes are made of medical quality polyvinyl chloride.

These tracheostomy tubes are made of thermo-sensitive, medical grade polymer, which develops its optimal material characteristics at body temperature.

We supply the Fahl® tracheostomy tubes in various sizes and lengths.

Fahl® tracheostomy tubes are reusable medical products intended for single-patient use.

Fahl® tracheostomy tubes are only to be used by the same patient and not by another patient.

With proper care and use, especially when complying with the guidance for cleaning, disinfection and care in this instruction for use, the service life of the tracheostomy tubes can last for up to 3 to 6 months.

Please see the size table in the appendix.

A particular feature of the Fahl® tracheostomy tube is the specially shaped neck flange.

The tube tip is rounded to prevent irritation of the mucus membranes in the trachea.

There are two lateral fastening eyelets on the neck flange for attaching a tube holder.

It may be advisable to use tubes of different lengths alternately so the tip of the tube does not always come in contact with the same spot in the trachea and cause irritation. This prevents pressure points or the formation of granulation tissue. Always discuss the exact procedure with your doctor/physician.

### MRT (Magnetic Resonance Tomography) information

These tracheostomy tubes can be worn by the patient during magnetic resonance tomography (MRT).

This only applies to tracheostomy tubes without metal parts.

In the case of our tracheostomy tubes with silver speaking valve, for example, the valve and its safety chain can be removed completely by detaching the inner cannula with speaking valve from the outer cannula.

## 1. Neck flange

A particular feature of the Fahl® tracheostomy tube is the specially shaped neck flange, which is designed to fit the anatomy of the neck. The sizes are printed on the neck flange.

On the neck flange of the tracheostomy tube you can find two lateral fastening eyelets for attaching a tube holder.

All Fahl® tracheostomy tubes with fastening eyelets include a disposable tube holder for single use. The tube holder secures the tracheostomy tube to the neck.

Please follow the instructions provided with the tube holder carefully when attaching and removing it from your tracheostomy tube.

Ensure that the Fahl® tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

## 2. Connectors and adapters

Connectors and adapters are used to connect compatible tube accessories.

The options available to a patient depend on the clinical condition, such as status post laryngectomy or tracheostomy.

Usually connectors/adapters are firmly connected to the inner cannula. These are standard 15 mm connectors, to which HMEs (Heat Moisture Exchangers) can be attached.

This specific connector is also available as a 15 mm rotatable variant. The rotatable connector is able to absorb torsion forces, which might occur with artificial ventilation tube systems. Therefore the cannula is released of any strain and also stabilized in its position, which helps to reduce irritations of the mucus membranes.

The 22 mm combi-adapter can be used to fasten compatible filter and valve systems with 22 mm attachments – such as, but not limited to the Humidophone® speaking valve with filter function (REF 46480), Humidotwin® Heat and Moisture Exchanger (HME) (REF 46460), Combiophon® speaking valve (REF 27131), Laryvox® HME filter cassettes (REF 49800).

## 3. Cannula tube

The cannula tube directly adjoins to the neck flange and the air flow is directed into the trachea.

The Easy-Lock hexagonal lock is used for fixing the inner cannula and compatible accessories securely.

The outer cannula features a lateral contrast strip which enables the exact localisation of the tracheostomy tube on x-rays.

## 4. Inner Cannula

By means of the hexagonal Easy-Lock connection, an inner cannula and any other compatible accessory can be securely fastened to the outer cannula.

Depending on the exact specification of each cannula type, the inner cannula is either equipped with a standard 15mm connector, firmly attached with specific adapters/connectors, or can be connected to detachable accessory parts, for example speaking valves.

The inner cannulas can be easily removed from the outer cannula, which enables an increased air supply in emergency cases such as respiratory distress.

In order to remove the inner cannula, it must first be detached from the outer cannula. By lifting the inner cannula slightly with thumb and index finger, the locking lines of the inner cannula are detached from the hexagonal Easy-Lock of the outer cannula.

An inner cannula may under no circumstances be used without the corresponding outer cannula.

## 4.1 Speaking Valves

Tracheostomy tubes with speaking valve (LINGO/PHON) are used for patients after a tracheostomy with an intact or partially intact larynx, thus enabling the user to speak.

Tracheostomy tubes with a silicone speaking valve have the valve directly attached at the inner cannula and it is included in the scope of delivery.

Furthermore, a Humidophone® or Combiophon® speaking valve is available, which can be used by tracheostomised patients who have a sieved tracheostomy tube with inner cannula and 22 mm combi-adapter.

**When using tracheostomy tubes with silver speaking valve, the valve can be detached by sliding it out from the inner cannula.**

**When using tracheostomy tubes with a silicone speaking valve, the valve can be detached by pulling off the valve from the inner cannula.**

## 5. Decannulation Plug

A decannulation plug is included in the scope of delivery of each tracheostomy tube with speaking valve. This decannulation plug may only be used by tracheostomised patients with an intact larynx and only under supervision of a doctor/physician. The decannulation plug allows a temporary sealing of the cannula and blocking of the air supply, which helps patients to relearn breathing through mouth and nose.

## ATTENTION

The decannulation plug is contra-indicated for laryngectomies and patients with COPD (chronic obstructive pulmonary disease). In this case the decannulation plug is not to be used under any circumstances! The decannulation plug may only be used with a sleeved outer cannula without inner cannula.

Under no circumstances is the decannulation plug to be used in combination with a tracheostomy tube with cuff! The decannulation plug may only be used with a sleeved outer cannula without inner cannula.

## CAUTION

In order to re-acustom a patient to normal breathing after a temporary tracheostomy, a decannulation plug is included with tracheostomy tubes with speaking valves. The plug can be inserted into the cannula for a short time so that the patient is encouraged to breathe in a natural manner, using nose and mouth. The decannulation may only be performed under the supervision of a doctor/physician! To avoid risk of suffocation, the decannulation plug may only be inserted in accordance with a doctor's/physician's instruction! Please also note the description of contra-indications for every tracheostomy tube variant!

## PRECAUTIONS

Make sure to read all accompanying product information, instructions for use, indications and contraindications. Discuss use of the product with your doctor/physician before first use.

## VIII. Directions for Tube Insertion and Removal

### For the Doctor/Physician

The correct tube must be selected by the physician or the trained medical professionals.

Select a tube that fits the patient's anatomy to optimise comfort and ventilation (breathing in and out).

The inner cannula can be detached from the outer cannula anytime for an enhanced air supply or for cleaning. Cleaning might become necessary e.g. if the cannula is clogged with secretion residues, which cannot be removed by coughing or by suction.

### For the Patient

## PRECAUTIONS

**Check the use-by or expiry date. Do not use after this date.**

Carefully examine the tube before first use to make sure that it is not damaged and that there are no loose parts.

Should you notice any anomaly or anything unusual, DO NOT use the tube. Return the tube to the manufacturer for inspection.

Tracheostomy tubes must be cleaned thoroughly prior to each insertion. This is also recommended prior to the first insertion unless the tube is sterile, then we do not recommend cleaning the tube before first time use.

Please note that the tube must always be cleaned and, if necessary, disinfected according to the instructions provided below before re-inserting.

If secretion collects in the lumen of the Fahl® tracheostomy tube and cannot be removed by coughing or aspiration, the tube should be removed and cleaned.

After cleaning and/or disinfection, carefully examine the Fahl® tracheostomy tube for sharp edges, cracks, or other signs of damage, since these may impair function and/or injure the mucus membranes in the airways.

Do not use damaged tubes.

If the tube material appears to be discoloured, discard the Fahl® tube immediately.

## WARNING

**Users must be trained by medical professionals in the safe use of the Fahl® tracheostomy tubes.**

### 1. Insertion of the tube

#### Step-by-step instructions to insert Fahl® tracheostomy tubes

Before application the users should clean their hands (see picture 3).

Remove tube from the package (see picture 4).

Next, push a tracheal compress, e.g. Sensotrach® DUO (REF 30608) or Sensotrach® 3-Plus (REF 30780), onto the tracheostomy tube.

In order to facilitate the insertion of the tracheostomy tube into the trachea, we recommend to lubricate the outer cannula with an Optiflu-id® Stoma Oil wipe (REF 31550), which ensures that the oil is applied evenly (see picture 4a and 4b).

If you are inserting the tube yourself, use a mirror to make insertion of the Fahl® tracheostomy tube easier.

When inserting the tube, hold it by the neck flange with one hand (see picture 5).

Pull the tracheostoma slightly apart with your free hand to allow the tip of the tube to fit into the tracheostoma more easily.

Medical aids to help dilating the tracheostoma are also available (e.g. tracheal dilator, REF 35500), which allows for an even and gentle dilating of the tracheostoma. This is recommendable especially in case of an emergency, e.g. a collapsing tracheostoma (see picture 6).

Please proceed with great caution when using a medical device for dilating the tracheostoma, to avoid damaging the cannula through friction.

Now carefully insert the tube into the tracheostoma during the inspiration phase (while breathing in) while tilting your head slightly back (see picture 7).

Advance the tube into the trachea.

Straighten your head once the tube has been inserted further into the trachea.

Tracheostomy tubes should always be fixed with a special tube holder. This tube holder stabilises the tracheostomy tube and ensures a secure fit within the tracheostoma (see picture 1).

## 2. Removing the tube

### PRECAUTIONS

**Accessories such as a tracheostoma valve or HME (Heat Moisture Exchanger) must be removed first before proceeding to remove the Fahl® tracheostomy tube.**

### ATTENTION

**In the case of an unstable tracheostoma or in emergencies, the tracheostoma may collapse after having removed the tracheostomy tube and may lead to an obstruction of the air supply. In this case, a new tracheostomy tube must be held ready for insertion. In order to secure the air supply temporarily, a tracheal dilator (REF 35500) may be inserted.**

Proceed very carefully to avoid injury to the mucus membranes.

The tube must always be cleaned and, if necessary, disinfected as follows before re-inserting according to the instructions provided below.

### Step-by-step instructions to remove the Fahl® tracheostomy tubes:

Remove the tracheostomy tube with the head tilted slightly back. Grip the tube at the side by the neck flange or the housing (see picture 7).

Carefully remove the tracheostomy tubes.

The Duravent® tracheostomy tube with silver speaking valve has a pre-formed guide bar in the valve adapter in which a flat Sterling silver valve plate with a circular lid is inserted. The inner cannula with the adjusted silver speaking valve can be removed from the outer cannula. The silver valve plate does not need to be removed first. A safety chain prevents the valve from being dropped.

In order to remove the inner cannula, it must first be detached from the outer cannula. By lifting the inner cannula slightly with thumb and index finger, the locking lines of the inner cannula are detached from the hexagonal Easy-Lock of the outer cannula.

For reinsertion the inner cannula is inserted in reverse order and attached to the hexagonal Easy-Lock of the outer cannula.

## IX. Cleaning and Disinfection

### PRECAUTIONS

**For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fahl® tracheostomy tube at least twice a day, more often during heavier secretion production.**

### CAUTION

**Do not use a dish washer, steamer, a microwave, a washing machine or similar for cleaning the tracheostomy tube!**

Keep in mind that your personal cleaning schedule, including additional disinfection if required, must always be determined in consultation with your doctor/physician and according to your individual needs.

Regular disinfection is generally only required, if this is indicated due to a medical prescription. The reason being that even upper airways of healthy patients are not germfree.

General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. We recommend chemical disinfection of the tubes according to the instructions provided below. Please consult your doctor/physician.

## CAUTION

**Residues on the tracheostomy tube caused by cleaning- and disinfection solution can lead to irritation of the mucus membranes or other adverse health effects.**

Tracheostomy tubes are to be viewed instruments with cavities, therefore during disinfection or cleaning it is important to ensure that the tube is fully moistened with the solution (see picture 8).

The cleaning and/or disinfection solution has to be prepared freshly each time.

### 1. Cleaning

Fahl® tracheostomy tubes must be cleaned/replaced regularly in accordance with the individual patient's needs.

Only use cleaning agents if the tube is outside the tracheostoma.

A mild, pH-neutral washing lotion can be used for cleaning the tube. We recommend using the special cannula-cleaning powder (REF 31100) as directed by the manufacturer.

Never clean the Fahl® tracheostomy tube with cleaning agents that are not approved by the manufacturer. Never use aggressive household cleaners, detergents, high-concentration alcohol or denture cleaners.

This may be very hazardous to health. The tube may also be destroyed or damaged.

Do not boil the tracheostomy tube. This could severely damage the tube.

### Cleaning Steps

Remove any additional items before cleaning.

First rinse the tube under lukewarm running water (see picture 9).

Use only lukewarm water to prepare the cleaning solution and observe the instructions for use of the cleaning agent.

We recommend the use of the cannula-cleaning tub with sieve insert (REF 31200) to facilitate the cleaning procedure.

Hold the sieve insert by the top edge to prevent contact with and contamination of the cleaning solution (see picture 10).

Place only one tube at a time in the sieve insert of the cannula-cleaning tube. If more than one tube is cleaned at same time, they may be compressed too strongly and damaged by excessive pressure.

Inner and outer cannula can be cleaned together.

Immerse the rinsing element with the cannula components into the cleaning agent.

After the immersion period (see instructions for cannula- cleaning powder), the tube is rinsed thoroughly several times under lukewarm clear water (see picture 9). There must be no residues of the cleaning agent left on the tube when it is inserted into the tracheostoma.

If necessary, for example, if persistent secretion residues cannot be removed by the cleaning bath, additional cleaning with a special cannula- cleaning brush (Optibrush®, REF 31850 or Optibrush® Plus with fibre top, REF 31855) may help. Only use the cleaning brush, if the tube is removed and already outside the tracheostoma.

Always insert the brush into the tube from the tube tip (see picture 11).

Use the brush as directed and proceed with great care to avoid damage to the soft tube material.

If you have a tracheostomy tube with speaking valve, always detach the speaking valve from the cannula prior to the cleaning process. The valve itself must never be cleaned with a cleaning brush as it could break off or could be otherwise damaged.

Thoroughly rinse the tube under lukewarm running water or sterile saline solution (0.9% NaCl solution).

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

Never use tubes with impaired functionality or damages such as sharp edges or cracks, as they may injure the mucus membranes in the airways. If damages are found, do not use the tube.

### 2. Chemical Disinfection Instructions

The Fahl® tracheostomy tube can be disinfected by cold disinfection with special chemical disinfectants.

Disinfection should always be done if ordered by a physician due to specific health concerns caused by disease, infection, or your specific situation.

Disinfecting is generally indicated to prevent cross-infection and in inpatient situations (e.g. hospitals, nursing homes, and/or other health care facilities) to limit infections.

## PRECAUTIONS

**Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).**

**Disinfectants that release chlorine or contain strong alkalis or phenol derivatives must never be used. They may seriously damage or even destroy the tube.**

#### Disinfection steps

Only use the Optic® tube disinfectant (REF 31180; not available in US) in accordance with the manufacturer's instructions. Alternatively we recommend a disinfectant on the basis of glutaraldehyd as active ingredient (also available in the US). Always observe the application range and spectrum of activity quoted by the manufacturer. Follow the instructions for the disinfectant.

After disinfection the tube must be rinsed inside and out very thoroughly with sterile saline solution (NaCl 0.9%) and then dried.

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

#### 3. Sterilization/Autoclave

These tracheostomy tubes can be sterilised with ethylene oxide (gas sterilisation) based on a validated process.

#### CAUTION

**Do not heat the tube over 65°C (149°F), boiling or steam sterilisation are not permitted as this will damage the tube.**

### X. Storage/Care and Maintenance

Currently not used and cleaned tubes should be stored in a clean plastic container, dry environment away from dust, sunlight and heat.

**After having cleaned, disinfected (if applicable) and dried the inner cannula, lubricate the exterior of the inner cannula on with Stoma Oil, e.g. Optifluid® Stoma Oil wipe in order to allow for easy insertion.**

**For lubrication, use only Stoma oil (REF 31500) or Optifluid® Stoma oil wipes (REF 31550).**

We highly recommend keeping at least two replacement tubes on hand to ensure continuous supply.

### XI. Service Life

These tracheostomy tubes are single-patient products.

With proper care and use, especially when complying with the guidance for cleaning, disinfection and care in this instruction for use, the service life of the tracheostomy tubes can last for up to 3 to 6 months.

The service life of a tube is influenced by many factors. The composition of the secretions, the quality of cleaning and other aspects, for example, are very important.

The maximum service life is not extended by using the tube intermittently (e.g. alternately with other tubes at normal cleaning intervals).

No later than after six months of use or no later than six months after opening the sterile packaging they must be replaced - or immediately if damaged - and must be disposed of correctly.

#### CAUTION

**Any changes, in specific, shortening and sieving, including repairs to the cannula, may only to be undertaken by the manufacturer or companies which are expressly authorised in writing by the manufacturer. Unprofessional work on tracheostomy tubes can lead to severe injuries.**

### XII. General Terms and Conditions

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, the Andreas Fahl Medizintechnik-Vertrieb GmbH will not take liability for damages, caused by alterations, mainly for shortening and sieving, through repairs, that were not undertaken by the manufacturer himself. This applies to damages caused to the tracheostomy tubes as well as for any caused consequential damage.

When using a tracheostomy tube for longer than the application period mentioned in Figure XI and/or using the tube contrary to the specifications, cleaning, care or storage directions as given in these instructions of use, the Andreas Fahl Medizintechnik-Vertrieb GmbH is released from any liability including the liability for defects – as far as legally permissible.

Sales and deliveries of all Andreas Fahl Medizintechnik-Vertrieb GmbH products are carried out exclusively according to our General Terms and Conditions of Business which are available directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

Duravent® is a trademark of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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
### Sizing information:

**O.D. Shield** = outer diameter behind the neck flange · **O.D. Tip** = outer diameter at the tube tip · **I.D.** = inner diameter at cannula tip · **O-tube** = outer cannula · **I-tube** = inner cannula · **Length** = length over the outer curve

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All measurements listed in the tables below are for standard cannulas.  
Customised cannulas are available upon request.

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\* available with or without inner cannula: UNI – with 15-mm connector; KOMBI – with 22 mm combi-adapter, size 5 and up, speaking cannula – depending on model – size 5 and up; multi-functional cannula, size 5 and up

\*\* available with or without inner cannula: UNI – with 15-mm connector; KOMBI – with 22 mm combi-adapter, size 5 and up, speaking cannula

\*\*\* available with or without inner cannula: KOMBI – with 22 mm combi-adapter, available only in 44-mm/62-mm lengths; LINGO – with sieve, available only in 44-mm/62-mm lengths

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