PORTEX™

0473

en EzPAP® Positive Airway Pressure Therapy System**fr Système de thérapie par pression expiratoire positive EzPAP®****de EzPAP® Atemwegs-Überdrucktherapiesystem****it Sistema terapeutico di ventilazione a pressione positiva EzPAP®****es Sistema de terapia EzPAP® de presión positiva en vías respiratorias****pt Sistema de Terapia de Pressão Positiva das Vias Respiratórias EzPAP®****nl EzPAP® positief luchtwegdruk-therapiesysteem****sv EzPAP®-system för PAP-behandling (Positive Airway Pressure)****da EzPAP® system til PAP-behandling****Smiths Medical ASD, Inc.**

Keene, NH 03431, USA

EC Smiths Medical International Ltd.

Hythe, Kent, CT21 6JL, UK

Australian Representative:

Smiths Medical Australia Pty. Ltd.

Brisbane, QLD 4113, Australia

www.smiths-medical.com

600-450-001 REV. 007 / 10

smiths medical

en English**1. DESCRIPTION:** The EzPAP® Positive Airway Pressure Therapy System is a non-invasive single patient use, hand held respiratory therapy device. EzPAP® Positive Airway Pressure Therapy System has the ability to provide supplemental oxygen when used with compressed oxygen. The device is also available with various mask sets.**2. INDICATIONS FOR USE:** The EzPAP® Positive Airway Pressure Therapy System is indicated for the treatment and prevention of apnea. The EzPAP® Positive Airway Pressure Therapy System is indicated for patients who are capable of following directions for Positive Airway Pressure (PAP) therapy.**3. CONTRAINDICATIONS:** There have been no contraindications to the use of PAP therapy have been reported, the following should be carefully evaluated before a decision is made to initiate therapy:

- Inability to tolerate increased breathing
- Hyperventilation syndrome
- Recent facial, oral or skull surgery or trauma
- Eosophagitis surgery
- Uncontrolled hypertension

- Known or suspected tamponade membrane rupture or other middle ear pathology
- Hemodynamic instability
- Deep sinusitis
- Epistaxis
- Active hemoptysis

4. WARNINGS: Use of this device or respiratory pressures may have adverse effects. Expert clinical judgement should be exercised in the selection of settings for each individual patient. If using compressed oxygen as the gas supply, do not exceed the clinically prescribed oxygen flowrate.**5. PRECAUTIONS:** It is important to use the EzPAP® Positive Airway Pressure Therapy System as directed by your physician. The EzPAP® Positive Airway Pressure Therapy System, as well as the intended nebulizer, are under product use. The selected nebulizer should be used according to the manufacturer's instruction for use in a positive pressure system.**6. INSTRUCTIONS FOR USE:** Please see the diagram and the product to become familiar with all the product features.6.1 Set up:
Please see the diagram and the product to become familiar with all the product features.

6.1.1 Connect one end of the tubing coming from the flow meter to the device's gas inlet port (on diagram).

6.1.2 Connect the patient end of the tubing coming from the flow meter to the device's gas inlet port (on diagram).

6.2 Pressure monitoring (recommended):

Connect the device pressure port (C) to a pressure gauge (not included). This is recommended for safety reasons to ensure a pressure transients to the patient. If pressure monitoring is not being done, do not use with after-the-flow pressure monitoring.

6.3 DO NOT OVERUSE AMBIENT AIR INLET: DO NOT EXCEED 15 LPM GAS INPUT OR 60 PSI COMPRESSED GAS INPUT.**6.4 Adverse reactions may include:** Decreased airway pressure may lead to hypotension and hypercapnia; increased cardiac output; hypoxia; and/or hypotension.**6.5 DO NOT DISASSEMBLE EXCEPT FOR MOUTHPIECE, MASK OR TUBING.****6.6 INSTRUCTIONS FOR NOSE:**

Please see the diagram and the product to become familiar with all the product features.

6.6.1 Insert Flow Meter Setting:

6.6.2 Insert the patient end of the tube while performing diaphragmatic breathing.

6.6.3 Set flow rate to 5 LPM in the wall or air oxygen flowmeter.

6.6.4 While performing airway pressure, slowly adjust flow meter until desired respiratory airway pressure is reached.

6.6.5 Maintain breathing through the EzPAP® Positive Airway Pressure Therapy System to accomplish goals.

6.6.6 Several repetitions may be needed to ensure that individual patient needs are met.

6.6.7 The greater the respiratory flow through the device - the higher the airway pressure.

6.6.8 The greater the respiratory flow through the device - the lower the airway pressure.

6.7 EzPAP® Procedure:

6.7.1 Insert the EzPAP® into the mouthpiece.

6.7.2 Have patient wear it with elbow resting comfortably on table.

6.7.3 Place mouthpiece lightly into mouth.

- Be sure mouthpiece is seated correctly on the mouthpiece during therapy.
- Disconnect the mouthpiece if a nose clip is necessary.
- If using a mask, apply mask tightly, but comfortably, over nose and mouth.

6.7.4 Perform breathing through the device as recommended.

6.7.5 EzPAP® Positive Airway Pressure Therapy System is intended to be used on an intermittent basis (not for prolonged periods of time).

6.7.6 Nose clip is 3.2 x 6.3 cm is prescribed.

7. FOR USE WITH NEBULIZER:**Suggested Instructions For Use With Small Nebulizer:**

It is important to use the EzPAP® Positive Airway Pressure Therapy System as directed by your physician. The EzPAP® Positive Airway Pressure Therapy System, as well as the intended nebulizer, are under product use. The selected nebulizer should be used according to the manufacturer's instruction for use in a positive pressure system.

Set up:

Review the diagrams contained with the device.

6.1 Connect the device pressure port (C) to a pressure gauge (not included). This is recommended for safety reasons to ensure a pressure transients to the patient. If pressure monitoring is not being done, do not use with after-the-flow pressure monitoring.

6.2 While performing airway pressure, slowly adjust flow meter until desired respiratory airway pressure is reached.

6.3 Maintain breathing through the EzPAP® Positive Airway Pressure Therapy System to accomplish goals.

6.4 Several repetitions may be needed to ensure that individual patient needs are met.

6.5 The greater the respiratory flow through the device - the higher the airway pressure.

6.6 The greater the respiratory flow through the device - the lower the airway pressure.

6.7 EzPAP® Procedure:

6.7.1 Insert the EzPAP® into the mouthpiece.

6.7.2 Have patient wear it with elbow resting comfortably on table.

6.7.3 Place mouthpiece lightly into mouth.

- Be sure mouthpiece is seated correctly on the mouthpiece during therapy.
- Disconnect the mouthpiece if a nose clip is necessary.
- If using a mask, apply mask tightly, but comfortably, over nose and mouth.

6.7.4 Perform breathing through the device as recommended.

6.7.5 EzPAP® Positive Airway Pressure Therapy System is intended to be used on an intermittent basis (not for prolonged periods of time).

6.7.6 Nose clip is 3.2 x 6.3 cm is prescribed.

8. CLEANING INSTRUCTIONS:

Clean as needed. Do not immerse in any liquid. Wash the parts by swabbing with warm soapy water. Rinse and shake dry. Air dry the device by laying the unit on its side with the mouthpiece removed until all liquid droplets are evaporated before reassembling use.

9. DISPOSAL:

Dispose of EzPAP® Positive Airway Pressure Therapy System in a safe manner according to Federal/State/Local regulations and guidelines for disposal of contaminated medical waste.

10. RETURNING INFORMATION:

EzPAP system includes EzPAP, one 7 ft. oxygen tube, one pressure port cap and one of the following:

• EzPAP® System with EzPAP (10 liter/minute) ... 23-0747

• EzPAP® System with medium mask (10 liter/minute) ... 23-0747

• EzPAP® System with Large Mask (10 liter/minute) ... 23-2747

• EzPAP® Kit (1 pro Koffer) ... 23-6000

The Smiths Medical and Portex design marks and EzPAP are trademarks of the Smiths Medical family of companies. The symbol " indicates that the trademark is registered in the U.S. Patent & Trademark Office & certain other countries.

11. WARRANTY:

EzPAP system is covered by a limited warranty.

• Apnea of respiratory pressure

• Hypoventilation syndrome

• Port of admission of the gas

• Port of control of the pressure

• Infrared sensor

• Tubing oxygen of 23 m, for recording the port of admission of gas (B) at circuit of alimentation mural

12. INDICATIONS:

The EzPAP® Positive Airway Pressure Therapy System (EzPAP®) is indicated for the treatment and prevention of the apnea. The EzPAP® Positive Airway Pressure Therapy System (EzPAP®) is intended for patients possed the following conditions:

• Tracheal obstruction or respiratory insufficiency due to the face, of the bouch or du crâne

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