Inhalog 1 plus Inhalog 2 plus

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Instructions for Use
Page 2

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Stationsgeräte
Medical Ward Equipment

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Contents

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Dräger Medical AG & Co. KGaA

For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observation of these instructions. The apparatus is only to be used for purposes specified here.

Maintenance

The apparatus must be inspected and serviced regularly by trained service personnel at one year intervals.
Repair and general overhaul of the apparatus may only be carried out by trained service personnel.
We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them.
Only authentic Dräger spare parts may be used for maintenance.
Observe chapter "Maintenance".

Not for use in areas of explosion hazard

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Liability for proper function or damage

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus is serviced or repaired by personnel not employed or authorized by DrägerService or if the apparatus is used in a manner not conforming to its intended use. Dräger cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Dräger are likewise not modified by the recommendations given above.

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Intended use

Inhalog 1 plus and Inhalog 2 plus

Apparatus for assisted intermittent positive pressure ventilation (IPPB) combined with periodic inhalation of aerosol medication, exclusively to assist patients with spontaneous breathing.

Caution

O2 cylinders must not be allowed to fall over and must not be exposed to heat.

Cylinder valves and pressure reducers on O2 cylinders must not be oiled or greased; do not touch with greasy hands – fire hazard!

Cylinder valves may only be opened / closed by hand, turning slowly and uniformly. Do not use tools.

Fires burn more fiercely in the presence of O2! No smoking, no naked lights!

Do not use mobile phones within a distance of 10 metres from the machine.

Mobile phones can cause interference to electrical and electronic medical appliances, thereby putting patients at risk.*

- * Dräger medical appliances comply with the interference immunity requirements of the specific standards for the products or EN 60601-1-2 (IEC 601-1-2). However, depending on the design of the mobile phone and situation of use, field strengths may occur in the immediate environment of a mobile phone that exceed the values of the standards quoted and therefore cause interference.

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
Preparing for use

Positioning Inhalog

- Place Inhalog 1 plus/Inhalog 2 plus on a level, horizontal surface, e.g. bedside table, trolley and table. Position the surface at the height of the patient's head.
- Do not position the apparatus near a radiator and avoid exposure to direct sunlight, as it may overheat. Ensure that air can circulate freely around the apparatus.

- 1** Do not cover the ventilation slits around the handle!

Inhalog 2 plus:

- 2** Air inlet must not become clogged.

Connecting to the power supply

- Check that the operating voltage specified on the rating plate matches the mains voltage.
- 3** Connect the power cord to the apparatus and
4 to a grounded socket outlet.

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


Connecting the medical gas supply

Only for Inhalog 1 plus

- 1 Screw medical air connecting hoses into sockets on rear of apparatus.

The medicament nebulizer of the Inhalog 1 plus cannot be used without medical air.

- 2 As an option, screw additional O₂ connecting hose to the connection at the back.

The medical gases must be dust and oil-free and dry. The gas pressure should be between 3 and 6 bar.

- Connect connecting hose to a central supply or to a medical gas cylinder with a pressure reducing valve.

Connecting the hose system

- Only use parts which have been disinfected and cleaned!

Dräger hose system

- 3 Plug breathing hose into the socket on the apparatus and into the »Inspiration« socket of the medicament nebulizer.
- 4 Plug nebulizer hose into middle socket on apparatus and into socket on medicament nebulizer.


In the case of recumbent patients or patients with a tracheal cannula, use short corrugated hose:

- 5 Fit corrugated hose between patient connection and cannula connector. If necessary, use adapter, cover and catheter connector, see »Accessories« Page 24.

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
Use mouthpiece on patients who are not subject to a tracheostomy:

- 1 Plug mouthpiece into patient connection

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
Disposable hose system

- 2 Plug breathing hose into socket on apparatus and into socket on nebulizer.
- 3 Plug nebulizer hose into middle socket on apparatus and into socket on medicament nebulizer.
- 4 Plug control hose for exhalation resistance into outer socket on apparatus and into socket of exhalation valve.

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
In the case of recumbent patients or patients with a tracheal cannula, use short corrugated hose:

- 5 Fit corrugated hose between patient connection and cannula connector

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Use mouthpiece on patients who are not subject to a tracheostomy:


- 6 Plug mouthpiece into patient connection.

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

Preparing the medicament nebulizer

Dräger medicament nebulizer:

When using the stenosis ring:

- 1 Open the medicament nebulizer by turning it anti-clockwise.
- 2 Fit the stenosis ring.
- 1 Close the medicament nebulizer clockwise.

- Unscrew nebulizer bottle and place on a horizontal surface.
- Fill with the prescribed amount of medication.
The filling levels are indicated on the nebulizer bottle:
3 mL – 6 mL – 10 mL – 20 mL
- Screw nebulizer bottle back onto medicament nebulizer.
- Hold medicament nebulizer vertically without tilting.

- Hook it into mount on apparatus.

Medicament nebulizer for disposable hose system:


- Unscrew nebulizer bottle and place on a horizontal surface.
- Fill with the prescribed volume of medication, noting the levels indicated on the bottle.
- Screw nebulizer bottle back on again.
- Hook it into mount on apparatus.

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
Settings

Only as directed by a doctor, otherwise therapy may prove unsuccessful

- Press power switch on rear of apparatus until it engages = ON.

Readiness for operation is indicated by the flashing symbol "*" on the display.


- Set »Trigger« knob for switch-on pressure.
Range of settings: -1 to -8 mbar
(-1 mbar = low application of pressure on inhalation starts inspirator, knob turned clockwise as far as possible.)
- The switch-on pressure appears as a mark on the left of the bar graph on the apparatus.
- Set »Pinsp.« knob for max. inspiratory pressure.
Range of settings: 5 to 40 mbar.
Inhalog switches over to expiration if the max. inspiratory pressure is reached during inspiration.
The set max. inspiratory pressure Pinsp. appears as a mark on the right of the bar graph.
- Set »Flow« knob (inspiratory flow).
Range of settings: 20 to 60 L/min.

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
When using the Dräger hose system with stenosis ring:

- 3 Set ring for expiration resistance:
Setting 0 = minimum resistance
Setting 3 = maximum resistance

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
Function check

- 1 Start inspiration = briefly pull »Trig-
ger« knob:

Inhalation commences, aerosol medica-
tion emerges from the mouthpiece.

- Seal mouthpiece with your hand and observe bar graph:

The indicated airway pressure (P_{aw}) increases and the apparatus switches to expiration when the set max. inspiratory pressure Pinsp. has been reached. The airway pressure drops to 0.

When using a disposable hose system:

- 2 The green LED »Rexp.« lights up continuously following the first inspiration.

Only for Inhalog 1 plus:

- 3 The green LED »O2« lights up continuously after the first inspiration when delivering O2.

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

Operation

- Ensure that mouthpiece is enclosed by the lips. Do not breathe through the nose; use a nose clamp if necessary.
- Breathe normally.

The patient's attempts to inhale generate a negative pressure which starts inspiration and medicament nebulization. Inhalog applies the aerosol medication to the lungs with the inspiratory flow set. When the max. inspiratory pressure P_{insp} is reached, the apparatus switches over to expiration, medicament nebulization is interrupted and the patient exhales – possibly against a preset resistance to prolong exhalation.

Inhalation can also be started manually:

- 1 Briefly pull »Trigger« knob.
- 2 The airway pressure (P_{aw}) is continuously displayed on the bar graph.

- 3 The patient's respiratory frequency f spont. bpm (breaths per minute) is measured by the apparatus and appears as a numerical value on the right of the bar graph.
- 4 The treatment time (symbol: ⊖, minutes) appears as a numerical value on the left of the bar graph.

Patients with tracheal cannulas and short corrugated hose

- The medicament nebulizer must be at the lowest point in the hose system to prevent condensation from reaching the patient.
- The nebulizer power is reduced by 1/3.
- Administer the treatment only in the presence and constant supervision of a nurse.

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Shutdown

- At the end of the prescribed treatment time
- Fully depress switch on rear of apparatus = OFF.
- Unscrew nebulizer housing, drain off any residual medication, fill with approx. 20 mL distilled water and screw nebulizer housing back into place.
- Press switch on rear until it engages = ON.
- Briefly pull »Trigger« knob = start inspiration. Medicament nebulizer sprays water into the air or onto a disposable tissue for approx. 10 seconds.
- Switch off apparatus and drain off remaining distilled water.
- The apparatus is now ready for the next therapy session with the same patient.

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
After use:

Dismantle parts

- Unplug power cable in order to avoid risk of electric shocks while stripping down the apparatus!

Disposable hose system:

- Remove and dispose of as directed by manufacturer.

Dräger hose system:

- Disconnect breathing hose and nebulizer hose from apparatus and breathing valve.
- Dismantle breathing valve and nebuli-zer:
 - 1 Pull off mouthpiece.
 - 2 Unscrew nebulizer bottle.
 - 3 Pull atomizer off nebulizer housing.
 - 4 Turn cap a quarter-turn anticlockwise and remove from nebulizer housing.
 - 5 Remove diaphragm, but do not dis-assemble further.
 - 6 Pull ring for expiration resistance off nebulizer housing (if used).

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Disinfection/cleaning

Only surface disinfectants based on

- aldehydes or
 - quaternary ammonium compounds
- should be used to disinfect the apparatus.

Disinfectants based on

- compounds containing alkyl amines,
 - compounds containing phenol,
 - compounds giving off halogens,
 - strong organic acids,
 - compounds giving off oxygen and
 - compounds containing alcohol
- should not be used.

Users in Germany are advised to use the disinfectants specified in the current DGHM* list.

This list also specifies the active agents on which the respective disinfectants are based. Users in those countries for which the DGHM list does not apply should use disinfectants based on the substances recommended above.

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Disinfection/cleaning in hospitals

Apparatus:

- Unplug the power plug.
- **Do not immerse the apparatus in liquids. Liquid must not be allowed to seep inside the apparatus, otherwise short-circuits may result when it is switched on!**
- Dirt can be wiped off the apparatus with a damp disposable cloth.

Disinfect by wiping

- e.g. with Buraton 10F (made by Schülke & Mayr, Norderstedt). Note the manufacturer's directions!

Parts of the dismantled Dräger hose system:

- Can be disinfected in the washing machine (93 °C / 10 minutes) **using detergent only!** or
- Disinfect in disinfectant solution, e.g. Buraton 10F (made by Schülke & Mayr, Norderstedt). Note the manufacturer's directions!

The various parts of the dismantled Dräger hose system can also be sterilized in hot steam at 134 °C.

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Cleaning/disinfection at home

Apparatus:

- **Do not immerse the apparatus in liquids. Liquid must not be allowed to seep inside the apparatus, otherwise short-circuits may result when it is switched on!**
- Dirt can be wiped off the apparatus with a damp disposable cloth.

Parts of the dismantled Dräger hose system:

- Clean in warm water, rinse thoroughly under running water, then:
- Immerse in a solution containing one part wine vinegar and two parts water. Leave for approx. 30 minutes.
- Rinse thoroughly under running water. Drain water out of hoses and allow all parts to dry completely.

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
Fault – Cause – Remedy

Fault	Cause	Remedy
Apparatus does not commence inspiration when triggered.	Leak or defect in hose connection or nebulizer. Medical air supply not connected to Inhalog 1 plus or supply pressure too low. Trigger faulty.	Plug hoses in securely, see page 6 et seq., replace defective hoses; check nebulizer, assemble correctly, replace defective parts. If used: check disposable patient system and replace if necessary. Connect medical gas hose, check supply pressure. Supply pressure must be 3.0 to 6.0 bar. Call DrägerService.
No aerosol from nebulizer.	Filling level in medicament nebulizer too low. Nebulizer soiled. Nebulizer hose not connected to apparatus or medicament nebulizer. Apparatus faulty.	Refill medication. Thoroughly clean nebulizer. Plug in nebulizer hose. Call DrägerService.
Inhalog 2 plus will not start.	No electric power. Compressor in Inhalog 2 plus faulty.	Check that mains connector is plugged in, check power supply. Check fuses and replace if necessary, see page 19. Call DrägerService.
Inhalog 1 plus Green LED »O ₂ « does not light up.	O ₂ supply not connected or supply pressure too low.	Connect O ₂ hose, check supply pressure.

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
Maintenance

Exchanging the fuse

- Unplug mains connector!
- 1 Press the two tabs on the fuse insert together simultaneously and
- 2 Pull out the insert at the same time.
- 3 Pull out the blown fuse and insert a new fuse.
Inhalog 1 plus: T 0,5 AL 250 V
IEC 127-2/III
Inhalog 2 plus: T 1 AL 250 V
IEC 127-2/III
- Press insert in until it engages.

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Exchanging the suction filter

- Only for Inhalog 2 plus
- once a year
 - or
 - in the event of high levels of contamination.
 - Unscrew the locking screws.
 - Remove protective sieve.
 - Exchange the suction filter behind it.
 - Fit protective sieve and secure with screws.

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
Disposing of the apparatus

- at the end of its service life.
- Inhalog 1 plus / Inhalog 2 plus can be disposed of in accordance with the applicable statutory regulations and in consultation with the relevant waste disposal companies.

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
What is what

Front view Inhalog 1 plus

- 1 »Trigger« knob for starting inspiration with negative pressure.
Can be started manually by briefly pulling the knob.
- 2 »Rexp.« knob for setting the expiration resistance.
- 3 »O₂« knob for setting a higher O₂ concentration.
- 4 »Flow« knob for setting the inspiratory flow.
- 5 »Pinsp.« knob for setting the maximum inspiratory pressure at which the apparatus switches over to expiration.
- 6 Numerical indication of the respiration rate in bpm = breaths per minute.
- 7 Bar graph for analog indication of the airway pressure (P_{aw}) with marks for trigger pressure and Pinsp.
- 8 Numerical indication of the desired therapy time.

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

**Front view Inhalog 2 plus**

- 1 »Trigger« knob for starting inspiration with negative pressure.
Can be started manually by briefly pulling the knob.
- 2 »Rexp.« knob for setting the expiration resistance.
- 3 »Flow« knob for setting the inspiratory flow.
- 4 »Pinsp.« knob for setting the maximum inspiratory pressure at which the apparatus switches over to expiration.
- 5 Numerical indication of the respiration rate in bpm = breaths per minute.
- 6 Bar graph for analog indication of the airway pressure (P_{aw}) with marks for trigger pressure and Pinsp.
- 7 Numerical indication of the desired therapy time.

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
Ambient conditions

During operation

Temperature	Inhalog 1 plus: 5 to 40 °C Inhalog 2 plus: 5 to 35 °C
Atmospheric pressure	600 to 1100 hPa
rel. humidity	10 to 95 %
During storage	
Temperature	-20 to 45 °C
Atmospheric pressure	600 to 1200 hPa
rel. humidity	0 to 95 %, without condensation

Performance characteristics

Control principle	pressure-controlled, assisted
Switch-on pressure	-1 mbar to -8 mbar*
Trigger volume	approx. 2 mL
Inspiratory Flow	20 to 60 L/min
Type of flow	decelerating
Nebulizer pressure (with disposable nebulizer type Lifeline 1616)	1 ± 0,5 bar
Nebulizer capacity	
$T_{\text{insp}} = 1.7 \text{ s}$, $f = 12 \text{ per min}$, Insp. flow 30 L/min and pressure 1,3 bar	2,5 mL/10 min

Spectrum of drops

Measured by the Fraunhofer Institute for Toxicology and Aerosol Research in Hannover.

The distribution of droplet diameter in an aerosol is described by a normal distribution.

The mean value equals 3.7 µm and the standard deviation equals 2.5 µm.

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Operating characteristics

Operating voltage	230 V ±10 %, 50 Hz
Protection class	II  ; Typ BF 
Electromagnetic compatibility	Tested to EN 60601-1-2
UMDNS-Code	11-634
Universal Medical Device Nomenclature System – Nomenklatur for medical products	
Classification as per EC Directive 93/42/EEC Annex IX	Class II a

Inhalog 1 plus

with O ₂ -and medical air from the central medical gas supply	
Power consumption	30 W
Dimensions (B x H x T)	260 x 210 x 260 mm
Weight	approx. 5 kg
Sound level	<60 dB (A)
Oxygen concentration	21, 30, 40, 50 Vol.% O ₂ ±10 Vol.-%
Air and O ₂ supply pressure from the central medical gas supply	3.0 to 6.0 bar*
O ₂ connection	NIST
Air connection	NIST

Inhalog 2 plus

with internal compressor	
Power consumption	160 W
Dimensions (B x H x T)	260 x 210 x 260 mm
Weight	approx. 7 kg
Sound level	<60 dB (A)

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Order List

Designation/description	Part No.
Inhalog 1 plus Scope of supply: basic unit, hose set plus, autoclavable patient set, carrying bag, power cable, Instructions for Use	2M 86 120
Inhalog 2 plus Scope of supply: basic unit, disposable patient set, carrying bag, power cable, Instructions for Use	2M 86 125
Accessories	
Patient set (nebulizer with mouthpiece, autoclavable)	2M 86 023
Hose set plus (corrugated hose, nebulizer hose, 3 hose clips, autoclavable)	2M 86 038
Disposable patient set (50 off)	MX 02 700
Trolley (feet and column)	2M 86 029
Table Inhalog/RespiCare	2M 86 030
Carrying bag	84 13 396
Accessories for tracheostomy application	
Adapter E-bracket	84 03 076
Cover	84 01 644
Corrugated hose	84 02 041
Catheter connector	M 23 840
Only for Inhalog 1 plus	
– required for connection to the central medical gas supply	
O2 connecting hose 1.5 m NIST (black)	M 34 401
O2 connecting hose 3.0 m NIST (black)	M 34 402
O2 connecting hose 5.0 m NIST (black)	M 34 403
Compressed air connecting hose 1.5 NIST (black)	M 34 407
Compressed air connecting hose 3.0 NIST (black)	M 34 408
Compressed air connecting hose 5.0 NIST (black)	M 34 409
Adapter NIST/DIN O ₂	
Adapter NIST/DIN Air	
Only for Inhalog 2 plus	
Suction filter	

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Index

Ambient	22
Disinfection/cleaning	15
Disposal	19
Expiration resistance	10
Function check	11
Fuse	19
Hose system	6
Intended	4
Maintenance	3
medical gas supply	6
Medicament nebulizer	8
Operating voltage	5
Operation	12
Order List	24
Positioning	5
Power supply	5
Preparing for use	5
Settings	9
Shutdown	13
Sterilisation	15
Stripping	14
Stripping down	14
Technical data	22
What is what	20

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© Dräger Medical AG & Co. KGaA