
# ISO-NEB®

## Filtered Nebulizer System

REF 41755

Rx ONLY

### ISO-NEB® Filtered Nebulizer System

Not made with natural rubber latex.

Contains DEHP

## Instructions for Use:

NOTE: These directions are general guidelines intended for use by qualified medical personnel.

- Remove the ISO-NEB from its package and be sure all components are intact.
- Rotate the nebulizer into upright position for right- or left-handed use.
- Verify that the arrows on one-way valves are oriented correctly: two on the inspiratory limb toward the mouthpiece and one on the expiratory limb toward the particle filter.
- Remove the nebulizer from the tee.
- Add prescribed medication to the nebulizer jar through the opening at the top of the cap (8 cc maximum).
- Reattach the nebulizer to the tee.
- Connect one end of the supply tubing to the fitting under nebulizer and the other end to the source gas outlet.
- Adjust the gas source to 6 to 8 LPM. When using an air compressor, adjust pressure for 23 to 25 psi.
- Proceed with treatment as prescribed. Exhale through the mouthpiece.
- Cover mouthpiece to divert flow toward filter when taking breaks during treatment.

## Warning Statements:

- Do not reuse, wash, soak or sterilize. Residuals from medication or cleaning agents retained within the device may result in contamination, dangerously increased breathing resistance and impaired valve function.
- Do not connect directly to an endotracheal or tracheostomy tube.
- Discard after single application.
- Re-processing of products intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only products may result in exposure to viral, bacterial, fungal, or prion pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these products. This product is not designed to be cleaned, disinfected, or re-sterilized.
- The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

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© 2006, 2009, 2011, 2013, 2016 Teleflex Incorporated