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1 Introduction

1.1 Intended use

The JOYCEone nasal masks are used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device.

1.2 Contraindications

In the case that any of the following symptoms are present the mask may not be used, or only with particular caution:

erosions and ulcerations, skin allergies, flush of the facial skin, pressure marks on the face, claustrophobia, anxiety, deformities of the face or nasopharynx, intake of medication which may induce vomiting, the necessity for immediate intubation. Please also observe the contraindications in the instructions for use of your therapy device.

1.3 Side effects

Nasal congestion, dry nose, morning xerostomia (dry mouth), feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks in the face, irritating noises when breathing

2 Safety

Read these instructions for use carefully. They are a constituent part of the mask and must be available at all times. In accordance with Directive 93/42/EEC, observe the points below.

Risk of injury from the supply of oxygen!

Oxygen can deposit in clothing, bed linen, and hair. In conjunction with smoking, naked flames, or electrical equipment, it can cause fires and explosions.

⇒ Do not smoke.

⇒ Avoid naked flames.

⇒ Do not use electrically conductive tubes (exception: breathing tubes which can be heated electrically).

Risk of injury from re-ingesting CO₂!

If the mask is used incorrectly, CO₂ may be re-ingested. ⇒ Do not close the openings of the mask.

⇒ Only apply the mask for a prolonged period if the therapy device is running.

⇒ Only use the mask within the quoted therapy pressure range.

⇒ Do not use the mask on patients who cannot take it off themselves.

⇒ Monitor patients with obstructive and restrictive lung diseases individually when using JOYCEone masks.

Risk of injury from slipping of the mask!

If the mask slips or falls off, the therapy is ineffective.

⇒ Monitor patients with restricted spontaneous respiration.

⇒ Activate low pressure/leakage alarms on the therapy device.

⇒ Monitor patients with obstructive and restrictive lung diseases individually.

Patient endangerment through application of mask during anesthesia!

If the nasal mask is used during anesthesia, the anesthetic gas supplied can escape through the mouth and endanger the patient.

⇒ Do not use nasal mask during anesthesia.

2.1 General information

- Use of third-party products may lead to functional failures and restricted fitness of purpose. Biocompatibility may also be compromised. Please note that in these cases, any claim under warranty and liability will be void if neither the accessories nor original replacement parts recommended in the instructions are used.
- To prevent infection, bacterial contamination, or functional impairments, follow the instructions in the section "Hygiene preparation".

3 Product description

3.1 Overview

Illustrations of the individual components can be found on the front cover.

Key

1. Headgear
2. Forehead cushion
3. Forehead support
4. Retaining ring
5. Elbow
6. Plug for pressure measurement port
7. Pressure measurement port
8. Rotating sleeve
9. Mask body
10. Headgear clip
11. Mask cushion

3.2 Compatible devices

The mask can be used with all therapy devices which are not for life-sustaining treatment of patients with sleep apnea or ventilatory insufficiency.

With some combinations of devices, the actual pressure in the mask will not correspond to the prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that the actual pressure in the mask corresponds to therapy pressure.

3.3 Exhalation system

The mask has an integrated exhalation system. The retaining ring and mask body are shaped in such a way that a gap forms between these parts. The exhaled air can escape through this gap.

3.4 Pressure measurement port

The mask has a pressure measurement port which can be used for measuring pressure or supplying oxygen. If you are not using this, close it with the plug for pressure measurement port to reach therapy pressure.

4 Operation

The following figures show you how to assemble, apply, remove, adjust, and dismantle the mask:

Figure 1: Assembling the mask

Figure 2: Applying the mask

Figure 3: Adjusting the mask

Figure 4: Removing the mask

Figure 5: Dismantling the mask

5 Hygienic preparation

⚠ WARNING

Risk of injury from inadequate cleaning!

Residues may congest the mask, impair the integrated exhalation system, and jeopardize therapy success.

- ⇒ For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.
- ⇒ Thoroughly clean the contact surfaces between the elbow and the rotating sleeve.

5.1 Cleaning the mask

1. Dismantle the mask (see figure 5).
2. Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask components with warm water and a mild detergent.	X	
When washing mask parts (exception: mask cushion and forehead cushion), clean thoroughly with a cloth or a soft brush.		X
Wash headgear by hand.		X

3. Rinse all parts with clear water.

NOTICE

Material damage due to abrasion!

Rubbing the mask cushion and forehead cushion can damage the coating.

⇒ Do not rub the mask cushion and forehead cushion dry.

4. Allow all parts to air-dry.
5. Perform a visual inspection.
6. Replace damaged parts if necessary.
7. Assemble the mask (see figure 1).

Discolorations of mask components do not impair the functionality of the mask.

5.2 Change of patient

You must subject the mask to a hygienic preparation in the event of a change of patient. For information on hygienic preparation in the event of a change of patient, please refer to our Internet brochure at www.weinmann.de. If you do not have Internet access, please contact Weinmann.

6 Disposal

You can dispose of all parts in domestic waste.

7 Faults

Fault	Cause	Remedy
Facial pressure pain.	Mask fits too tightly.	Loosen headgear.
Draft in the eyes.	Mask fits too loosely.	Tighten headgear.
	Mask does not fit.	Contact your specialist retailer.
Therapy pressure is not reached.	Mask is not adjusted correctly.	Re-adjust mask (see figure 3).
	Mask cushion is damaged.	Replace mask cushion.
	Tube system is leaking.	Check connection and fit of tubes.
	Air is escaping from pressure measurement port.	Close pressure measurement port with plug (included in scope of supply).

8 Technical data

	Nasal mask
Product class according to 93/42/EEC	IIa
Dimensions (W x H x D)	96 mm x 115 mm x 83 mm
Weight	87 g
Therapy pressure range	4 hPa - 30 hPa
Hose connection: Cone according to EN ISO 5356-1	Ø 22 mm (male)
Pressure measurement port	Ø 4 mm
Width of headgear clips	24 mm
Temperature range: Operation Storage	+5°C to + 40°C -20°C to +70°C
Flow resistance at 50 l/min at 100 l/min	0.05 hPa 0.28 hPa
Quoted two-figure noise emission value according to ISO 4871: - Sound pressure level - Sound power level uncertainty factor	15 dB(A) 23 dB(A) 3 dB(A)
Service life	Up to 12 months ¹
Standards applied	EN ISO 16510-2:2009

¹ The materials of the mask deteriorate if exposed to e.g. aggressive detergents. In individual cases, therefore, it may be necessary to replace mask parts sooner.

CE 0197

The right to make design modifications is reserved.

9 Warranty

Weinmann gives the customer a limited manufacturer warranty on new genuine Weinmann products and any replacement part fitted by Weinmann in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from www.weinmann.de on the Internet. We can also send you the warranty conditions on request.

In the event of a claim under warranty, please contact your specialist dealer.

Product	Period of guarantee
Masks including accessories	6 months

10 Declaration of conformity

Weinmann Geräte für Medizin GmbH + Co. KG hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The complete text of the declaration of conformity can be found at: www.weinmann.de

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