JOYCE vented

JOYCE vented 40 hPa

JOYCE non vented

JOYCE non vented 40 hPa

JOYCE GEL vented

JOYCE GEL non vented

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1 Device description

Legend

(Description of individual parts: see enclosed leaflet)

- 1** Headgear
- 2** Forehead cushion
- 3** Fine adjustment component of forehead support
- 4** Coarse adjustment component of forehead support
- 5** Forehead support
- 6** Headgear clips
- 7** Pressure measurement ports
- 8** Port cap
- 9** Rotating sleeve (non vented: blue coloured)
- 10** Elbow (non vented: blue coloured)
- 11** Retaining ring (non vented: blue coloured)
- 12** Mask frame
- 13** Mask cushion

Intended use

The nasal mask is an accessory and acts as a seal between the patient and the therapy device. It is used for the treatment of sleep apnoea and for the non-invasive ventilation of patients with ventilatory insufficiency.

Before use by the patient

Patients must be instructed in the use of the mask by qualified personnel before therapy begins.

2 Safety instructions

Safety regulations

Following the requirements of Richtlinie 93/42 EWG please observe the following points:

- Please read these instructions for use very carefully. They are a part of the product and should be available at all times.

- Use the mask only in the therapy pressure ranges required for the individual variants:

- 4-30 hPa: JOYCE vented
JOYCE non vented
- 4-40 hPa: JOYCE vented 40 hPa
JOYCE non vented 40 hPa
JOYCE GEL vented
JOYCE GEL non vented

- If third-party items are used, functional failures may occur and fitness for use may be restricted. Biocompatibility requirements may also not be met. Please note that in such cases, any claim under warranty and liability will be voided if neither the accessories nor genuine replacement parts recommended in the instructions for use are used.
- Patients with limited spontaneous respiration must be monitored at all times. The appropriate negative pressure/leak alarm on the ventilation device must be activated.
- Patients with obstructive and restrictive lung diseases should be examined on an individual basis during use of the JOYCE masks (for example with regard to re-inhalation of CO₂).
- Do not close the openings in this mask except in situations described in these instructions for use. Otherwise there is a risk of re-inhaling CO₂ and of suffocation.
- Only put the nasal mask over your face for a prolonged period when the therapy device is switched on. Otherwise there is a risk of re-inhaling CO₂.
- Do not use a breathing hose that is electrically conductive.
- Follow the instructions in section 4 "Hygiene treatment" to avoid infection, bacterial contamination and poor functionality of the mask.

Side-effects

The following side-effects may occur:

nasal congestion, dry nose, dry mouth in the morning, sinus pressure, irritation of the conjunctiva, irritating noises during exhalation, pressure sores on the face, reddening of the skin.

3 Operation

Fig. A-E: see enclosed leaflet.

3.1 Adjust forehead support (Fig. A)

Fine adjustment

Squeeze the upper side tabs of the coarse adjustment component inwards. Fit the fine adjustment component into one of the three notches.

To change the height of the forehead cushion, i.e. the fine adjustment (by approx. 1 cm), turn the fine adjustment component through 180°. Then position it as described above.

Coarse adjustment

Squeeze the lower side tabs of the coarse adjustment component inwards. Place it in one of the slots in the mask frame.

3.2 Putting on the mask

1. Unhook one of the lower headgear clips from the mask.
2. Guide the headgear over your head and hook the headgear clips back onto the strap.
3. Adjust the upper and lower headgear straps so that the mask sits firmly, but not too tightly, on your face.

HEADstrap only:

The HEADstrap headgear provides the following options for you to do this (see also "Overview" on the attached leaflet):

- coarse adjustments: slot "a" or "b" of the headgear.
 - fine adjustments: four hook-and-loop straps."
4. If necessary, re-adjust the forehead support.

3.3 Exhaled air system (Fig. E)

vented variants

vented variants (all mask components are transparent) have a so-called **integrated exhaled air system**: the retaining ring and

mask frame are shaped in such a way that a gap is formed between these parts. Through this gap the exhaled air can escape. The special geometry of the parts guarantees that when putting on the retaining ring the gap cannot be closed. **Do not** connect an external exhaled air system, otherwise the therapy pressure will not be reached.

non vented variants

non vented variants (retaining ring, elbow and rotating sleeve in blue) have no exhaled air system. It is essential to connect an external exhaled air system. Otherwise there is a risk of reinhaling CO₂ and of asphyxia.

Follow the instructions for use of your exhaled air system.

3.4 Connecting the mask to the therapy device (Fig. E)

Push the hose of your therapy device onto the rotating sleeve (vented variants) or onto your external exhaled air system (non vented variants). Ensure that the connection to the hose is secure and airtight. The connection must not be loosened by the therapy pressure.

If you want to detach the mask temporarily from the hose, pull the rotating sleeve off the elbow.

3.5 Gel-filled mask components

On JOYCE GEL vented and JOYCE GEL non vented, the forehead cushion and the mask cushion are gel-filled as standard. Gel forehead cushions and gel mask cushions are also available as replacement parts for other JOYCE variants.

Caution!

Gel-filled mask components are highly sensitive components. Always handle them with care to prevent deformation and damage.

3.6 Therapy pressure ranges

The individual JOYCE mask variants are each available for a particular therapy pressure range (see "7. Scope of supply" on page 17);

- 4-30 hPa: JOYCE vented
JOYCE non vented
- 4-40 hPa: JOYCE vented 40 hPa
JOYCE non vented 40 hPa
JOYCE GEL vented
JOYCE GEL non vented

The JOYCE masks for the therapy pressure range 4 - 30 hPa have a soft mask cushion, the JOYCE masks for the therapy range 4 - 40 hPa have a harder mask cushion or a gel mask cushion (gel-filled variants).

The harder mask cushion for the therapy pressure range 4 - 40 hPa is marked with the number "40".

The soft (up to 30 hPa) and harder (up to 40 hPa) mask cushion and the gel-filled mask cushion (gel variants) are also available as individual replacement parts to enable you to adapt your mask to changed therapy requirements.

3.7 Ports for pressure measurement (Fig. C)

JOYCE masks are fitted with two ports that can be used to measure pressure or to feed in oxygen.

Caution!

- Oxygen therapy can have side-effects if incorrectly administered or used in the wrong dosage. Only carry out this sort of treatment if prescribed by a physician.
- When carrying out oxygen therapy, strictly follow the instructions for use of your O₂ system and your therapy device.

When the pressure measurement ports are not in use, always close them with the cap provided. Otherwise air will escape from the ports and the therapy pressure will not be reached.

3.8 Compatible devices

Therapy devices

JOYCE masks can be combined with all devices designed for carrying out the therapies described under "Intended Use".

Note that if unsuitable device combinations are used, the actual pressure in the mask may not correspond to the therapy pressure you have been prescribed.

This may also happen even if the correct pressure is indicated on the therapy device. Have your device combination adjusted by a physician or distributor so that the pressure in the mask corresponds to your therapy pressure.

Air humidifiers

JOYCE masks can be used with either cold-air or warm-air humidifiers.

Never set the humidification level so high that condensation forms in the air hose. Follow the instructions for use of your humidifier and therapy device.

3.9 After use

1. Unhook one of the lower headgear clips from the mask and remove the mask.
2. Loosen the headgear clips and remove the headgear from the mask.
3. Proceed as described in 4 "Hygiene treatment".

3.10 Disassembling/reassembling the mask

Disassembling the mask (Fig. D - A)

1. Detach the mask cushion from the mask frame.
2. Remove the cap from the pressure measurement ports.
3. Rotate the retaining ring counter-clockwise. Now pull the retaining ring and the elbow away from the mask frame.
4. Press the elbow through and out of the retaining ring and pull the rotating sleeve off the elbow.

5. Detach the coarse adjustment component from the mask frame.
6. Detach the fine adjustment component from the coarse adjustment component. Then detach the forehead cushion from the fine adjustment component.

disinfecting the contact surfaces between the mask frame and the retaining ring. Residues on these surfaces can reduce the success of the therapy and the functioning of the integrated exhaled air system (vented variants).

Reassembling the mask (Fig. A - D)

1. Attach the forehead cushion to the fine adjustment component. Then attach the fine adjustment component to the coarse adjustment component.
2. Attach the coarse adjustment component to the mask frame.
3. Push the rotating sleeve onto the elbow until you hear it click into place.
4. Press the elbow into the retaining ring until you hear it click into place.
5. Attach the retaining ring and elbow onto the mask frame. Turn the retaining ring clockwise until you hear it click into place.
6. Push the cap onto the pressure measurement ports.
7. Attach the mask cushion to the mask frame.
8. Fasten the headgear to the mask with the clips.

Weekly

Disassemble the mask and subject the parts of the mask to a hygiene treatment (see "4.2. Cleaning, disinfection, sterilisation" on page 15).

4.2 Cleaning, disinfection, sterilisation

Permissible procedures

Procedure	Dishwasher	Hand wash	Disinfection*	Sterilisation (alternative to disinfection)
Components				
Mask cushion	•	•	•	•
Mask cushion GEL	•	•	•	
Forehead cushion	•	•	•	
Forehead cushion GEL		•	•	
Port cap	•	•	•	•
Mask frame	•	•	•	
Coarse adjustment component	•	•	•	
Fine adjustment component	•	•	•	
Elbow	•	•	•	
Rotating sleeve	•	•	•	
Retaining ring	•	•	•	
Headgear clips	•	•	•	
Headgear HEADstrap		•		
Headgear JOYCEstrap		•		

Caution!

Always clean the components of the mask thoroughly before disinfecting or sterilising them as indicated in the above Table.

4 Hygiene treatment

4.1 Schedule

Daily

Clean the mask in warm water to remove all residues. You can disassemble the mask for cleaning. Use a mild cleaning agent (e.g. dishwasher detergent). After cleaning, rinse all parts carefully with clean water.

Caution!

- For patients with a weakened immune system or a special medical history, **daily** disinfecting or sterilization of the mask components may be necessary. In this case, please consult your doctor.
- Please take great care when cleaning and

Gel-filled mask components are highly sensitive components. Always handle them with care to prevent deformation and damage.

Description of procedures

Dishwasher	Rinse cycle at 65 °C. Then rinse all parts carefully with clean water.
Hand wash	Hand wash in warm water. Use a mild detergent. Then rinse all parts in clean water and allow to air-dry.
Disinfection*	Clean mechanically, e.g. with a brush, in dilute solution. After disinfection rinse all parts thoroughly with distilled water and dry thoroughly.
Sterilisation (alternative to disinfection)	Steam sterilisation in equipment rated according to EN 285. Temperature 134 °C, minimum sterilisation time 3 minutes.

*We recommend the disinfectants GIGASEPT FF, CIDEX OPA and ANIOSYME DD1. Evidence of tolerance to these disinfectants has been provided. It is essential to follow the manufacturer's instructions.

Notes on the hygiene treatment and on the hygiene treatment for a change of patient can be found in a brochure from our information portal at www.weinmann.de. If you do not have Internet access, contact Weinmann.

Note

Wear suitable gloves (e.g. disposable gloves) during the hygiene treatment.

Wash the HEADstrap headgear before using for the first time, as the colour might run.

The headgear may only be washed by hand!

Do not iron the headgear, otherwise the hook-and-loop closures will no longer work!

4.3 Following treatment

Thoroughly dry all parts following the hygiene treatment. Check the components for residual contamination or damage. Repeat the hygiene treatment if required. Replace damaged parts. Reassemble the mask.

4.4 Multiple patient use

If the mask is used by a different patient, it must be subjected to a hygiene treatment first.

Disassemble the mask. Disinfect or sterilize all parts as described in section 4.2 "Cleaning, disinfection, sterilisation". Use a new headgear.

5 Functional checks

Carry out a visual check after every hygiene treatment.

If the mask cushion is damaged and is no longer airtight, replace it with a new one.

Check the functionality and air-tightness of the ball-and-socket joint on the elbow. If it no longer moves easily or is no longer airtight, replace the elbow and retaining ring.

Parts that are no longer usable can be disposed of as normal household waste.

The forehead cushion and the mask cushion may become discolored, but this does not impair function.

6 Troubleshooting

Problem	Cause	Solution
Painful pressure on the face.	Mask is too tight.	Loosen headgear slightly. Adjust position of forehead support.
Draught in the eyes.	Mask is too loose.	Tighten headgear slightly. Adjust position of forehead support.
	Mask does not fit.	Use a different size of mask; contact your distributor.
	High therapy pressure.	Use firmer mask cushion (for pressures up to 40 hPa). Available from your distributor.
	Mask not adjusted correctly.	Adjust headgear. Adjust position of forehead support.
Therapy pressure is not reached.	Mask cushion is damaged.	Replace mask cushion.
	Hose system is leaking.	Check push-fit connector and fit of hose sleeves.
	Air is escaping from pressure measurement ports.	Cover pressure measurement ports with cap (included in package).

7 Scope of supply

Note

All masks are supplied fully assembled and with headgear.

vented variants

Article	Article no. with HEADstrap	Article no. with JOYCEstrap
JOYCE vented, Size S	WM 26110	WM 26110/1
JOYCE vented, Size M	WM 26120	WM 26120/1
JOYCE vented, Size L	WM 26130	WM 26130/1
JOYCE vented, Size XL	WM 26140	WM 26140/1
JOYCE vented 40 hPa, Size S	WM 26111	WM 26111/1
JOYCE vented 40 hPa, Size M	WM 26121	WM 26121/1
JOYCE vented 40 hPa, Size L	WM 26131	WM 26131/1
JOYCE vented 40 hPa, Size XL	WM 26141	WM 26141/1
JOYCE GEL vented, Size S	WM 26112	WM 26112/1
JOYCE GEL vented, Size M	WM 26122	WM 26122/1
JOYCE GEL vented, Size L	WM 26132	WM 26132/1

non vented variants

Article	Article no. with HEADstrap	Article no. with JOYCEstrap
JOYCE non vented, Size S	WM 26160	WM 26160/1
JOYCE non vented, Size M	WM 26170	WM 26170/1

Article	Article no. with HEADstrap	Article no. with JOYCEstrap
JOYCE non vented, Size L	WM 26180	WM 26180/1
JOYCE non vented, Size XL	WM 26190	WM 26190/1
JOYCE non vented 40 hPa, Size S	WM 26161	WM 26161/1
JOYCE non vented 40 hPa, Size M	WM 26171	WM 26171/1
JOYCE non vented 40 hPa, Size L	WM 26181	WM 26181/1
JOYCE non vented 40 hPa, Size XL	WM 26191	WM 26191/1
JOYCE GEL non vented, Size S	WM 26162	WM 26162/1
JOYCE GEL non vented, Size M	WM 26172	WM 26172/1
JOYCE GEL non vented, Size L	WM 26182	WM 26182/1

Replacement parts

Article	Article no.
Ball-and-socket joint set: retaining ring, elbow, rotating sleeve	WM 26206
Ball-and-socket joint set non vented: retaining ring, elbow, rotating sleeve	WM 26207
Rotating sleeve vented	WM 26254
Forehead support: Forehead cushion, forehead support coarse adjustment component, forehead support fine adjustment component	WM 26201
Forehead support GEL: Forehead cushion GEL, forehead support coarse adjustment component, forehead support fine adjustment component	WM 26215

Article	Article no.
Forehead cushion	WM 26200
Forehead cushion GEL	WM 26209
Mask cushion 30 hPa, Size S	WM 26210
Mask cushion 30 hPa, Size M	WM 26220
Mask cushion 30 hPa, Size L	WM 26230
Mask cushion 30 hPa, Size XL	WM 26240
Mask cushion 40 hPa, Size S	WM 26211
Mask cushion 40 hPa, Size M	WM 26221
Mask cushion 40 hPa, Size L	WM 26231
Mask cushion 40 hPa, Size XL	WM 26241
Mask cushion GEL, Size S	WM 26212
Mask cushion GEL, Size M	WM 26222
Mask cushion GEL, Size L	WM 26232
Headgear HEADstrap	WM 26360
JOYCEstrap headgear, blue, incl. headgear clips	WM 26290
Headgear clips (4 pcs.)	WM 26205
Port cap	WM 26208
Instructions for use DE, EN, FR, NL, IT, DA, NO, SV, ES, PT, TR RU (pdf, Internet: www.weinmann.de)	WM 16950 WM 16814

Accessories

Article	Article no.
Headgear, blue: JOYCEstrap EXTRA, 3-point, Size M, incl. headgear clips	WM 26291
Exhaled air system: Silentflow 2 Noise reduction system	WM 23600 WM 23685

8 Specifications

Technical data

JOYCE	Value
Device class according to Richtlinie 93/42 EWG:	II a
Dimensions (WxHxD):	80 x 130 x 90 mm ¹
Weight:	c. 85 grams ¹
Therapy pressure range: mask cushion soft: mask cushion hard: mask cushion GEL:	4 to 30 hPa 4 to 40 hPa 4 to 40 hPa
Hose connector: cone complying with EN ISO 5356-1 vented: non vented:	Ø 22 mm (male) Ø 22 mm (female)
Pressure measurement ports:	Ø 4 mm
Width of headgear clips:	max. 20 mm
Temperature range Operation: Storage:	+ 5 °C to + 40 °C - 20 °C to + 70 °C
Flow resistance vented: at 50 l/min at 100 l/min non vented: at 50 l/min at 100 l/min	0.02 hPa 0.11 hPa 0.09 hPa 0.37 hPa
Quoted two-figure noise emission value to ISO 4871: – sound pressure level – sound power level (uncertainty factor: 3 dB(A))	15 dB(A) 23 dB(A)
Working life	up to 12 months ²

¹ depending on size and design

² If the materials used to produce masks are exposed to aggressive cleaning agents, for example, they age. In individual cases,

therefore, you may need to replace your mask sooner (see "5. Functional checks" on page 16).

CE 0197

Subject to design changes


Materials

Article	Material
Mask cushion, port cap	Silicone
Mask cushion GEL	Silicone, silicone gel
Forehead cushion, retaining ring	TPE, Polypropylene
Forehead cushion GEL	Silicone, silicone gel, polycarbonate
Mask frame, elbow	Polycarbonate
Headgear JOYCEstrap	Cotton, polyamide, elastane
Headgear HEADstrap	Nylon/ Spandex Polyurethane Polyester/Lycra
Rotating sleeve, forehead support coarse adjustment component, forehead support fine adjustment component	Polypropylene
Headgear clips	POM

All mask components are latex-free.

Pressure/flow characteristic curve

The pressure/flow characteristic curve shows outlet flow as a function of therapy pressure.

9 Warranty

- Weinmann guarantees that the product is free from defects, if used as intended, for a period of 6 months from date of purchase.
- Warranty claims must be accompanied by proof of purchase showing merchant and date of purchase.
- We acknowledge no liability under warranty in the event of:
 - Non-compliance with the instructions for use
 - Incorrect operation
 - Improper use or improper treatment
 - Repairs to the device attempted by unauthorised persons
 - Force majeure
 - Transport damage due to inadequate packing of returned items
 - Deterioration in normal use and customary wear
 - Failure to use original replacement parts.
- Weinmann does not accept liability for damage consequential upon defects unless due to willful acts or gross negligence or in the event of harm to life or limb due to minor negligence.
- Weinmann reserves the right to choose whether to rectify the defect, to supply an item free of defects or to make an appropriate reduction in the purchase price.

- If a warranty claim is refused, we will not accept the transport cost to or from our address.
- This does not affect your statutory rights.

10 Declaration of conformity

Weinmann Geräte für Medizin GmbH + Co. KG declares herewith that the product complies fully with the respective regulations of the Medical Device Directive 93/42/EEC.

The unabridged text of the Declaration of Conformity can be found on our website at www.weinmann.de

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Geräte für Medizin GmbH+Co.KG