0104-0545  
0104-0544

0104-0545  
0104-0544

**KimVent**  
*Closed Suction System for Neonates / Pediatrics*

BALLARD<sup>®</sup> TRACH CARE<sup>®</sup> Technology

*Directions for Use*


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6

CE 0086

Distributed in the U.S. by Kimberly-Clark Global Sales, LLC, Roswell, GA 30076 USA  
 In USA, 1-800-KCHELPS • www.kchealthcare.com  
 Kimberly-Clark, Roswell, GA 30076 USA  
 Kimberly-Clark N.V., Belgicastraat 13, 1930 Zaventem, Belgium  
 Sponsored in Australia by Kimberly-Clark Australia Pty Limited; 52 Alfred Street, Milsons Point, NSW 2061 • 1-800-101-021

**GENERAL DIRECTIONS FOR USE**

These instructions apply to the following families of KimVent® products:

- Neonates / Pediatrics Y
- Neonates / Elbow

**Warning:**

1. Do not trim or cut the endotracheal tube (not supplied) while the KimVent® Closed Suction System is attached, otherwise the KimVent® catheter may also be cut and that portion of the catheter may be aspirated into the lower respiratory tract of the patient and may cause death or serious injury.
2. Do not reuse, reprocess, or resterilize this medical device. Reuse, reprocessing, or resterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.

**Cautions:**

1. Inspect KimVent® catheter package before opening. Do not use product if packaging has been compromised. Non-sterile contents may cause infection.
2. Excess fluid in heat and moisture exchanger (HME) may increase gas flow resistance. When introducing fluid into airway connection, ensure that fluid does not enter HME.
3. Stop withdrawal when black mark on the tip of the catheter is visible within the dome (Fig. 2). Withdrawal of the black mark past the dome, may inflate the protective sleeve of the suction catheter and may cause oxygen desaturation.
4. Always place the thumb valve in the locked position when not in use to prevent inadvertent activation.
5. Single patient use only.
6. Do not use for more than 24 hours.
7. Inspect Sodium Chloride vial prior to opening. Do not use product if vial has been compromised. Compromised contents may cause infection.
8. Rx Only.
9. Select the appropriate size KimVent® catheter. Most experts suggest that the catheter selected should occupy no more than one half of the internal diameter of the artificial airway.
10. Do not leave the catheter within the airway. Always withdraw until black mark on the tip of the catheter is visible within the dome. Any catheter left extended into the airway will cause increased airway resistance.
11. Use appropriate regulated vacuum levels. Though experts in the neonatal community cannot agree on a set regulated level of suction, some experts recommend you should use no more than -100 mm Hg.
12. Use appropriate suction technique. Though experts in the neonatal community cannot agree on duration of the suction procedure, some experts recommend the entire suction procedure should last no longer than 5 to 10 seconds. Actual duration of negative pressure should be no more than 5 seconds per episode.
13. Always use caution and good clinical judgement no matter what ventilator mode is in use. If the clinician notes any signs of suction intolerance such as oxygen desaturation, negative ventilator system pressures, patient stress or excessive discomfort, adjustments to the ventilator settings may need to be made. These adjustments (please refer to the ventilator's instructions for use) may include manipulation of the inspiratory trigger sensitivity, inspiratory

volume or flow rate require the use of precautions may include:  
 14. This medical device classified in the EC based on data for that exposure to D in humans. A risk indicated patient risk, was performed when used as directed

**Setup for Manifold**

1. Select appropriate
2. Remove original ei (mm) Kimberly-Cl
3. Reconnect the ven
4. Depress and hold t regulator to the di
5. Release thumb cor circuit and ET adaj

**Setup for Neonat**

1. Select appropriate tube).
2. Attach KimVent® cat
3. Attach thumb cont
4. Hold thumb contrc regulator to the di
5. Remove the origin ventilator circuit.

**Suggested Suctio**

1. Stabilize KimVent® catheter into the e opposite hand (Fi
2. Release catheter a
3. Depress thumb cor withdrawal when dome (Fig. 2).
4. Repeat steps 1-3 a

**Patient Lavage Pr**

1. Advance the cathe
2. Instill desired amo
3. Advance catheter t as above.

**Catheter Irrigatio**

- (A) Neonatal / Pedia  
 Use same port for  
 (B) Manifold and EI  
 Use separate port I (nearest the patie mode is in use. If the clinician notes any signs of suction intolerance such as oxygen desaturation, negative ventilator system pressures, patient stress or excessive discomfort, adjustments to the ventilator settings may need to be made. These adjustments (please refer to the ventilator's instructions for use) may include manipulation of the inspiratory trigger sensitivity, inspiratory
1. Withdraw the cath should be comple
  2. Introduce fluid slor thumb control val
  3. Continue irrigatio

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