### Instruction Sheet

English • Deutsch • Français • Italiano • Español • Português •  
Nederlands • Svenska • Dansk • Norsk • Suomi • Ελληνικά • Polski •  
Magyar • Русский • Türkçe • Česky • 简体中文 • 繁體中文

- Refer to your positive airway pressure device User Guide for details on setting up the Leak Valve.
- The Leak Valve can be connected to components with a 22 mm female ISO connector on the patient side (C) and device side (D) and to components with a 15 mm male ISO connector on the patient side (C) only.
- Position the Leak Valve so that air from the vent holes does not blow directly onto the patient's chest.

**Note:** Connecting accessories with a small internal diameter to the device side of the Leak Valve may affect valve function.

### Cleaning

- Wash the Leak Valve daily/after each use (or more often if required) in warm water using mild detergent or by wiping all accessible surfaces with a 70% v/v isopropyl alcohol swab.
- Rinse thoroughly and allow to air dry out of direct sunlight.
- Inspect the Leak Valve on a regular basis. Discard and replace after 7 days of use or if visibly soiled after cleaning.

**Note:** Do not try to disassemble the Leak Valve.

### Disposal

The Leak Valve does not contain any hazardous substances and may be disposed of with normal household refuse.

### Technical Specifications

Pressure range:	2 to 40 cm H <sub>2</sub> O (bilevel), 4 to 20 cm H <sub>2</sub> O (CPAP) (see also Pressure-Flow diagram)
Operating temperature:	0°C to +40°C
Operating humidity:	0 to 100%
Storage and transport temperature:	-20°C to +60°C
Storage and transport humidity:	0 to 95% non condensing
Dimensions:	81 mm (length), 40 mm (diameter)
Weight:	18 g
Pneumatic impedance:	≤0.5 cmH <sub>2</sub> O at 50 L/min and ≤0.3 cmH <sub>2</sub> O at 20 L/min
Physical deadspace:	≤30 mL

### Symbols


- Manufacturer; REF Catalogue number; Humidity limitation;
- Temperature limitation; CE CE labelling in accordance with EC directive 93/42/EEC, class II a; Disposable: For single patient use only - do not reuse; Does not contain natural rubber latex;
- Consult accompanying documents; LOT Batch code;
- Handle with care; Keep dry; This way up

For more information please refer to [www.resmed.com](http://www.resmed.com).

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
### English

### Intended Use

The ResMed Leak Valve (24988) is an accessory for single patient use (1 to 7 days), intended to provide leak ventilation during therapy with a positive airway pressure device. It is designed with an anti-asphyxia valve (A) to allow for spontaneous breathing and prevent CO<sub>2</sub> rebreathing in the event of ventilator failure. It also allows continuous airflow through the vent holes (B).

### ⚠️ WARNINGS

- Read this Instruction Sheet and your device User Guide before use.
- Ensure that the vent holes and the anti-asphyxia openings are not covered.

### Setup

### ⚠️ WARNING

Always set up the Leak Valve in the breathing circuit with the arrows and the symbol pointing in the direction of air flow from the device to the patient.

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