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


Accessory


Legend

- | | | | |
|----------|----------------------------|-----------|-----------------------|
| 1 | Hose holder | 7 | Plug |
| 2 | Corrguated hose | 8 | O ₂ outlet |
| 3 | O ₂ supply hose | 9 | LED |
| 4 | Oxygen valve | 10 | Connecting cable |
| 5 | Velcro strips | 11 | O ₂ inlet |
| 6 | Socket | | |

Accessory

- | | |
|-----------|----------------------------|
| 12 | Double hose |
| 13 | Corrguated hose |
| 14 | Sleeve for corrugated hose |
| 15 | Adapter (transparent) |
| 16 | Sleeve for adapter |
| 17 | Adapter (grey) |

Special marks on the device

Unit	
□	Unit inlet
□	Unit outlet
Rating plate	
☒	Do not dispose of the unit in domestic waste!
做人	Protection class B
□	Protection class II, protective insulation
年	Year of manufacture
□	O ₂ -pressure max 1.5 bar = 1500 hPa

Safety information in these instructions

Safety information is marked as follows in these instructions for use:

Warning!

Warns of a risk of injury and potential material damage.

Caution!

Warns of material damage and potentially false therapy results.

Note

Contains useful tips.

Contents

1.	Description	33
1.1	Purpose	33
1.2	Functional description	34
2.	Safety instructions	35
2.1	Safety regulations	35
3.	Setting up the device	36
4.	Operation	38
4.1	Switching the device on and off	38
4.2	Using the device	38
4.3	After use	38
5.	Hygienic preparation	39
5.1	Intervals	39
5.2	Cleaning and disinfection	39
5.3	Changing patients	41
6.	Performance check	42
6.1	Intervals	42
6.2	Method	42
7.	Troubleshooting	43
8.	Maintenance	44
8.1	Intervals	44
8.2	Changing the double hose	44
8.3	Waste disposal	45
9.	Supply Schedule	46
9.1	Standard supply schedule	46
9.2	Spare parts	46
9.3	Accessory	46
10.	Technical specifications	47
11.	Warranty	48
12.	Declaration of Conformity	49
13.	Index	50

1. Description

1.1 Purpose

Important:

Like the therapy devices, the oxygen valve does not have an alarm that warns in the event of any faults. It must therefore be used only for patients who can endure the absence of oxygen therapy over several hours, so that an unnoticed failure of the oxygen supply will not represent a risk.

During therapy the oxygen valve can be used to safely supply up to 4 l/min oxygen into the mask. It can be used together with an oxygen concentrator, an oxygen cylinder with a pressure reducing valve or liquid oxygen.

You can use the valve with the following devices:

CPAP devices

- SOMNOsmart 2
- SOMNOset
- SOMNObalance (e)
- SOMNOsoft 2 (e)

nCPAP devices

- SOMNOcomfort
- SOMNOcomfort 2 (e)
- CPAP 20(e)
- SOMNOsoft +

BiLevel devices

- SOMNOvent S
- SOMNOvent ST
- SOMNOvent auto-S
- SOMNOvent auto-ST

CS devices

- SOMNOvent CR

The valve opens when the therapy device is switched on and is working properly. In the event of a fault (e.g. power loss) the supply of oxygen is interrupted to prevent an accumulation of oxygen in the device and thus to prevent the risk of fire that may result with oxygen.

1.2 Functional description

The valve is electrically connected with the therapy device and scans the latter constantly to ensure it is working correctly. If the therapy device is switched on, the internal solenoid valve opens.


This is indicated by a green LED. If the therapy device is switched off or does not work properly, the valve closes. The control lamp goes out. Once the device is working properly again, the valve opens again.

2. Safety instructions


2.1 Safety regulations

- Please read this instruction manual carefully. It is part of the device and must always be available.
- Please also follow the instruction manuals included with your therapy device, your oxygen system and all accessories.
- Use the oxygen valve for the described purpose only (see "1. Description" on page 33).
- During therapy with oxygen, your clothes, bedclothes or hair can come into contact with oxygen-enriched air. This means easy flammability. Therefore, do not smoke nor use candles or fire nearby. This also applies after therapy until oxygen has been removed by thorough ventilation.
- Do not exceed the maximum supply quantity of 4 l/min. Higher rates can result in oxygen accumulation in the therapy device.
- If used or dosed improperly, oxygen therapy can have side-effects. Therefore, use the device only if it has been prescribed by a doctor. Always ensure that the flow prescribed by the physician is applied (maximum 4 l/min.).
- Always follow the safety instructions in the instruction manual supplied with your oxygen system.
- Have repairs carried out only by the manufacturer Weinmann or refer to your dealers for repairs.
- Please follow section "5. Hygienic preparation" on page 39 to prevent infection or bacterial contamination.
- Do not allow liquids to get into the therapy device or into the housing of the oxygen valve. The electronic components could be damaged.
- If third-party items are used, functional failures may occur and fitness for use may be restricted. Biocompatibility requirements may also not be met. Please note that in such cases, any claim under warranty and liability will be voided if neither the accessories nor genuine replacement parts recommended in the instructions for use are used.
- As there may be problems with communication in the course of setting therapy parameters (e.g., Softstart time, humidity level, therapy pressure), pull out the plug of the connecting cable to the therapy unit.



3. Setting up the device



1. Set up your therapy device and oxygen system as described in the relevant instruction manuals. Secure the O₂ supply hose to your hose system with the hose holders.
2. Put the oxygen valve in a position so that you can connect the connecting cable and the O₂ supply hose. If you intend to use the Velcro fastening on the therapy device, proceed as follows:

- The surfaces for adhesive must be dry and free of grease and dust, etc.
 - Use adhesive to attach the two soft strips of Velcro to each side of the label underneath the oxygen valve.
 - Attach the hard strips to them.
 - Remove the protective film from the back of the hard strips (1.) and press the valve firmly on to the desired position (2.) on the therapy device, e.g., as shown at the back on the right.
 - Leave the valve there for at least one hour before opening the Velcro fastening to allow the adhesive time to bind with the surfaces.
3. Insert the adapter **17** of the hose system into the output of the device or output of the humidifier.

4. Use the connecting cable to connect the oxygen valve to the socket **6** of your therapy device. The control lamp may light up briefly as the connection is being made.
5. Attach the device-side end of the O₂ supply hose to the outlet side of the valve (small socket on the left) **8**.

6. Insert the hose from your oxygen system onto the inlet side of the valve (large socket on the right) **11**.
 7. Assemble the exhalation system and the mask according to the relevant instruction manuals.
 8. Remove the seal from the port on the mask and attach the O₂ supply hose there.
-
9. To humidify respiratory air, use the humidifier obtainable for your therapy device. Additional humification of the oxygen with a bubble humidifier is allowed, but is not necessary.

The use of an atomizer is not permitted.

A hose system with an internal O₂ supply hose is available as an accessory.

You can use the hose system with the following devices:

CPAP devices

- SOMNOsmart 2
- SOMNOset

nCPAP devices

- SOMNOcomfort
- SOMNOsoft +

BiLevel devices

- SOMNOvent S
- SOMNOvent ST
- SOMNOvent auto-S
- SOMNOvent auto-ST

CS devices

- SOMNOvent CR

4. Operation

4.1 Switching the device on and off

The oxygen valve is ready for use when it has been connected as described above. You simply have to switch on your therapy device and your oxygen system (with the prescribed flow). The valve automatically switches on the oxygen supply to the therapy device and switches it off again.

4.2 Using the device

No action is required when the device is in operation.

4.3 After use

Please switch off the devices in the following sequence.

1. Switch off the oxygen supply by switching off the oxygen concentrator or the cylinder valve (bottled oxygen) or the flow controller (liquid oxygen).
2. Switch off your therapy device.
3. Clean the mask and the exhalation system as described in the relevant instruction manuals.
4. To clean the hose system, close or remove the O₂ supply hose.


5. Hygienic preparation

5.1 Intervals

Clean the hose system whenever it is dirty or at least once a month. The O₂ supply hose has to be changed at least every six months or earlier if it's soiled. (If a hose system is used, the internal O₂ supply hose has to be replaced according to this schedule (see "8.2 Changing the double hose" on page 44).

5.2 Cleaning and disinfection

Cleaning the housing


Caution! Before cleaning the device, remove the connecting cable from the socket **6** on the therapy device.

Remove it by pressing slightly on the small lever underneath the plug. Do not allow any water or disinfectant to penetrate the housing. Never place the valve in water.

Use a soft damp cloth to wipe the housing including connecting sockets and the connecting cable. The valve must be completely dry before use.

Clean the tube system

The tube system WM 23737 may be cleaned in warm water with a temperature of up to 70°C. Do not allow any water to enter the double tube. Proceed as follows:

1. Detach the tube system from the device and from the exhalation system.
2. Pull the double tube a little way out of the creased tube. Seal the opening with the sealing plug supplied.
3. Close the small opening of the grey adapter and the O₂ supply tube with the sealing plugs supplied.
4. Clean the creased tube thoroughly to remove all residues, using warm water containing a little washing-up liquid. Thoroughly wash out the inside of the tube.

5. Rinse the creased tube thoroughly inside and out with clear warm water.
6. **Thoroughly shake out the tube system.**
7. Hang up the tube system and allow it to drip dry, in order to prevent moisture entering your therapy device.
8. Use the therapy device to dry the tube system completely. For details see the instructions for use for your therapy device.

Disinfecting the housing

The housing, the connecting sockets and the connecting cable can be cleaned by simply wiping them with disinfectant. We recommend TERRALIN.


Disinfect the tube system

As a disinfectant we recommend GIGASEPT FF. You are recommended to wear suitable gloves (e.g. household or disposable gloves) during disinfection procedures. Observe the instructions for the disinfectant used.

No disinfectant may penetrate the double tube.


To disinfect, proceed as follows:

Remove double tube:

1. Release both sleeves from the adapters.
2. Pull the double tube out of the creased tube.
3. Take the short piece of the double tube off the gray adapter. Then take the tube out of the transparent adapter.
4. Release the long part of the double tube from the transparent adapter.

Disinfect double tube



1. Seal all four openings of the double tube using the sealing plugs supplied.
2. Clean the double tube with a little disinfectant in hot water so that there are no residues.
3. After disinfecting, rinse the double tube thoroughly with distilled water.
4. Leave the parts to dry completely.

Disinfect creased tube and adapter

1. Proceed as described in "Clean the tube system".
2. After disinfecting, rinse thoroughly with distilled water.

Fit double tube

1. Push the short part of the double tube through the transparent adapter and then push it as far as possible onto the gray adapter.
2. Push the sleeve back onto the gray adapter.
3. Push the long part of the double tube into the groove of the transparent adapter. The two ends of the tube must lie smoothly next to one another.
4. Insert the free end of the double tube in the creased tube.
5. Push the sleeve of the creased tube onto the transparent adapter.

After disinfecting

Dry the assembled tube system completely with the aid of the therapy device. The procedure can be found in the instructions for use for your therapy device.

5.3 Changing patients

To prepare the device for use by a new patient, replace the O₂ supply hose or the hose system and clean the outside surfaces of the O₂ valve by wiping them with disinfectant.

6. Performance check

6.1 Intervals

Check daily that the oxygen valve works properly. If you detect a fault, do not continue to use it until the fault has been remedied.

6.2 Method

1. Assemble the valve so that it works with the therapy device. The control lamp **9** on the top of the valve housing does not light up.
2. Close the opening of the mask, for example, using your thumb or hand.
3. Switch on your therapy device and then your oxygen supply. When the therapy device is switched on, the valve opens with a slight "click"; this is indicated by the green indicator lamp **9**. You can now set the prescribed flow on the flowmeter of your oxygen system. If this is not possible, first check the function of your oxygen system. (For example: Is the cylinder empty? Are any of the hoses kinked?)
4. Switch the device off again. The valve closes audibly, the control lamp **9** goes out and the flow displayed for the oxygen system drops to "0".

If correct function is not possible, contact your dealer or the manufacturer for repairs to the equipment.

7. Troubleshooting

Fault	Cause	Remedy
Display does not light up or only sometimes although the therapy device is working properly. The prescribed flow of max. 4 l/min is not reached	No power	Check that connecting cable 10 is securely attached.
	Mask leaks	Adjust the head cover so that the mask is sealed. Replace mask if defective.
The prescribed flow is not reached	Hoses are bent	Check hose connections.
	Excessive resistance in the hose system, output pressure of the concentrator is too low.	Shorten oxygen supply hoses. If necessary, use other concentrator (e.g., OXYMAT)

In the event of any faults that cannot be remedied immediately, contact the manufacturer Weinmann or your dealer immediately to have the valve repaired.

Do not continue to use it because a fault may increase the risk of fire.

8. Maintenance



8.1 Intervals

The oxygen valve requires no maintenance.

For hygienic reasons, we recommend replacing the O₂ supply hose or the double hose after six months, or earlier if it is soiled.

8.2 Changing the double hose

The double tube consists of a short part which is used for pressure measurement during therapy, and a long part that is used for O₂ supply. To replace the double tube, proceed as follows:

1. Detach both sleeves from the adapters.
2. Pull the double tube out of the creased tube.
3. Pull the shorter part of the double tube off the grey adapter. Then pull the tube out of the transparent adapter.
4. Remove the long part of the double tube from the transparent adapter.
5. Dispose of the old double tube.
6. Push the short part of the new double tube through the transparent adapter, then push it as far as possible onto the grey adapter.
7. Push the sleeve onto the grey adapter again.
8. Press the long part of the double tube into the groove in the transparent adapter. The two tube ends must lie smoothly side by side.
9. Insert the free end of the double tube into the creased tube.
10. Push the sleeve of the creased tube onto the transparent adapter.

8.3 Waste disposal

To dispose of the unit properly, please contact a licensed and certified electronic waste disposal merchant. Names and addresses can be obtained from your Environmental Officer or municipal authorities. The unit packaging (cardboard box and inserts) can be disposed of as waste paper.

9. Supply Schedule

9.1 Standard supply schedule

Oxygen valve, packed

WM 24042

Part	Order No.
Oxygen valve, Basic unit	WM 24102
O ₂ supply hose	WM 24015
6 hose holders	WM 24268
Velcro strips for fastening	WM 24687
Instruction manual	WM 16792

9.2 Spare parts

Part	Order No.
O ₂ supply hose	WM 24015
Velcro strips for fastening	WM 24687

9.3 Accessory

Part	Order No.
Hose system for oxygen valve, consisting of: – Corrugated hose 1830 mm – Double hose – Adapter (transparent) – Adapter (gray) – Pressure connector	WM 23737 WM 23146 WM 24626 WM 23901 WM 24123 WM 24129

Hose system for oxygen supply valve WM 23737 is not suitable for the following devices:
SOMNOcomfort 2 (e), SOMNObalance (e), CPAP 20(e), SOMNOsoft 2 (e).

10. Technical specifications

Product class in compliance with 93/42/EEG	IIa
Dimensions (W x H x D) in cm	4.5 x 8.5 x 3.0
Weight	approx. 123 g
Temperature range	
– Operation	+5 °C to +40 °C
– Storage	-20 °C to +70 °C
Air pressure range	700 – 1060 mbar
Permissible humidity during operation and storage	≤95 % rF (no condensation)
Electrical connection	7–15 V DC
Current consumption	
– In operation	96.0 mA
– Standby	8.5 mA
Classification in compliance with EN 60601-1	
– Protection type for electrical shock	Protection class II
– Protection type for electrical shock	Type B
Electromagnetic compatibility (EMC) in compliance with EN 60601-1-2	
– Radio interference suppression	EN 55011
– Radio interference immunity	IEC 1000-4 Part 2-6 and 11
Max. permissible flow	<4 l/min O ₂
Min. O ₂ advance pressure at the valve inlet for 4 l/min O ₂	250 hPa
Max. permissible ram pressure of the O ₂ system	1500 hPa
Permissible therapy pressures	4 bis 30 hPa (1hPa = 1mbar)

Subject to modification

11. Warranty

- Weinmann offers a warranty that the product, when used in accordance with requirements, will remain free from defects for a period of two years from date of purchase. For products whose durability is clearly indicated as less than two years, the warranty expires on the expiration date indicated on the packaging or in the user's manual.
- Claims against the warranty can be made only when accompanied by the sales receipt, which must show salesperson and date of purchase.
- We offer no warranty in the case of:
 - Disregard of usage instructions
 - Operating errors
 - Improper use or improper handling
 - Third-party intervention by non-authorized persons for the purpose of device repair
 - Acts of God, e.g., lightning strikes, etc.
 - Transport damage as a result of improper packaging of returned items
 - Lack of maintenance
 - Operational and normal wear and tear, which includes, for example, the following components:
 - Filters
 - Batteries / rechargeable batteries
 - Articles for one-time usage, etc.
 - Failure to use original spare parts
- Weinmann is not liable for consequential harm caused by a defect if it is not based on intention or gross negligence. Weinmann is also not liable for minor physical injury to life or limb resulting from negligence.
- Weinmann reserves the right to decide whether to eliminate defects, to deliver a defect-free item or to reduce the purchase price by a reasonable amount.
- If Weinmann rejects a claim against the warranty, it assumes no expense for transport between customer and manufacturer.
- Implied warranty claims remain unaffected by these changes.

12. Declaration of Conformity




Weinmann Geräte für Medizin GmbH + Co. KG declares herewith that the product complies fully with the respective regulations of the Medical Device Directive 93/42/EEC. The unabridged text of the Declaration of Conformity can be found on our website at www.weinmann.de

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Geräte für Medizin GmbH+Co.KG**