**S<sub>P</sub>O<sub>2</sub> Accuracy**

70 - 100%: A<sub>rms</sub> < 3.0%

0 - 69%: Unspecified

Accuracy test methods available from the manufacturer upon request

**LED Specifications**

LED	Wavelength	Radiant Power (typical)	Brand
Red	660 nm	3.5 mW	All below
Infrared	885 nm	3.5 mW	Nellcor, CSI
	890 nm	3.5 mW	HP / Philips
	905 nm	3.5 mW	Datex, BCI, Minolta, Bluepoint Medical
	910 nm	3.5 mW	Nonin
	935 nm	3.5 mW	Nihon-Kohden
	940 nm	3.5 mW	Ohmeda, Novamatrix

**DIRECTIONS FOR USE**

**Intended Use:** Bluepoint Medical SoftCap® oximetry sensors are indicated for use in continuous non-invasive monitoring of arterial oxygen saturation and pulse rate. They are original or fully compatible replacement sensors intended for use with major brands of pulse oximeters and patient monitors. Select the sensor model number that matches the brand of oximeter or patient monitor that you are using.


**Precautions:** For further instructions, warnings and precautions refer to the operating instructions provided with your oximeter. Do not autoclave.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**Warnings:** Match the sensor model with your brand of oximeter. A mismatch may cause patient injury, inaccurate measurements or equipment damage. Check the sensor every 8 hours and reposition to an alternate location every 24 hours. Avoid application of the sensors to oedematous or fragile tissue. Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure. Excessive patient motion, excessive ambient light, electromagnetic interference, dysfunctional haemoglobin, low perfusion, intravascular dyes, finger nail polish and long or artificial finger nails may affect the sensor performance and the accuracy of the measurement. Note: The materials used in the manufacture of the sensor contain no natural latex protein. The materials have undergone extensive bio-compatibility testing. Further information is available upon request.

**Application of the Sensor**

1. Connect the sensor cable to the Bluepoint Medical adaptor cable that fits your oximeter or directly to the patient monitor. Then switch on the oximeter and verify proper operation.
2. Insert the patient's digit into the sensor. The preferred digit to use for adults is the index finger. The digit is correctly inserted when:
  - The tip of the digit touches the rear guide posts.
  - The sensor cable extends along the top of the patient's hand.
3. Visually monitor the sensor site to ensure the integrity of the skin.
4. Secure the lead with medical tape.

**Cleaning the Sensor**

To clean the sensor between uses or between patients, wipe the sensor surfaces with alcohol.

**Environmental Requirements**

- Operating temperature range: 5° - 40° C (41° -104° F)
- Storage and transportation temperature range -20° - 70° C (-4° -158° F)
- Relative humidity: 15 - 95% non-condensing

**Disposal**

Follow local regulations governing disposable recycling plans regarding disposal or recycling of device components. A list of component materials is available by contacting Bluepoint Medical.

All registered trademarks recognized.

Bluepoint Model - Adult	Bluepoint Model - Paediatric	Oximeter Brand	Cable Length
SC4001	SCP4001	Nellcor	0.90 m
SC4002	SCP4002	Nellcor	2.00 m
SC4006	SCP4006	HP/Philips	2.95 m
SC4007	SCP4007	HP/Philips	2.00 m
SC4008	SCP4008	HP/Philips	2.95 m
SC4011	SCP4011	Datex	0.90 m
SC4012	SCP4012	Datex	2.95 m
SC4013	SCP4013	Datex S/5	2.95 m
SC4016	SCP4016	Ohmeda	2.95 m
SC4017	SCP4017	Ohmeda	1.00 m
SC4021	SCP4021	BCI	1.20 m
SC4026	SCP4026	Nonin	0.90 m
SC4031	SCP4031	Minolta	0.33 m
SC4036	SCP4036	CSI	0.90 m
SC4037	SCP4037	CSI	2.95 m
SC4041	SCP4041	Novamatrix	2.95 m
SC4046	SCP4046	Nihon-Kohden	2.95 m
SC4047	SCP4047	Nihon-Kohden	1.60 m
SC6500	SCP6500	Bluepoint Medical	1.20 m
SC6501	SCP6501	Bluepoint Medical	0.30 m

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#### Bluepoint Medical GmbH & Co. KG