
# DEVILBISS®

## 9055 Series *DeVilbiss*® RPM™ Bilevel CPAP System Instruction Guide

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**TABLE OF CONTENTS**

Symbols	4
Important Safeguards	6
Introduction	7
Important Parts	7
Accessories/Replacement Items	7
Unpacking Your Bilevel CPAP	8
Other Required Components	8
Optional Accessories	8
Assembly of Your Bilevel CPAP	8
Normal Operation of Your Bilevel CPAP	8
Additional Operating Modes and Features	9
Epworth Sleepiness Scale Questionnaire	10
Travel Information	11
DC Operation	11
Supplemental Oxygen	12
Cleaning	12
Cabinet	12
Air-Inlet Filter	12
Mask, Headgear and Tubing	12
Recommended Maintenance Interval for Equipment Owner	13
Troubleshooting	13
Specifications	14
Warranty	15
Declaration of Conformity	15

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
**Symbols**

-  ON / OFF– turns blower on and off
-  Select Next Item on LCD display
-  Increase value shown on LCD display
-  Decrease value shown on LCD display
-  Delay– starts the delay cycle, if prescribed
-  Attention, consult instruction guide
-  Type BF equipment– applied part
-  Class II electrical protection– double insulated
-  Data port
- ~ 100 – 240 Volts AC, 50/60 Hz; 120 Volts AC, 400 Hz

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


## IMPORTANT SAFEGUARDS

When using electrical products, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

**▲ DANGER** Urgent safety information for hazards that will cause serious injury or death.

**▲ WARNING** Important safety information for hazards that might cause serious injury.

**CAUTION** Information for preventing damage to the product.

**NOTE** Information to which you should pay special attention.

**READ ALL INSTRUCTIONS BEFORE USING.**

### **▲ DANGER**

- Electric Shock Hazard – do not use while bathing.
- Electric Shock Hazard – do not place or drop into water or any other liquid.
- Electric Shock Hazard – do not attempt to open or remove cabinet; there are no user-serviceable internal components. If service is required, return the product to a qualified Sunrise Medical provider or service center. Opening or tampering with the product will void the warranty.

### **▲ WARNING**

- This bilevel CPAP unit should be used only with masks recommended by Sunrise Medical, your physician, or respiratory therapist. To avoid rebreathing of exhaled air 1) a mask should not be used unless the bilevel CPAP unit is turned on and operating properly, and 2) the vent hole(s) associated with the mask should never be blocked. When the bilevel CPAP unit is turned on and functioning properly, new air from the unit flushes the exhaled air out through the mask vent hole(s). However, when the bilevel CPAP unit is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation. This warning applies to most bilevel and CPAP units.
- This is not a life support device and is intended to be used on spontaneously breathing adults.
- To avoid electric shock, always unplug power cord from wall outlet power source when performing cleaning.

### **CAUTION**

- The circular DIN connector located on the back of the Model 9055 bilevel CPAP is used to attach external accessories to the product. The DIN connector must only be used with Sunrise Medical approved accessories that comply with IEC 60601-1 or an appropriate safety standard.
- Do not place this bilevel CPAP unit near other equipment or devices that create or attract electromagnetic fields. Placing the unit in electromagnetic fields greater than 10 V/m can affect its operation. Examples of such equipment are defibrillators, diathermy equipment, cellular telephones, CB radios, radio-controlled toys, microwave ovens, etc.
- Never rinse or place the unit in water. Never allow liquids to get into or around any of the ports, switches, or air filter; doing so will result in product damage. If this occurs, do not use the unit. Contact your provider for service.
- Do not place the bilevel CPAP where it can be bumped or where the power cord may trip someone.
- Only DeVilbiss humidifiers are compatible for use with the DeVilbiss Model 9055 bilevel CPAP. Do not use any other humidifier systems.

### **▲ OXYGEN USE WARNING**

- Oxygen supports combustion. To avoid possible physical injury, do not smoke while using your bilevel CPAP unit with supplemental oxygen. Do not use near hot, sparking objects or open sources of flame.
- Always turn the bilevel CPAP unit ON prior to turning on the oxygen flow and turn off the oxygen before turning the bilevel CPAP OFF.
- The flow setting of the oxygen source shall be specified by a physician.
- The bilevel CPAP is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

### **CAUTION**

- Oxygen is a prescription gas and should only be administered under the supervision of a physician.

**SAVE THESE INSTRUCTIONS.**

## INTRODUCTION

### Indications

The 9055 series DeVilbiss RPM bilevel CPAP is intended for use in treating obstructive sleep apnea in adult patients.

### What is Obstructive Sleep Apnea (OSA)?

Obstructive Sleep Apnea (OSA) is a condition where the upper airway relaxes during sleep causing the airway to obstruct. "Apnea" is the Latin word meaning "without breath." People diagnosed with sleep apnea stop breathing repeatedly during sleep, resulting in extreme exhaustion. High blood pressure, cardiac enlargement, and cardiac arrhythmia are common, as are shortness of breath and loud snoring. In a 1993 report, the National Commission on Sleep Disorders Research stated that as many as 18 million Americans may have signs and symptoms of sleep apnea.

### How Does The Bilevel CPAP System Work?

The bilevel CPAP instrument provides non-invasive continuous air pressure to your upper airway, acting as an air "splint" to keep the airway passages open during sleep. This reduces and eliminates the obstructive events, allowing uninterrupted sleep. The DeVilbiss 9055 Series is designed to provide this therapy to you quietly and comfortably so that your sleep is restorative and refreshing. The bilevel mode provides two levels of pressure established by the physician. There is an inspiratory positive airway pressure (IPAP) and an expiratory positive airway pressure (EPAP) setting. These two pressures allow for a more comfortable breathing pattern when treating your obstructive sleep apnea. In turn, this should allow you to better comply with the therapy so you receive the most benefit from this treatment.

### IMPORTANT PARTS (Figures 1, 2, 3)

- A. Keypad/LCD Display (refer to **Figure 3** for LCD symbol definition)
- B. Outlet Port
- C. Carrying Case (not shown)
- D. Air-Inlet Filter
- E. AC Power Connector
- F. DIN connector
- G. Power Cord (assure proper configuration for outlet)


Figure 1


Figure 2

### ACCESSORIES/REPLACEMENT ITEMS (Figure 4)

- |                                  |             |
|----------------------------------|-------------|
| 1. CPAP Humidifier               | 9100D       |
| 2. Humidifier Heater             | 9200 Series |
| 3. DeVilbiss® FlexAire™ Mask     | 9353D       |
| DeVilbiss FlexAire Shallow Mask  | 9353S       |
| DeVilbiss Serenity™ Mask         | 9352D       |
| DeVilbiss Serenity Shallow Mask  | 9352S       |
| DeVilbiss Serenity Gel Mask      | 9352G       |
| DeVilbiss Serenity Gel Shallow   | 9352GS      |
| 4. Air-Inlet Filter (4/pk)       | 8000D-602   |
| 5. Optional Fine Particle Filter | 8000D-603   |
| 6. Oxygen Adapter                | 7353D-601   |
| 7. 6' 22 mm Tubing               | 7351D-616   |


Figure 3

**NOTE**—The standard Air-Inlet Filter is designed to filter particles larger than 3 microns in size. The bilevel CPAP must have the standard filter in place during operation. It is recommended that the filter be replaced every 6 months or if it is damaged.

**NOTE**—The optional Fine Particle Filter is designed to filter particles larger than 0.3 microns in size. The fine particle filter is placed in the filter area of the bilevel CPAP air intake with the standard filter placed over the fine particle filter. The optional fine particle filter should be replaced every 30 days or if it appears dirty or damaged.

Ask your Sunrise Medical provider for information about other DeVilbiss mask styles and sizes.


Figure 4

## UNPACKING YOUR 9055 BILEVEL CPAP

Check to be sure that you have the following items:

- Model 9055 bilevel CPAP
- 6-foot long, 22 mm diameter air tubing
- Instruction guide
- Power cord appropriate to your wall outlet
- Carrying bag

### Other Required Components

- You will need to have a nasal CPAP mask from Sunrise Medical or as recommended by your physician or clinical specialist.

### Optional Accessories (Figure 4)

- Epworth Sleepiness Scale (ESS) Guide
- DeVilbiss Model 9100D Cool Passover Humidifier
- DeVilbiss Model 9200 Series Humidifier Heater

## ASSEMBLY OF YOUR BILEVEL CPAP

**NOTE**—Refer to *Important Parts* for part identification.

1. Place the bilevel CPAP on a stable surface.

**NOTE**—Ensure that the unit is away from heating vents, radiators, drapes and bedspreads so that the air intake is not restricted.

2. Plug the power cord into the connector on the back of the unit **Figure 2**.
3. Plug the other end of the cord into the wall outlet.
4. Prepare your mask for use, in accordance with the mask instruction guide.
5. Connect one end of the air tubing to the mask and the other end to the outlet port of the bilevel CPAP **Figure 5**.


Figure 5

**NOTE**—Some tubing is prone to kinking, which can prevent flow from reaching the mask. Use only 6' 22 mm diameter, smooth-bore, kink-resistant tubing.

6. If you are using a humidifier, assemble the system in accordance with the instructions provided with the DeVilbiss humidifier.

**NOTE**—Only the DeVilbiss 9100D humidifier and 9200 series heaters are compatible for use with the 9055 series DeVilbiss bilevel CPAP. Do not use any other humidifier systems.

## NORMAL OPERATION OF YOUR BILEVEL CPAP

1. Put your mask on and press the ON/OFF button on the keypad to start the blower **Figure 3**. The pressure in your mask will come to prescription pressure. If you do not turn the unit ON using the ON/OFF button, it will automatically start after you have taken a few breaths with the mask on. At this time you may see one of four settings on the display:

**Actual Pressure Fig. 6:** This is the actual pressure displayed while breathing on the device. **NOTE**—Due to fluctuations while breathing, the actual pressure will vary from breath to breath.

**Epworth (optional) Fig. 7:** This represents the Epworth Sleepiness Scale (ESS). This is a measurement of sleepiness and will display periodically. Follow the instructions in the Epworth section of this Instruction Guide on page 10.

**Mask Check (optional) Fig. 8:** This is a reminder that masks should be cleaned regularly and changed periodically. Press any button to acknowledge the notification.

**Filter Check (optional) Fig. 9:** This is a reminder that filters should be cleaned regularly and changed periodically. Press any button to acknowledge the notification.

2. Position and adjust mask on face according to the mask instructions **Figure 10**.

**NOTE**—The bilevel CPAP will operate normally and increase to the prescribed CPAP pressure even though the “Mask Check” (if applicable) and/or “Filter Check” (if applicable) have not been acknowledged, or “Epworth” (if applicable) has not been completed. Once the “Mask Check,” and/or “Filter Check” have been acknowledged, or “Epworth” completed, the display will show “Act Pres.”

**NOTE**—The 9055 series bilevel CPAP will turn off if power is turned off or disconnected. When power is restored, the bilevel CPAP will return to the previous on or off state within 2 seconds.

3. If you have difficulty falling asleep with the pressure at the prescription level, press the Delay button on the keypad and the pressure will stay at a low level for a preset time interval. You may repeat this delay time by pushing the button again if you have not fallen asleep during the first delay period.

**NOTE**—Your provider has set a Delay Time in the unit. Delay time can be from 0 to 45 minutes in 5-minute increments. To change: scroll to Delay using the SELECT button and use the UP/DOWN buttons to select the delay

A	c	t	P	r	e	s
X	X	.	X	.	c	m

Figure 6

E	p	w	o	r	t	h
S	t	a	r	t	>	>

Figure 7

	M	a	s	k		
	C	h	e	c	k	

Figure 8

	F	i	l	t	e	r
	C	h	e	c	k	

Figure 9


Figure 10

time. If you are unable to change the delay time, your provider has locked this setting.

4. If you take your mask off during the night, your bilevel CPAP will automatically stop providing pressure. If the Mask Alert is enabled, an audible alert will be present for 20 seconds. Pressing any keypad button will cancel this audible alert. Pressing any keypad button a second time will remove the visible alert from the display. When you put your mask back on and start breathing again, the bilevel CPAP will start operating. You can use the Delay feature again at this time to get back to sleep.

**NOTE**—A large system leak will result in a Mask Alert. If the cause of the alert is corrected, the visible and audible alert will cancel and be reset.

5. In the morning, take your mask off and press the ON/OFF button on the keypad to shut the blower OFF. The automatic shutoff feature will shut the blower OFF if you do not push the ON/OFF button. There will be a 20-second audible alert during the shutdown period if the Mask Alert is enabled. Pushing any keypad button will cancel this audible alert.

**ADDITIONAL OPERATING MODES AND FEATURES**

The Display Menu—With the unit ON, press the SELECT button on the keypad to scroll through the items on the display. Not all of these will display on your unit. Some are optional, depending on what your physician and/or home medical equipment provider has selected. The table below shows the items you may see and describes the feature associated with each item.

Item	Displayed Text Examples	Description	Range	Comments
Actual Pressure	Act Pres 20.0 cmH	Actual mask pressure calculated from CPAP pressure and flow	3-25 cmH <sub>2</sub> O	The value displayed may vary slightly from your prescription based on operating conditions.
Compliance Bar Graph	Compl % 70	The bar graph is displayed to give feedback on device usage. Graph covers usage over the last 30 days of operation (rolling 30-day basis).	0-100%	The value displayed is your actual use time versus the time you are expected to use the bilevel CPAP.
CPAP Pressure	CPAP 20.0 cmH	Prescription setting – this will be seen only if your bilevel CPAP is to be used in the CPAP mode.	Delay Pressure to 25 cmH <sub>2</sub> O, in 0.5 cm increments	Your provider will set the CPAP prescription for your bilevel CPAP.
IPAP Pressure	IPAP 20.0 cmH	Prescription setting – this will be your inspiration pressure; it is seen if your bilevel CPAP is to be used in the bilevel mode.	EPAP to 25 cmH <sub>2</sub> O, in 0.5 cm increments	Your provider will set the IPAP prescription for your bilevel CPAP.
EPAP Pressure	EPAP 10.0 cmH	Prescription setting – this will be your exhalation pressure; it is seen if your bilevel CPAP is to be used in the bilevel mode.	Delay pressure to IPAP, in 0.5 cm increments	Your provider will set the EPAP prescription for your bilevel CPAP.
Delay Time	Delay 10 min	Selectable delay time setting	0 to 45 minutes in 5 minute increments	You can select the delay you desire using the UP/DOWN button, if the provider has unlocked this setting.
Delay Pressure	Delay Pr 4.0 cmH	Pressure during the delay	3 cmH <sub>2</sub> O to EPAP (or CPAP in CPAP mode), in 0.5 cm increments	Your provider has set the delay pressure value.
IPAP Rounding	IPAP Rnd 1	IPAP rounding waveform profile	0 to 5, in increments of 1	Your provider has set the IPAP rounding value.
EPAP Rounding	EPAP Rnd 1	EPAP rounding waveform profile	0 to 5, in increments of 1	Your provider has set the EPAP rounding value.
Enable Mask Alert	Msk Alrt Enabled	Indicates whether the 20 second audible mask alert feature is enabled or disabled.	Enabled or Disabled	The mask alert can be enabled or disabled by using the UP/DOWN button, unless the provider has locked the alert.
Humidifier Present	9100 Hum Attached	Indicates that a DeVilbiss humidifier is connected to the CPAP.	Attached or Removed	The indication of the attachment or removal of the humidifier can be changed using the UP/DOWN button, unless the provider has locked this setting.
Compliance Meter	ComplHrs 123456.7	The total number of hours the bilevel CPAP has sensed breathing.	N/A	
Hour Meter	Unit Hrs 123456.7	The total number of hours that the bilevel CPAP blower has been ON.	N/A	
Firmware Version	Firmware M5 V0.17	Shows version of installed firmware.	N/A	
Serial Number	Serial # 95000001	Shows the serial number of the bilevel CPAP unit.	N/A	
Language Selection	Language English	Shows the selected LCD display language.	Up to 15 selections	The display language can be changed using the UP/DOWN button, unless the provider has locked the language.

The display indications below are not in the scroll menu, but they are normal operating display items:

Item	Displayed Text	Description	Range	Comments
Epworth Sleepiness Scale (only seen on a periodic basis)	Epworth Start >>	Patient Epworth Sleepiness Scale – see the description of how the Epworth Sleepiness Scale Questionnaire is used in “Epworth Sleepiness Scale” on page 10.	N/A	Your provider may enable this feature, and it may be the first item displayed on a periodic basis.
Filter Check Notification (only seen on a periodic basis)	Filter Check	This display reminds you to check the inlet filter. See “Cleaning” section on page 12.	N/A	Your provider may enable this feature, and it may be the first item displayed on a periodic basis.
Mask Check Notification (only seen on a periodic basis)	Mask Check	This display reminds you to check the condition of your mask. See the instructions that came with your mask.	N/A	Your provider may enable this feature, and it may be the first item displayed on a periodic basis.
Blower OFF Indication	OFF	Indication that AC power is connected but blower is turned OFF	N/A	
Delay ON Indication	Delay ON	Informs the user that the delay is running (active)	N/A	

The bilevel CPAP unit will operate (provide bilevel or constant-pressure CPAP therapy) while the Epworth, or mask or filter alerts are displayed. The blower ON/OFF and DELAY keypad buttons will continue to function regardless of these alert prompts.

#### EPWORTH SLEEPINESS SCALE®

**NOTE**—Please read all instructions before proceeding.

The Epworth Sleepiness Scale Questionnaire is a tool used to assess your sleepiness and your benefit from using the bilevel CPAP regularly. The unit will record and display the results of a patient-completed Epworth Sleepiness Scale (ESS) at various intervals. The interval is determined by your home medical equipment provider and/or physician. The ESS feature can be turned on or off by the health care provider at device set-up.

The Epworth Sleepiness Scale (ESS) is an eight-question sleepiness assessment that has four possible numerical answers for each question (0, 1, 2, or 3).

The questionnaire can be temporarily skipped by turning the blower OFF and then back ON. The unit must be off for two hours before ESS is displayed again. For all eight (8) ESS questions, the keypad UP/DOWN buttons will be used to scroll through the possible answers (0, 1, 2, or 3), and the SELECT button will be used to enter your answer and advance sequentially through the process. The eight (8) ESS questions that are to be answered and the scale for your answers of 0 to 3 are provided on the next page and in a separate booklet. The questions are to be referenced as you proceed through the ESS survey in steps 1 to 4 below. The question number is referenced on the bilevel CPAP LCD display.

If the ESS is enabled and due for completion, an ESS prompt will be shown on the LCD display when the blower is turned ON **Figure 11**. In addition, the LCD back light will turn on and several short beeps will be sounded from the audible indicator.

E	p	w	o	r	t	h	
S	t	a	r	t	>	>	

Figure 11

1. With the display as shown in **Figure 11**, press the SELECT button to proceed to ESS Q1. The display shown in **Figure 12** references Question number 1, “Sitting and reading.” See questions (situations) on the next page.

E	S	S		Q	1		
			*				

Figure 12

2. Use the UP/DOWN arrow buttons to display the most appropriate number as it relates to you for the situation.

3. Press the SELECT button to enter your answer and advance to the next question.

4. Repeat steps 2 and 3 until you have answered all eight (8) questions of the ESS. The total will be displayed once all of the questions are answered **Figure 13**. If “??” is displayed, then one or more of the questions has not been answered. **NOTE**—Record the date and total into the chart on page 2 of the ESS booklet. If instructed to do so by your homecare provider, complete the ESS Reply Card and mail it to the address printed on the front of the card.

E	p	w	o	r	t	h	
T	o	t	a	l	=	X	X

Figure 13

5. Press the SELECT button to scroll past the Epworth Total. An Epworth “Finished” prompt will be displayed **Figure 14**.

F	i	n	i	s	h	e	d
	Y	e	s				

Figure 14

6. Use the UP/DOWN buttons to display YES or NO. If it is necessary for you to change or review your answers, press the SELECT button with NO on the display and return to the questionnaire. If you are finished, select YES on the display, press the SELECT button and the display will return to the “Act Pres” (Actual Pressure) reading.

**NOTE**—The bilevel CPAP will operate throughout the ESS assessment. The blower ON/OFF and DELAY keypad buttons will continue to function regardless of the ESS prompts. If you do not answer the ESS, it will remain on the display for 15 minutes before timing out and being removed. Any keypad activity will reset the 15-minute time-out period. Turning the blower OFF and then back ON will temporarily skip the questionnaire; the unit must be off for two hours before ESS will be displayed again.

**Epworth Sleepiness Scale® (ESS)**

**The Epworth Sleepiness Scale® (ESS)**

How likely are you to doze off or fall asleep in the following situations, in contrast to feeling just tired? This refers to your usual way of life in recent times.

Even if you haven't done some of these things recently try to work out how they would have affected you. Use the following scale to choose the most appropriate number for each situation:

- 0 = would never doze
- 1 = slight chance of dozing
- 2 = moderate chance of dozing
- 3 = high chance of dozing

It is important that you answer each question as best you can.

SITUATION	CHANCE OF DOZING (0-3)
1. Sitting and reading	_____
2. Watching TV	_____
3. Sitting, inactive in a public place (e.g. a theater or a meeting)	_____
4. As a passenger in a car for an hour without a break	_____
5. Lying down to rest in the afternoon when circumstances permit	_____
6. Sitting and talking to someone	_____
7. Sitting quietly after a lunch without alcohol	_____
8. In a car, while stopped for a few minutes in the traffic	_____

THANK YOU FOR YOUR COOPERATION

© M. W. Johns 1997

Used with the permission of Dr. Murray W. Johns

**TRAVEL INFORMATION**

This unit is equipped with a universal power supply capable of accepting line voltages of 100-240V~ 50/60Hz, and 120V~ 400Hz. No adjustments are required to the AC power connector or the line fuses to accept this entire range of voltages. If traveling to another country where the line voltage is different than the domestic supply voltage, simply obtain a line cord configured for that particular country and substitute it for the line cord supplied with your bilevel CPAP. Operate unit as normal.

**NOTE**—Your home medical equipment provider has set your prescription pressure based on your prescription. The Model 9055 bilevel CPAP has onboard automatic altitude compensation. Therefore, no adjustments for altitude are required.

**DC Operation of the 9055**

The 9055 Series DeVilbiss bilevel CPAP can be operated from a 12V DC power source in a motor home, recreational vehicle or boat provided that the appropriate power inverter has been installed. Connect the DeVilbiss bilevel CPAP to the inverter. The appropriate power rating of the inverter should be at least 120 watts. There are no switches or fuses to change for DC operation. Do not use the bilevel CPAP in a moving vehicle.

## SUPPLEMENTAL OXYGEN

Supplemental oxygen may be prescribed with nasal CPAP for treatment of obstructive sleep apnea.

**CAUTION**—Oxygen is a prescription gas and should only be administered under the supervision of a physician.

### ▲ WARNING

Oxygen supports combustion. To avoid possible physical injury, do not smoke while using your bilevel CPAP unit with supplemental oxygen. Do not use near hot, sparking objects or open sources of flame.

### ▲ WARNING

Always turn the bilevel CPAP unit on prior to turning on the oxygen flow and turn off the oxygen before turning the unit off.

1. Press the ON/OFF button on the keypad to start the blower. Connect the optional oxygen adapter (part #7353D-601) to the outlet port of the unit. Always turn CPAP unit on before the oxygen source and turn the oxygen source off before turning off the bilevel CPAP unit.
2. Connect one end of the 22 mm diameter tubing to the open end of the adapter and the other end to the CPAP mask.
3. Oxygen may be injected into the air supply by connecting one end of standard oxygen tubing to the oxygen adapter and the other end to the oxygen source (**Figure 15**).
4. If the bilevel CPAP is to be left for an extended period of time, shut off the oxygen flow. Do not leave the oxygen source unattended.


Figure 15

### ▲ WARNING

The flow setting of the oxygen source shall be specified by a physician.

### ▲ WARNING

Equipment is not suitable in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

**NOTE**—Equipment is suitable for use with oxygen if prescribed by your physician.

**NOTE**—At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on the pressure setting, patient breathing pattern, mask selection and leak rate. This applies to most types of bilevel and CPAP machines.

## CLEANING

### ▲ WARNING

To avoid electric shock, always unplug power cord from wall outlet power source.

**CAUTION**—Never rinse or place the unit in water. Never allow liquids to get into or around any of the ports, switches, or air filter; doing so will result in product damage. If this occurs, do not use the unit. Contact your provider for service.

## Cabinet

1. Unplug the power cord from the wall outlet.
2. Wipe the cabinet with a clean, damp cloth and mild dishwashing detergent every few days to keep dust free. Do not use alcohol or any alcohol-based solutions to clean the unit. These cleaning solutions may damage the plastic enclosure. **Allow the unit to dry completely before plugging the unit into a power source.**

## Air-Inlet Filter

1. The air-inlet filter should be checked daily and cleaned as necessary. The filter is located in the rear of the unit (**Figure 16**).
2. To clean, remove the dirty foam filter from the frame.
3. Wash the filter in a solution of warm water and dishwashing detergent and rinse with tap water. Allow the filter to dry completely. If the filter (part #8000D-602) is damaged, contact your home medical equipment provider for a replacement.
4. Insert clean filter into housing.


Figure 16

## Mask, Headgear, and Tubing

Clean according to the manufacturer's instructions.

**RECOMMENDED MAINTENANCE INTERVAL FOR THE EQUIPMENT OWNER**

**▲ DANGER**

Electric shock hazard. Do not attempt to open or remove cabinet; there are no user-serviceable internal components. If service is required, return unit to a qualified Sunrise Medical provider or authorized service center. Opening or tampering with the unit will void warranty.

Pressure– Check outlet pressure once per year or as needed with a calibrated gauge, calibrate unit if the difference between the set and actual pressures is  $> \pm 0.5$  cmH<sub>2</sub>O

Standard Filter– Check filter daily, clean if necessary. It is recommended that the filter be replaced every 6 months or if it is damaged. The standard filter is designed to filter particles larger than 3.0 microns in size. The bilevel CPAP must have the standard filter in place during operation.

Optional Fine Particle Filter– Check filter daily, replace if dirty. Replace every 30 days or if damaged. The optional fine particle filter is designed to filter particles larger than 0.3 microns in size. The fine particle filter is placed in the filter area of the bilevel CPAP air intake with the standard filter placed over the fine particle filter.

**TROUBLESHOOTING**

**NOTE**–DeVilbiss 9055 Series contains no user-serviceable parts. If you believe your unit is not working properly, BEFORE YOU RETURN IT TO THE HOME MEDICAL EQUIPMENT PROVIDER WHERE YOU PURCHASED IT OR TO SUNRISE MEDICAL, please take a few moments to check for these possible causes:

BASE UNIT		
Symptoms	Possible Causes	Remedies
Unit does not turn on.	<ol style="list-style-type: none"> <li>Unit not connected to a live power source.</li> <li>Power cord not attached properly.</li> </ol>	<ol style="list-style-type: none"> <li>Plug unit into a live receptacle. Check if receptacle operated by a switch.</li> <li>Firmly attach power cord.</li> </ol>
Unit display shows E05 error.	<ol style="list-style-type: none"> <li>Unit operated with outlet and/or mask exhalation port blocked (no leak flow).</li> <li>Unit calibration incorrect.</li> </ol>	<ol style="list-style-type: none"> <li>Remove occlusion, then turn unit off, unplug it, plug it back in, and turn it on.</li> <li>Contact your Sunrise Medical provider or authorized service center.</li> </ol>
Unit display shows any other error code.	<ol style="list-style-type: none"> <li>Unit malfunction.</li> </ol>	<ol style="list-style-type: none"> <li>Contact your Sunrise Medical provider or authorized service center.</li> </ol>
No air flow emitted from flow generator when power is turned on, motor runs.	<ol style="list-style-type: none"> <li>Air filter may be occluded by dirt or bed clothing.</li> <li>Mask clogged.</li> </ol>	<ol style="list-style-type: none"> <li>1a. Clean or replace the air filter.</li> <li>1b. Move unit away from loose material.</li> <li>2. Check mask for obstructions.</li> </ol>
Symptoms of OSA have recurred.	<ol style="list-style-type: none"> <li>Air filter may be occluded.</li> <li>Incorrect air pressure setting.</li> <li>Sleep apnea condition has changed.</li> </ol>	<ol style="list-style-type: none"> <li>Clean or replace the air filter and relocate unit away from drapes, etc.</li> <li>Contact your Sunrise Medical provider or authorized service center.</li> <li>Call your physician; you may need to be retested.</li> </ol>
Skin becomes irritated where mask touches skin.	<ol style="list-style-type: none"> <li>Headgear is too tight.</li> <li>Incorrect mask size.</li> <li>Allergic reaction to the mask material.</li> <li>Worn mask.</li> </ol>	<ol style="list-style-type: none"> <li>Loosen the headgear.</li> <li>Contact your Sunrise Medical provider.</li> <li>Call your physician, or contact your Sunrise Medical provider.</li> <li>Replace mask.</li> </ol>
Dryness of throat or nose.	<ol style="list-style-type: none"> <li>Inadequate room humidity.</li> <li>Room temperature too high.</li> </ol>	<ol style="list-style-type: none"> <li>Increase room humidity or add a DeVilbiss CPAP humidifier; see page 7.</li> <li>Decrease room temperature.</li> </ol>
Air from flow generator seems warm.	<ol style="list-style-type: none"> <li>Air filter may be occluded.</li> <li>Air filter and/or outlet port may be obstructed.</li> <li>Unit located too close to a heat register or radiator.</li> </ol>	<ol style="list-style-type: none"> <li>Clean or replace the air filter.</li> <li>Move the unit away from obstruction.</li> <li>Move the unit away from heat source.</li> </ol>
Nasal, sinus, ear pain, or runny nose.	<ol style="list-style-type: none"> <li>Reaction to air flow pressure, sinus infection, or nasal congestion.</li> </ol>	<ol style="list-style-type: none"> <li>Discontinue use. Contact your physician.</li> </ol>

**SPECIFICATIONS**

Size .....	4.0"H x 7.5"W x 10.8"D (10.2 cm x 19.6 cm x 26.2 cm)
Weight.....	3.8 lbs. (1.7 kg)
Electrical Requirements.....	100-240V~, 50/60 Hz or 120V~, 400 Hz
Maximum Power Consumption.....	65 watts max from AC power source
Bilevel CPAP Pressure Range .....	3-25 cmH <sub>2</sub> O
Pressure Delay Time .....	0 to 45 minutes, in 5 minute increments
Delay Pressure Range.....	3 cmH <sub>2</sub> O to EPAP or CPAP prescription pressure
Operating Temperature Range.....	41°F to 93°F (5°C to 34°C)
Operating Humidity Range .....	0% to 95% RH non-condensing
Operating Atmospheric Conditions .....	Sea level to 9,000 feet
Storage & Transportation Temperature Range. ....	-40°F to 158°F (-40°C to +70°C)
Storage & Transportation Humidity Range.....	0% to 95% RH non-condensing
Maximum Limited Pressure.....	30 cmH <sub>2</sub> O under normal use
Sound Level .....	< 30 dBA
Filter specs.....	Standard Filter > 3.0 micron particles Optional Fine Particle Filter > 0.3 micron particles

## Max flow rates

83 L/m .....	@ 8 cmH <sub>2</sub> O (1/3 max pressure)
134 L/m.....	@ 17 cmH <sub>2</sub> O (2/3 max pressure)
163 L/m.....	@ 25 cmH <sub>2</sub> O (max pressure)
Pressure Accuracy.....	± 0.5 cmH <sub>2</sub> O

Pressure Swings (cmH<sub>2</sub>O peak to peak at 500mL tidal volume)

Breaths per minute	<u>10</u>	<u>15</u>	<u>20</u>
@ 8 cmH <sub>2</sub> O	0.5	0.5	0.5
@ 17 cmH <sub>2</sub> O	1.0	1.0	1.0
@ 25 cmH <sub>2</sub> O	1.0	1.0	1.0

**NOTE**–For results of extreme testing of static and dynamic pressure regulation refer to the Sunrise Medical web site ([www.sunrisemedical.com](http://www.sunrisemedical.com)).

Warranty.....	Two-year limited
Equipment classification with respect to protection from electric shock.....	Class II <input type="checkbox"/>
Degree of protection from electric shock.....	Type BF Applied Part
Degree of protection against ingress of liquids.....	None
Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.	
Mode of operation .....	Continuous

**Clinical Specifications for Digital Outputs**

Mask Pressure.....	0 to 25.5 cmH <sub>2</sub> O ± 0.5
Patient Flow.....	-127 to +127 L/min ± 10 L/min
Leak Flow.....	0 to +127 L/min
Estimated Tidal Volume .....	0 to 1023 mL

**NOTE**–Leak and Tidal Volume are estimates. They are provided for trending purposes only.

**WARRANTY**

The 9055 series DeVilbiss RPM bilevel CPAP system is warranted to be free from defective workmanship and material for a period of two years from date of purchase. Any defective part(s) will be repaired or replaced at Sunrise Medical's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your Sunrise Medical provider or Sunrise Medical Service Department for instructions:

U.S.A. 800-333-4000, 814-443-4881

Canada 905-660-2459

Europe 44-138-444-6688.

**NOTE**—This warranty does not cover providing a loaner unit, compensating for costs incurred in rental while said unit is under repair, or costs for labor incurred in repairing or replacing defective part(s).

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

**DECLARATION OF CONFORMITY**

**Manufacturer:** Sunrise Medical  
**Address:** Respiratory Products Division  
100 DeVilbiss Drive  
Somerset, Pennsylvania 15501-2125 USA  
**Product Designation:** Bilevel CPAP  
**Type, Model:** DeVilbiss® 9055 series

We herewith declare that the above-mentioned product complies with the requirements of EC Directive 93/42/EEC and the following:

**Quality System Standards Applied:** ISO9001/EN46001/ISO13485

**Notified Body RWTÜV**

**MDD Annex II Applied**

**Safety Standards Applied:** IEC 601-1/EN60601-1

**EMC Compliance to:**  
EMC Directive 89/336/EEC  
IEC 801-2 through 5  
CISPR 11/EN55011

**European Contact:**  
Sunrise Medical Ltd.  
Sunrise Business Park  
High Street  
Wollaston, West Midlands DY8 4PS  
ENGLAND  
44-138-444-6688

®Registered U.S. Patent and Trademark Office and other countries.

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®Déposé auprès du U.S. Patent and Trademark Office et dans d'autres pays.
®Eingetragen beim U.S. Patent and Trademark Office sowie in anderen Ländern.
®Gedeponeerd bij U. S. Patent and Trademark Office en in andere landen.
**Sunrise Medical Canada, Inc.**
**Sunrise Medical Pty. Limited**