USER GUIDE • MANUEL D'UTILISATION • ISTRUZIONI D'USO • GEBRAUCHSANLEITUNG • GEBRUIKSAANWIJZING
BRUKSANVISNING • MANUAL DE USO • INSTRUÇÕES DE USO • INSTRUKCJA OBSŁUGI • ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

Respireo Primo N

NASAL MASK • MASQUE NASAL • MASCHERA NASALE • NASENMASKE • NEUSMASKER
NESEMASKE • MASCARILLA NASAL • MÁSCARA NASAL • MASKA NOSOWA • ΠΙΝΙΚΗ ΜΑΣΚΑ


- ENGLISH
- FRANÇAIS
- ITALIANO
- DEUTSCH
- NEDERLANDS
- NORSK
- ESPAÑOL
- PORTUGUÊS
- POLSKI
- ΕΛΛΗΝΙΚΑ


A


B

1


2

3


4


B

5.1


5.2


B

6

7

Respireo Primo N

NASAL MASK

Vented

for use by a single patient

Instructions for the patient

1 - APPLICATION

Respireo Primo N is a nasal mask with calibrated vent ports for use by a single patient. Respireo Primo N is designed for use by patients who have been prescribed continuous positive airway pressure (CPAP) or Bi-level positive airway pressure (BiPAP) treatment at home or a hospital. Respireo Primo N is suitable for adult patients (>30 kg).

2 - IMPORTANT NOTES

- 2.1 The mask must only be used with CPAP or Bi-level flow generators recommended by a doctor or respiratory therapist.
- 2.2 The operating pressure range is 4-24 cmH₂O.
- 2.3 The mask must only be used with the appliance on and operating at a pressure of at least 4 cmH₂O. With low pressures using a CPAP appliance, the flow through the vent ports may be insufficient to eliminate all the gas exhaled from the tube, which could lead to partial re-inhalation.
- 2.4 The vents in the mask must never be obstructed.
- 2.5 If you experience discomfort or irritation, or suffer an allergic reaction to any of the components of the mask, please report this to your doctor or respiratory therapist.
- 2.6 If oxygen is used, the supply must be turned off when the appliance is not in operation, to prevent oxygen accumulating inside, which is a fire risk.
- 2.7 Do not smoke or light naked flames when using oxygen.
- 2.8 Do not use the mask in the presence of nausea or vomiting.
- 2.9 The mask must be cleaned and/or disinfected. For further details refer to the section entitled "Cleaning and disinfecting the mask".
- 2.10 Check the mask before use. If there is any visible damage, report it to the dealer.
- 2.11 The mask must only be used in accordance with these instructions. No other instructions are applicable.
- 2.12 Some of the components of the mask are small enough to be swallowed by children and must never be left unattended. An adult must always be present.

3 - COMPONENTS OF THE MASK (Fig. A)

- a) cushion; b) main body; c) forehead support; d) adjusting knob; e) adjusting screw; f) cap; g) curved fitting; h) swivel fitting; i) headgear.

4 - HOW TO PUT ON THE MASK (Fig. B)

The mask comes fully assembled ready for use. Proceed as follows to put on the mask.

Turn the adjusting knob (d) until the forehead support (c) is fully extended (Fig. 1).

Hold the mask firmly in one hand and extend the headgear (i) with the other enough to pull it over your head (Fig. 2).

Place the mask gently on your face so that the cushion (a) fits snugly round the nose and position the headgear (i) on your head (Fig. 3).

Connect the two lower ends of the headgear (i) to the main body (b) (Fig. 4).

Tighten the upper straps first, then the lower straps, until the mask fits on your face comfortably, without overtightening (Fig. 5).

Connect the mask to a CPAP or a Bi-level appliance at the prescribed pressure (Fig. 6)

After lying down, turn the adjusting knob (d) until the mask feels comfortable but still provides a good seal (Fig. 7).

Readjust the lower straps if necessary to eliminate any leaks, but make sure the mask is not too tight.

5 - HOW TO REMOVE THE MASK

Release one of the lower straps. Remove the mask from your face and lift the headgear over your head. The other straps do not need to be loosened or disconnected.

6 - HOW TO DISASSEMBLE THE MASK

Before disassembling the mask, disconnect the patient circuit from the CPAP or Bi-level appliance.

Make a note of the position of the adjusting knob (d) for reference next time.

Undo the four straps and detach the headgear (i) from the main body (b).

Pull the curved fitting (g) off the main body (b).

Pull the swivel fitting (h) off the curved fitting (g).

Detach the cushion (a) from the main body (b).

Remove the cap (f) from the main body (b).

7 - CLEANING AND DISINFECTING THE MASK

The mask must be cleaned and/or disinfected each time before use.

N.B. The headgear should be cleaned once a week.

CLEANING THE MASK

First wash your hands thoroughly.

Disassemble all the components (see 6 - HOW TO DISASSEMBLE THE MASK) and wash carefully in water at about 30°C using mild soap.

Rinse thoroughly with water and leave in the air to dry.

The mask can be cleaned in a washing machine at 30°C (programme for delicate fabrics) but must not be spin- or tumble-dried. The mask can be washed in this way up to 20 times.

DISINFECTING THE MASK

The mask must be taken apart for disinfection (see 6 HOW TO DISASSEMBLE THE MASK). Immerse the components in a cold solution of disinfectant (e.g. 5.25% sodium hypochlorite solution), following the manufacturer's instructions.

IMPORTANT. Check that the disinfectant solution is compatible with the materials of which the mask is made (see 12 DISPOSAL INSTRUCTIONS).

N.B. The headgear (i) must not be disinfected.

8 - PRECAUTIONS

Do not clean the mask with solutions containing bleach, chlorine, alcohol, aromatic compounds, humidifiers, antibacterial agents or aromatic oils. Such solutions would damage the product and reduce its life.

Prolonged exposure to direct sunlight or heat may cause damage to the mask.

If any visible damage is noted (cracks, tears, etc.), the damaged part must be replaced.

Check carefully to make sure there are no soap deposits or foreign bodies inside the mask.

Check regularly that the vents are not obstructed. Clean with a soft brush if necessary.

9 - HOW TO REASSEMBLE THE MASK

Insert the curved fitting (g) in the main body (b), pressing gently.

Insert the swivel fitting (h) in the curved fitting (g), pressing gently.

Fit the grooved side of the cushion (a) in the frame of the main body (b). Insert the top, then the sides of the mask, pressing the cushion until the whole edge fits snugly in the main body.

Fit on the cap (f).

Attach the four straps of the headgear (i) to the main body (b).

Turn the adjusting knob (d) to the reference position you noted previously.

10 - PRESSURE/OXYGEN PORT

If additional oxygen or a pressure reading is required, use the port in the main body (b).

Open the cap (f) by lifting upwards, without removing it.

Push the pressure/oxygen tube firmly into the port.

11 - TROUBLE-SHOOTING CHART

PROBLEM	POSSIBLE CAUSES	SOLUTIONS
<ul style="list-style-type: none">No air or insufficient air comes out of the mask.	<ul style="list-style-type: none">Flow generator off.Flow generator regulated wrongly.Patient circuit tube not inserted.	<ul style="list-style-type: none">Switch on the flow generator.Regulate the flow generator.Connect the tube to the flow generator and mask.
<ul style="list-style-type: none">The mask is noisy.	<ul style="list-style-type: none">Cushion mounted incorrectly and air leak.	<ul style="list-style-type: none">Mount the cushion correctly.
<ul style="list-style-type: none">The mask does not fit the face snugly.	<ul style="list-style-type: none">Mask applied incorrectly.	<ul style="list-style-type: none">Reposition the headgear.
<ul style="list-style-type: none">Exhaled air does not flow out and difficulty inhaling.	<ul style="list-style-type: none">Vent ports obstructed.	<ul style="list-style-type: none">Remove the obstruction.

12 - DISPOSAL INSTRUCTIONS

The components of the mask are made of the following materials and must be disposed of in accordance with the applicable laws.

COMPONENT	MATERIAL
Main body	Polycarbonate
Forehead support	Polycarbonate
Adjusting knob	Polypropylene
Adjusting screw	Nylon
Curved fitting	Polycarbonate
Swivel fitting	Polycarbonate
Forehead rest	Silicone
Cushion	Silicone
Cap	Silicone
Headgear	Nylon/polyurethane

13 - TECHNICAL DATA

MINIMUM CONTINUOUS FLOW

Mask Pressure (cmH ₂ O)	4	8	12	16	20	24
Flow Rate (lpm) (± 3 lpm)	18	27	34	40	45	49

N.B. Due to manufacturing differences, the exhalation flow rate may differ from the values shown in the table.

Normative references	CE marking in accordance with European Directive 93/42/EEC, class IIa EN ISO 14971, ISO 5356-1, ISO 10993-1, EN ISO 17510-2
Resistance	at 50 lpm: 0.1 cmH ₂ O - at 100 lpm: 0.5 cmH ₂ O
Operating pressure	4-24 cm H ₂ O
Temperature range	5°-40°C
Dead space volume inside the mask	114 ml

WARRANTY

Refer to the general warranty conditions. The Manufacturer is liable for the safety, reliability and performance of the apparatus, provided that it is used in accordance with the instructions and for the intended purpose only, and that any repairs are made by the Manufacturer or an approved service centre.

The Manufacturer declines all liability and the warranty will be invalidated if the user fails to follow the instructions and the important notes contained herein.

COPYRIGHT

The information contained in this booklet may not be used for anything other than the intended purpose. This booklet is the property of Air Liquide Medical Systems S.p.A and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

TECHNICAL UPDATES

Air Liquide Medical Systems periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated to ensure that they continue to be consistent with the characteristics of the devices marketed.

If the booklet accompanying the apparatus is mislaid, a replacement copy can be obtained from the Manufacturer by quoting the data shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.p.A.