**SPIRO**

**Instructions for use** en

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**Bruksanvisning** sv

**Brugsanvisning** da

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**RÜSCH** CARE

en

## **Speech valve and humidifier**

**SPIRO Speech valve,  
Ref. 85 89 11, 85 89 13,  
85 89 14**

Daily check of function and possible damage. When damaged or not functioning, replace the valve. Never take out the membrane. Cleaning of humidifier housing daily.

Sharp instruments must not be used. Use your fingers or plastic tweezers. Humidifier housing and cap are cleaned with water and air-dried.

Detergents, solvents or other fluids must not be used.

Always check function after cleaning. The air-conditioning medium needs to be changed at least twice daily, if necessary more often. It is situated at the bottom of the humidifier housing, parallel to the phonation valve.

Air preheater is used when needed. The opening for inhalation is turned downwards.

The nipple cover should be fixed to the speech valve and only be taken away when oxygen supply is attached.

The speech valve should be replaced after approximately two weeks, more often if needed.

All products are to be stored in a dark, dry place.

The medium for air-conditioning may show a change of colour from white to light yellow. This is normal.

**Instructions and warnings relating to the use of speech valves in combination with tracheostomy tubes**

A speech valve may only be used in connection with a tracheostomy tube if the upper airways are unobstructed.

If a speech valve is used, tracheostomy tubes with a phonation window should be given preference.

If a speech valve is used together with a tracheostomy tubes with cuff, the cuff must be completely deflated!

Before each use, the speech valve must be checked for correct functioning, and cleaned or replaced, as applicable.

The speech valve may be only used when the patient is conscious. The patients have to be trained to use the speech valve by qualified staff.

In the case of an incipient dyspnoea, the speech valve must be removed immediately from the tracheostomy tube and appropriate measures must be taken at once to correct the dyspnoea.

**Humidifier housing  
(with speech valve)**

**Ref. 85 89 11**

**Humidifier housing  
(without speech valve)**

**Ref. 85 89 13**

**Medium for air-conditioning,  
plastic bag of 100**

**Ref. 85 89 14**

**Complete Set**

**Ref. 12 18 00**

**Deflector**

**Ref. 12 18 04**

**Note:**

The speech valve must never  
be used without medium for  
air-conditioning!!

Do not use the speech valve  
while sleeping.

Then use SPIRO humidifier  
without phonation valve.

### **Presentation and Storage:**

The products Spiro humidifier housing with speech valve and Spiro humidifier housing, Ref. 858911 and 858913, can be reprocessed for single patient use.

The filter medium, Ref. 858914, can not be reprocessed and is intended for single use only.

Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only medical devices may result in exposure to viral, bacterial, fungal or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these medical devices. This product is not designed to be cleaned, disinfected or sterilized.

The products are supplied non sterile, packaged in individual pouches.





Do not use the product when expiry date is reached.

Discharge the product!

This product should be stored in a dry place, preferably away from direct and indirect sources of light and heat.

Store in a cool and dry place and away from light and ozone (no UV-light)! Avoid high temperatures during storage.

**Attention-Phonation valve**

Do not use without conditioning medium.  
Do not use when tracheostomy cannula is blocked.  
Do not use when patient's unconscious or asleep.  
The patient's physician is responsible for correct application.

**Assembly**

- 1 Tube
- 2 Connector to silver tube
- 3 Humidifier housing
- 4 Medium for air conditioning
- 5 Cap with or without speech valve
- 6 Air preheater

**Change of the air-conditioning medium**  
**Cleaning****Air preheating**

Inhaled air is preheated under the clothes

**Extra humidifier**

The condensate from the exhaled air is reused

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## SWITZERLAND TELEFLEX MEDICAL AG
Willy Rüsch GmbH,
Henkel GmbH, Druckerei, 70499 Stuttgart