RÜSCH / RÜSCH CARE

TracheoFlex Tracheostomy Tubes
Instructions for Use

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tr

Instructions for Use
TracheoFlex Tracheostomy Tubes

en

Note:

The enclosed instructions for use are valid for the following products:
TracheoFlex tracheostomy tubes with adjustable fixation flange with or without cuff.

Warning

DEHP Contains DEHP, DBP

DBP The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction.

Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical devices containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Set components:

Component	With Cuff	Without Cuff
Fixed standard connector	X	X
Tightly fitting high-pressure cuff, valve for Luer and Luer-Lock syringe tip	X	
Mounted adjustable fixation flange	X	X
Velcro strap to fix the fixation flange to the neck	X	X
Insertion aid	X	X

This tracheostomy tube set is a sterile, reusable product
for single-patient use only.

These instructions for use are primarily directed at the tracheostomised patient himself as well as at the nursing staff. They should guide them in how to safely use the tube set, clean it and above all correctly resterilise it. For this reason, the instructions for use on hand should always be read by the user, the nursing staff as well as the patient himself.

Caution

This product is delivered in a sterile condition by the manufacturer. However, sterility is only guaranteed if the packaging is intact. Do not expose to temperatures above 50 °C.

Indication

Cannulation of tracheotomised patients, particularly in the case of a narrow tracheostoma and narrow-lumen tracheas. When using cannulas with cuff, for short-term intermittent inflation.

Diagnosis only by doctor!

Contraindications

- Patients in whom the cuff needs to be inflated for a longer period of time.
- For patients during radiation therapy and magnetic resonance imaging.

Caution

- During surgery on the trachea, use of electrocautery may cause the tube to burn.
- Protect the tracheostomy cannula from direct laser bombardment (risk of tube burning)!
- **Caution:**
In the event of Radiation Therapy where the target organ (cancer) is adjacent to the location of the cuff on the tube, the cuff material can be compromised by Radiation and can thus lose its resilient properties.

In this instance, a visual inspection and functional performance testing should be conducted to verify that the integrity of the cuff has not been compromised **during radiation therapy**.

As the duration and dosage of this therapy can influence the material capabilities of Poly Vinyl Chloride (PVC), and can therefore increase the risk of cuff damage, we recommend to check the tube and its function thoroughly after each applied dose fraction, and to exchange the tube in the event that performance is affected and / or damage to the material is observed.

Description

- **Shaft**

This product consists of a highly flexible, armoured tracheostomy cannula made of polyvinylchloride (medical grade PVC).

- **With cuff / without cuff**

When using cannulas with cuff a high filling pressure must be used to inflate the cuff as it fits very tightly. This pressure does not give any real information about the actual pressure of the cuff on the tracheal mucosa. The inflation status of the cuff is indicated by the pilot balloon.

- **Fixation Flange**

The fixation flange may be moved along the tube shaft and can thus be individually adapted to the patient's individual anatomy and can be secured to the tube in optimum position.

Warning!

The graduation on the tube shaft can be used for orientation. It is, however, not an indication of measurements!

The neck band (velcro band) ensures patient-friendly, gentle fixation of the flange to the patient's neck and can be adapted to the thickness of the neck.

Caution:

The tube neckband can not be cleaned and reused (single use only)!

Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only medical devices may result in exposure to viral, bacterial, fungal or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these medical devices. This product is not designed to be cleaned, disinfected or sterilized.

Delivery form and storage recommendations

This tracheostomy tube set is a sterile reusable product **to be used with one patient only**.

Each set is supplied in an individual sterile pouch.

For the sterilisation method used see product label.

Sterility of the set and its components can only be guaranteed as long as the packaging has not been opened or damaged.

Do not resterilise this product and do not use beyond its expiry date, but dispose!

This product should be stored in dry surroundings, preferably away from direct or indirect sources of heat and light.

Store in cool surroundings, protected from light and ozone (if applicable, no ultraviolet light)! Avoid high temperatures.

A) Suggested Directions for Use

1. Selection of the correct tube size

The correct choice of tube size depends on the size and shape of the individual tracheostoma, the diameter of the trachea and the patient's underlying disease and is primarily the responsibility of the attending physician. In general, the largest possible tube that will fit both the trachea and the tracheostoma on the basis of the patient's anatomy should be selected to ensure optimum ventilation of the lungs and, at the same time, the lowest possible airway resistance.

2. Preparation of the tracheostomy tube

- Prior to use, check that the set is complete and functions correctly.
- When using cannulas with cuff, the patency of the tube lumen and cuff filling system should be checked before intubation.

Prior to use, check that the cuff, the pilot balloon and the valve are intact and have no leaks.

In order to check that the system is leaktight, fill the cuff with air using a syringe (see table for inflation volume of leakage test), then remove the syringe and leave for approx. 5 minutes. Leaktightness is demonstrated if there is no pressure or volume loss. In order to speed up this test procedure, the entire tube can also be temporarily immersed in sterile water. The absence of ascending bubbles will definitely indicate leaktightness.

Caution!

To prevent the cuff membrane from overstretching (herniation), the maximum filling volume should not be exceeded (see table).

In order to prevent damage to the cuff, avoid using pointed objects or cutting instruments during the leakage test.

After this procedure, completely evacuate the air from the cuff.

- Please also ensure that the fixation flange may be adjusted on the shaft and that secure fixation to the tube shaft can be achieved by arresting the fixaton screw. The area of the tube shaft where the fixation flange is to be secured must not be wetted with lubricants, grease or oil. If necessary, clean briefly with ethyl alcohol.
- Insert the introducer into the tube until it fits tightly in the tube and the olive tip projects beyond the end of the tube.

Caution

Transport damage or incorrect storage may cause damage to the cuff and other tube components. Never use damaged tubes and leaky cuffs!

3. Inserting the tracheostomy tube

- Intubate the patient following currently accepted medical techniques noting the specific "Warnings and precautions" in these instructions for use.
- The introducer must be held in fixed position throughout the entire procedure. After introducing the tube into the trachea and positioning it correctly, the introducer should be removed. When removing the introducer, the tracheostomy tube must be held in position.

Note

The intubation process can be facilitated by prior wetting of the insertion aid and tube tips with a small amount of a water-soluble lubricant.

Caution

When using cannulas with cuff please ensure that the cuff is not inadvertently damaged by sharp-edged objects during intubation.

To avoid subsequent injury to cartilage, the tube should be in the tracheal opening without tension.

Too deep an intubation should be avoided as otherwise there is the danger of one-sided ventilation of the lung!


4. Adapting and securing the fixation flange to the tube shaft

The information given below with regard to the fixation screw is only valid from tube size I.D. 7 mm for the new fixation flange system. Smaller sizes are supplied with threaded ring.


Without changing the tracheal position of the tube, the fixation flange can now be adjusted in the longitudinal axis to adapt it to the anatomical conditions of the neck surface.

To secure the flange, the fixation screw must be tightened clockwise by approx. 90° (a quarter of a complete turn) until the fixation flange is tightly connected to the shaft.

After correct fixation the marking on the fixation screw points upwards. (see ill. 1)

ill. 1

ill. 2

Warning!

The graduation on the tube shaft can be used for orientation. It is, however, not an indication of measurements!

Warning

The intubation depth can be set according to individual requirements via the fixation screw on the adjustable neck plate. Ideally, this is done after previous fibroscopic inspection of the trachea. The fixation screw needs to be screwed tightly to safely prevent the tube from sliding out of the neck plate.

Caution (cuff related!)

When using cannulas with cuff, incorrect securing of the fixation screw can lead to the air supply channel being squeezed together somewhere on the shaft of the tube, thus rendering inflation or deflation of the cuff impossible.

In this case the fixation screw must be loosened and the position of the fixation flange corrected. After correct fixation the marking on the fixation screw points upwards. (see ill. 2)

Caution

If the fixation flange cannot be secured due to contamination by lubricants, it is recommended to briefly clean the shaft.

Caution

Do not immerse in alcohol!

5. Fixing the fixation flange to the patient's neck

After fixing the fixation flange on the tube shaft, it must be secured to the patient's neck with the velcro band in order to prevent the tube from slipping out. To achieve this, thread the ends of the neck band through the slots in the fixation flange and secure at the sides with the velcro fasteners.

To avoid irritation of the skin in the area of the fixation flange, especially in the case of a wet tracheostoma, it is obligatory to place a commercially available slitted compress between the skin and the fixation flange.

When using cannulas with cuff:

For safety reasons, the tracheostomy tube should be fixed before the cuff is inflated, if it is introduced by the patient himself or by auxiliary medical staff.

Caution:

- Avoid exerting high tensile forces on the trachea through the ventilation tubing and the connector on the tube shaft.
Warning! Danger of extubation!
- Avoid all alcoholic disinfectants when looking after the tracheostoma.
- Never wet the slitted compress with disinfectants containing alcohol!

6. Inflating the cuff (Cuff related!)

The cuff of this tracheostomy tube set may only be inflated intermittently or for a short time (e.g. during eating) in order to avoid damage to the tracheal mucosa. To inflate the cuff, use a syringe to inject air into the filling system through the valve. The amount of air should be just sufficient to seal the trachea (minimum inflation volume). If ventilation gases still escape between the cuff and tracheal wall during active ventilation, the minimum inflation volume of the cuff has not yet been achieved. In this case the user or patient helper can hear a hissing noise using a stethoscope or with the ear at the patient's mouth or stoma.

Air is then carefully added to the filling system until the hissing noise can no longer be heard during the inhalation phase.

It should be noted that the cuff should not be filled with more than the minimum inflation volume of air. Otherwise the proportion of cuff pressure which is passed on to the mucosa is too high, mucosal circulation at the contact point with the cuff can be impaired and major sequelae can be expected in the long term.

Caution:

If the trachea cannot be sealed off despite reaching the maximum filling volume, the entire filling volume should be removed from the cuff and a renewed filling and inflation attempt should be made. If this is again not successful, the next largest tracheostomy cannula should be selected.

Warning:

- Anaesthetic gases administered via the tracheostomy tube may lead to gas diffusion with a pressure increase in the cuff system. As the actual pressure on the mucosa exerted by the inflated cuff cannot be established by direct pressure monitoring, the cuff should only be inflated intermittently or for a short time. (The intervals between inflations should be determined by the physician.)
- At all costs, avoid closing the machine end of the tube with a phonation valve or sealing plug or similar when the cuff is inflated. Even with an empty cuff, connecting a phonation valve may sometimes lead to considerable difficulty in respiration (dyspnoea). Closing the tube (with an empty cuff), e.g. with a plug to wean the patient off the tube, may only be performed under medical supervision and responsibility.

If the desired inflation has been achieved, the syringe should be removed from the valve.

7. Deflating the air (Cuff related!)

Firstly, the tracheal section above the cuff should be cleaned by draining accumulated secretion and mucus. If this is not possible in a conscious patient with maintained reflexes, it is recommended to advance a suction catheter into the trachea through the tracheostomy tube and, while suctioning, to completely remove the air out of the cuff using a syringe. Thus escaping secretion can be effortlessly suctioned up and does not cause coughing. When the air has been completely removed from the cuff, the cannula can be safely removed through the stoma.

B) Cleaning and disinfecting the tube set

If the tube is heavily soiled with encrusted secretion, blood or other body fluids which cannot be removed completely when the cannula is in place, the tracheostomy tube should be removed and cleaned outside the body. Simple cleaning (see below) should always be carried out, followed by disinfection if required. These tube sets must not be resterilized.

1. Simple cleaning

The lumen of the removed tube may be cleaned with a soft brush and neutral soap solution and then rinsed out with clear, hand-hot water.

Caution:

Temperatures should not exceed 50 °C!

2. Disinfection and cleaning

When using disinfectants (e.g. disinfectants based on active oxygen) ensure that their exposure times and concentrations are exactly observed. The disinfectant solution must be freshly prepared before each new application.

Thorough rinsing is indispensable to remove any possible residues!

Important!

We recommend the usage of disinfectants based on active oxygen (e.g. Sekusept classic, Sekusept aktiv). Cleaning with these agents has been tested and approved by us. This way you make sure that the tube set is not damaged through the use of disinfectants. Unsuitable disinfectants (especially those containing alcohol) can damage the tube and lead to dangerous malfunction.

Basically all disinfectants can penetrate the tube material and lead to irritation of the mucous membrane later when used again.

Thorough rinsing to remove the residues is therefore essential!

Warning

- When using cannulas with cuff: The use of brushes in the cuff region should be avoided to prevent injury to the cuff.
- Cleaning in an alcoholic bath and in solvents must be avoided at all costs due to the danger of material hardening and shrinking.
- The prescribed exposure times and concentrations of disinfecting solutions must not be exceeded (please observe manufacturer's instructions).
- Only use freshly prepared disinfecting solutions!

3. Resterilisation

Resterilisation is not allowed!

4. Frequency of cleaning

The cannula and its components should be cleaned as required. The underlying disease, mucus production, state of the tracheostoma and

the subjective hygiene feeling of the patient greatly determine how often cleaning is performed. As a rule of thumb we recommend cleaning the tracheostomy cannula at least twice a week, disinfecting it if necessary. Carefully dry the individual components after preparing the tube set.

C) Storage until reuse

Long-term tracheotomized patients often have several complete tracheostomy sets available for personal use. Thus the cleaned and dried tube sets must be stored in dry, closed containers which are protected from light and ozone. The containers should prevent soiling and contamination of the tubes and thus should be stored in a cool and dry place.

Caution (cuff related!)

The cuff must not be inflated during storage.

D) Safety testing before reuse

The function of the tube set must be tested before every renewed use in order to ensure that the tube set has not been damaged due to being reused too many times or being prepared incorrectly for reuse.

The following must be tested in particular:

- Kinking
- Working connection system at the tube
- Intact edges at the tip of the tube (no sharp edges or chipped material)
- Intact cuff and filling system (see preparation of cannula set)
- Intact fixation flange and tube shaft
- The inner lumen must be patent (no encrustation stuck to the walls)

Warning

On no account should you continue to use set components

if you observe the following:

- Damage, especially to the outer or inner surface of the shaft and cuff system (e.g. surface cracks, leaking cuff, parts of the armouring penetrating the tube wall, etc.).
- Brittleness or hardening of cuff or other set components.
- Separation of the standard connector from shaft.
- Damage of the fixation screw and the neck plate.

Safe use can no longer be expected under these circumstances!

Table: Important parameters of the TracheoFlex sets according to DIN EN ISO 5366-1

Valid for all sizes:

Residual cuff volume (ml air without excess pressure) = 0,

Cuff inflation volumen at 20 mbar excess pressure = 0 and angle $\Phi = 90^\circ$

Nominal cannula size (ID)	O.D. max. approx. (mm)	Leakage test; see inflation volume with air (ml)	Maximum inflation volume with air (ml)	Middle line length (CLL) approx. (mm)	Size A approx. (mm)	Size B approx. (mm)	Size C approx. (mm)
7.0	10.3	8	10	52 – 62	35	0 – 10	17
8.0	11.3	8	10	60 – 85	40	0 – 25	20
9.0	12.4	10	15	72 – 112	45	0 – 40	27
10.0	13.5	10	15	76 – 111	50	0 – 35	26
11.0	14.8	10	15	80 – 110	45	0 – 30	35
12.0	15.9	10	15	80 – 110	45	0 – 30	35

The inflation volumes given in the table were determined during experiments and thus should only be used as a guideline. The actual filling volume required in each case depends on the anatomical situation (see also par. A. 6, "Inflating the cuff").

Exceeding the above-mentioned maximum filling volumes may lead to non-reversible overstretching (herniation) of the cuff (in extreme cases, the cuff will burst!).

All information is in accordance with the state of the art when going to press. Subject to technical alterations.

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SPAIN TELEFLEX MEDICAL IBERIA S.A.
SWITZERLAND TELEFLEX MEDICAL AG
**Willy Rüscher GmbH**
Druckerei Täuber GmbH, Waiblingen