RÜSCH/RÜSCH CARE TRACHEOFLEX S

Tracheostomy tubes and accessories en

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Important : Please, read the entire instructions prior to use.

Note:

These instructions for use apply to the following products:

- TracheoFlex S Tracheostomietuben tracheostomy tubes with or without low-pressure cuff,
- Inner cannula for TracheoFlex S tracheostomy tubes,
- Accessories for TracheoFlex S tracheostomy tubes such as cough cap, phonation valve, sealing cap, connectors, neckbands and introducers.

Indications

Cannulation of tracheostomised patients through an existing tracheostoma.
Diagnosis only by a doctor!

Contraindications

For patients during radiation therapy and magnetic resonance imaging (MRI).

Warning

 Contains DEHP, DBP
The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical devices containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Form of supply and storage:

This tracheostomy tube set is a sterile reusable product and **only intended for single patient use.**

Each set is supplied in individual sterile packaging.

Sterilisation procedure see product label!

The sterility of the set and its components can only be guaranteed as long as the packaging is undamaged and has not been opened.

Do not resterilise the set! Do not use the set later than the expiry date indicated on the label (dispose of product!).

This product should be stored in dry surroundings, preferably protected from sources of direct and indirect light or heat.

Store in cool and dry surroundings, protected from light and ozone (no ultraviolet light)!
Avoid high temperatures during storage.

Set components

The composition of the TracheoFlex S set varies in quantity and components **depending on the model.**

Depending on set configuration, it can contain the following:

- Tracheostomy tube, armoured, with flexible, non-displacing neck plate
 - with / without phonation window,
 - with / without low-pressure cuff
- Inner cannula with low-profile connector)
- 15-mm standard connector
- Sealing cap
- Cough cap
- Phonation valve
- Introducer
- Velcro® neckband

Note:

The neckband can be ordered separately.

Additional accessories:

The following accessories can also be ordered separately:

- Tracheostomy wound dressing
- Tracheostomy neckbands
- Artificial noses, filters, etc.
- Cleaning brushes for tracheostomy tubes
- Endotest filling and pressure monitoring device

Please refer to our current catalogues for an overview of the available accessories. To order these, please contact your Rüscher representative.

Description

Tube and shaft

Rüscher TracheoFlex S tracheostomy tubes made of thermosensitive polyvinyl chloride (PVC) are anatomically shaped and can be supplied, depending on the model, as tubes with or without cuff, with or without a phonation opening or a combination of either of these.

The tubes are reinforced with highly flexible armoring. The tube is radiopaque due to the armoring.

Cuff

The filling system for the low-pressure cuff consists of a filling tube, a pilot balloon and a self-closing valve with a Luer connection.

Note:

The notes mentioned in connection with an existing cuff do not apply to tracheostomy tubes without cuff!

Neck plate

The neck plate is anatomically shaped and, firmly fixed to the tube shaft. The neck plate cannot be adjusted, but can still adapt individually to the soft tissues of each patient's neck thanks to its special design.

Neckband

The non-stretch and non-irritating tube neckband (Velcro band) provides patient-friendly and gentle fixation of the neck plate on the patient's neck and can be individually adapted to the neck's circumference.

The tube neckband contained in the set as well as other tracheostomy neckbands can also be ordered separately.

Inner cannula

The tracheostomy tube can be used optionally with or without an inner cannula.

The advantage of the inner cannula is that in patients who produce increased secretions because of their underlying disease, the inner lumen of the tube can be made patent again immediately without changing the tube if it is narrowed or blocked by secretions.

In this case, the inner cannula is removed and replaced by a new one. After insertion of the inner cannula with the low-profile connector into the tube, it is locked at the tube with the 15 mm connector or the cough cap as accompanying set components.

A) Suggested directions for use:

1. Select the correct tube size

The correct choice of tube size depends on the size and shape of the individual tracheostoma, the diameter of the trachea and the patient's underlying disease and is primarily the responsibility of the attending physician. In general, the largest possible tube that will fit both the trachea and the tracheostoma on the basis of the patient's anatomy should be selected to ensure optimum ventilation of the lungs and, at the same time, the lowest possible airway resistance.

2. Preparation of the tracheostomy tube

- Prior to use, check that the set with its components is complete and functions correctly.

- **When using tubes with cuff:** The patency of the tube lumen and cuff filling system should be checked before intubation. Prior to use, check that the cuff, the pilot balloon and the valve are intact and have no leaks. Insert a Luer syringe into the valve which belongs to the filling system and fill with enough air to arrive at a full cuff that shows no creases.
- After this procedure, completely empty the cuff again. Do not overinflate the cuff!
- Insert the introducer into the tube until it sits tightly in the tube and the olive tip projects beyond the end of the tube.

Caution!

Transport damage or incorrect storage can lead to damage to the cuffs or other components of the tube.

Never use damaged tubes, leaky cuffs or other damaged set components!

3. Introducing the tracheostomy tube

- Intubate the patient following currently accepted medical techniques noting the specific „Warnings and precautions“ in these instructions for use.

Note:

The intubation procedure can be facilitated by previously wetting the introducer and the tube, each in the area of the tip, with a small amount of water-soluble lubricant.

When using tubes with cuff, a small amount of water-soluble lubricant can also be applied to the cuff in order to ease introduction.

Caution!

Do not use aerosols based on lidocaine (e.g. Xylocain Spray!).

- The introducer must be held in fixed position throughout the entire procedure. After introducing the tube into the trachea and positioning it correctly, the introducer should be removed.
- When removing the introducer, the tracheostomy tube must be held in position.
- **When using tubes with cuff:** Following intubation fix the neck plate of the tube with the neckband as described under 4. After that, fill the cuff with air using a Luer syringe or a suitable cuff pressure gauge (e.g. Rüscher Endotest) until the trachea is sufficiently sealed. The use of minimal occluding volume, leak-minimising techniques and regular monitoring of the cuff pressure help to reduce the incidence of many of the adverse reactions associated with the use of cuffed tracheostomy tubes (See also „Warnings, precautions, cuff-related“).

Caution!

Please ensure that, when using tubes with cuff, the cuff is not accidentally damaged by sharp-edged objects during intubation.

Caution!

In order to avoid possible late-effect cartilage damage, the tube should be positioned without tension in the tracheal opening.

4. Fixing the neck plate on the patient's neck

After removing the introducer, the entire system must be secured to the patient's neck in order to prevent the cannula from slipping out. To do this, the ends of the neckband are passed through the corresponding slots in the neck plate of the tube and fastened at the side with the Velcro closure.

To avoid irritation of the skin where it is in contact with the plate – particularly in the case of a moist tracheostoma – a commercially available perforated pad (e.g. Rüscher wound dressing for tracheostoma) must be placed between the skin and the neck plate.

Caution!

When using tubes with cuff, the tracheostomy tube should, for safety reasons, be fixed before the cuff is inflated, if it is introduced by the patient himself or only by auxiliary medical staff!

Attention!

- Avoid high tractive forces, which act on the trachea via the insufflation tubes and the connector screwed on the tube shaft.
Attention! Danger of extubation!
- On principle, avoid alcoholic disinfectants during care of the tracheostoma.
- Never wet the slit compress with disinfectants containing alcohol!

5. Introduction and fixation of the inner cannula

When the inner cannula with the low-profile connector is used, it is introduced carefully into the tracheostomy tube up to stop position and then fixed by means of the 15 mm connector or the cough cap to the tube's connection system in the manner described above.

Note:

For safety reasons, the 15 mm connector is secured by means of a separate clip mechanism. Please, make sure that the clip snaps in completely. Due to manufacturing tolerances, it is possible that, after the clip mechanism at the 15 mm standard connector has snapped in, an artificial positive pressure ventilation cannot yet achieve complete sealing of the system. In this case, continue turning the 15 mm standard connector after the clip mechanism has snapped in, until complete sealing of the system is achieved. Under these circumstances, the clip mechanism prevents accidental loosening and unintentional disconnection of the 15 mm connector through possible torque on the connection system by the ventilation tubes connected to it.

6. Cough cap

Depending on the set configuration, a cough cap can be included. The cough cap can be used either with or without the low-profile inner cannula. It is secured as described previously by means of the screw cap.

The cough cap provides for even more security for the patient, since secretions which may be coughed up are collected by the cap and then drained off.

Note: The cough cap should be attached so that the opening faces downwards.

7. Sealing cap

The sealing cap is used for weaning from the tracheostomy tube. It must be used only under medical/nursing supervision. It is secured as described previously by means of the screw cap.

Warning:

If the sealing cap is employed, ensure that a tracheostomy tube with phonation window and an inner cannula with phonation window are used.

Therefore, before application of the sealing cap, always make sure that the low-profile inner cannula is being removed from the TracheoFlex S tracheostomy tube, since it does not have a phonation window!

In the case of cuffed tubes, the cuff must be deflated before using the sealing cap!

If there is not adequate ventilation through the natural airways after the cap has been attached, the cap must be removed again immediately.

Caution!

The sealing cap must under no circumstances be used together with tubes with an inflated cuff, as complete closure of the airway can result (risk of suffocation)!

8. Phonation valve

The phonation valve can be used in patients with preserved speech. During inhalation, the air passes only through the tracheostomy tube. During exhalation, the valve closes automatically so that the exhaled air can escape through the natural airway.

In this connection it is vital that the tracheostomy tube is also equipped with a phonation window!

Before application always make sure that the low-profile inner cannula is being removed from the tracheostomy tube, since it does not have a phonation window!

During application, always ensure that the phonation openings of the tracheostomy tube are not obstructed or blocked des by secretion. If this is the case, the tube should be cleaned or replaced by a new one.

9. Use in association with techniques for percutaneous dilative tracheostomy (PDT), (Rüsch PercuQuick-Set or PercuTwist-Set)

The TracheoFflex S tracheostomy tubes can be used in principle with most percutaneous dilative tracheostomy techniques. When introducing the tube immediately after dilation of the tracheostoma, it must be borne in mind that this can sometimes be more difficult than in a surgical tracheostomy.

When using cuffed tubes:

When introducing a TracheoFflex S tracheostomy tube with cuff during a percutaneous dilative tracheostomy procedure, ensure in particular that the cuff is completely deflated, i.e. empty of air. Introduction of the tube is facilitated by pulling the cuff carefully upwards after deflating it and additionally wetting it with a water-soluble lubricant if necessary.

B) Warnings and precautions

1.) Warnings and precautions (cuff-related):

- Ensure that the filling system has no leaks. The cuff pressure of the low-pressure cuff should be measured regularly. During use, you should make sure at regular intervals that the system is free from defects. Any deviation from the selected sealing pressure should be pursued and promptly corrected.
- Empty the cuff prior to extubation by attaching a syringe to the cuff filling system and removing all of the air in the cuff until a certain vacuum is obtained within the syringe and the pilot balloon has collapsed completely.
- Inflation of the cuff by "feel" alone or by using a predefined volume of air is strongly advised against, as resistance alone is an unreliable guide during inflation. The diffusion of nitrous oxide, oxygen or air may either increase or decrease the cuff volume and cuff pressure.
- The use of a cuff pressure gauge (e.g. Rüsch ENDOTEST) will help to monitor and adjust the cuff pressure.
- Do not overinflate the low-pressure cuff. The cuff pressure should normally not exceed 25 cm H₂O. Overinflation can result in tracheal injury, damage to the cuff with subsequent deflation or cuff distortion, which may lead to airway obstruction.
- The use of lidocaine-based aerosols (e.g. Xylocain spray) has been associated with the development of holes in PVC cuffs (Jayasuiya KD, Watson WF: "P.V.C. cuffs and lidocaine-based aerosol"; Br J Anesth. 1981 Dec; 53 (12): 1368).
The use of a water-soluble lubricant gel is recommended.
- Various anatomical structures within the intubation route or any intubation aid with rough surfaces can also damage the integrity of the cuff. Particular care is required to avoid damaging the thin-walled cuff during intubation. If the cuff is damaged, the tube must not be used.
- **Caution:**
In the event of Radiation Therapy where the target organ (cancer) is adjacent to the location of the cuff on the tube, the cuff material can be compromised by Radiation and can thus lose its resilient properties. In this instance, a visual inspection and functional performance testing should be conducted to verify that the integrity of the cuff has not been compromised **during radiation therapy**.
As the duration and dosage of this therapy can influence the material capabilities of Poly Vinyl Chloride (PVC), and can therefore increase the risk of cuff damage, we recommend to check the tube and its function thoroughly after each applied dose fraction, and to exchange the tube in the event that performance is affected and / or damage to the material is observed.
- Deflate the cuff prior to repositioning the tube. Movement of the tube with the cuff inflated could result in damage to the cuff or in patient injury possibly requiring medical intervention.
- Syringes, stopcocks or other devices should not be left permanently connected to the cuff inflation system for extended periods of time.

2.) Warnings and precautions (general):

- Contains DEHP, DBP
The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical devices containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

- The use of tracheostomy tubes during procedures that involve the use of LASERS or electrocautery, especially in the presence of oxygen-enriched mixtures or gas mixtures containing nitrous oxide, can result in the ignition of the inflammable mixture with resulting emission of cauterising and toxic products, including hydrochloric acid (HCL) and result in burns in the patient.
- Avoid strong lateral and axial traction on the tube and neck plate. Particularly in ventilated patients high tensile forces can act on the connector and tube via the ventilation tubes, with the resulting danger of the tube being pulled out of the trachea (risk of extubation)!
- If the position of a patient is altered following an intubation, the tube must always be checked again to ensure that it is in correct position.
- Expert clinical judgement should be exercised in the selection of the appropriate size of tracheostomy tube for each individual patient.
- Intubation and extubation should be performed in accordance with currently accepted medical techniques.
- If the tracheostomy tube is lubricated prior to intubation, it is essential to ensure that the lubricant does not get into and occlude the tube lumen or cuff filling system, thereby preventing ventilation or causing damage to the cuff.
- Lubricant residues on the connection system of the tube must always be removed in order to avoid accidental disconnection or loss of fixation of the tube with the set components employed (e.g. phonation valve, etc.).
- The device should be disposed of in accordance with applicable national regulations.

3.) Warnings and precautions: Phonation valve

- The phonation valve may be used together with a tracheostomy tube only if the patient's upper airway is unobstructed.
- The phonation valve should only be used in combination with tracheostomy tubes that have a phonation window.
- Before application, please make sure that the low-profile inner cannula is removed from the tracheostomy tube, since it does not have a phonation window!
- If a cuffed tracheostomy tube is used, the **cuff must previously be completely deflated!** The phonation valve may only be used if the cuff is deflated!
- The function of the phonation valve must be checked prior to every use and cleaned or replaced by a new phonation valve if necessary.
- The phonation valve may only be used when the patient is conscious. Patients must be trained by qualified staff on how to use the phonation valve.
- In the event of incipient dyspnoea, the phonation valve must be removed immediately from the tracheostomy tube and appropriate measures to correct the dyspnoea must be instituted promptly.

C) Summary of important dimensions TracheoFlex S Sets (in accordance with DIN EN ISO 5366-1:2009)

Size / ID Inner diameter in mm	OD Nominal external diameter in mm	Nominal centre-line length Dimension A+B+C in mm	Dimension B for tubes <u>without</u> phonation opening in mm	Dimension B for tubes <u>with</u> phonation opening in mm	Angle in degrees°	Cuff resting diameter in mm	Inner cannula Nominal inner diameter in mm
6	9.4	49	5	8	110-120	18	4.0
7	10.4	67	11	15	110-120	22	5.0
8	11.4	77	16	15	110-120	24	6.0
9	12.4	77	15	11	110-120	25	7.0
10	13.4	88	18	11	110-120	26	8.0
11	14.4	88	20	15	110-120	28	9.0

D) Cleaning and disinfection of the tracheostomy tube set

**Cleaning of the tracheostomy tube and the set components
(only for products which are expressly suited for it – single patient use)**

If the tube is heavily soiled with encrusted secretion, blood or other body fluids, which cannot be removed completely with the tube in situ, the tracheostomy tube should be removed and cleaned outside the body. Simple cleaning should always be carried out, followed by disinfection if required.

Caution!

This tube set must not be resterilised!

1. Simple cleaning

The lumen of the removed tube can be cleaned with a soft brush (e.g. Rüschi cleaning brushes for tracheostomy tubes) and neutral domestic cleaning agents and then rinsed with sterile clear (lukewarm) water.

Warning! No temperatures above 50° C!

2. Disinfection and cleaning

When using disinfectants (e.g. disinfectants based on active oxygen) ensure that their exposure times and concentrations are exactly observed.

The disinfectant solution must be freshly prepared before each new application.

Thorough rinsing is indispensable to remove any possible residues!

Important!

We recommend the usage of disinfectants based on active oxygen (e.g. Sekusept classic, Sekusept aktiv). Cleaning with these agents has been tested and approved by us. This way you make sure that the tube set is not damaged through the use of disinfectants. Unsuitable disinfectants (especially those containing alcohol) can damage the tube and lead to dangerous malfunction.

Basically all disinfectants can penetrate the tube material and lead to irritation of the mucous membrane later when used again.

Thorough rinsing to remove the residues is therefore essential!

Warning!

- When using tubes with cuff: The use of brushes should be avoided in the cuff area, in order not to damage it.
- Cleaning in an alcohol bath or in solvents must be avoided in all circumstances because of the risk of material hardening and material shrinking.
- The stated exposure time and concentration of disinfectant solutions must not be exceeded (please adhere to the manufacturer's instructions!).
- Only use freshly prepared disinfectant solutions!

3. Resterilisation

Resterilisation is **not possible!**

4. Simplified home care cleaning procedure

In the home care of tracheostomy patients, routine disinfection of the tube is often unnecessary (to be decided only by the attending physician). In these cases, cleaning can be limited to a simple procedure.

In the hospital setting, disinfection is always recommended in addition to simple cleaning.

5. Frequency of cleaning

The tracheostomy tube and its components should be cleaned as required. The underlying disease, mucus production, the state of the tracheostoma and the patient's subjective sense of hygiene greatly determine how often cleaning needs to be performed. As a **rule of thumb**, we recommend cleaning the tracheostomy tube at least twice a week, disinfecting if necessary.

Important:

Carefully dry the individual components after preparing the suitable set components of the tube set for re-use.

6. Restrictions concerning the cycles of cleaning and reuse

(only for single patient use components!)

The design of the tube and the choice of material allow for repeated cleaning and re-use for one individual patient.

If the above-mentioned precautions are observed, approximately 20 cleaning cycles are normally possible, without impairing the efficiency of the tube. However, in particular cases the lifespan of a tracheostomy tube can be considerably shorter due to extraordinary pressure during wrong handling or due to improper cleaning. In such cases the tube must be replaced (see E, „Safety check before re-use“, page 9).

Uninterrupted indwelling time should not exceed approx. 14 days. The maximum total service life including cleaning cycles should not exceed 28 days.

Caution:

The tube neckband can not be cleaned and reused (single use only)!

- Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only medical devices may result in exposure to viral, bacterial, fungal or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these medical devices. This product is not designed to be cleaned, disinfected or sterilized..

7. Storage until re-use

Long-term tracheostomised patients, especially in the home care area, often have several complete tracheostomy tube sets available for personal use. The cleaned and dried tracheostomy tube sets must therefore be stored temporarily in dry closed containers, protected from light and ozone. The containers should prevent tube soiling and contamination and should be stored in a cool and dry place protected from dust.

Note:

The components of individual sets should be stored separately and should not be jumbled up.

Caution (with respect to cuffs)!

During storage the cuff must not be blocked.

E) Safety check before re-use

The function of the tube set must be checked before every use in order to ensure that the tube has not been damaged during cleaning. In particular, the following must be checked:

- Firm and intact connection of neck plate and tube shaft
- No kinks
- Tube's connecting system functioning correctly
- Inner lumen and phonation window, if present, must be patent (no adherent incrustations)
- Undamaged edges at the tip of the tube (no sharp edges or chipped material)
- When using tubes with cuff: Undamaged low-pressure cuff and filling system (see 2 „Preparation of the tracheostomy tube, page 3)

Warning:

On account should you continue to use set components if you observe the following:

- Damage, especially to the outer and inner surface of the shaft and cuff system (e.g. surface cracks, leaking cuff, parts of the armouring protruding through the tube wall, etc.)
- Brittleness or hardening of the tube, the cuff or other set components
- Kinks
- Leaks and/or damage to the cuff system
- Damage to the connection system or the neck plate

Safe use can no longer be expected under these circumstances!

All information is in accordance with the state of the art when going to press. Subject to technical alterations.

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SPAIN TELEFLEX MEDICAL IBERIA S.A.
SWITZERLAND TELEFLEX MEDICAL AG
**Willy Rüsch GmbH**
Druckerei Täuber GmbH, Waiblingen