FAHL

## TRACHEOTEC®

---

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EN

## I. FOREWORD

These instructions for use are valid for Fah!® tracheostomy tubes. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling of Fah!® tracheostomy tubes.

**Please read the instructions for use carefully before using the product for the first time!**

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the tracheostomy tube. It contains important information on the product!

The Fah!® tracheostomy tubes are medical devices as defined by the MPG (German medical devices act) and conform to, among others, the "Basic Requirements" as per Annex IX of Directive 93/42 EEC.

## II. INTENDED USE

Fah!® tracheostomy tubes are intended to stabilise the tracheostoma following laryngectomy or tracheostomy.

The tracheostomy tube is designed to keep the tracheostoma open.

Tracheostomy tubes with low pressure cuff are always indicated after tracheostomy has been performed for any cause if the cannula needs to be sealed against the wall of the trachea.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

**Users must be trained by medical professionals on how to use the Fah!® tracheostomy tubes safely and appropriately.**

Fah!® LINGO tube variants are intended for tracheotomised patients who retain their larynx or for laryngectomees who use a shunt valve (voice prosthesis users).

## III. WARNINGS

**Patients must be trained by medical professionals on how to use the Fah!® tracheostomy tubes safely.**

**Fah!® tracheostomy tubes must never be blocked, e.g. by mucus or encrustation. Blockage can cause asphyxiation!**

**Mucus in the trachea can be aspirated (removed by suction) through the tracheostomy tube with a tracheal suction catheter.**

**Damaged tracheostomy tubes must not be used and must be discarded immediately. Use of a damaged tube may result in airway compromise.**

**When the tracheostomy tubes are inserted or removed, this can lead to irritations, coughing or slight bleeding. If bleeding persists, consult your doctor immediately!**

**Tracheostomy tubes are single-patient products and intended for use by one patient only.**

**Do not re-use with another patient. Do not reprocess for re-use by another patient.**

**Do not use tracheostomy tubes during laser treatment (laser therapy) or with electrosurgical devices! Contact with a laser beam may damage the tube.**

**In the presence of oxygen-rich mixtures the tube may ignite, which may cause injury from heat and exposure to decomposing and toxic combustion products.**

**CAUTION!**

**Tracheostomy tubes with speech function are only recommended for tracheotomised patients with normal secretion and normal mucous tissue.**

**CAUTION!**

**If there is strong secretion, a tendency to form granulation tissue, if the patient is undergoing radiation therapy or if scabs have formed, a perforated cannula version is only recommended under regular surveillance by a doctor and provided that shorter replacement intervals are observed (as a rule once a week), because the perforations in the outer tube can increase the formation of granulation tissue.**

#### IV. COMPLICATIONS

The following complications could emerge when using this device:

Contamination of the stoma may make it necessary to remove the tube. Contaminations can furthermore cause infections that necessitate the use of antibiotics.

Unintentional Inhalation of a tube/button that was not sized correctly may need to be removed by a doctor/physician.

If mucus blocks the tube it should be removed and cleaned.

#### V. CONTRAINDICATIONS

**Do not use if the patient is allergic to the material.**

##### **CAUTION!**

**Never under any circumstances use cannula variants without low-pressure cuff for patients on mechanical respiration!**

##### **CAUTION!**

**Only use perforated/fenestrated cannula variants during artificial respiration after consultation with the treating doctor/physician.**

##### **CAUTION!**

**Tracheostomy tubes with speaking valves must not be used on any account by laryngectomised patients, as this could lead to serious complications, even suffocation!**

#### VI. PRECAUTIONS

The correct tracheostomy tube size should be selected by the treating physician or by trained medical professionals.

We highly recommend always keeping at least two replacement tubes on hand to ensure continuous supply.

Do not perform any type of adjustment, repair or alteration on a tracheostomy tube or compatible accessories. If a product is damaged, it must be discarded immediately.

##### **CAUTION!**

**United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).**

#### VII. PRODUCT DESCRIPTION

The TRACHEOTEC® tracheostomy tubes are made of medical-grade synthetic materials (polyvinyl chloride and others).

We supply the Fahl® tracheostomy tubes in various sizes and lengths.

Fahl® tracheostomy tubes are reusable medical products intended for single-patient use.

Fahl® tracheostomy tubes are only to be used by the same patient and not by another patient.

The maximum period of use is 29 days (reckoned from the date on which the package was opened). If a pathogenic microbe (e.g. MRSA) is present, the period of use/indwelling time is reduced accordingly.

Please see the size table in the appendix.

A particular feature of the Fahl® tracheostomy tubes is the specially shaped neck flange.

The tube tip is rounded to prevent irritation of the mucus membranes in the trachea.

There are two lateral fastening eyelets on the neck flange for attaching a tube holder.

It may be advisable to use tubes of different lengths alternately so the tip of the tube does not always touch the same spot in the trachea and cause irritation. This prevents pressure points or the formation of granulation tissue. Always discuss the exact procedure with your doctor/physician.

## MRT (Magnetic Resonance Tomography) information

### CAUTION!

As the tracheostomy tubes with cuff have a small metal spring in the non-return valve of the pilot balloon with inflation tube, the tracheostomy tube with cuff should not be used while an MRT (magnetic resonance tomography) examination is being performed. MRT is a diagnostic imaging technique in which inner organs, tissues and joints are depicted with the help of magnetic fields and radio waves. Metallic objects can be drawn into the magnetic field and can induce changes by their acceleration. Although the metal spring is extremely small and lightweight, interactions that could give rise to health impairments or to malfunction or damage of the technical equipment used or of the tracheostomy tube itself nevertheless cannot be ruled out in this context. If wearing a tracheostomy tube is indicated for keeping the tracheostoma open, we recommend use of a metal-free tracheostomy tube instead of the tracheostomy tube with cuff while an MRT is being performed.

### CAUTION!

**Do not use a tracheostomy tube with low-pressure cuff while undergoing a nuclear magnetic resonance tomography / MRT!**

#### 1. Neck flange

A particular feature of the Fahl® tracheostomy tube is the specially shaped neck flange, which is designed to fit the anatomy of the neck.

The size specifications are printed on the neck flange.

The neck flange of the tracheostomy tubes has two side eyelets for inserting a tube holder.

All Fahl® tracheostomy tubes with fastening eyelets include a disposable tube holder for single use. The tube holder secures the tracheostomy tube to the neck.

The insertion aid (obturator) included in delivery facilitates insertion of the tracheostomy tube.

Please follow the instructions provided with the tube holder carefully when attaching and removing it from your tracheostomy tube.

Ensure that the Fahl® tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

#### 2. Connectors and adapters

Connectors and adapters are used to connect compatible tube accessories.

The options available to a patient depend on the clinical condition, such as status post laryngectomy or tracheostomy.

A 15 mm standard connector enables secure connection to required accessories.

#### 3. Cannula tube

The cannula tube borders directly on the neck flange and directs the airflow into the trachea.

The X-ray contrast strip in the tube that runs along the side allows the tube to be depicted radiographically and the position of the tube to be checked.

In the TRACHEOTEC® Suction tracheostomy tube versions, the X-ray contrast strip runs along the middle of the cannula tube.

#### 3.1. Low-pressure cuff

In the product versions with low-pressure cuff, the very thin-walled, high-volume low-pressure cuff adapts well to the trachea and ensures reliable sealing if inflated correctly.

The low-pressure cuff can be inflated like a balloon. The small pilot balloon on the inflation tube indicates whether the tracheostomy tube is in sealed (inflated) or non-sealed condition.

The low pressure cuff is inflated via a pilot line with one-way valve and pilot balloon.

##### 3.1.1 Leakproofness test of the cannula and the low-pressure cuff (if present)

The cannula and the low pressure cuff must be tested for leakages directly before and after every insertion and thereafter at regular intervals. For this purpose, inflate the low-pressure cuff to 22 to 24 mm Hg (1 mm Hg corresponds to 1.35951 cm H<sub>2</sub>O) and then observe whether a spontaneous drop in pressure occurs (recommended for inflating and testing: Cuff pressure gauge control inflator, REF 19500). There should be no substantial drop in pressure during the observation period. This leakproofness test must also be performed prior to every renewed insertion (e.g. after cleaning the tracheostomy tube) (see picture 7c).

The following are possible signs indicating leakages in the cuff (balloon):

- Externally visible damage to the balloon (holes, fissures etc.)
- Audible and/or tactile perception of air escaping from the balloon
- Water in inlet tubes of the tracheostomy tube (after cleaning!)
- Water in the cuff (after cleaning!)
- Water in the pilot balloon (after cleaning!)
- No cough stimulus when pressure is applied to the pilot balloon

#### **CAUTION!**

**When testing the balloon or when inserting, removing or cleaning the tracheostomy tube, never under any circumstances use sharp or pointed objects, such as for instance forceps or clamps, since these can damage or even destroy the balloon. Should one of the above-named signs of leakage be detectable, the tracheostomy tube must under no circumstances be used any more since it is no longer functioning properly!**

### **3.2 Obturator**

Please check whether the obturator can easily be removed from the cannula before inserting the tracheostomy tube!

After checking the obturator for ease of movement, reinsert the obturator into the cannula for inserting the tracheostomy tube.

### **3.3 Suction opening (only in Suction tracheostomy tube versions)**

Secretions that have accumulated above the inflated cuff can be removed via the suction opening in the outer tube of the Suction tracheostomy tubes.

The suction tube leading to the outside can be connected to a syringe or a suction device. The user must decide in consultation with the treating physician and after appropriate risk analysis, which of these suction variants to use. The individual disease condition of the patient should be considered here.

Suction using a suction device may in any case only be performed if the suction device is equipped with vacuum adjustment. The suction pressure must never exceed – 0.2 bar.

When using the suction variants of the tracheostomy tubes, the secretions that have accumulated above the cuff must be suctioned off via the suction opening on the tracheostomy tube directly before deflating the cuff in order to prevent aspiration of secretions. In this way, additional suctioning with a suction catheter while deflating the cuff may be rendered unnecessary, which facilitates handling significantly.

#### **CAUTION!**

**When using any of the tracheostomy tube versions with a suctioning device („SUCTION“), take special care during the suctioning procedure to generate a vacuum for as short a time as possible; drying out of the subglottic space as a result of this must be avoided.**

#### **CAUTION!**

**Contraindicated in patients with increased bleeding tendency (e.g. during anticoagulant therapy). In such cases, the Suction tracheostomy tube with suction opening may not be used, since suctioning entails an increased risk.**

Tracheostomy tubes with speaking valve (LINGO/PHON) are used after tracheotomy with a complete or partially retained larynx and enable the user to speak.

#### **CAUTION!**

**Make sure to read all accompanying product information, instructions for use, indications and contraindications. Discuss use of the product with your doctor/physician before first use.**

## **VIII. DIRECTIONS FOR TUBE INSERTION AND REMOVAL**

### **For the Doctor/Physician**

The appropriate tracheostomy tube must be selected by a doctor/physician or trained medical professionals.

Select a tube that fits the patient's anatomy to optimise comfort and ventilation (breathing in and out).

## For the Patient

### CAUTION!

**Only insert the tracheostomy tube with the cuff in fully deflated condition (see picture 7a)!**

**Carefully examine the packaging to ensure that it has not been tampered with or damaged. Do not use the product if the packaging has been damaged.**

It is advisable to use sterile disposable gloves.

Carefully examine the tube before first use to make sure that it is not damaged and that there are no loose parts.

Should you notice any anomaly or anything unusual, DO NOT use the tube. Return the tube to the manufacturer for inspection.

The tube must always be cleaned and, if necessary, disinfected as follows before re-inserting according to the instructions provided below.

If secretions collect in the lumen of the Fahl® tracheostomy tube and these cannot be removed by coughing or suctioning, the tube or stoma button should be removed and cleaned.

After cleaning and/or disinfection, carefully examine the Fahl® tracheostomy tube for sharp edges, cracks, or other signs of damage, since these may impair function and/or injure the mucus membranes in the airways.

Never under any circumstances continue using damaged tracheostomy tubes.

If the tube material has become discoloured, the tracheostomy tube should be discarded immediately.

### WARNING

**Patients must be briefed by trained medical professionals in the safe use of the Fahl® tracheostomy tubes.**

## 1. Insertion of the tube

### Step-by-step instructions to insert Fahl® tracheostomy tubes.

Before application, users should clean their hands (see picture 3).

Remove tube from the package (see picture 4).

If an obturator is to be used, this must first be fully inserted into the cannula tube so that the collar on the gripping piece of the obturator comes into contact with the outer edge of the 15 mm connector and the tip of the olive projects beyond the tip of the cannula (proximal end of cannula). The obturator must be held in this position during the entire procedure.

When using tracheostomy tubes with low-pressure cuff, pay special attention to the following points:

Before inserting the tracheostomy tube, check the cuff (balloon) as well - it must be free of damages of any kind and must be leakproof so as to ensure tight sealing as required. We therefore recommend to perform a leakproofness test prior to every insertion (see Section VII, No. 3.1.1). The balloon must be emptied completely prior to insertion of the tracheostomy tube (see picture 7b)! When using an aid for dilating the tracheostoma, take care to ensure that the tracheostomy tube, and especially the cuff, is not damaged by friction.

Next, push a tracheal compress, e.g. SENSOTRACH® DUO (REF 30608) or SENSOTRACH® 3-Plus (REF 30780), onto the tracheostomy tube.

To facilitate insertion of the tracheostomy tube, it is advisable to lubricate the outer tube by wiping it with an OPTIFLUID® stoma oil wipe (REF 31550) which allows the stoma oil to be applied evenly to the whole surface of the tube (see picture 4a and 4b).

If you are inserting the tube yourself, use a mirror to make insertion of the Fahl® tracheostomy tube easier.

When inserting the Fahl® tracheostomy tube, hold it by the neck flange with one hand (see picture 5).

Pull the tracheostoma slightly apart with your free hand to allow the tip of the tube to fit into the tracheostoma more easily.

Special aids for dilating the tracheostoma are also available (tracheal dilator, REF 35500). These allow the tracheostoma to be dilated gently and evenly, for instance also in emergency situations with collapsing tracheostoma (see picture 6).

Make sure that the tube is not damaged by friction when using an instrument for assistance.

Now carefully insert the tube into the tracheostoma during the inspiration phase (while breathing in) while tilting your head slightly back (see picture 7).

Advance the tube into the trachea.

Straighten your head once the tube has been inserted further into the trachea.

Once the tube has been brought into the correct position in the trachea, remove the obturator and tighten the locking screw.

If an obturator is used, this must then immediately be removed from the tracheostomy tube.

The tracheostomy tubes should always be attached with a special tube holder. This stabilises the tube and thus ensures that the tracheostomy tube is securely seated in the tracheostoma (see picture 1).

### **1.9 Inflating the low-pressure cuff (if present)**

The low-pressure cuff is inflated by applying a defined pressure to the cuff via the Luer connection (standardised conical connection) of the inflation tube by means of a cuff pressure gauge (e.g. MUCOPROTECT®, REF 19500). If not instructed otherwise by the doctor/physician, we recommend a cuff pressure of at least 15 mm Hg (20 cm H<sub>2</sub>O) to 22 mm Hg (30 cm H<sub>2</sub>O). The cuff pressure should never under any circumstances exceed 22 mm Hg (approx. 30 cm H<sub>2</sub>O).

**Inflate the low pressure cuff at most to this target pressure and check to make sure that sufficient air is supplied via the tracheostomy tube.**

**Always make sure that the low-pressure cuff is undamaged and in perfect working order.**

If the desired sealing is not achieved even after trying repeatedly with the specified limit volume, a tracheostomy tube with larger diameter may be indicated.

The correct cuff pressure must be checked regularly, i.e. at least every 2 hours.

#### **CAUTION!**

All instruments used for inflating the cuff must be clean and free of foreign particles! Detach the instruments from the Luer connection of the inflation tube as soon as the cuff has been inflated and close the connection with the cap.

#### **CAUTION!**

If the maximum pressure is exceeded for longer periods of time, the blood circulation in the mucus membrane can be impaired (risk of ischaemic necrosis, pressure ulcers, tracheomalacia, tracheal stenosis, pneumothorax). In patients undergoing artificial respiration, the cuff pressure should not be allowed to drop below the cuff pressure value specified by the doctor/physician in order to prevent unnoticed aspiration. Hissing noises in the region of the balloon, especially during expiration, indicate that the trachea is insufficiently sealed by the balloon. If the trachea cannot be sealed with the pressure values specified by the doctor/physician, the entire air should be withdrawn again from the balloon and the sealing process should be repeated. If this does not lead to success, we recommend to use the next larger tracheostomy tube with balloon. Due to the permeability of the balloon wall for gases, it is normal for the pressure in the balloon to decline slightly over time, but it can on the other hand also rise unintentionally during gas anaesthesia. Regular pressure monitoring is therefore urgently recommended.

The cuff must never under any circumstances be inflated with excessive amounts of air, since this can lead to damage of the tracheal wall, tears in the low pressure cuff with subsequent deflation, or distortion of the cuff, in which case airway obstruction cannot be ruled out.

#### **CAUTION!**

During anaesthesia, the cuff pressure can rise/fall due to nitrous oxide (laughing gas).

## 2. Removing the tube

### CAUTION!

EN Accessories such as a tracheostoma valve or HME (Heat Moisture Exchanger) must be removed first before proceeding to remove the Fah!® tracheostomy tube.

### CAUTION!

If the tracheostoma is unstable, or in emergency situations (puncture/dilation tracheostomy), the tracheostoma can collapse after withdrawal of the tracheostomy tube, thereby impairing air supply. A fresh tracheostomy tube must be kept ready for use in such cases and must be quickly inserted if necessary. A tracheal dilator (REF 35500) can be used for temporarily securing the air supply.

The cuff must be emptied before removing the tracheostomy tube. The head should be tilted back slightly for removal of the tube.

### CAUTION!

Never use a cuff pressure gauge to empty the low-pressure cuff. Always use a syringe for this.

Before the air is removed from the balloon by means of a syringe and the tracheostomy tube is withdrawn, the region of the trachea above the balloon must first be cleaned by suctioning off secretions and mucus. If the patient is responsive and reflexes are intact, it is recommended that the patient be suctioned while at the same time unblocking the tracheostomy tube. Suctioning is performed by inserting a suction catheter through the cannula tube into the trachea. In this way, suctioning can be performed without any problems and gently for the patient and cough stimulus and the risk of aspiration are minimised.

**Next, deflate the low pressure cuff while suctioning off at same time.**

If secretions are present, these are now taken up by the suction tube and can no longer be aspirated. Please note that the tracheostomy tube must in every case be cleaned, if necessary disinfected, and lubricated with stoma oil as specified below prior to reinsertion.

Proceed very carefully to avoid injury to the mucus membranes.

The tube must always be cleaned and, if necessary, disinfected as follows before reinserting according to the instructions provided below.

### Step-by-step instructions to remove the Fah!® tracheostomy tubes:

The tracheostomy tubes should be removed with the head slightly tilted back. Grip the tube or button at the side by the neck flange or the housing (see picture 7).

Carefully remove the tracheostomy tubes.

## IX. CLEANING AND DISINFECTION

### CAUTION!

For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fah!® tracheostomy tube at least twice a day, more often during heavier secretion production.

### CAUTION!

For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fah!® tracheostomy tube at least twice a day, more often during heavier secretion production.

Keep this in mind, your personal cleaning schedule, including additional disinfection if required, must always be determined in consultation with your doctor/physician and according to your individual needs.

Regular disinfection is only required if this is indicated from the medical point of view as ordered by the doctor/physician. The reason for this is that the upper airways are not free of microbes in healthy patients either.

General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. We recommend chemical disinfection of the tubes according to the instructions provided below. Please consult your doctor/physician.

**CAUTION!**

**Cleaning agent and disinfectant residues on the tracheostomy tube can lead to irritations of the mucus membranes or other health impairments.**

Tracheostomy tubes are considered as belonging to the category of instruments with hollow cavities. Special care must therefore be taken when performing cleaning or disinfection that the tube is fully wetted with the solution used and is unobstructed (see picture 8).

The solutions used for cleaning and disinfection must be freshly prepared prior to each use.

**1. Cleaning**

Fahl® tracheostomy tubes must be cleaned / replaced regularly in accordance with the individual patient's needs.

Only use cleaning agents if the tube is outside the tracheostoma.

A mild, pH-neutral washing lotion can be used for cleaning the tube. We recommend using the special cannula-cleaning powder (REF 31100) as directed by the manufacturer.

Never clean Fahl® tracheostomy tubes with cleaning agents that are not approved by the manufacturer. Never use aggressive household cleaners, detergents, high-concentration alcohol or denture cleaners.

This may be very hazardous to health. The tube may also be destroyed or damaged.

Never under any circumstances expose the tubes to boiling water. This could seriously damage the tracheostomy tube.

Heating to 65°C can cause slight clouding of the material; this normalises again after some time though.

**Cleaning Steps****Cleaning the tracheostomy tubes without low-pressure cuff**

Remove any additional items before cleaning.

First, rinse the tube thoroughly under running water (see picture 9).

Use only lukewarm water to prepare the cleaning solution and observe the instructions for use of the cleaning agent.

We recommend the use of the cannula-cleaning tub with sieve insert (REF 31200) to facilitate the cleaning procedure.

To do this, hold the sieve insert by the upper edge to prevent contact with and contamination of the cleaning solution (see picture 10).

Place only one tracheostomy tube at a time in the sieve insert of the cannula-cleaning tub. If more than one tube is cleaned at a time, they may be compressed too strongly and damaged by excessive pressure.

Immerse the sieve insert with the tracheostomy tube components in the prepared cleaning solution.

After the immersion time has expired (see instructions for use of the cannula cleaning powder), wash the tracheostomy tube thoroughly several times with lukewarm, clear water (see picture 9). There must be no residues of the cleaning agent on the tube when it is inserted into the tracheostoma.

If necessary, for example, if persistent secretion residues cannot be removed by the cleaning bath, additional cleaning with a special cannula-cleaning brush (OPTIBRUSH®, REF 31850 or Optibrush® Plus with fibre top, REF 31855) may help. Only use the cleaning brush, if the tube is removed and already outside the tracheostoma.

Always insert the cannula cleaning brush into the cannula from the cannula tip (see picture 11).

Use the brush as directed and proceed with great care to avoid damage to the soft tube material.

The outer cannula must only be cleaned with the cuff inflated, and the safety balloon must not come into contact with the cleaning solution in order to prevent intrusion of cleaning solution into the balloon, since this would lead to serious functional impairment and health risks for the user.

Thoroughly rinse the tube or button under lukewarm running water or sterile saline solution (0.9% NaCl solution).

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

Never use tracheostomy tubes with impaired functionality or with damages such as sharp edges or cracks, because this may lead to injuries of the mucus membranes in the trachea. If damages are found, the tube must not be used under any circumstances.

The obturator can be cleaned in the same way as the tracheostomy tube.

### **Cleaning the tracheostomy tubes with low-pressure cuff**

The indwelling time and thus the replacement intervals of the tracheostomy tube must be determined individually in consultation with the treating doctor/physician. The outer cannula must be cleaned or the tracheostomy tube replaced after 1 week at the latest though, since after that the risk of formation of granulation tissue or of tracheomalacia, etc. can increase significantly, depending on the patient's disease condition.

The outer cannula with cuff can be cleaned/rinsed with sterile saline solution.

A brush must not be used for tracheostomy tubes with cuff in order to avoid damaging the balloon!

Damage to the balloon and hence a defect of the tracheostomy tube can only be avoided by careful and gentle handling.

## **2. Chemical Disinfection Instructions**

### **2.1 Disinfecting the inner cannula / cleaning the outer cannula without low-pressure cuff**

The Fahl® tracheostomy tube can be disinfected by cold disinfection with special chemical disinfectants.

Disinfection should always be done if ordered by a doctor/physician due to specific health concerns caused by disease, infection, or your specific situation.

Disinfecting is generally indicated to prevent cross-infection and in inpatient situations (e.g. hospitals, nursing homes, and/or other health care facilities) to limit infections.

#### **CAUTION**

**Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).**

**Disinfectants that release chlorine or that contain strong alkalis or phenol derivatives must under no circumstances be used. This could badly damage or even destroy the tracheostomy tube.**

### **2.2 Disinfecting the outer cannula with low-pressure cuff**

Tracheostomy tubes with low-pressure cuff may only be disinfected provided that this is performed and checked with utmost care. The balloon must always be inflated beforehand.

Since even minor carelessness during disinfection of outer cannulas with low-pressure cuff can lead to damage and functional impairment of the same, we recommend replacement of the cannula instead of disinfection.

#### **Disinfection steps**

For this purpose, OPTICIT® tube disinfectant (REF 31180; not available in the USA) should be used in accordance with the manufacturer's instructions.

As alternative, we recommend a disinfectant based on glutaric aldehyde as active ingredient (also available in the USA). Always observe the area of application and spectrum of activity specified for the disinfectant by the manufacturer.

Follow the instructions for the disinfectant.

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

## **3. Sterilisation / Autoclaving**

**Resterilisation is not allowed.**

#### **CAUTION!**

**Heating to over 65°C, disinfection with boiling water or steam sterilisation are not allowed and lead to damage of the tracheostomy tube.**

## X. STORAGE/CARE

Currently not used and cleaned tubes should be stored in a clean plastic container, dry environment away from dust, sunlight and heat.

Still packed replacement tubes should be stored in a dry environment away from dust, sunlight and heat.

After cleaning tracheostomy tubes, make sure that the balloon has been deflated for storage.

**Only use stoma oil (REF 31500) as lubricant, or OPTIFLUID® stoma oil wipes (REF 31550).**

We highly recommend keeping at least two replacement tubes on hand to ensure continuous supply.

## XI. SERVICE LIFE

These tracheostomy tubes are products intended for single-patient use.

The maximum period of use is 29 days (reckoned from the date on which the package was opened). If a pathogenic microbe (e.g. MRSA) is present, the period of use/indwelling time is reduced accordingly.

The service life of a tube is influenced by many factors. The composition of the secretions, the quality of cleaning and other aspects, for example, are very important.

The maximum service life is not extended by using the tube intermittently (e.g. alternately with other tubes at normal cleaning intervals).

The tracheostomy tube must be replaced at the latest 29 days after opening the package.

This applies regardless of whether the tube was used every day, only on some days or not at all after the package was opened. After expiry of the 29-day period of use, the tubes must not be used anymore but must be disposed of correctly. Damaged tubes must be replaced immediately.

### CAUTION!

**Modifications to the tracheostomy tube of any kind, in particular perforations as well as repairs to the tube, may only be carried out by the manufacturer himself or by companies expressly authorised for this in writing by the manufacturer! Inexpertly performed work on tracheostomy tubes can lead to serious injuries.**

## XII. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the tracheostomy tube, especially due to shortening of the same or perforations, if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to the tubes themselves caused thereby and to any consequential damages caused thereby.

If the tracheostomy tube is used beyond the period of use specified under Section XI, and/or if the tracheostomy tube is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

TRACHEOTEC® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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
### Sizing information:

**O.D. Shield** = outer diameter behind the neck flange · **O.D. Tip** = outer diameter at the tube tip · **I.D.** = inner diameter at cannula tip · **O-tube** = outer cannula · **I-tube** = inner cannula · **Length** = length over the outer curve · **Length at mid radius** = length of centre line from neck flange to tip

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All measurements listed in the tables below are for standard cannulas.  
Customised cannulas are available upon request.

**GRÖSSENTABELLE TRACHEOTEC® TRACHEALKANÜLE /  
TRACHEOTEC® TRACHEOSTOMY TUBE – SIZE TABLE**

Größe / Size	A.D. Schild mm / OD of flange mm	A.D. Spitze mm / OD of tip mm	I.D. mm / ID mm	Länge mm / Length mm	Länge mittl. Radius mm / Length at mid radius	TRACHEOTEC®	TRACHEOTEC® UNI	TRACHEOTEC® SUCTION	TRACHEOTEC® LINGO	TRACHEOTEC® LINGO CUFF
3,0	4,7	4,7	3,0	47,2	43,7	x	-	-		
3,5	5,3	5,3	3,5	50,6	46,7	x	-	-		
4,0	6,0	6,0	4,0	57,1	52,7	x	-	-		
4,5	6,7	6,7	4,5	57,7	52,7	x	-	-		
5,0	7,3	7,3	5,0	61,1	55,7	x	x	x		
5,5	8,0	8,0	5,5	63,6	57,7	x	x	x		
6,0	8,7	8,7	6,0	67,1	60,7	x	x	x		
6,5	9,3	9,3	6,5	73,5	66,7	x	x	x		
7,0	10,0	10,0	7,0	77,0	69,6	x	x	x	x	x
7,5	10,7	10,7	7,5	81,5	73,6	x	x	x	x	x
8,0	11,3	11,3	8,0	87,9	79,5	x	x	x	x	x
8,5	12,0	12,0	8,5	90,4	81,5	x	x	x	x	x
9,0	12,7	12,7	9,0	94,9	85,5	x	x	x	x	x
10,0	14,0	14,0	10,0	105,3	95,0	x	x	-	x	x

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UNI with cuff and 15 mm swivel connector · SUCTION with cuff and suction line  
LINGO with sieve · LINGO CUFF with cuff and sieve

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