RÜSCH / RÜSCH CARE
ULTRA TRACHEOFLEX, EXTRA TRACHEOFLEX

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pl

Instructions for Use for TracheoFlex Tracheostomy Tubes:

ULTRA TracheoFlex, EXTRA TracheoFlex with low-pressure cuff

Important:

Please, read the entire instructions prior to use.

Sterile, as long as packaging has not been opened or damaged.

Discard after use, for single patient use only.

The tube can be cleaned and re-used for this patient
(see par. B "Cleaning and disinfecting the tube set").

Note:

Products which have been labelled differently (use once) may only be used according to what is stated on the label!

Sterilisation method, see packaging.

Exposure to elevated temperatures, UV light and excessive moisture should be avoided during storage.

Warning

 DEHP Contains DEHP, DBP DBP The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical devices containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Note:

The enclosed instructions for use are valid for the following products:

TracheoFlex tracheostomy tubes with adjustable fixation flange with low-pressure cuff.

This tracheostomy tube set is a sterile, reusable product

for single-patient use only. Resterilisation is not allowed!

These instructions for use are primarily directed at the tracheostomised patient himself as well as at the nursing staff. They should guide them in how to safely use the tube set, clean it and above all correctly resterilise it. For this reason, the instructions for use on hand should always be read by the user, the nursing staff as well as the patient himself.

Caution!

This product is delivered in sterile condition by the manufacturer.

However, sterility is only guaranteed if the packaging is intact.

Do not expose to temperatures above 50 °C.

Indication

Cannulation of tracheotomised patients through an existing tracheostoma.
Diagnosis only by a doctor!

Contraindications

For patients during radiation therapy and magnetic resonance imaging.

Set components:

- Fixed standard connector
- Low-pressure cuff, valve for Luer and Luer-Lock syringe tip
- Mounted adjustable fixation flange
- Velcro strap to fix the fixation flange to the neck
- Insertion aid

Description:

- Shaft

This product consists of a highly flexible, armoured tracheostomy cannula made of polyvinylchloride (medical grade PVC).

- **Cuff**

The cuff filling system of the low-pressure cuff consists of a filling tube, a pilot balloon, a self-closing valve with a Luer connection.

- **Fixation flange**

The fixation flange may be moved along the tube shaft and can thus be individually adapted to the patient's individual anatomy and can be secured to the tube in optimum position.

Warning!

The graduation on the tube shaft can be used for orientation.

It is, however, not an indication of measurements!

The neck band (velcro strap) ensures patient-friendly, gentle fixation of the flange to the patient's neck and can be adapted to the thickness of the neck.

Caution:

The tube neckband can not be cleaned and reused (single use only)!

Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only medical devices may result in exposure to viral, bacterial, fungal or prionic pathogens.

Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these medical devices. This product is not designed to be cleaned, disinfected or sterilized.

Delivery form and storage recommendations

This tracheostomy tube set is a sterile reusable product
to be used with one patient only.

Each set is supplied in an individual sterile pouch.

For the sterilisation method used see product label.

Sterility of the set and its components can only be guaranteed as long as the packaging has not been opened or damaged.

Do not resterilise this product and do not use beyond its expiry date, but dispose!

This product should be stored in dry surroundings, preferably away from direct or indirect sources of heat and light.

Store in cool surroundings, protected from light and ozone (if applicable, no ultraviolet light)! Avoid high temperatures.

A) Suggested Directions for Use

1. Selection of the correct tube size

The correct choice of tube size depends on the size and shape of the individual tracheostoma, the diameter of the trachea and the patient's underlying disease and is primarily the responsibility of the attending physician. In general, the largest possible tube that will fit both the trachea and the tracheostoma on the basis of the patient's anatomy should be selected to ensure optimum ventilation of the lungs and, at the same time, the lowest possible airway resistance.

2. Preparation of the tracheostomy tube

- Prior to use, check that the set is complete and functions correctly.
- The patency of the tube lumen and cuff filling system should be checked before intubation.

Prior to use, check that the cuff, the pilot balloon and the valve are intact and have no leaks.


Insert a Luer syringe into the adapter belonging to the filling system and inject enough air to inflate the cuff fully (until it shows no creases).

- After this procedure, completely evacuate the air from the cuff.
- Please also ensure that the fixation flange may be adjusted on the shaft and that secure fixation to the tube shaft can be achieved by arresting the fixaton screw. The area of the tube shaft where the fixation flange is to be secured must not be wetted with lubricants, grease or oil. If necessary, clean briefly with ethyl alcohol.


Caution!

Transport damage or incorrect storage may cause damage to the cuff and other tube components. Never use damaged tubes and leaky cuffs!


3. Preparation of the new fixation flange

ill.1

ill.2

ill.3

The slightly curved wings of the fixation flange must always point upwards, the lettering must not be upside down, the fixation flange must be aligned at a right angle to the cannula shaft (ill.1).

The marker on the swivel lock must be in open position (9 o'clock, ill. 2). Before use, check whether the fixation flange can be freely moved along the tube shaft (marker on swivel lock must be in open position, 9 o'clock, ill.3).

4. Inserting the tracheostomy tube

- Insert the introducer into the tube until it sits tightly in the tube and the olive tip projects beyond the end of the tube.
- Intubate the patient following currently accepted medical techniques noting the specific "Warnings and precautions" in these instructions for use.
- The introducer must be held in fixed position throughout the entire procedure. After introducing the tube into the trachea and positioning it correctly, the introducer should be removed. When removing the introducer, the tracheostomy tube must be held in position.

Note:

The intubation process can be facilitated by prior wetting of the insertion aid and tube tips with a small amount of a water-soluble lubricant.

Caution!

Please ensure that the cuff of the tube is not inadvertently damaged by sharp-edged objects during intubation.

To avoid subsequent injury to cartilage, the tube should be in the tracheal opening without tension.

Too deep an intubation should be avoided as otherwise there is the danger of one-sided ventilation of the lung!


5. Adapting and securing the fixation flange to the tube shaft

Without changing the tracheal position of the tube, the fixation flange can now be adjusted in the longitudinal axis of the tube to adapt it to the anatomical conditions of the neck surface.


After correct alignment of the fixation flange, the white swivel lock is turned clockwise a quarter of a revolution, i.e. by approx. 90°, from its open position (9 o'clock, ill.4a) into closed position (12 o'clock) in order to secure the fixation flange.

This is indicated by the marker on the swivel lock pointing upwards and the fixation flange being firmly connected with the shaft (ill.4b).


After correct fixation the marking on the fixation screw points upwards (see ill.5).

ill.4a

ill.4b

ill.5

Warning!

The graduation on the tube shaft can be used for orientation.
It is, however, not an indication of measurements!

Warning!

The intubation depth can be set according to individual requirements via the fixation screw on the adjustable neck plate. Ideally, this is done after previous fibroscopic inspection of the trachea. The fixation screw needs to be screwed tightly to safely prevent the tube from sliding out of the neck plate.


Once the patient is intubated and the fixation screw correctly fixated, use a Luer syringe or a suitable cuff pressure gauge (e.g. Rüsch ENDOTEST) to inflate the cuff with just enough air to provide an effective seal of the trachea. The use of minimal occluding volume, leak-minimising techniques and continuous monitoring of the cuff pressure can help to reduce the incidence of many of the adverse reactions associated with the use of cuffed tracheal tubes (see also "Warnings, precautions").

Caution!

Incorrect securing of the fixation screw can lead to the air supply channel of the cuff system being squeezed together somewhere on the shaft of the tube, thus rendering inflation or deflation of the cuff impossible.

How to identify a squeezed cuff inflation channel:

Due to incorrect positioning of the fixation flange or the swivel lock, the control balloon and cuff can only be filled / deflated with difficulty or not at all!

Necessary counteractive measures:

In this case the fixation screw must be loosened and the position of the fixation flange corrected. After correct fixation the marking on the fixation screw points upwards (see ill. 5).

Caution!

If the fixation flange cannot be secured due to contamination by lubricants, it is recommended to briefly clean the shaft.

Caution!

Do not immerse in alcohol!

6. Fixing the fixation flange to the patient's neck

After fixing the fixation flange on the tube shaft, it must be secured to the patient's neck with the velcro band in order to prevent the tube from slipping out. To achieve this, thread the ends of the neck band through the slots in the fixation flange and secure at the sides with the velcro fasteners.

To avoid irritation of the skin in the area of the fixation flange, especially in the case of a wet tracheostoma, it is obligatory to place a commercially available slitted compress between the skin and the fixation flange.

Caution!

For safety reasons, the tracheostomy tube should be fixed before the cuff is inflated, if it is introduced by the patient himself or only by auxiliary medical staff.

Caution!

- Avoid exerting high tensile forces on the trachea through the ventilation tubing and the connector on the tube shaft.
Warning! Danger of extubation!
- Do not sharply kink the tube shaft on the outside of the secured fixation flange. If, for example, the fixation flange is placed close to the connector end of the tube – in the case of very deep intubation – lateral kinking of the tube shaft can lead to a strong leverage effect on the tube shaft, which in turn can cause initial or total tearing of the shaft.
- Avoid all alcoholic disinfectants when looking after the tracheostoma.
- Never wet the slitted compress with disinfectants containing alcohol!

7. Warnings and precautions

Warnings and precautions (cuff related):

- Ensure that the filling system has no leaks. The cuff pressure of the low-pressure cuff should be measured regularly or permanently. During use, you should make sure at regular intervals that the system is free from defects. Even the least deviation from the selected sealing pressure should be pursued and promptly corrected.
- Empty the cuff prior to extubation by attaching a syringe to the cuff filling system and removing all of the air in the cuff until a certain vacuum is obtained within the syringe and the pilot balloon has collapsed completely.
- Inflation of the cuff by “feel” alone or by using a predefined volume of air is strongly advised against, as resistance alone is an unreliable guide during inflation. The diffusion of nitrous oxide, oxygen or air may either increase or decrease the cuff volume and cuff pressure.
- The use of a cuff pressure gauge (Rüsch ENDOTEST Ref. 11 27 00) will help to monitor and adjust the cuff pressure.
- Do not overinflate the low-pressure cuff. The cuff pressure should normally not exceed 25 cm H₂O. Overinflation can result in tracheal injury, damage to the cuff with subsequent deflation or cuff distortion, which may lead to airway obstruction.

- The use of lidocaine-based aerosols has been associated with the development of holes in PVC cuffs (Jayasuiya KD, Watson WF: "P.V.C. cuffs and lidocaine-base aerosol"; Br J Anesth. 1981 Dec; 53 (12): 1368). Expert clinical judgement is required when using this substance in order to prevent cuff leaks.
 - Various anatomical structures within the intubation route or any intubation aid with rough surfaces can also damage the integrity of the cuff. Particular care is required to avoid damaging the thin-walled cuff during intubation. If the cuff is damaged, the tube must no longer be used.
- **Caution:**
In the event of Radiation Therapy where the target organ (cancer) is adjacent to the location of the cuff on the tube, the cuff material can be compromised by Radiation and can thus lose its resilient properties. In this instance, a visual inspection and functional performance testing should be conducted to verify that the integrity of the cuff has not been compromised **during radiation therapy**.
As the duration and dosage of this therapy can influence the material capabilities of Poly Vinyl Chloride (PVC), and can therefore increase the risk of cuff damage, we recommend to check the tube and its function thoroughly after each applied dose fraction, and to exchange the tube in the event that performance is affected and / or damage to the material is observed.
- Deflate the cuff prior to repositioning the tube. Movement of the tube with the cuff inflated could result in damage to the cuff or in patient injury possibly requiring medical intervention.
 - Syringes, stopcocks or other devices should not be left inserted in the cuff inflation system for extended periods of time.

Warnings and precautions (general):

- Contains DEHP, DBP
The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical devices containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.
- The use of tracheostomy tubes during procedures that involve the use of lasers or electrocautery, especially in the presence of oxygen-enriched gas mixtures or mixtures containing nitrous oxide, can result in the tube igniting with resulting burns and with emission of cauterising and toxic products, including hydrochloric acid (HCL).
- Avoid strong lateral and axial traction on the cannula and neck plate. Particularly in ventilated patients, high tensile forces can act on the connector and tube via the ventilation tubes, with the resulting danger of the tracheostomy tube being pulled out of the trachea (risk of extubation)!
- If the position of a patient's tube is altered following an intubation, the tube must always be checked to ensure that it is in correct position.
- Expert clinical judgement should be exercised in the selection of the appropriate size of tracheal tube for each individual patient.
- Intubation and extubation should be performed in accordance with currently accepted medical techniques.

- If the tracheostomy tube is lubricated prior to intubation, it is essential to ensure that the lubricant does not get into and occlude the tube lumen or cuff filling system, thereby preventing ventilation or causing damage to the cuff.
- Lubricant residues on the connection system, on the fixation of the neck plate and on the area of the tube shaft within which the adjustable neck plate can be moved about, must always be removed in order to avoid accidental disconnection or loss of fixation of the tube in the neck plate.
- The device should be disposed of in accordance with applicable national regulations.

8. Summary of important dimensions in accordance with DIN EN ISO 5366-1: 2009

| ULTRA TracheoFlex | | | | | | | | |
|--------------------------|---------------|---------------|--|----------------------|----------------------|----------------------|--------------------------------|--------------|
| Order size | I.D. in mm | O.D. in mm | Overall length (CLL=A+B+C) in mm | Dimension A in mm | Dimension B in mm | Dimension C in mm | Angle in degrees° (approx.) | Cuff Ø in mm |
| 7.0 | 7.0 | 10.4 | 63.0 | 32.0 | 3.0 | 28.0 | 90 | 22 |
| 8.0 | 8.0 | 11.4 | 88.0 | 32.0 | 21.0 | 35.0 | 90 | 24 |
| 9.0 | 9.0 | 12.4 | 117.0 | 42.0 | 40.0 | 35.0 | 90 | 25 |
| 10.0 | 10.0 | 13.4 | 117.0 | 46.0 | 36.0 | 35.0 | 90 | 26 |
| 11.0 | 11.0 | 14.4 | 116.0 | 43.0 | 34.0 | 39.0 | 90 | 28 |

ULTRA TracheoFlex length variants:

Intra tracheal Length: +20 mm

Extra tracheal Length: +30 mm

| Order size | I.D. in mm | O.D. in mm | Overall length (CLL=A+B+C) in mm | Dimension A in mm | Dimension B in mm | Dimension C in mm | Angle in degrees° (approx.) | Cuff Ø in mm |
|------------|---------------|---------------|--|----------------------|----------------------|----------------------|--------------------------------|--------------|
| 7.0 | 7.0 | 10.4 | 113.0 | 52.0 | 33.0 | 28.0 | 90 | 22 |
| 8.0 | 8.0 | 11.4 | 138.0 | 52.0 | 51.0 | 35.0 | 90 | 24 |
| 9.0 | 9.0 | 12.4 | 167.0 | 62.0 | 70.0 | 35.0 | 90 | 25 |
| 10.0 | 10.0 | 13.4 | 167.0 | 66.0 | 66.0 | 35.0 | 90 | 26 |
| 11.0 | 11.0 | 14.4 | 166.0 | 63.0 | 64.0 | 39.0 | 90 | 28 |

ULTRA TracheoFlex length variant:

Extra tracheal Length: +40 mm

| Order size | I.D. in mm | O.D. in mm | Overall length (CLL=A+B+C) in mm | Dimension A in mm | Dimension B in mm | Dimension C in mm | Angle in degrees° (approx.) | Cuff Ø in mm |
|------------|---------------|---------------|--|----------------------|----------------------|----------------------|--------------------------------|--------------|
| 7.0 | 7.0 | 10.4 | 103.0 | 32.0 | 43.0 | 28.0 | 90 | 22 |
| 8.0 | 8.0 | 11.4 | 128.0 | 32.0 | 61.0 | 35.0 | 90 | 24 |
| 9.0 | 9.0 | 12.4 | 157.0 | 42.0 | 80.0 | 35.0 | 90 | 25 |
| 10.0 | 10.0 | 13.4 | 157.0 | 46.0 | 76.0 | 35.0 | 90 | 26 |
| 11.0 | 11.0 | 14.4 | 156.0 | 43.0 | 74.0 | 39.0 | 90 | 28 |

EXTRA TracheoFlex

| Order size | I.D. in mm | O.D. in mm | Overall length (CLL=A+B+C) in mm | Cuff Ø in mm |
|------------|---------------|---------------|--|--------------------|
| 7.5 | 7.5 | 10.3 | 155.0 | 22 |
| 8.5 | 8.5 | 11.3 | 155.0 | 24 |
| 9.5 | 9.5 | 12.3 | 155.0 | 25 |
| 10.5 | 10.5 | 13.3 | 155.0 | 26 |

B) Cleaning and disinfecting the tube set

Cleaning the tracheostomy tube (only with products which are expressly suitable for cleaning – single patient use)

If the tube is heavily soiled with encrusted secretion, blood or other body fluids which cannot be removed completely when the cannula is in place, the tracheostomy tube should be removed and cleaned outside the body. Simple cleaning (see below) should always be carried out, followed by disinfection if required. These tube sets must not be resterilized.

1. Simple cleaning

The lumen of the removed tube may be cleaned with a soft brush and neutral soap solution and then rinsed out with clear, hand-hot water.

Caution!

Temperatures should not exceed 50 °C!

2. Disinfection and cleaning

When using disinfectants (e.g. disinfectants based on active oxygen) ensure that their exposure times and concentrations are exactly observed. The disinfectant solution must be freshly prepared before each new application.

Thorough rinsing is indispensable to remove any possible residues!

Important!

We recommend the usage of disinfectants based on active oxygen (e.g. Sekusept classic, Sekusept aktiv). Cleaning with these agents has been tested and approved by us. This way you make sure that the tube set is not damaged through the use of disinfectants. Unsuitable disinfectants (especially those containing alcohol) can damage the tube and lead to dangerous malfunction.

Basically all disinfectants can penetrate the tube material and lead to irritation of the mucous membrane later when used again.

Thorough rinsing to remove the residues is therefore essential!

Warning!

- The use of brushes in the cuff region should be avoided to prevent injury to the cuff.
- Cleaning in an alcoholic bath and in solvents must be avoided at all costs due to the danger of material hardening and shrinking.
- The prescribed exposure times and concentrations of disinfecting solutions must not be exceeded (please observe manufacturer's instructions).
- Only use freshly prepared disinfecting solutions.

Warning!

Disinfectants can penetrate the tube material and lead to mucosal irritation when reused. Thorough rinsing is absolutely necessary to remove residues.

3. Resterilisation

Resterilisation is not allowed!

4. Frequency of cleaning

The cannula and its components should be cleaned as required. The underlying disease, mucus production, state of the tracheostoma and the subjective hygiene feeling of the patient greatly determine how often cleaning is performed. As a rule of thumb we recommend cleaning the tracheostomy cannula at least twice a week, disinfecting it if necessary. Carefully dry the individual components after preparing the tube set.

C) Storage until reuse

Long-term tracheotomized patients often have several complete tracheostomy sets available for personal use. Thus the cleaned and dried tube sets must be stored in dry, closed containers which are protected from light and ozone. The containers should prevent soiling and contamination of the tubes and thus should be stored in a cool and dry place, protected from dust.

Caution!

The cuff must not be inflated during storage.

D) Safety testing before reuse

The proper function of the tube set must be tested before every renewed use in order to ensure that the tube set has not been damaged due to being reused too many times or being prepared incorrectly for reuse.

The following must be tested in particular:

- Kinking
- Working connection system at the tube
- Intact edges at the tip of the tube (no sharp edges or chipped material)
- Intact low-pressure cuff and filling system
(see preparation of cannula set)
- Intact fixation flange and tube shaft
- The inner lumen must be patent (no encrustation stuck to the walls)

Warning!

On no account should you continue to use set components if you observe the following:

- Damage, especially to the outer or inner surface of the shaft and cuff system (e.g. surface cracks, leaking cuff, parts of the armouring penetrating the tube wall, etc.)
- Brittleness or hardening of cuff or other set components
- Kinking
- Leaks and/or damage to the cuff system
- Separation of the standard connector from the shaft
- Damage of the fixation screw and the neck plate.

Safe use can no longer be expected under these circumstances!

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ULTRA TracheoFlex

| Misura | D.I. in mm | D.E. in mm | Lunghezza totale (CLL=A+B+C) in mm | Dimen- sione A in mm | Dimen- sione B in mm | Dimen- sione C in mm | Angolo in gradi ° (ca.) | Cuffia Ø in mm |
|--------|---------------|---------------|---|----------------------------|----------------------------|----------------------------|-------------------------------|-------------------------|
| 7,0 | 7,0 | 10,4 | 63,0 | 32,0 | 3,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 88,0 | 32,0 | 21,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 117,0 | 42,0 | 40,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 117,0 | 46,0 | 36,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 116,0 | 43,0 | 34,0 | 39,0 | 90 | 28 |

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| Misura | D.I. in mm | D.E. in mm | Lunghezza totale (CLL=A+B+C) in mm | Dimen- sione A in mm | Dimen- sione B in mm | Dimen- sione C in mm | Angolo in gradi ° (ca.) | Cuffia Ø in mm |
|--------|---------------|---------------|---|----------------------------|----------------------------|----------------------------|-------------------------------|-------------------------|
| 7,0 | 7,0 | 10,4 | 113,0 | 52,0 | 33,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 138,0 | 52,0 | 51,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 167,0 | 62,0 | 70,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 167,0 | 66,0 | 66,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 166,0 | 63,0 | 64,0 | 39,0 | 90 | 28 |

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| Misura | D.I. in mm | D.E. in mm | Lunghezza totale (CLL=A+B+C) in mm | Dimen- sione A in mm | Dimen- sione B in mm | Dimen- sione C in mm | Angolo in gradi ° (ca.) | Cuffia Ø in mm |
|--------|---------------|---------------|---|----------------------------|----------------------------|----------------------------|-------------------------------|-------------------------|
| 7,0 | 7,0 | 10,4 | 103,0 | 32,0 | 43,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 128,0 | 32,0 | 61,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 157,0 | 42,0 | 80,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 157,0 | 46,0 | 76,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 156,0 | 43,0 | 74,0 | 39,0 | 90 | 28 |

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8. Resumo dos tamanhos importantes

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| ULTRA TracheoFlex | | | | | | | | |
|--------------------------|---------------|---------------|--|--------------------------|--------------------------|--------------------------|---------------------------------|---------------------|
| Ta- manho | D.I. em mm | D.E. em mm | Comprimento total (CLL=A+B+C) em mm | Dimen- são A em mm | Dimen- são B em mm | Dimen- são C em mm | Ángulo em graus° (aprox.) | Balão Ø em mm |
| 7,0 | 7,0 | 10,4 | 63,0 | 32,0 | 3,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 88,0 | 32,0 | 21,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 117,0 | 42,0 | 40,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 117,0 | 46,0 | 36,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 116,0 | 43,0 | 34,0 | 39,0 | 90 | 28 |

ULTRA TracheoFlex variantes de comprimento:

Comprimento intratraqueal: +20 mm

Comprimento extra traqueal: +30 mm

| Ta- manho | D.I. em mm | D.E. em mm | Comprimento total (CLL=A+B+C) em mm | Dimen- são A em mm | Dimen- são B em mm | Dimen- são C em mm | Ángulo em graus° (aprox.) | Balão Ø em mm |
|--------------|---------------|---------------|--|--------------------------|--------------------------|--------------------------|---------------------------------|---------------------|
| 7,0 | 7,0 | 10,4 | 113,0 | 52,0 | 33,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 138,0 | 52,0 | 51,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 167,0 | 62,0 | 70,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 167,0 | 66,0 | 66,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 166,0 | 63,0 | 64,0 | 39,0 | 90 | 28 |

ULTRA TracheoFlex variante de comprimento:

Comprimento extra traqueal: +40 mm

| Ta- manho | D.I. em mm | D.E. em mm | Comprimento total (CLL=A+B+C) em mm | Dimen- são A em mm | Dimen- são B em mm | Dimen- são C em mm | Ángulo em graus° (aprox.) | Balão Ø em mm |
|--------------|---------------|---------------|--|--------------------------|--------------------------|--------------------------|---------------------------------|---------------------|
| 7,0 | 7,0 | 10,4 | 103,0 | 32,0 | 43,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 128,0 | 32,0 | 61,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 157,0 | 42,0 | 80,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 157,0 | 46,0 | 76,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 156,0 | 43,0 | 74,0 | 39,0 | 90 | 28 |

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ULTRA TracheoFlex

| Bestelmaat | I.D. in mm | O.D. in mm | Gehele lengte (CLL=A+B+C) in mm | Afmeting A in mm | Afmeting B in mm | Afmeting C in mm | Hoek in graden° (ongeveer) | Cuff ∅ in mm |
|------------|---------------|---------------|---------------------------------------|---------------------|---------------------|---------------------|----------------------------------|--------------------|
| 7,0 | 7,0 | 10,4 | 63,0 | 32,0 | 3,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 88,0 | 32,0 | 21,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 117,0 | 42,0 | 40,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 117,0 | 46,0 | 36,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 116,0 | 43,0 | 34,0 | 39,0 | 90 | 28 |

Ultra TracheoFlex lengte variant:

Intratracheale lengte: +20 mm

Extratracheale lengte: +30 mm

| Bestelmaat | I.D. in mm | O.D. in mm | Gehele lengte (CLL=A+B+C) in mm | Afmeting A in mm | Afmeting B in mm | Afmeting C in mm | Hoek in graden° (ongeveer) | Cuff ∅ in mm |
|------------|---------------|---------------|---------------------------------------|---------------------|---------------------|---------------------|----------------------------------|--------------------|
| 7,0 | 7,0 | 10,4 | 113,0 | 52,0 | 33,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 138,0 | 52,0 | 51,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 167,0 | 62,0 | 70,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 167,0 | 66,0 | 66,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 166,0 | 63,0 | 64,0 | 39,0 | 90 | 28 |

Ultra TracheoFlex lengte variant:

Extratracheale lengte: +40 mm

| Bestelmaat | I.D. in mm | O.D. in mm | Gehele lengte (CLL=A+B+C) in mm | Afmeting A in mm | Afmeting B in mm | Afmeting C in mm | Hoek in graden° (ongeveer) | Cuff ∅ in mm |
|------------|---------------|---------------|---------------------------------------|---------------------|---------------------|---------------------|----------------------------------|--------------------|
| 7,0 | 7,0 | 10,4 | 103,0 | 32,0 | 43,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 128,0 | 32,0 | 61,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 157,0 | 42,0 | 80,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 157,0 | 46,0 | 76,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 156,0 | 43,0 | 74,0 | 39,0 | 90 | 28 |

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ULTRA TracheoFlex

| Bestillings størrelse | I.D. mm. | U.D. mm. | Total længde (CLL=A+B+C) i mm. | Dimen- sion A i mm. | Dimen- sion B i mm. | Dimen- sion C i mm. | Vinkel ° i grader (ca.) | Cuff Ø i mm |
|--------------------------|-------------|-------------|--------------------------------------|---------------------------|---------------------------|---------------------------|-------------------------------|-------------------|
| 7,0 | 7,0 | 10,4 | 63,0 | 32,0 | 3,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 88,0 | 32,0 | 21,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 117,0 | 42,0 | 40,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 117,0 | 46,0 | 36,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 116,0 | 43,0 | 34,0 | 39,0 | 90 | 28 |

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| Bestillings størrelse | I.D. mm. | U.D. mm. | Total længde (CLL=A+B+C) i mm. | Dimen- sion A i mm. | Dimen- sion B i mm. | Dimen- sion C i mm. | Vinkel ° i grader (ca.) | Cuff Ø i mm |
|--------------------------|-------------|-------------|--------------------------------------|---------------------------|---------------------------|---------------------------|-------------------------------|-------------------|
| 7,0 | 7,0 | 10,4 | 113,0 | 52,0 | 33,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 138,0 | 52,0 | 51,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 167,0 | 62,0 | 70,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 167,0 | 66,0 | 66,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 166,0 | 63,0 | 64,0 | 39,0 | 90 | 28 |

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| Bestillings størrelse | I.D. mm. | U.D. mm. | Total længde (CLL=A+B+C) i mm. | Dimen- sion A i mm. | Dimen- sion B i mm. | Dimen- sion C i mm. | Vinkel ° i grader (ca.) | Cuff Ø i mm |
|--------------------------|-------------|-------------|--------------------------------------|---------------------------|---------------------------|---------------------------|-------------------------------|-------------------|
| 7,0 | 7,0 | 10,4 | 103,0 | 32,0 | 43,0 | 28,0 | 90 | 22 |
| 8,0 | 8,0 | 11,4 | 128,0 | 32,0 | 61,0 | 35,0 | 90 | 24 |
| 9,0 | 9,0 | 12,4 | 157,0 | 42,0 | 80,0 | 35,0 | 90 | 25 |
| 10,0 | 10,0 | 13,4 | 157,0 | 46,0 | 76,0 | 35,0 | 90 | 26 |
| 11,0 | 11,0 | 14,4 | 156,0 | 43,0 | 74,0 | 39,0 | 90 | 28 |

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SPAIN TELEFLEX MEDICAL IBERIA S.A.
SWITZERLAND TELEFLEX MEDICAL AG
**Willy Rüsch GmbH**
Druckerei Täuber GmbH, Waiblingen