Respironics Inc.
1001 Murry Ridge Lane
Murrysville, PA 15668 USA

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1 Safety Information

Intended Use

The Wisp Pediatric Nasal Mask and SE elbow accessory are intended to provide an interface for application of non-invasive ventilation, including CPAP and bi-level therapy, to patients. The mask is for single patient use in the home and multipatient use in the hospital/institutional environment. The mask is to be used on patients greater than 10 kg/22 lbs (>10 kg/22 lbs) for whom CPAP, bi-level therapy, or non-invasive ventilation has been prescribed. An exhalation port is not built into the elbow and a separate exhalation device must be used with this elbow.

Note: This mask is not made with natural rubber latex, DEHP, or BPA.

Contraindications

This mask is not intended for patients who are dependent on mechanical ventilation for their life support.


⚠ Warnings

- When using the Wisp Pediatric Nasal Mask with the SE elbow, it is the responsibility of the healthcare professional to ensure proper device set-up in accordance with patient condition, including the use of proper alarms.
- The mask must be used under qualified supervision for patients who are unable to remove the mask by themselves.
- This elbow requires a separate exhalation device.
- Use of this mask while the system is not on and operating may cause rebreathing of exhaled air. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation.
- Attaching an exhalation device requires therapy pressure level adjustment to compensate for increased leak.

Note: For a complete list of warnings associated with the Wisp Pediatric Nasal Mask, refer to the instructions for use that came with your mask.

2 Accessory Parts

Features


3 Before Use/Accessory Care

- Hand wash the accessory before use and daily (refer to the Wisp Pediatric Nasal Mask instructions for use for additional care information).
- Inspect the mask. Replace any parts that are broken or worn.

4 Assembly/Disassembly

Attaching the Elbow

Push the elbow into the cushion until the quick release tabs click into place **1****2**.

If the removable adapter is required to connect to the circuit, attach the adapter to the standard elbow as shown in the Accessory Parts section.

Removing the Elbow

Squeeze the quick release tabs and pull the elbow away from the cushion **2****1**. Detach the adapter from the elbow.

5 Health Care Provider & Clinician Information

Philips Respironics System One Resistance Control

This symbol may appear on your Philips Respironics therapy device or mask. The System One setting is not applicable to this mask with the SE elbow accessory, or any mask requiring a separate/additional exhalation device.

Multi-Patient Use


Refer to the Disinfection Guide for Professional Users to reprocess between patients in a clinical setting. Access the latest version of the Disinfection Guide at www.healthcare.philips.com or by contacting customer service at 1-800-345-6443 (USA or Canada only) or +1-724-387-4000.

Specifications

The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your therapy device, and that alerts and alarms function properly.

⚠ Warning: Use outside of these specifications may result in ineffective therapy. The pressure flow curve shown below is an approximation of expected performance. Exact measurements may vary.

Pressure Flow Curve

Resistance

Drop in Pressure at:	
50 SLPM	100 SLPM
All Sizes 0.3 cm H ₂ O	1.3 cm H ₂ O

Disposal

Dispose of in accordance with local regulations.

Deadspace

SCS	28.6 mL
SCM	34.1 mL
SCL	45.6 mL

Storage Conditions

Temperature: -4° F to 140° F (-20° C to 60° C)
Relative Humidity: 15% to 95%, non-condensing

Sound Levels

A-weighted Sound Power Level	15.0 dBA
A-weighted Sound Pressure Level at 1 m	8.1 dBA

6 Symbols

Manufacturer

Humidity Limits

Warning or Caution

Temperature Limits

Not Made with Natural Rubber Latex

7 Limited Warranty

Respironics, Inc. warrants that its mask systems (including mask frame and cushion) (the "Product") shall be free from defects of workmanship and materials for a period of ninety (90) days from the date of purchase (the "Warranty Period"). If the Product fails under normal conditions of use during the Warranty Period and the Product is returned to Respironics within the Warranty Period, Respironics will replace the Product. This warranty is nontransferable and only applies to the original owner of the Product. The foregoing replacement remedy will be the sole remedy for breach of the foregoing warranty.

This warranty does not cover damage caused by accident, misuse, abuse, negligence, alteration, failure to use or maintain the Product under conditions of normal use and in accordance with the terms of the product literature, and other defects not related to materials or workmanship.

This warranty does not apply to any Product that may have been repaired or altered by anyone other than Respironics. Respironics disclaims all liability for economic loss, loss of profits, overhead, or indirect, consequential, special or incidental damages which may be claimed to arise from any sale or use of the Product. Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

THIS WARRANTY IS GIVEN IN LIEU OF ALL OTHER EXPRESS WARRANTIES. IN ADDITION, ANY IMPLIED WARRANTY, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, IS SPECIFICALLY DISCLAIMED. SOME JURISDICTIONS DO NOT ALLOW DISCLAIMERS OF IMPLIED WARRANTIES, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS UNDER LAWS APPLICABLE IN YOUR SPECIFIC JURISDICTION.

To exercise your rights under this limited warranty, contact your local authorized Respironics, Inc. dealer or Respironics, Inc. at 1001 Murry Ridge Lane, Murrysburg, Pennsylvania 15668, USA.

Limited Warranty for Australia and New Zealand only

Note: For Australian and New Zealand customers this warranty replaces the warranty contained above.

Respironics, Inc., a Philips Healthcare company warrants that the Products shall be free from defects of workmanship and materials and will perform in accordance with the Product specifications. 2. This warranty is valid for a period of 30 days from the date of purchase from an authorised Respironics dealer. 3. If the Product is found to contain a defect of workmanship or materials or fails to perform in accordance with the Product specifications during the applicable warranty period, Respironics will repair or replace, at its option, the defective material or part. 4. The customer is responsible for returning the product to an authorised Philips Respironics dealer, and collecting the product from the authorised Philips Respironics dealer after repair or replacement, at its own cost. Philips Respironics is responsible only for the freight cost of transporting the product between the authorised Philips Respironics dealer and Respironics. Respironics reserves the right to charge an evaluation and postage fee for any returned Product as to which no problem is found following investigation. 5. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to materials or workmanship. 6. The warranty provided by Respironics herein is not transferrable by the Buyer in the event of any sale or transfer of Products purchased by the Buyer from an authorised Respironics dealer. 7. To the extent permitted by law, where the Buyer has the benefit of an implied guarantee under the Australian Consumer Law, but the Product is not of a kind ordinarily acquired for personal, domestic or household use or consumption Respironics' liability shall be limited, at the option of Respironics, to the replacement or repair of the Product or the supply of an equivalent Product. 8. To exercise your rights under this warranty, contact your local authorised Philips Respironics dealer. A list of all authorised dealers is available at the following link: http://www.healthcare.philips.com/au_en/homehealth/distributors_index.wpd.

Alternatively, you can make a claim under this warranty by contacting Respironics directly at: Philips Electronics Australia Limited, 65 Epping Road, North Ryde NSW 2113, Australia. Tel: 1300 766 488, Email: prcontact@philips.com. 9. The following statement is provided to a Buyer who is a "consumer" under the Australian Consumer Law: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the good repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. 10. The following statement is provided to a Buyer who is a "consumer" under the Consumer Guarantees Act 1993, New Zealand: Our goods come with guarantees that cannot be excluded under the Consumer Guarantees Act 1993. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

AUSTRALIAN SPONSOR DETAILS:

Philips Electronics Australia Ltd.
65 Epping Road, North Ryde, NSW 2113
Australia

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©2016 Koninklijke Philips N.V. All rights reserved.