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


Artificial nose for infants

Indication

For the efficient filtering, warming and humidifying of inhaled air for patients with spontaneous breathing through an endotracheal tube or tracheal cannula, in hospital or at home. Diagnosis by a medical specialist is a prerequisite of use.

Description

The sanabelle HME-A2 comprises a plastic housing with an insert of hydroscopic foam material. The HME has a 15 mm female connector, which can be fitted to the standardised 15 mm connector of the tracheal cannula. The HME is latex-free, DEHP-free and electrically non-conductive.

Usage

Unpack the sanabelle HME-A2 only before immediate use carefully. Easy to insert and remove thanks to the patented **clip-to-fit** function. Lightly squeeze the case on the side with the foam ① and slide the HME back to the catch on the 15 mm attachment on the cannula or the tube ②. The normal humidifying effect sets in after a few breaths. To easily remove the HME, squeeze the top of the case and remove it from the 15 mm attachment ②.

The HME must not be used together with a humidifier operating with heated or nebulised water or in combination with a nebuliser.

IMPORTANT: The instructions for use of the sanabelle HME-A2 must be carefully read prior to use, to minimise the potential risk of dyspnoea, labored breathing, hypercapnia and hypoxia.

Contraindications

- The sanabelle HME-A2 must not be used outside the recommended tidal volume range, in order to prevent an increased concentration of CO_2 in case of a low tidal volume, or insufficient humidification, if the tidal volume is too high. Please observe the doctor's recommendation. Please consult a specialist in case of questions.
- The device is not to be used on patients with strong exsiccosis or excessive secretions in the lungs and airways or in case of leaking gas flow (expired gas does not flow through the humidifier).
- Do not connect the HME to an artificial ventilation system. The HME does not have a machine connection.

Precautions:

- Always ensure the patency of the HME.
- Replace the sanabelle HME-A2 after a maximum of 24 hours to counteract any possible bacterial or fungicidal growth. The HME must be replaced more frequently depending on the needs, for example, if the breathing resistance is in danger of increasing due to accumulations of secretion.
- The HME must be slid back to the catch on the 15mm connection because otherwise the air can escape unfiltered through the slots on the bottom.
- Monitor patients for potential fatigue! Upon indication of difficulty breathing, remove the HME immediately! Reason: the HME's resistance requires increased breathing effort.
- The sanabelle HME-A2 is for single use only! Cleaning, disinfecting or re-sterilising it for reuse may impair the function and increase the risk of transmitting infections.
- Dispose of used HMEs immediately in a hygienically proper manner!

Technical Data

Tidal volume (V_T), recommended range	7,5 - 50 ml
Internal volume	3 ml
Humidity loss	$V_T = 50$ ml: 8,7 mg/l
Pressure drop inspiration	5 l/min: 0,2 hPa (~ cm H_2O)
Weight	2,5 g
Length	30 mm
Connection	Endotracheal or tracheostomy tubes and tracheal cannulae with a 15 mm universal attachment
Environmental impact	Hydrocarbon, carbon dioxide, water and carbon oxide are produced when the product is burned.

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