# FRANÇAIS

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## ENGLISH

**Prior to first use** as well as after each use clean or disinfect (sterilize) the products.

### Cleaning and Disinfection

For cleaning soap-suds or mild alkaline cleaning agents can be used. Hard brushes or other materials which might damage the silicone surface must not be used. For disinfection use commercially available disinfectants for instruments recommended for use with silicone and plastic materials. The immersion time and concentration stated by the manufacturer have to be followed strictly.

Disinfectants containing compounds similar to phenol or alkylamines (glucoprotamine) are unsuitable. Eliminate all residues of the cleaning agents and disinfectants used by thoroughly rinsing with water after each cleaning and disinfecting procedure.

### For treatment of our reusable products from silicone and polysulfone, we recommend the following agents of the company ECOLAB:

#### 1. For manual treatment:

Sekusept® powder classic 2% concentration  
Sekusept® active 2% concentration

#### 2. For mechanical treatment:

Sekumatic® FR 0.5% measured out by machine (Miele Vario programme)  
Sekumatic® FRE 0.5% measured out by machine (Miele Vario programme)  
Sekumatic® FNZ 0.3% neutralisation (Miele Vario programme)

The agents mentioned above have to be used following the instructions for use of ECOLAB (ECOLAB Germany - Anwendungstechnik - Phone: +49 (211) 98 93-0 - www.ecolab.com).

### Sterilization

For sterilization a validated steam sterilization procedure according to EN 554 at 134°C can be applied.

#### Attention:

To connect or disconnect the catheter mount, hold it at the connector or cuff in order to avoid cracking. During steam sterilization the connectors do not have to be removed from the silicone tubing, however they must not be connected to other standard connectors. Furthermore no hard and sharp-edged objects must come in contact with the product.

Please note that each processing accelerates natural ageing of the material. Furthermore product lifetime and reusability of the sterilized devices depend importantly on a careful pretreatment (cleaning/disinfection). Therefore we recommend not to reuse the devices more than 100 times.

**Prior to each use** the products need to be checked carefully for damage (e.g. material cracks). The connectors must be securely attached to the silicone tubing. **Damaged devices must not be reused!** Swivelling connectors have to be tested for easy rotation. For **safe disconnection** (e.g. for swivel connectors), VBM offers the **uncoupling wedge** (REF 60-60-999).