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
## 1 General instructions

These operation instructions are valid for all medical pressure regulators of the type "Exact Ergo". Medical pressure regulators are manufactured with great care. Improper handling, e.g. the use of force or dropping the device can lead to malfunctions of the pressure display or at the flow adjustment.

### Particular care is required when used in conjunction with oxygen!

For this reason, use of this device requires a thorough understanding of and strict adherence to the operating instructions.

Notify our customer service team of any malfunctions that may occur that you are unable to solve on your own.

(see: "Eliminating possible malfunctions".)

## 2 Intended use

The pressure regulator is intended to administer medical gases from gas cylinders during the treatment and care of patients, independently of a fixed gas supply system.

It is connected to the outlet connection of a cylinder valve and reduces the cylinder pressure to a low, constant value. This constant low pressure (back pressure) enables the provision of a defined flow that is largely independent of the cylinder pressure.

The Exact Ergo pressure regulator is suitable for use by emergency services, in clinics and hospitals as well as in the area of home care - **not** for use in aseptic, sterile areas!

## 3 Product description

This medical pressure regulator is designed for a cylinder pressure of up to 300 bar (see type plate (7)). It is equipped with an inlet connection (1) that is specific to the type of gas, via which the connection to the cylinder valve is created. The current cylinder pressure can be read from the pressure gauge (2) as soon as the cylinder valve is opened. The pressure gauge rotates and can therefore be adjusted for optimum readability depending on the situation.

The outlet provides a fixed flow (e.g.: 5 l/min). It is equipped with a fixed hose nipple (4a) or a threaded connection (4b).

The pressure regulator is **not** suitable for use in driving surgical tools!

Each pressure regulator is equipped with a pressure relief valve (safety valve) (3). This is pre-adjusted at the factory and must not be altered. Manipulation can affect the flow rate accuracy or lead to the uncontrolled escape of gas.

**The pressure regulator, as well as all equipment coming into contact with the gas, must be kept absolutely free of oil and grease!**

**There is a risk of fire or explosion in the case of oxygen, oxygen mixtures and nitrous oxide!**

## 4 Technical data

Gases (depending on version – see type plate)	Oxygen 200 bar, oxygen 300 bar, nitrous oxide, medical air, all medical gases and mixtures
Inlet connection in accordance with ISO 5145 or national standards	Country- and gas-specific
Nominal Inlet pressure (cylinder pressure) P1 at 15 °C	Up to 300 bar (gas-specific)
Outlet connection	Fixed hose nipple or 9/16“ UNF (DISS) + hose nipple adaptor
Nominal outlet pressure	1.5 bar +/-10%
Flow Accuracy	See marking near the outlet connection +/- 10%
Pressure relief valve opening pressure	approx. 7.5 bar (first regulator stage)
Temperature range – operation/storage	-20 °C to +60 °C / -20 °C to +70 °C
Service life	15 years, subject to the recommended inspection and maintenance periods
Maintenance	Technical inspection after 5 years Maintenance after 10 years (see inspection plate)
Checkup	Every 6 months
Weight	Approx.670 g
Metallic materials	Housing: Brass, matt-chromed Pistons: Brass Regulating springs: Stainless steel No aluminium and no stainless steel in gas wetted areas.
Plastic parts	Polyamide, Pressure gauge protection: Rubber Contains no halogenated plastic parts
Elastomeres	EPDM, Silicone
Filter	Sintered bronze

Complies with / approved according to:

- DIN EN ISO 10524-1 (and referenced standards)
- ASTM G175-3

Complies with the essential requirements of the European Medical Device Directive

Equipment class IIb  
93/42/EEC

**CE 0482**

## 5 Identification

Every pressure regulator displays a type plate (7) with the following information:

Manufacturer: VTI  
Type designation: Exact Ergo  
Service pressure: e.g. 300 bar (depending on gas)  
Gas type: O<sub>2</sub> (or other gas)  
Serial number (6-digit): XXXXXX  
EC conformity mark: CE 0482  
Flow outlets are labelled with the litre capacity in l/min.

## 6 Operation

**Attention: Before connecting the gas cylinder, you must wash your hands.  
Do not use hand cream or fatty ointments!**

### • Connecting the pressure regulator to the cylinder valve

Secure the cylinder so that it cannot fall or tip over. Unscrew the protective cap on the gas cylinder. Remove the pressure regulator from the packaging. Check first of all:

- Is the seal present on the connection (1) of the pressure regulator?
- Is the connection of the pressure regulator and its counterpart on the cylinder valve free of dirt and contamination?

The connections on the cylinder valve and on the pressure regulator are standardised and specific to each gas type. As a result, they cannot be confused. Nevertheless, you must check:

- Does the gas cylinder actually contain the correct gas (label/sticker on cylinder shoulder)?
- Is the pressure regulator designed for this gas (see type plate)?

By means of the screw connection (1), screw the pressure regulator by hand clockwise onto the thread of the cylinder valve. Screwing on by hand is sufficient to create a gas-tight connection between the pressure regulator and the cylinder valve.

**Do not use any tools!**

Pressure gauges (2) must not be used as levers! Never use the pressure gauge to tighten/loosen the pressure regulator from the cylinder!

### • Commissioning

Make certain that the equipment to be connected is suitable for medical purposes, for use with medical gases and for the pressure range. Make sure that hose fittings and screwed connections are securely seated and pressure-tight (is the seal present and intact?).

Thin hose lines may be bent if they are not routed carefully. Make sure that they are lying free and that they cannot be subjected to mechanical stress or crushed!

Open the cylinder valve **slowly** with at least one turn.

The pressure gauge (2) displays the available pressure in the cylinder in bar.

If the pointer is in the red zone, the supply can only be guaranteed for a limited period. In this case, the cylinder should be replaced as soon as possible. For information, refer to section , "Duration of Gas Cylinders".

If the cylinder valve is opened quickly, this can lead to a brief pressure surge which is bled off via the pressure relief valve (3). The noise created as a result does not indicate a fault. The pressure regulator is operating normally.

**The uncontrolled discharge of oxygen in confined spaces can lead to a risk of fire!  
Spontaneous combustion may be caused if high concentrations of gas accumulate in  
bedspreads, items of clothing etc!**

**Make sure that all threaded and hose connections are securely and tightly seated and that  
respirator masks, nasal canulas etc. are correctly fitted.**

**If possible, make patients aware of the potential hazards!**

**All medical gases – including oxygen – are medicines!  
For this reason, medical prescriptions and the instructions of nursing staff and paramedics must be observed!**

#### • Turning off the pressure regulator and connecting it to a new gas cylinder

As soon as the pointer on the pressure gauge (2) enters the red zone and continued supply is required, you should prepare to change the cylinder. Always leave some residual pressure in the cylinder. This helps prevent possible contamination of the cylinder.

Also read the section "Duration of Gas Cylinders" for information.


Close the valve on the gas cylinder by hand. Allow the remaining gas to escape from the pressure regulator via the outlet until the pressure gauge (2) indicates "0". The pressure regulator can now be released from the valve by turning the screw connection (1) anti-clockwise and, if necessary, attached to a new cylinder.

Each time the cylinder is changed, check that the pressure regulator is functioning (read the section "Checkup")!

**Do not set up gas cylinders in the immediate vicinity of sources of heat (e.g. heaters, ovens)!**

**Naked flames and smoking are strictly forbidden in the vicinity of valves carrying gas!**

**Do not use external components, e.g. adhesive plaster, to seal the inlets and outlets!**

**Always secure cylinders against falling over and uncontrolled movement!**

**Attention: Liquid gas cylinders (e.g. nitrous oxide) must only be used in the upright position!**

### 7 Cleaning

The device must be cleaned regularly. Cleaning agents must be used in moderation and must **not** enter the inlets or outlets.

**Under no circumstances use agents containing oil, grease or wax! Use a clean, dry or damp cloth.**

### 8 Checkup

Even devices not in use must be checked at least every six months as follows:

- Are there any signs of external damage?
- Is the O-ring for the inlet connection present and intact?
- Does the pressure gauge display a realistic pressure when the pressure regulator is connected and the cylinder valve open?

- Does gas flow out of the flow outlet?

If necessary, inform the customer service of any faults that occur!

### 9 Maintenance

The pressure regulator must be inspected and maintained at regular intervals by trained personnel. We recommend an inspection every five years and a maintenance 10 years after manufacture – or sooner in case of visible defects.

The interval specified on the inspection plate (6) must also be observed to devices in storage!

Only genuine VTI spare parts and lubricating agents specified by VTI may be used!

### 10 Storage and transport

The pressure regulator must always be stored in a clean and dry environment!

Take particular care to ensure that the inlet connection is protected against contamination and damage. Suitable transport protection, e.g. the original packaging, must be used when the devices are transported!

## 11 Disposal

Metal and plastic parts left over following maintenance and repair work must be properly disposed off. If the pressure regulator is no longer used, it is disassembled by the technical department responsible for carrying out maintenance or by the manufacturer and sent for appropriate disposal.

## 12 Warranty

For the pressure regulator supplied, we offer a warranty for a period in accordance with our General Conditions of Sale and Delivery, provided that the operating instructions are adhered to and that the device is not subject to improper external interference. The warranty is limited to error-free functionality, impermeability and flawless material. All other claims are excluded, unless otherwise stipulated by mandatory legal provisions.

Also excluded are defects that are not attributable to our product, but to external influences, improper handling or the failure to adhere to our operating instructions.

## 13 Eliminating possible malfunctions

Fault	Possible cause	Remedy
Pressure gauge (2) does not display a reading	Cylinder valve not opened	Open the valve
	Cylinder is empty	Connect a new cylinder
	Pressure gauge faulty	Inform customer service
Safety valve (3) is activated (allows gas to escape continuously)	Malfunction in regulator system	Inform customer service
Leak at the screw connection between the pressure regulator and valve	Seal missing or damaged	Inform customer service
Screw connection to the cylinder valve cannot be loosened by hand	Residual pressure between valve and pressure regulator	See the section "Turning off the pressure regulator and connecting it to a new gas cylinder"

## 14 Spare parts


Only use genuine spare parts!

O-rings for high-pressure connection (connection to the cylinder valve)

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Your specialist retailer/customer service

## 16 Duration of gas cylinders (hours:minutes)

Does not apply to liquid bottled gases, e.g. carbon dioxide, nitrous oxide!

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**Zertifikat – Nr.: QS – 1332**

We declare under sole responsibility, that  
the a.m. products are in accordance with  
the requirements of directive 93/42/EWG.

The management system is certified by the  
notified body

**MEDCERT GmbH**  
**Pilatuspool 2**  
**20355 Hamburg**

according to the **DIN EN ISO 13485:2003**  
and **93/42/EWG annex II**.

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This certificate is valid for all products with part numbers and revision status of the product  
documentation named above, manufactured until June, 2019.

Director Business Unit MEDcontrol

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VTI Ventil Technik GmbH
HRB 5075 AG Arnsberg
Commerzbank AG, Filiale Menden