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Produkt „Made in Germany“

These instructions must be read carefully before every use of the products. They must be provided for informational purposes to anyone entrusted with the care of patients.

Note!

Intended for use with one single patient only.

These instructions apply to the following product models:

All of the product variants are available in the models Basic, Standard, and Plus. The Basic model does not include inner cannulas; the Standard and Plus models are each supplied with 2 inner cannulas. Additional accessory items are supplied with the Plus model (see the item Specifications of the Product Models).

	Basic REF	Standard REF	Plus REF
Optima® Basic	102213-102217	102013-102017	102113-102117
Optima® Basic - Cuff	102003-102007	102023-102027	102123-102127
Optima® Voice	102233-102237	102033-102037	102133-102137
Optima® Voice - Cuff	102243-102247	102043-102047	102143-102147
Optima® Suction	102253-102257	102053-102057	102153-102157
Optima® Suction Voice	102263-102267	102063-102067	102163-102167

Product Description

The *Optima*® tracheostomy tubes without/with low-pressure cuff are classified in Class IIb in accordance with Annex IX of the European Directive regarding medical devices 93/42 EEC Medical Devices and are always delivered in sterile (EO) condition. Tracheostomy tubes without low-pressure cuff can be delivered in non-sterile condition on request of specific customers.

The characteristic feature for all models of the *Optima*® tracheostomy tubes without/with low-pressure cuff is the transparent PVC outer tube which does not contain DEHP.

The 15-mm universal connector is permanently attached to the cannula outer tube shield. Translucent inner cannulas click into the 15-mm universal connector and can easily be grasped for removal by two slots on the universal connector. The geometry and material of the shield allow it to adapt to the patient's anatomical features and fit snugly against the contours of the neck, assuring good movability and a high degree of comfort when wearing the device. The cannula shield has lateral eyes for fixing the tracheostomy tube in place with a tube holder.

The tubes for inflating the low-pressure cuff (*Optima*® tracheostomy tubes with low-pressure cuff) as well as the tubes for the models with suction are embedded in a molded depression in the outer tube so that they can be positioned securely. They emerge laterally through the shield on the side away from the patient. The very thin walls and large volume of the low-pressure cuff in the models

with low-pressure cuffs ensure dependable sealing when correctly inflated (blocking). The low-pressure cuff seals the tube against the walls of the trachea and can be inflated with air like a balloon. As only a low level of pressure is required to ensure a good seal, the low-pressure cuff serves as a low-pressure cuff and the pressure imposed on the tracheal mucous membrane can be kept to a minimum after inflation. When the air is removed, the balloon collapses closely to the outer tube thanks to its thin walls. There is less irritation of the mucous membrane and the cannulas can be inserted and removed more easily.

A control balloon with a check valve is connected directly to the low-pressure cuff via a tube embedded in an inflation channel of the outer tube. This control balloon is used to inflate/deflate the low-pressure cuff and to check the applied low-pressure cuff pressure.

The Standard and Plus models of the *Optima*® tracheostomy tubes are always provided with 2 inner cannulas which can easily be removed as required, e.g., for cleaning or in the event of respiratory distress. The fenestrated inner cannula in the cannula models with a voice function (Voice) is colored red. Since the inner cannula is visible on the side away from the patient even when it has been inserted, the user and/or nursing personnel are given an immediately visible signal that a fenestrated inner cannula is in use.

The standardized 15-mm universal connector assures the secure connection of accessories with a 15-mm connection such as HME, speaking valve, or ventilation machines to all *Optima*® tracheostomy tubes. The inner cannulas are mutually interchangeable (depending on size). All of the tracheostomy tubes include an obturator of the appropriate size and an atraumatic tip.

The tracheostomy tubes are of high-quality manufacture. All of the materials used are free of DEHP and latex. The careful selection of materials assures optimal characteristics without any risk of collapse when in use by adjustment of the cannula tube at body temperature. These cannulas feature low weight and simple use. All of the materials used in their manufacture are biocompatible and meet applicable legal standards.


The patient's anatomy and particular illness must be carefully considered when selecting the cannula to ensure optimal fit and the best possible air flow. The treating physician should select the product. Special product information, indications, and contraindications in the instructions must be taken into account.

Sizes/Lengths

The sizes of the tracheostomy tubes correspond to the European standard DIN EN ISO 5366-1 (equivalent ISO 5366-1:2000).

The designated sizes refer to the inside diameter in mm of the outer cannula at the tip of the cannula.

The necessary information about a specific product can be viewed on the packaging and on the cannula shield.

LB — length measured along the outer arc
L — length measured along the center line

Size	ID Outer cannula	OD	IDI Inner cannula	L	LB	Low-pressure Cuff
6	6.0	9.2	4.4	64.5	73	20
7	7.0	10.5	5.4	70.0	79	22
8	8.0	11.6	6.4	75.5	87	26
9	9.0	12.6	7.4	81.0	93	30
10	10.0	13.6	8.4	87.5	100	30

All measurements shown in mm

Specifications of the Product Models

1. Tracheostomy tubes without Voice Function:

Optima® Basic

Comprising an outer cannula with a 15-mm universal connector, without low-pressure cuff, 2 inner cannulas (one radiopaque),¹ obturator, neck strap, shower cap² and cough cap².

Size	Basic Model REF	Standard Model REF	Plus Model REF
6	102213	102013	102113
7	102214	102014	102114
8	102215	102015	102115
9	102216	102016	102116
10	102217	102017	102117

Optima® Basic Cuff

Comprising an outer cannula with a 15-mm universal connector, with low-pressure cuff, 2 inner cannulas (one radiopaque),¹ obturator, neck strap, shower cap² and cough cap².

Size	Basic Model REF	Standard Model REF	Plus Model REF
6	102003	102023	102123
7	102004	102024	102124
8	102005	102025	102125
9	102006	102026	102126
10	102007	102027	102127

¹ Included only with the Standard and Plus models

² Included only with Plus model

Optima® Suction

Comprising an outer cannula with a 15-mm universal connector, with low-pressure cuff, 2 inner cannulas (one radiopaque),¹ obturator, neck strap, shower cap,² cough cap,² integrated suction line and separate suction funnel.

Size	Basic Model REF	Standard Model REF	Plus Model REF
6	102253	102053	102153
7	102254	102054	102154
8	102255	102055	102155
9	102256	102056	102156
10	102257	102057	102157

2. Tracheostomy tubes with Voice Function:

Tracheostomized patients fitted with these cannulas are able to speak.

The tracheostomy tube with voice function must not be used for laryngectomized patients.

Optima® Voice

Comprising a fenestrated outer cannula with 15-mm universal connector, without low-pressure cuff, 2 inner cannulas (one radiopaque and one fenestrated),¹ obturator, neck strap, speaking valve, shower cap,² cough cap² and decannulation cap.

Size	Basic Model REF	Standard Model REF	Plus Model REF
6	102233	102033	102133
7	102234	102034	102134
8	102235	102035	102135
9	102236	102036	102136
10	102237	102037	102137

¹ Included only with the Standard and Plus models

² Included only with Plus model

Optima® Voice Cuff

Comprising a fenestrated outer cannula with 15-mm universal connector, with low-pressure cuff, 2 inner cannulas (one radiopaque and one fenestrated),¹ obturator, neck strap, speaking valve, shower cap,² cough cap² and decannulation cap.

Size	Basic Model REF	Standard Model REF	Plus Model REF
6	102243	102043	102143
7	102244	102044	102144
8	102245	102045	102145
9	102246	102046	102146
10	102247	102047	102147

Optima® Suction Voice

Comprising a fenestrated outer cannula with 15-mm universal connector, with low-pressure cuff, 2 inner cannulas (one radiopaque and one fenestrated),¹ obturator, neck strap, speaking valve, shower cap,² cough cap,² decannulation cap, integrated suction line and separate suction funnel

Size	Basic Model REF	Standard Model REF	Plus Model REF
6	102263	102063	102163
7	102264	102064	102164
8	102265	102065	102165
9	102266	102066	102166
10	102267	102067	102167

NOTE!

The shower cap is not suitable for bathing and swimming!

Accessories

Inner Cannula Size 6 transparent	REF 102018
Inner Cannula Size 7 transparent	REF 102019
Inner Cannula Size 8 transparent	REF 102020
Inner Cannula Size 9 transparent	REF 102021
Inner Cannula Size 10 transparent	REF 102022

¹ Included only with the Standard and Plus models

² Included only with Plus model

Inner Cannula Size 6 red, fenestrated	REF 102028
Inner Cannula Size 7 red, fenestrated	REF 102029
Inner Cannula Size 8 red, fenestrated	REF 102030
Inner Cannula Size 9 red, fenestrated	REF 102031
Inner Cannula Size 10 red, fenestrated	REF 102032

Inner Cannula Size 6 white, radiopaque	REF 102038
Inner Cannula Size 7 white, radiopaque	REF 102039
Inner Cannula Size 8 white, radiopaque	REF 102040
Inner Cannula Size 9 white, radiopaque	REF 102041
Inner Cannula Size 10 white, radiopaque	REF 102042

All inner cannulas are manufactured free of any DEHP!

Shower cap	REF 102080
Cough cap	REF 102081
Decannulation cap	REF 102082
PRIMEDI speaking valve	REF 20489

Indication

The treating physician is always responsible for the selection of the cannulas best suited to the patient.

Optima® tracheostomy tubes can be used following any tracheostomies, regardless of the cause, and serve to stabilize/hold open the tracheostoma. They are inserted through the tracheostoma into the trachea and held in place on the throat by a tube holder.

Optima® tracheostomy tubes with voice function have been designed for tracheostomized patients whose larynx is intact.

Optima® Voice, *Optima*® Voice Cuff and *Optima*® Suction Voice are used after tracheostomies when the larynx remains completely or partially intact and enable the user to speak.

NOTE!

Observe the instructions regarding cannulas with voice function (Voice)!

The tracheostomy tubes *Optima*® Voice, *Optima*® Voice Cuff, and *Optima*® Suction Voice can be outfitted with a speaking valve with a matching connection via the integrated, standardized 15-mm universal connector on the outer cannula (**NOTE: follow manufacturer's instructions!**). These cannula models are provided with fenestration/screen (speech field) as standard. The speech field of the outer cannula has a number of small holes, the sieving. One of the two included inner cannulas has an oval opening (fenestration) reaching over the speech field of the outer cannula so that air can flow freely for speaking when the patient breaths out. The Basic model of the Voice tubes does not have an inner cannula of this type.

The use of the *Optima*® tracheostomy tubes with low-pressure Cuff following tracheostomies is always indicated whenever a seal between the tracheal wall and cannula is required.

These tubes must not be used during MRI/MRT procedures because of the check valve in the control balloon of the low-pressure cuff inflation system.

Do not attempt to reposition a tracheostomy tube with inflated low-pressure cuff while it is being worn! Risk of injury!

NOTE!

Optima® Voice, *Optima*® Voice Cuff, and *Optima*® Suction Voice are recommended only for tracheostomized patients with normal secretions and normal mucous membrane. These cannula models are recommended for patients with heavy secretions, a tendency to granulation tissue, and clogging only if they are regularly examined by a physician and the intervals between replacement are shorter (weekly, as a rule) because the screen in the outer tube can intensify the formation of granulation tissue.

Follow carefully the indication instructions for the specific cannula models.

Solely the unfenestrated inner cannula included with every cannula model with low-pressure cuff (except the Basic model) must be used for patients who are being artificially ventilated.

NOTE!

In the case of mechanical ventilation, fenestrated outer cannulas must be used only together with unfenestrated inner cannulas.

NOTE!

Optima® Voice, *Optima*® Voice Cuff, and *Optima*® Suction Voice must never be used in patients with recurrent aspiration pneumonia because of the increased risk of aspiration.

Optima® Voice, *Optima*® Voice Cuff, and *Optima*® Suction Voice are recommended only for tracheostomized patients with normal secretions and normal mucous membrane.

Only the manufacturer is permitted to shorten lengths, add screens, and otherwise modify the cannulas as patients may otherwise be put at risk.

Any work on tracheostomy tubes which is not performed properly can lead to serious injury! Any and all liability claims are excluded if modifications are carried out by unauthorized persons.

A decannulation cap is included with the voice cannulas for use in preparing a possible decannulation after a temporary tracheostomy. This temporarily stops the air flow through the cannula so that the patient becomes re-accustomed to breathing through mouth/nose. The decannulation and the discontinuation of breathing through the tracheostomy tube when it is being worn must be done solely under a physician's supervision and in accordance with the physician's instructions. **RISK OF SUFFOCATION!**

Tracheostomy tubes with suction function (*Optima*® Suction, *Optima*® Suction Voice)

These models with low-pressure cuff allow suctioning of the tracheostomy tube while being worn without requiring additional suction tubes or suction catheters. The suction takes place immediately above the blocked low-pressure cuff using the suction system permanently fixed to the cannula. The suction funnel enclosed separately is connected to the cannula using a Luer connection before suction and can then be connected to a suction device outfitted with a vacuum regulator (**NOTE: follow manufacturer's instructions!**). The suction funnel is connected with the Luer connector on the tube connector/fingertip of the device suction tube.

The suction pressure must not exceed -0.2 bar (-20 kPa). The suction can also be performed by attaching a syringe to the Luer connector. The procedure must be determined in consultation with the treating physician so that the circumstances of the specific illness are given consideration and risks can be avoided. A procedure which is as optimal and gentle as possible must be chosen.

The design of the Suction model also allows the conduct of a sub-glottal flushing (e.g., with saline solution, NaCl 0.9%, **exact dosage according to physician's instructions**). This allows a reduction of mucous membrane irritation or the softening of clogging. During the flushing procedure, the usual suction process via the integrated suction system above the low-pressure cuff is reversed and the solution is introduced into the trachea. At the conclusion of the procedure, the solution must be removed using a syringe or suction device as described above.

NOTE!

The tightness of the seal and the adequacy of the block must be checked before the flushing procedure is carried out. Be sure to prevent any aspiration of the solution. Suction off any secretions above the low-pressure cuff before beginning the flushing procedure.

The design of the *Optima*® Suction Voice model combines the advantages of a voice model combined as a unit with the integrated suction function.

NOTE!

The cannula model *Optima*® Suction Voice may be used solely for patients who have only a slight tendency to the formation of secretions. Only the Suction model may be used with patients who have a tendency to form heavy secretions.

Note:

Cannulas with Suction and Suction Voice function may be used in patients with an increased risk of bleeding (e.g., when taking anti-coagulants) solely in consultation with the physician because of the greater risks during suction procedures.

Special Models

Special models taking into account anatomical circumstances are technically not possible and are not offered.

Contraindications

The *Optima*® Basic, *Optima*® Basic Cuff, and *Optima*® Suction models of the tracheostomy tubes have been designed for use exclusively in patients with an intact larynx.

Make sure that the cannula is fixed firmly in place in patients who are being artificially ventilated because a disconnection puts them at risk of death!

NOTE!

The decannulation cap may be used solely with tracheostomized patients with an intact larynx. It may be used only under a physician's supervision.

NOTE!

Do not use the tracheostomy tube when conducting medical laser therapy. There is a significant risk of injury. There is also the risk of tube burn when electrocautery is used during operations on the trachea. Contact can lead to emission of toxic gases, burns, and/or harm to health.

Observe carefully the description of contraindications in the relevant product descriptions/specifications!

NOTE!

Do not insert tracheostomy tubes in patients with abnormal anatomy or pathology of the upper airways because there is a risk of partial or total blockage of the airways.

Scope of Delivery

A cannula pass is enclosed with all of the tracheostomy tubes in addition to the contents of the packaging described above. All of the necessary data of the tracheostomy tube are shown on the self-adhesive label and are provided for insertion in the patient pass.

Cannula Pass

The cannula pass is enclosed in every package containing a tracheostomy tube. The pass, when the provided label is glued into it, contains information about:

- Identification number
- Size
- Length
- Lot number

Data about the patient and the treating physician can be entered on the reverse side.

The cannula pass should be readily available so that the required information about the cannula can be easily obtained for orders or subsequent questions.

Insertion of a Tracheostomy tube

Preparation

Exercise great diligence and care before and during the insertion of the tracheostomy tubes. When removing a cannula from the sterile packaging, check the date for use shown on the package and the integrity of the sterile packaging to ensure sterility of the cannula. The use of sterile disposable gloves is recommended.

The treating physician selects the size and model of the tracheostomy tube which is to be used. The correctness of the cannula should be checked by comparing the information on the label and the cannula shield.

Prepare the tracheostomy tube before insertion. This includes the inspection of the entire cannula for external damage and loose parts so that the possibility of injury is excluded and the required sealing is assured. Never use the cannula if there is anything unusual about it.

Checking the Airtightness of the Seal in the *Optima*® Basic Cuff, *Optima*® Voice Cuff, *Optima*® Suction and *Optima*® Suction Voice Tracheostomy tube models

These cannulas with a low-pressure cuff must first be tested for airtightness as described below.

The low-pressure cuff must be completely deflated before insertion of the cannula!

The airtightness of the cannula and low-pressure cuff should be checked before and after every insertion.

Inflate the low-pressure cuff with air up to a pressure of 22 to 24 mm Hg (~ 30 to 33 cm H₂O or 29 to 32 hPa). If cannula and low-pressure cuff are airtight, there will not be any spontaneous, significant loss of pressure. A low-pressure cuff pressure measuring device (**NOTE: follow manufacturer's instructions!**) is recommended for the test.

The following signs during/after cleaning or in inserted cannulas may appear as further indications that cannula and/or low-pressure cuff are not airtight:

- Water in the input tubes to the cannula
- Water in the low-pressure cuff
- Water in the control balloon
- No inducing of coughing when pressure is applied to the control balloon

NOTE!

If one or more of the above-mentioned signs appear indicating that the cannula is not airtight, use of this cannula must not be continued under any circumstances. The proper functioning of this cannula is no longer assured!

Tracheostomy tubes with low-pressure cuff which display the above-mentioned defects or anomalies must be returned to the manufacturer for inspection.

Insertion of the Tracheostomy tube

If the inspection determines that the cannula is in good condition, the tracheostomy tube can be inserted. We recommend drawing a cannula protective bib (compress) over the cannula shaft to protect the tracheostoma and to keep the skin and clothing dry.

We recommend applying stoma oil to the outer cannula to reduce the friction on the tracheostomy tube and simplify its insertion into the trachea. This will also reduce the patient's urge to cough. When using water-soluble lubricants, avoid any clogging/blockage of the cannula (**NOTE: follow manufacturer's instructions!**).

We recommend tilting the head back slightly when inserting the tracheostomy tubes. Hold the cannula by the shield. Use the other hand to spread open the tracheostoma slightly if needed, then carefully insert the end of the cannula in to the stoma opening, axially at first, during the inspiration phase (when the patient is breathing in). Special tools (trachea spreaders) are available to spread open the tracheostoma (**NOTE: follow manufacturer's instructions!**) which enable an even and gentle spreading of the tracheostoma, e.g., including emergency cases when the tracheostoma collapses. When using any kind of tool, be sure that the cannula, especially the low-pressure cuff, is not damaged by rubbing against the tool or sharp edges. While pushing the cannula further into the opening, the patient simultaneously moves his/her head into its normal position and pushes the cannula, following the curve of the cannula, into the trachea. If the patient is inserting the cannula him-/herself, he/she should use a mirror to simplify the insertion of the tracheostomy tube. He/she can observe every action.

NOTE!

Never insert a cannula unless the low-pressure cuff is completely deflated! The cannula shield must not be pushed into the tracheostoma. Make sure that the cannula shield is always outside of the tracheostoma.

When inserting and using the models with low-pressure cuff and suction, make sure that no damage occurs to the low-pressure cuff and the corresponding tube or control balloon and that no damage can occur to the suction system.

Tracheostomy tubes should always be held firmly in place with a special tube holder (one tube holder is always included). The holder stabilizes the cannula and ensures its secure fit in the tracheostoma. It is attached to the lateral eyes on the cannula shield using Velcro fasteners, is adjustable, and is comfortable to wear. When cannulas with low-pressure cuff and/or suction are being used, we recommend fixing the inflation tube for the low-pressure cuff inflation and, if present, the suction tube in place so that the connections are not subject to tension. This will maintain the patient's mobility.

If needed and according to model, commercially available filter and valve systems can be adapted using the 15-mm universal connector on the outer cannula.

NOTE!

If a ventilation system is attached to the models with low-pressure cuff using the 15-mm universal connector, make sure that any significant tension or twisting of this connection is avoided. There is a risk of disconnection or obstruction of the tracheostomy tube. The patient

is at risk of death if there is a disconnection! Never use screened/fenestrated tracheostomy tubes with low-pressure cuff in patients who are being artificially ventilated.

Optima® Basic Cuff, Optima® Voice Cuff, Optima® Suction, and Optima® Suction Voice Tracheostomy tube Models

Inflating the Low-pressure Cuff

Once the cannula has been inserted and fixed in place, a low-pressure cuff pressure meter is used to inflate the low-pressure cuff with air. The inflation pressure must be adjusted for the specific requirements and is determined by the treating physician. Unless otherwise instructed by the physician, we recommend a low-pressure cuff pressure of no less than 16 to 18 mmHg / 22 to 24 cmH₂O/21 to 24 hPa.

The pressure should be checked at regular intervals, approx. every two hours.

NOTE!

The pressure used to inflate the low-pressure cuff should never be greater than 25 mm Hg (approx. 34 cm H₂O or 33 hPa) to exclude the possibility of injury and/or blockage of the airways.

The low-pressure cuff must be undamaged and function absolutely correctly. If the desired seal is not achieved with the maximum inflation pressure given above even after repeated attempts, it may be necessary to use a cannula with a larger diameter. The correct low-pressure cuff pressure must be checked regularly, i.e., at least every two hours.

Removal of a Tracheostomy tube

NOTE!

The low-pressure cuff model may not be removed from the tracheostoma until the low-pressure cuff has been completely deflated; otherwise, the mucous membranes of the trachea may be injured.

Before deflating the low-pressure cuff using a syringe, any secretions which have collected above the low-pressure cuff should be suctioned off to prevent choking/aspiration. After removing the inner cannula, insert a suction catheter/suction tube through the cannula tube to below the cannula tip while the tracheostomy tube with low-pressure cuff is still in place. Any secretions which are present can be suctioned off by simultaneously applying suction and deflating the low-pressure cuff. This prevents any aspiration of the secretions. When the Suction model of the cannula is used, the secretion is suctioned off using the suction system integrated in the cannula. After the suction has been completed, the suction tube in the cannula should be closed using the provided cap so that the risk of infection is minimized.

The patient's head should be tilted back slightly when the cannula is removed. Extreme care must be taken so that the mucous membrane is not injured. Grasp the cannula on the sides of the shield for removal.

NOTE!

If the tracheostoma is unstable or in emergency situations (puncture, dilatation tracheostoma), the tracheostoma may collapse when the cannula is removed, causing difficulty in breathing. In this case, a new cannula must be available for immediate use. A trachea spreader (NOTE: follow manufacturer's instructions!) is useful for temporarily assuring the air flow.

The inner cannula can become clogged/closed by liquids/secretions. In this case, the inner cannula must be replaced and/or cleaned as necessary. After suctioning off any liquid/secretions which have collected above the low-pressure cuff, open the lock between inner and outer cannula by lifting it slightly so that the inner cannula is released. Before removing the inner cannula of cannulas with voice function, any speaking valve or filter system being used must be disconnected. After it has been cleaned, the inner cannula is returned by reversing the procedure. Note that the inner cannula must always be examined for damage and cleaned and/or disinfected according to the following instructions before being re-inserted. The lubrication of the inner cannula with stoma oil ensures low friction between inner and outer cannula. Stoma oil is absolutely harmless.

NOTE!

Do not twist the inner cannula along its axis when removing it as it may otherwise be damaged. Do not use force when removing inner cannulas. If the inner cannula cannot be released, remove the entire cannula.

Usually the outer cannula remains in the trachea when an inner cannula is changed.

The inner cannula is returned by reversing the procedure described above.

Cleaning

Optima® tracheostomy tubes must be cleaned regularly for hygienic reasons and to ward off the risk of infection. We recommend cleaning the outer cannulas at least once a day and cleaning the inner cannulas thoroughly at least twice a day. The cannulas can also be cleaned several times daily if there has been heavy formation of secretions or for other hygienic reasons.

Always follow the instructions below exactly!

Cleaning the Inner Cannulas

Use warm water and a mild, pH-neutral washing lotion to clean the cannula. It is important to rinse the inner cannula thoroughly with clear water several times after washing.

NOTE!

If special cannula cleaning agents for tracheostomy tube are used, follow the instructions of the agent manufacturer. The agent must be suitable and approved for the treatment of plastic cannulas.

The use of a cannula soaking container with a screen insert simplifies handling.

The screen insert holding the inner cannulas (outer and inner cannulas next to each other for cannulas without a low-pressure cuff) is immersed in the prepared cleaning solution. Hold the screen insert by the upper edge so that there is no contact with the cleaning solution. After expiration of the soaking period (**see instructions for using the cleaning agent**), the inner cannula must be rinsed thoroughly a number of times with clear, lukewarm water.

If necessary (e.g., if there are stubborn and tough secretion residues which cannot be removed by the cleaning bath), additional cleaning using a special cannula cleaning brush is possible. This brush must also be used in accordance with instructions. When using a brush, take special care so that the cannula material is not damaged. Always insert the cannula cleaning brush into the cannula from the tip of the cannula. Always remove speaking valves beforehand and clean them according to manufacturer's instructions because they can otherwise be damaged.

The cannulas must subsequently be rinsed thoroughly again under running water.

The cannulas can be thermally disinfected at a minimum temperature of 65° C (149° F) using clean water and commonly available household agents without the addition of any chemical agents. After the cannulas have been placed in the water bath, the temperature must be monitored and maintained (temperature control using commonly available household appliances with adequate regulator). Precautions must be taken to ensure that the temperature does not go significantly higher than the aforementioned level and that the water does not begin to boil. This can cause substantial damage to the tracheostomy tubes.

NOTE!

Do not use a dishwasher, steam cooker, or microwave oven to clean the inner cannulas!

After washing the cannula, dry it thoroughly with a clean, lint-free cloth.

NOTE!

Do not use any cleaning agents which have not been approved by the cannula manufacturer for cleaning tracheostomy tube. The use of caustic household cleaners, high-proof alcohol, hydrogen peroxide, or denture cleaning agents is absolutely prohibited. This would cause an acute health risk! Moreover, the cannula might be destroyed or damaged.

If special cannula cleaning agents for tracheostomy tube are used, follow the instructions of the agent manufacturer. The agent must be suitable and approved for the treatment of cannulas made of medical-quality PVC.

Solutions must not be allowed to overheat/come to a boil during the cleaning. Use only lukewarm water for the preparation of the cleaning solution and follow the instructions for the cleaning agent.

Never use a dishwasher, steam cooker, or microwave oven to clean the tracheostomy tubes!

The following instructions for disinfecting procedures and for care of the tracheostomy tube must be followed exactly to prevent any health risks and to maintain the usability of the cannula.

Disinfection

Special chemical disinfectants may be used for cold disinfection of *Optima*® tracheostomy tubes. However, such a procedure should be used only if the treating physician has instructed its use because of the specific illness or if it is indicated by the particular care situation. As a rule, disinfection is appropriate to avoid cross-infections and during use in in-patient areas (clinic, nursing home) as a means of reducing the risk of infection.

Tracheostomy tubes must be disinfected after completion of the cleaning procedure described above by placing them in a disinfectant (e.g., PRIMEDICLEAN®) suitable for use with the material and the cannulas. Manufacturer's instructions for the use of the disinfectant must be followed exactly. After disinfection, the cannulas must be rinsed very thoroughly with boiled water which has cooled completely or with a sterile saline solution, then dried.

NOTE!

Never use disinfectants which release chlorine or contain strong bases or phenol derivatives. This can cause major damage to, or even destroy, the cannula tube.

After cleaning and any disinfection, the tracheostomy tube can be re-inserted into the tracheostoma. To prevent any risks to the patient's health, all residue of any kind from the cleaning agent/disinfectant must be removed from the cannulas before they are re-inserted into the tracheostoma.

Care

After cleaning, any disinfection, and drying of the cannula, the outer surface of the outer cannulas, with the exception of outer cannulas with low-pressure cuff, should be lubricated by the application of stoma oil. This ensures safe sliding during insertion, and the suppleness of the cannula material of the outer cannula is maintained by regular care.

Sterilization

The *Optima*® tracheostomy tubes are provided in sterile packaging and are single-patient products. Sterility is guaranteed only as long as the unopened package is undamaged and the written shelf life has not been expired. Resterilization is prohibited.

NOTE!

Never use cannulas whose functional performance has been impaired or which have visible damage such as sharp edges or cracks; they can cause injury to the mucous membranes in the trachea.

Damaged cannulas must not be used and must be disposed of.

Regular checks prevent risks!

Useful Life, Storage

Optima® tracheostomy tubes are sterile single-patient products.

The useful life of the *Optima*® tracheostomy tubes is influenced by many different factors. The composition of secretions, the conscientiousness of the cleaning, and other aspects can be of decisive importance.

We recommend replacement of the cannulas after 29 days (from the opening of the packaging).

This wearing period must be observed without exception for *Optima*® tracheostomy tubes with low-pressure cuff.

If *Optima*® tracheostomy tubes without a low-pressure cuff are maintained very carefully, the expected lifetime is extended and the cannulas can be used for more than 29 days. The

Optima® tracheostomy tube without low-pressure cuff must be replaced with a new *Optima*® tracheostomy tube no later than 6 months after its first use. Replaced tracheostomy tubes must be properly disposed of.

We strongly recommend always having at least two replacement cannulas on hand to ensure availability of replacements at all times.

Unused cannulas should be stored in a dry environment and protected from sunlight and/or heat.

Warranty

We warrant the freedom of the products from defects within the scope of our general terms and conditions of business (TOB). You can view a copy of our TOB on our business documents or on our Web site on the Internet (www.primed-halberstadt.de).

Repairs and other work on the tracheostomy tubes may be carried out exclusively by Primed Halberstadt Medizintechnik GmbH. Failure to observe these terms will cause the warranty to expire immediately!

Reuse/Single-patient Product

NOTE!

Tracheostomy tubes are single-patient products and as such are intended exclusively for use with a single patient. Any reuse with other patients and any preparation for reuse are prohibited.

Catalog/Accessories

Primed Halberstadt Medizintechnik GmbH specializes in the manufacture and supply of aids and devices for tracheostomized and laryngectomized patients. A list of the products currently available can be seen in our product catalog, which we will be happy to send to you free of charge, or can be viewed at www.primed-halberstadt.de.

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Products are subject to manufacturer's modifications at any time.

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Optima® is a registered trademark of Primed Halberstadt Medizintechnik GmbH.