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
PRIMEDI-PHON VENT® - Instructions for use

PRIMEDI-PHON VENT®

Speaking valve with respiration connection, rotatable

REF 200320

CE 0482

Content

Read the following details for warnings, precautionary measures and instructions before use!

The PRIMEDI-PHON VENT® speaking valve is individually packaged as a single patient product. The package at hand contains a PRIMEDI-PHON VENT® speaking valve (REF 200320) as well as warning labels to attach to the inlet tube of the control balloon of a tracheal cannula with cuff, and instructions for use.

The content is free of latex and not sterile.

The instruction manual provides information for doctors, health care professionals, patients and their relatives for the ensurance of correct use.

Indications for selecting a product are determined by the attending physician, in consultation with the patient!

Please consult the attending physician if there are questions or complications.

⚠ Caution:

**US Federal law restricts the sale of this product to doctors, or at the direction of a physician.
(Applies only to United States.)**

Description

PRIMEDI-PHON VENT® is a speaking valve with valve diaphragm, that has been developed for patients wearing a tracheostomy tube to enable finger-free, normal speaking, whilst breathing.

It consists of a plastic housing, a lid with valve function, as well as a rotating, 15mm connection on the patient side for connecting to a ventilator. The valve diaphragm is centrally fixed on the lid and, in contrast to an open one-way speaking valve, this membrane closes the valve in the default position. Only by inhalation (inspiration), does the valve open so that air flows into the tracheostomy tube and the lungs. After the inspiration, the diaphragm closes again. Additionally, the increase in pressure at exhalation (expiration) supports a closing of the valve. During the exhalation process, the air around the tracheostomy tube passes through the voice box and throat (larynx and pharynx) to the mouth and nasal cavity. When passing through the vocal cords, speech (phonation) is made possible.

PRIMEDI-PHON VENT® is suitable for use in adults, children and newly-born and can be used in conjunction with the following types of tracheostomy tube:

- Tracheostomy tube without fenestration or sieving
- Tracheostomy tube with sieving in the outer cannula, without fenestration in the inner cannula
- Tracheostomy tube with fenestration or sieving **(only on a doctor's orders)**
- Tracheostomy tube without cuff
- Tracheostomy tube with cuff

① Warning:

when using tracheostomy tube with cuff the cuff must be completely deflated, otherwise the patient cannot breathe!

① Warning:

In every case check beforehand, especially in children and newly born, if sufficient air passageway surrounds the cannula.

Fenestrated/ sieved tracheostomy tube:

PRIMEDI- PHON VENT® can be used with sieved outer cannulas, in which case a non-fenestrated inner cannula is to be used.

△ Warning:

when using PRIMEDI- PHON VENT® in conjunction with a tracheostomy tube, which has an inner cannula with circlip, or similar, which protrudes beyond the 15 mm-connector, it should be examined, as to whether a firm seating of the valve is assured. If not, the inner cannula must be removed prior to the use of PRIMEDI-PHON VENT®.

The plastic housing of the speaking valve has a centre opening with 15 mm connector on the patient's side that enables a connection to the 15 mm connector of a tracheostomy tube (e.g. Optima® and Ultrasoft® tracheostomy tube). A second connector on the patient side, with a diameter of 22 mm, allows also the connecting of tracheostomy tubes with a 22 mm fitting (e.g. appropriate models of a Primedistom® tracheostomy tube).

Through the closed default position of the PRIMEDI- PHON VENT® a column of air forms in the tracheostomy tube, which prevents an ingress of secretions into the cannula and a blocking of the speaking valve. Air and secretions are led up the windpipe (trachea) to mouth and/or nose.

This provides for the restoration of positive subglottic pressure, which eases swallowing and allows effective coughing/ expectoration and may cause a reduction of aspiration. Also the evaporation of oral secretions is favoured as a result since the air passes through the upper respiratory tract during exhalation. The need for removal by suction may thus be reduced.

Because PRIMEDI- PHON VENT® reinstates a normal closed respiratory system, by which a positive airway pressure, with no manual shutter, can be created, a clear voice with normal diction is permitted and which, especially in children, supports normal speaking development.

The sense of smell can be improved through the use of a PRIMEDI- PHON VENT® speaking valve, because the airflow while exhaling again passes through the mouth and nasal cavity. An improved sense of smell can lead to an improved sense of taste, a bigger appetite and a higher calorie intake.

After a little practice, the subglottic pressure receptors are innervated and may cause spontaneous swallowing. As a result swallowing therapy can be effectively supported.

Because a manual or finger closure of the tracheostomy tube is superfluous when using this speaking valve the risk of infection is reduced and tracheal hygiene thereby supported. Also, the valve functions as a filter to prevent the ingress of particles into the trachea.

△ Warning:

While PRIMEDI- PHON VENT® can ease swallowing and in some patients reduce aspiration, the presence or risk of aspiration should be carefully judged individually with each patient, to determine appropriate use of this speech valve with regard to the swallowing function.

⚠ Warnung:

the use of PRIMEDI- PHON VENT® speaking valves must be assessed carefully in patients with severe secretions that may cause blockages of the air passageways.

This speaking valve can be used also as an additional aid in the weaning of patients from artificial respiration. During exhalation, the tight fitting ensures the restoration of the physiological PEEP pressure (positive end expiratory pressure), which can improve oxygen uptake. Increasing familiarisation by exhaling through the upper airway gives the patient more confidence and supports breathing muscle training.

Also the PRIMEDI- PHON VENT® speaking valve can represent an alternative to plugging the tracheostomy tube, particularly in patients who cannot tolerate plugging because of physiological and psychological reasons. If the patient tolerates plugging for only short periods of time, the speaking valve can serve as an aid between the plugging during the transition from an open tracheostomy tube to tracheal plugging. Because it is possible for a patient who has the use of a PRIMEDI- PHON VENT® speaking valve to exhale regularly through the upper respiratory tract he will gain confidence and the doctor will be able to determine the consistency of the air passage.


Fig. Air flows when using a PRIMEDI- PHON VENT® speaking valve

- Inspiration through speaking valve
- Expiration bypassing the tracheostomy tube through the rhino pharynx

Indication

PRIMEDI- PHON VENT® speaking valves are intended for use with tracheostomised and ventilated patients. They are suitable for awake, alert tracheostomised adults, children and newly born, if not otherwise contraindicated.

It is possible to use with ventilated patients.

The airflow around the tracheostomy tube and through the upper airway during exhalation must be sufficient.

⚠ Warning:

for use on one patient only. This product is to be sold and used only for the purpose specified.

⚠ Warning:

patients with a PRIMEDI-PHON VENT® speaking valve must be examined and/or monitored according to the directions of the physician.

⚠ Warning:

not for use with severe blockage of the lung passageways (e.g. tracheal and/or laryngeal stenosis). Take care with end-stage pulmonary disease. Do not use when there is a severe mucus secretion from the lungs.

⚠ Warning:

PRIMEDI- PHON VENT® is not suitable for laryngectomised patients (larynx extirpated).

Indications for the use of a PRIMEDI- PHON VENT® speaking valve (including, but not limited to):

- Quadriplegia
- Chronic obturation pulmonary disease
- Slight tracheal and laryngeal stenosis
- Double-sided paralysis of the vocal cords without significant airway obstruction (recurrence paresis)
- Non-obturate laryngeal tumours (including patients with vocal chord function after resection of the tumour)
- Patients with sleep apnoea, which instead of a plug in the waking state receive a tracheostomy
- Patients who mentally or physically cannot tolerate a tracheal plugging
- Neuromuscular disease
- Head injury
- Pulmonary disease
- Mild tracheomalacia

Contraindications

- Unconscious, sleeping and / or comatose patients
- Laryngectomised (voice box removed) patients
- Inflated tracheostomy tube cuff
- Tracheostomy tubes with foam-filled cuffs
- Severe airway obstruction, which may prevent adequate exhalation (e.g. subglottic stenosis)
- Severe and increased mucus secretion
- Significantly reduced lung elasticity, which can cause an air pocket
- Difficult breathing
- Not suitable for use with endotracheostomy tube.

Use

The speaking valve PRIMEDI- PHON VENT® can be used in patients with portable or conventional ventilators whereby it is connected to the rotating 15 mm-connection.

The following application guidelines should be observed together with the advice of the physician: If the tracheostomy tube has been replaced, the insertion of the PRIMEDI- PHON VENT® speaking valve may need to be delayed 48-72 hours in case a tracheal swelling and/or bronchial spasm have been caused. Universal precautions should be observed.

1. Check before inserting

The following assessment criteria must be observed before the insertion of a PRIMEDI- PHON VENT® speaking valve:

Cognitive condition:

The patient must be awake and responsive, and he must try to communicate. The speaking valve must not be inserted while the patient is sleeping.

⚠ Warning:
do not use while sleeping.

Medical/pulmonary status:

Patient must have the appropriate lung mechanisms to exhale around the tracheostomy tube and through the nose and mouth. The following criteria should be considered when assessing the patient:

- Signs of life
- Patient response
- Air passageway consistency
- The correct position of the patient and the tracheostomy tube
- Oxygen saturation
- Respiratory function
- Breathing sounds
- Psychological condition and motivation of the patient

Deflation of the cuff (if fitted):

When using PRIMEDI- PHON VENT® the air must be let out from the cuff completely, so that the exhaled air can flow around the tracheostomy tube and through the mouth, nose, and throat.

Secretion management:

The use of a PRIMEDI- PHON VENT® speaking valve can facilitate the flow and oral expectoration of secretions by the patient. Excessive amounts, viscosity and an acute infection can affect a secretion treatment. To what extent increased quantities and/or the different viscosity of secretions can be treated depends on the individual patient. It may be that the use of the valve must be temporarily restricted or suspended, until a satisfactory treatment of the mucous secretion is possible.

⚠ Warning:

the use of PRIMEDI- PHON VENT® speaking valves must be assessed carefully in patients with severe mucous secretions, which may cause blockages of the air passageways.

Swallowing:

The patient's risk of aspiration should be evaluated, as it may affect the amount, the thickness and the ability to treat secretions. The presence of a strong aspiration plays an important role under certain circumstances in determining whether the patient is suited for letting out the air from the cuff and the insertion of a PRIMEDI- PHON VENT® speaking valve. The safety and efficacy of the swallowing process can be reduced by the use of a tracheostomy tube. While some tracheostomised individuals exhibit no swallowing difficulties, others suffer from dysphagia and aspiration, although their primary diagnosis would not typically indicate swallowing difficulties. Use of a PRIMEDI-PHON VENT® speaking valve can improve the safety and efficiency of swallowing and may reduce aspiration. The closed default position restores the normal, closed system of the patient. This eases the swallowing process, as the pharyngeal / laryngeal sensation is increased and the positive subglottic air pressure restored.

⚠ Warning:

While PRIMEDI- PHON VENT® can ease swallowing and in some patients reduce aspiration, the presence or risk of aspiration should be carefully and individually assessed with each patient, to determine appropriate use of this speaking valve with regard to the swallowing function.

Air passageway consistency:

The patient should only use a PRIMEDI- PHON VENT® speaking valve, if he is able to exhale effectively around the tracheostomy tube upwards through the larynx and pharynx and out of the nasal and oral cavity.

The diagnosis must be checked to ensure that no known air passageway restrictions are present (e.g., tumours, stenosis, granulation tissue).

The patient's ability to exhale effectively depends on the size of the tracheostomy tube used. Therefore, a size should be selected, which ensures sufficient airflow around the tracheostomy tube and thus eases speaking and the use of a PRIMEDI- PHON VENT® speaking valve.

The tracheostomy tube cuff can also cause an occlusion, even if the air has been let out. It should therefore be included in the assessment of air passageway consistency. A patient who uses a tracheostomy tube with cuff, should be considered for one without a cuff if this is medically appropriate, in order to avoid the necessity of having to let the air out in order to use the speaking valve.

⚠ Warning:

before inserting a PRIMEDI- PHON VENT® speaking valve the air must be completely let out of the cuff of a cuffed tracheostomy tube. Otherwise, the patient cannot breathe.

Do not use tracheostomy tube with a foam-filled cuff.

Patients must always be examined to ensure that the breathing passage is sufficient.

Assessment of the respiratory tract consistency at the bedside:

1. Let the air completely out of the tracheostomy tube cuff (if any).
2. Instruct the patient to breathe through the tracheostomy tube.
3. Close the tracheostomy tube manually with a gloved finger and at the same time instruct the patient to exhale through mouth and nose to ensure an adequate exhalation. This can be checked in that the patient blows onto a handkerchief, a mirror or a feather. The patient should be asked to vocalise (e.g. say „Ah“) to assess the presence and quality of speech. Some patients can perhaps breathe out adequately, but are initially not able to speak and require a voice assessment and/or new speaking exercises.
4. In some patients the steps 1-3 must be performed many more times, until they grow accustomed to exhaling through the upper respiratory tract. If it has been determined that the patient can adequately breathe out and speak, and all other criteria have also been met, the use of PRIMEDI- PHON VENT® speaking valve can be considered.

Compliance of the lungs:

The lungs of seriously ill patients and those with chronic lung disease have an altered compliance (compliance). Therefore the usage of PRIMEDI- PHON VENT® speaking valve must be confined to short-term intervals during the day, and it must be closely monitored. Serious lung diseases cause a loss of lung elasticity and a poor natural stiffness (recoil). As a result, exhalation is prolonged. A careful assessment of the suitability of a PRIMEDI- PHON VENT® speaking valve is required to avoid potential complications due to a case of air traps caused by non-elastic lungs. For these patients, the choice of the correct size of tracheostomy tube for the use of such a speaking valve is especially important because this allows the flow of exhaled air to be encouraged.

Patient care:

The PRIMEDI- PHON VENT® speaking valve can be used at all levels of health care. Whether a patient is suitable for the use of a PRIMEDI- PHON VENT® speaking valve, can already be assessed 48 to 72 hours after a tracheotomy. The valve, on the order of the physician, can be used depending on the severity of existing tracheal oedema and secretions, as soon as the patient's health condition has stabilized and first communication attempts are being made. PRIMEDI- PHON VENT® can be used with infants from two weeks old, provided that the assessment criteria are met.

2. Installing / removing without ventilation

After the assessment criteria to be observed before inserting have been satisfied, the speaking valve can be installed under observation of the medical instructions and in accordance with the following instructions.

⚠ Warning:

Under no circumstances is the PRIMEDI- PHON VENT® speaking valve to be installed directly into the tracheostoma without cannula!

⚠ Warning:

do not use with endotracheostomy tube.

⚠ Warning:

do not use while sleeping.

⚠ Warning:

patients with a PRIMEDI- PHON VENT® speaking valve must be examined and / or monitored according to the directions of the physician.

⚠ Warning:

the PRIMEDI-PHON VENT® is to be removed before the start of medicated nebuliser treatments. If the speaking valve is inadvertently used during the medical nebuliser, the valve must be removed immediately and rinsed thoroughly to remove medication residue, as this can affect the function of the valve membrane.

1. Clarification

For the prevention/reduction of anxiety in patients and to ensure a successful application, the patient, family members, and supervising medical staff (day and night) should be adequately informed and taught about the use of PRIMEDI- PHON VENT®.

2. Patient assessment

The patient should be examined before and after the installation of the speaking valve on the following criteria:

- Signs of life (e.g. heart rate, respiratory rate, oxygen saturation)
- Breathing sounds
- Change in the colour and responsiveness of the patient
- Respiratory function
- Tracheal and oral secretions

3. Suction

Both oral and tracheal suctioning is recommended as needed.

4. Deflation of the cuff

Should a tracheostomy tube with cuff be used, the air should be let out slowly and completely. After letting out the air the patient may have to be suctioned again to remove secretions on and/or above the cuff.

For patients with a tracheostomy tube with cuff, the use of one without a cuff should be considered, if this is medically justifiable, in order to dispense with the need to let the air out when using the PRIMEDI- PHON VENT®.

5. Size of the tracheostomy tube, tracheostomy tube with or without cuff

At the direction of the physician a smaller size and/or a cannula without cuff must occasionally be used to ensure a sufficient flow of exhaled air for the operation of the PRIMEDI- PHON VENT®.

6. Warning label

The warning label supplied is attached to the supply tube of the control balloon of a cuffed tracheostomy tube cuff to emphasise to medical staff the proper use of the PRIMEDI- PHON VENT® speaking valve.

7. Installing the PRIMEDI- PHON VENT® speaking valve

Remove the PRIMEDI- PHON VENT® speaking valve from the packaging.

To use, set the casing of the PRIMEDI PHON VENT® speaking valve on the 15 mm connector or in the 22 mm connector of the tracheostomy tube being used. When mounting on the 22 mm connection on the patient side the casing of the PRIMEDI- PHON VENT® valve clicks noticeably into place.

⚠ To ensure unobstructed breathing it is important to make sure that the valve opening is not blocked by clothing or other things.

Use the PRIMEDI- PHON VENT® only with a secure connection and check the seating!

8. Patient monitoring

The patient should be kept under observation to ensure that the valve diaphragm opens during inspiration and remains closed during exhalation. Patients with a speaking valve should be examined to see that there is sufficient airflow. If the patient has difficulty breathing the PRIMEDI- PHON VENT® should be removed immediately and the consistency of the air passages checked.

⚠ Warning:

If the patient indicates difficulty in the use of PRIMEDI- PHON VENT®, this may be because of an airway obstruction due to a stenosis, tissue mass, tracheomalacia, granulation, vocal cord paralysis in the midline position, secretions, or that the tracheostomy tube is too big for the patient's trachea.

After the blockage has been corrected, the patient should be examined again with regard to the insertion of a PRIMEDI- PHON VENT®.

⚠ Warning:

Be careful when using PRIMEDI- PHON VENT® in conjunction with a heat and moisture exchanger (HME) or hygroscopic condenser humidifier (HCH). These aids/devices absorb humidity from the exhaled air. Because with a PRIMEDI- PHON VENT® speaking valve the air does not flow through the tracheostomy tube when breathing out the performance of the HME or HCH may be compromised and so an additional humidification may be required.

9. Removing the PRIMEDI- PHON VENT® speaking valve

To remove PRIMEDI- PHON VENT®, hold the tracheostomy tube with one hand and with the other hand gently twist out the speaking valve. If a tracheostomy tube with a rotary fixing is used it may be necessary to remove it with a to and fro tipping movement instead of simply twisting. Then remove the ventilation connection to the valve.

10. Patient familiarisation

Although usually there is an instant and easy adaptation to the PRIMEDI- PHON VENT® speaking valve, a gradual familiarisation may be necessary in some cases. If the patient has not vocalised for a long time, regular breathing and, where necessary, speaking, may have to be re-learned. Here, a speech therapist and respiratory therapist can help.

Patients will perceive more normal respiratory sensations (e.g. air flow through the mouth and nasal cavities) and the impact of increased respiratory muscle activity.

Some patients will initially cough more, because the subglottic pressure and normal exhaled airflow in the oral and nasal cavities has been restored. This facilitates secretion, which serves to support pulmonary hygiene.

⚠ If the patient has to cough excessively, the speaking valve should be removed and consistency in the air passages checked.

Anxiety of the patient:

Tracheostomised patients may suffer anxiety with the first use of the PRIMEDI- PHON VENT® speaking valve. Therefore patients should have it explained prior to the use of a PRIMEDI- PHON VENT® speaking valve that they will feel the airflow through the upper respiratory tracts during exhalation and that in the beginning secretions may be secreted through air passages and mouth. These details can help to reduce the anxiety of the patient. Also distractions should be catered for (phone calls, visits, doctor visits, etc.) to facilitate exhalation and/or speech. Optical support techniques, such as simple spirometry or the use of mirrors, cotton wool, feathers, whistles or bubbles are other ways.

3. Use in ventilation

1. Examine ventilator

Examining the settings of the ventilation unit before, during, and after the insertion of the PRIMEDI- PHON VENT® comprises the following checks (including, but not limited to):

- Mode
- Tidal volume (VT)
- Rate
- Proportion of inhaled oxygen
- Oxygen saturation
- Positive end expiratory pressure (PEEP)
- Peak inspiratory pressure (PIP)
- Sensitivity
- Alarm settings

⚠ Adjustments to the ventilator are allowed only on doctor's orders.

2. Deflating the cuff

During the deflation of the cuff (when using a cannula with cuff) adjustments to the ventilator may be required to compensate for leaks around the cannula and to address the needs of the patient.

If the PIP pressure significantly decreases after letting out the air, it may be that air is escaping through the upper respiratory tract. Regulation of the tidal volume may be required to compensate. Increase the tidal volume in small increments, (to avoid over-compensation), until the PIP pressure corresponds to the values from before the deflation of the balloon.

⚠ Warning:

the PIP values from before the deflation of the balloon must not be exceeded.

⚠ Warning:

Before inserting the PRIMEDI- PHON VENT® the tracheostomy tube must be completely deflated.

3. Attaching the PRIMEDI- PHON VENT® when ventilating

Connect the patient side of the PRIMEDI- PHON VENT® speaking valve directly to the 15 or 22mm connection of the tracheostomy tube by inserting until a secure fitting is achieved. The connection is implemented via a friction connection.

A ventilation tube with a 15mm adapter can then be connected to the rotating 15 mm connector of the ventilator side.

⚠ Warning:

Connect the PRIMEDI- PHON VENT® speaking valve directly to the tracheostomy tube to minimize dead space in the system and to prevent a blockage of the speaking valve due to the formation of condensation in the ventilator tube.

⚠ Warning:

If the patient has difficulties with the use of the PRIMEDI- PHON VENT® speaking valve this may be due to an obstruction of the air passageway caused by a stenosis, tissue mass, tracheomalacia, granulation, vocal cord paralysis in the midline position, secretions, or a tracheostomy tube too big for the patient's trachea. After the blockage has been corrected the patient should be examined again regarding the use of the PRIMEDI- PHON VENT® speaking valve.

4. Air passageway pressures

Air passageway pressures can rise with the use of the PRIMEDI- PHON VENT® speaking valve as the patient breathes out through the oronasopharynx and thus creates an almost natural physiological PEEP-pressure. This is part of the normal physiology that has been restored with the closed respiratory system, which is made through the form fitting „sealed“ speaking valve design. Mechanical PEEP- requirements can thus be reduced. Furthermore, (normal) turbulent air currents in the cannula are increased and a higher pressure created. However, increased air passageway pressure should nevertheless remain within the allowable limits for the patient.

If the peak pressure is above the permissible limit, the PRIMEDI- PHON VENT® speaking valve must be removed immediately and the consistency of the air passageway checked. Because some patients with a speaking valve experience a slightly higher air passageway pressure, the pressure values for a separation, as with every change to the ventilation tube, must be re-checked.

5. Alarm settings of the ventilator

Before, during, and after the usage of the PRIMEDI- PHON VENT® speaking valve all alarm settings on the ventilator must be re-checked and, if necessary, changed.

⚠ Warning:

If the ventilator alarm settings are not checked and adjusted patient safety may be impaired.

If the PRIMEDI- PHON VENT® speaking valve is integrated into the ventilation system, the patient no longer breathes out into the ventilation hose. For this reason, ventilators of acute patients must have the upper and lower values for tidal volume and minute volume and apnoea-alarm re-checked and adjusted. The upper and lower pressure alarm settings on the ventilator must be re-evaluated at this time and re-set (more sensitive to separation and blocking) to ensure the safety of the patient. The manufacturer's recommendations of self-testing the ventilator are to be followed (e.g. quicker, extended self-test tube changes for some ventilators of acute patients). Some manufacturers make a speaking valve mode available, which should be used with the PRIMEDI- PHON VENT®.

6. Patient monitoring and speaking valve removal

Examine patients with a PRIMEDI- PHON VENT® speaking valve to ensure there is a sufficient airflow around the tracheostomy tube. If the patient has difficulty breathing, remove the speaking valve immediately and check the consistency of the air passageway. To remove the PRIMEDI- PHON VENT® speaking valve take this out of the ventilation tube and restore the initial state. Re-set the ventilator to all its previous values, before the tracheostomy tube cuff can be re-inflated.

⚠ Warning:

Remove the PRIMEDI- PHON VENT® speaking valve and before re-inflating the cuff of the tracheostomy tube adjust the settings on the ventilator to the initial values.

7. Patient familiarisation

Excessive air discharge from nose and mouth: If the patient feels discomfort due to a continuous stream of air through the mouth and nose, observe the following recommendations:

- Due to reduced control of the glottis, vocal cord and breathing training may be required
- In this case, a compensation on the ventilator can help to control the flow of air through the upper respiratory tracts. Settings on the flow rate, the tidal volume and/or mode can increase the well being of the patient;

4. Troubleshooting

If the patient cannot exhale properly through the upper respiratory tract, the following criteria may have to be reconsidered under certain circumstances:

Cuff assessment:

Check that the air has been completely let out of the cuff of the tracheostomy tube. A tracheostomy tube without a cuff, although not required, can, however, ensure an optimum consistency of the air passages and should be seriously considered when using a PRIMEDI- PHON VENT® speaking valve, if the patient is a suitable candidate.

Tracheostomy tube assessment:

The size of the tracheostomy tube should be evaluated carefully to determine whether a smaller tracheostomy tube would be more suitable in case the size of the tracheostomy tube or that of the deflated cuff hinders adequate exhalation.

Closure of the air passageway:

A medical examination (e.g., bronchoscopy) should be considered to ensure that no undetected respiratory tract closure exists (tissue mass, stenosis, granulation, vocal cord paralysis, etc.).

Review of the situation:

The optimal situation of patient and tracheostomy tube should be examined again.

Useful life

The PRIMEDI- PHON VENT® can be reinstalled after proper cleaning. The cleaning instructions are to be complied with.

If handled and cleaned correctly and all precautionary measures are taken into account and observed, the durability/useful life is at least two months

In the case of incorrect cleaning or use the normal service life cannot be guaranteed. If the speaking valve, through operation and care not recommended by the manufacturer, begins to stick, get noisy or vibrate in less than two months, it must be replaced.

PRIMEDI- PHON VENT® can be used longer if it does not begin to stick or get noisy or vibrate, exhibit increased resistance by inspiration, or show other problems.

Unused speaking valves should be stored in a dry environment protected from sunlight and/or heat.

⚠ Warning:

A substituted PRIMEDI- PHON VENT® must be correctly disposed of promptly, according to the instructions of the physician.

⚠ Warning:

PRIMEDI- PHON VENT® speaking valves should be stored protected from sunlight, heat and dust.

Disclaimer:

Primed® Halberstadt Medizintechnik GmbH (hereinafter „Primed®“) declares that reasonable care has been used in the manufacture of this product. This warranty is the only valid warranty and is in lieu of any other, whether direct or indirect, whether written or oral, including but not limited to any indirect warranties of merchantability or fitness for a particular purpose. No product can be completely effective under all circumstances because of biological differences in humans. For this reason and due to the fact that we have control neither over the conditions under which the product is used, nor over the diagnosis of patients, the methods of administration or handling once the product has left our possession, Primed® guarantees not the prevention of adverse effects by use of the product. Primed® shall not be liable for any incidental or consequential damages or costs, directly or indirectly resulting from the use of the product. Primed® will replace any product which has been proven faulty at the time of delivery. No representative of Primed® is entitled to modify the foregoing information relating to this product, or to alter or accept any additional liability or responsibility.

Cleaning instructions

PRIMEDI- PHON VENT® is individually wrapped. Ideally, the patient should have an additional speaking valve as a reserve so that one PRIMEDI- PHON VENT® speaking valve can be worn, while the other is being cleaned.

Always clean the PRIMEDI- PHON VENT® speaking valve separately. Before cleaning it is therefore necessary to remove the valve from the 15mm/22mm connection of the tracheostomy tube.

The cleaning frequency of the PRIMEDI- PHON VENT® speaking valve is dependent on the secretion tendencies of the patient. It should, however, be cleaned at least once a day.

A mild pH neutral washing lotion or warm soap solution can be used for daily cleaning. The water should not be hot (above 65°C). The speaking valve should then be thoroughly rinsed with lots of warm water. We recommend air-drying the PRIMEDI- PHON VENT® speaking valve.

⚠ Warning:

Do not dry the speaking valve with heat!

⚠ Warning:

Do not use a brush or cotton swab for cleaning a PRIMEDI- PHON VENT® speaking valve as this may damage the valve.

⚠ Warning:

The PRIMEDI- PHON VENT® speaking valve should not be cleaned with hot water above 65°C, nor with peroxide, bleach, vinegar or alcohol.

Under no circumstances should cleaning agents be used that release chlorine, or contain strong lye or phenol derivatives. The PRIMEDI- PHON VENT® speaking valve could be significantly damaged or even destroyed.

⚠ Warning:

PRIMEDI- PHON VENT® is not autoclavable!

Technical Data

<u>Indication:</u>	tracheotomised and ventilated patients
<u>Weight:</u>	ca. 6 g
<u>Length/height:</u>	34 mm
<u>Connections:</u>	<ul style="list-style-type: none">- 15 mm connection (ISO- inner cone)- 22 mm connection (mounting flange)- 15 mm connection (ISO-outer cone), rotatable
<u>Environmental impact:</u>	Hydrocarbons, carbon dioxide, water, and carbon oxides are released during the combustion of the product.

MRT / CT

Due to the materials used, there are no restrictions in examinations by MRI or CT for patients using a PRIMEDI-PHON VENT® speaking valve.

Reuse / single-patient product

The PRIMEDI-PHON VENT® speaking valve is a single-patient product and is intended exclusively for use with only one patient. A re-use with another person is not allowed.

Warranty

We accept the guarantee for the flawlessness of the delivered item, in accordance with our Terms and Conditions (GTC).

You will find a print-out of our terms and conditions and more on our business stationery. Failure to comply with the use according to the instructions voids any warranty claims.

Any unauthorized modification of the product leads to the exclusion of any warranty or liability on the part of the manufacturer.

The PRIMEDI-PHON VENT® speaking valve REF 200320 is a class IIa medical product pursuant to the European Medical Product Directive 93/42 EEC.

Primed® Halberstadt Medizintechnik GmbH specialises in the supply of aids for patients who have had a tracheotomy and/or laryngectomy. You can see a current product range in our product catalogue, which you can request for free or download from www.primed-halberstadt.de.

PRIMEDI-PHON VENT® is a registered trademark of Primed® Halberstadt Medizintechnik GmbH.

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PRIMEDI-PHON VENT® es una marca registrada de Primed® Halberstadt Medizintechnik GmbH.