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# Instructions for use for PRIMED® tracheostomy tubes for neo- nates and children

<b>NEO</b>	Neonatal tracheostomy tubes without cuff
<b>PED</b>	Pediatric tracheostomy tubes without cuff
<b>PEDL</b>	Pediatric tracheostomy tubes without cuff, long
<b>PEDC</b>	Pediatric tracheostomy tubes with cuff
<b>PEDCL</b>	Pediatric tracheostomy tubes with cuff, long

**These instructions for use must be read through carefully each time before the products are used. They must be provided to all persons entrusted with care, so that they can inform themselves and follow the instructions contained.**

The medical professionals must instruct parents and carers in detail in the correct care of the stoma and the suction technique. The changing of the tubes must be practised in the hospital under all circumstances.

## **Product description**

PRIMED® tracheostomy tubes for neonates and children are a Class IIb medical device in accordance with Annex IX of the European Medical Device directive 93/42/EEC.

All product versions are delivered sterile (EO).

PRIMED® tracheostomy tubes for neonates and children are exclusively produced from high quality, DEHP-free polyvinyl chloride. All the materials employed are free from latex, DEHP and bisphenol A. The tube shafts are radiopaque.

The ergonomically designed neck flange has an integrated rotatable 15 mm connector for connection to standard devices, accessories and/or hose systems (flexible extension adapter for 15 mm connectors REF 200090). All tubes are delivered individually packaged with an obturator, neck strap and disconnecting wedge. The tips of the tubes and obturators areatraumatically rounded.

Cuffed tube models have a cuff with a very thin wall, which adapts to the natural shape of the trachea safely and protectively when inflated in the low pressure range and seals off the airway securely. The thinness of the cuff wall allows it to sit tightly against the shaft of the tube (tight to shaft) when deflated. The pilot balloon is used to display and verify the filling of the cuff.

**PRIMED® tracheostomy tubes for neonates and children** are exclusively produced without inner cannulae.

To guarantee optimal positioning and the best air supply possible, the attending doctor should select a tube based on the anatomy of the child and its condition. In case of queries or complications, the attending doctor should be consulted immediately.

#### **CAUTION!**

**PRIMED® tracheostomy tubes for neonates and children are intended for use on one patient only**

**PRIMED® tracheostomy tubes for neonates and children must only be used by persons who have been properly trained and/or have experience in the use of the procedures in paediatric tracheotomy.**

#### **Indication**

**PRIMED® tracheostomy tubes for neonates and children** are indicated following tracheotomies, comparable with the indication for adults depending on problems caused by obstruction, secretion and ventilation, and permit access to the airways and their monitoring. Cuffed tubes are indicated in cases where it is necessary to create a seal between the wall of the trachea and the tube.

#### **Contraindication**

**PRIMED® tracheostomy tubes for neonates and children should not be used if the patient is allergic to the materials used.**

**PRIMED® tracheostomy tubes for neonates and children should not be used during treatments involving lasers (laser therapy).** If the laser were to come into contact with the tube, it could cause damage. If oxygen-rich mixtures are present, the tube could catch fire, with it not being possible to rule out injuries resulting from the influence of the heat and the release of corrosive and toxic combustion products. There is a considerable risk of injury.

**In the case of surgical interventions on the trachea, the use of an electrocautery also poses the risk of a tube fire.**

**Pediatric tracheostomy tubes with cuff** are exclusively designed for patients with a preserved larynx.

**Due to the one-way valve in the filling system's pilot balloon, the cuffed tubes are not suitable for use in MRI scans. (Artefacts in imaging).**

**We recommend consulting the doctor and using a suitable tracheostomy tube for the period of the investigation.**

#### **Warnings**

**In the case of ventilated patients, special attention should be paid to correct fastening, as disconnection could endanger the patient's life. The connectors should be checked at regular intervals after connection.**

We recommend setting a cuff pressure of max. 25 cm H<sub>2</sub>O (approx. 18.8 mm Hg or 25 hPa) when filling the cuffs of paediatric tubes. Exceeding the maximum pressure for extended periods of time can affect the blood supply to the mucosa. There is a risk of tracheal injuries and/or impaired breathing.

In the case of ventilated patients, the cuff pressure must not fall below the minimum set by the doctor, in order to prevent silent aspiration. Hissing noises in the area of the cuff, especially during exhalation, are a sign that the cuff is not sealing off the trachea properly. Should it prove impossible to seal off the trachea with these pressure values, the air should be removed from the cuff again completely and the inflation procedure repeated. If the repeated procedure is also unsuccessful, the next largest pediatric tracheostomy tube with cuff should be selected.

The correct cuff pressure must be checked regularly, i.e., approx. every 2 hours.

All the devices and aids used to fill the cuff must be clean and free from foreign particles. The pilot balloon should be sealed using the attached LL plug immediately after inflation of the cuff.

Never reposition a tracheostomy tube with the cuff inflated! Risk of injury! Pediatric tracheostomy tubes with cuff (PEDC and PEDCL) should never be sealed, as under certain circumstances the patient might not have a sufficient airway available to breathe despite the cuff being completely deflated. When using lubricants, it must be ensured that the tube lumen does not become blocked, as this can impede breathing.

The lubricants used must be suitable for tracheostomy.

### **Leakage test**

Cuffed tubes should always be checked for leaks before use.

The tube and cuff should be checked for leaks before and after every use.

This is done by filling the cuff with air at a pressure of 30 to 33 cm H<sub>2</sub>O (22-24 mm Hg or 29 to 32 hPa). If the tube and cuff are free from leaks, there will be no significant drop in pressure. It is recommended to use a low-pressure cuff pressure measuring device (**OBSERVE the manufacturer's instructions for use!**).

### **Caution!**

**The tube must not be used under any circumstances if it displays leaks.**

**Following the leakage test, the cuff must be completely deflated again before the tube is used!**

### **Use of a PRIMED® tracheostomy tubes for neonates and children**

Extreme care and caution should be exercised when using the PRIMED® tracheostomy tubes. When removing the tube from the sterile packaging, ensure that the

specified expiry date has not elapsed and the sterile packaging does not display any signs of damage. The wearing of sterile disposable gloves is recommended. Check the tube for external damage and loose parts. If you identify any problems, the tube must not be used any further. It should be returned to the manufacturer for checking/repairs.

The size should be checked against the size specifications on the label and the neck flange to make sure it matches.

Tracheostomy tubes must be cleaned thoroughly and checked for damage before every subsequent use (see Cleaning). Following confirmation that the tube is in proper working order, the obturator can be inserted into the tube.

To improve the lubricity of the PRIMED® tracheostomy tube, facilitate introduction into the trachea and minimise the urge to cough which develops, it is recommended to rub the external tube and the tip of the obturator with suitable lubricant (e.g., stoma oil). Once prepared in this way, the tracheostomy tube can be inserted.

We recommend placing a protective tube cover (compress, e.g., REF 200234 paediatric foam compress made of PU foam) on the tube shaft.

When placing the tube, it is important to ensure that the soft tube material / cuff is not damaged by friction or sharp edges and that the cuff is completely deflated.

As soon as the tracheostomy tube is in the right position, the obturator should be removed immediately and the tube fixed in place using the neck strap.

**Caution!**

**Make sure that the neck flange is always outside of the tracheostoma.**

**Inflating the cuff**

Once the tube has been inserted and fixed in place, the cuff is inflated with a low-pressure cuff pressure measuring device to a pressure of max. 25 cm H<sub>2</sub>O (approx. 18.8 mm Hg or 25 hPa) with air, unless the doctor has prescribed a different pressure.

**Removing a PRIMED® tracheostomy tube**

It is recommended to change tracheostomy tubes for neonates and children including the accessories frequently and routinely. The time and interval of the change should be prescribed and monitored by the attending doctor.

Before removing/changing a cuffed tube, any secretions that have collected above the cuff must be removed, unless suctioning is not advisable from a medical perspective.

The cuff is deflated by removing the air slowly using a syringe via the pilot balloon with one-way valve.

Grasp the PRIMED® tracheostomy tube by the connector when removing. Following removal of the tracheostomy tube, it should be cleaned in accordance with the instructions (see Cleaning).

**Caution!**

Cuffed tracheostomy tubes must only be removed from the tracheostoma once the cuff has been completely deflated, as the tracheal mucosa could otherwise be injured.

We recommend keeping at least one replacement tube at the patient's bedside for emergencies at all times.

In case of queries or complications, the attending doctor should be consulted immediately.

Whenever PRIMED® tracheostomy tubes are used, there is a risk that candida, bacteria, etc., may colonise the products.

It must therefore be ensured that the tube is always cleaned and, if necessary, disinfected before reinsertion. The instructions for cleaning and disinfection must always be followed closely.

**Caution!**

When disconnecting accessory parts such as an artificial nose (heat and moisture exchanger) or a ventilation tube, avoid pulling, twisting and pressing on the tracheostomy tube excessively, as this could result in disconnection or blockage of the tracheostomy tube or even damage it.

Each PRIMED® paediatric tracheostomy tube comes with a disconnecting wedge. Pushing the disconnecting wedge between the 15 mm-connector of the tracheostomy tube and the connector of the accessory or ventilation tube allows disconnection of the items without excessive force.

Once fully applied, the disconnecting wedges centres itself on the 15 mm-connector of the tracheostomy tube and, following disconnection and removal of the accessory / ventilation tube, can be removed from the patient in an axial direction.

**Cleaning**

Tracheal cannulas and accessories must always be cleaned as appropriate to the needs of the specific patient.

**The following information must be complied with!**

– Remove the rotatable 15 mm connector from the neck flange using the disconnecting wedge supplied prior to cleaning.

To do so, push the tip of the wedge under the lower edge of the 15 mm connector. Continuing pushing, press one side of the top collar of the neck flange on the side furthest away from the patient inwards towards the connector until the connector releases.

– It is best to preclean the PRIMED® tracheostomy tube, the 15 mm connector, the obturator and the disconnecting wedge under lukewarm running water. When cleaning, be especially careful not to damage the soft material of the tube.

- We recommend a pH-neutral liquid soap for the subsequent fine cleaning. However, it is also possible to use a suitable tube cleaning agent (**OBSERVE the manufacturer's instructions for use!**). In this case, the external cannula, 15 mm connector, obturator and, if applicable, disconnecting wedge should be immersed alongside each other in the prepared cleaning solution.

Use of a tube cleaning container with tray insert facilitates handling.

#### **Caution!**

##### **Clean cuffed tracheostomy tubes inflated.**

**Under no circumstances should the tracheostomy tube be cleaned with agents which have not been approved by the tube manufacturer.**

**Under no circumstances should you use aggressive domestic cleaning agents, high-percentage alcohol, hydrogen peroxide or agents for cleaning dentures. This presents an acute health risk! In addition, there is a risk of the tube being destroyed or damaged.**

**The agent must be suitable for the treatment of tubes made of medical quality PVC and be approved.**

**Only one tracheostomy tube should be placed in the tray insert of the tube cleaning container at a time.**

**Under no circumstances should the tracheostomy tubes be cleaned in a dishwasher, steamer or microwave oven.**

#### **Disinfection**

**We recommend always using cold disinfection. The effect of thermal disinfection (see below) is disputed in professional circles, and if the specified temperatures are exceeded, it can even result in damage to the tube.**

PRIMED® tracheostomy tubes can be cold disinfected with suitable disinfectants (e.g., PRIMEDICLEAN® REF 70810 for combined cleaning and disinfection). (**OBSERVE the manufacturer's instructions for use!**)

However, it should only be performed if the attending doctor has prescribed it due to a specific condition or it is indicated by the respective care situation.

Following the cleaning described above, the tracheostomy tubes must then be immersed in a disinfectant suitable for the tubes for disinfection. The manufacturer's specifications must be observed closely.

Alternatively, the tubes can also be thermally disinfected. To do so, disinfect the tube thermally in clean water without the addition of chemical additives at a temperature of min. 65°C but max. 68°C under domestic conditions. The temperature must be maintained in the water bath for **max.** 15 minutes after immersion of the tubes and not exceeded. It must be ensured that the temperature of 68°C is not exceeded and boiling out does not occur. This can cause serious damage to the tracheostomy tube.

### **Disinfection of cuffed tracheostomy tubes**

#### **CAUTION!**

**As even minor inobservances during the disinfection of cuffed tubes can damage them or impede their function, we recommend replacing the tube for a new PRIMED® cuffed tracheostomy tube instead of disinfection.**

Cold disinfection of a cuffed tube should always only be performed if deemed necessary according to the ruling of the attending doctor and exercising extreme caution after performance of the cleaning described. The disinfectant PRIMEDI-

CLEAN® should only be used in accordance with the manufacturer's specifications for disinfection.

#### **Disinfection procedure:**

- The cuff should always be inflated at a pressure of 30 to 33 cm H<sub>2</sub>O (22-24 mm Hg or 29 to 32 hPa) prior to possible disinfection.
- The pilot balloon must not be immersed in the disinfection solution together with the inflated cuff, as disinfectant solution could enter the filling line and the cuff via the valve in the pilot balloon. Should this occur nevertheless, the tube cannot be used any further and must be disposed of properly.

#### **CAUTION!**

**Under no circumstances should disinfectants be used which release chlorine or contain strong alkalis or phenol derivatives or hydrogen peroxide. The use of such could seriously damage or even destroy the tubes.**

**If disinfection of a cuffed tube is performed, extreme caution should be exercised and the instructions followed closely. There is otherwise a risk of considerable damage to the tracheostomy tube which would present a serious risk to the user's health!**

**Following every cleaning/disinfection procedure, the tubes (inside and outside), obturator, etc., must be rinsed through with cold, boiled water and then dried.**

**To rule out health risks, there should be no more cleaning agent or disinfectant residues on the tube when it is inserted in the tracheostoma.**

Following cleaning and/or disinfection, the rotatable 15 mm connector should be connected to the neck flange again. To do so, coat the upper edge of the flange on the side facing away from the patient minimally with a lubricant suitable for tracheostomy and push the 15 mm connector onto the flange with the largest diameter first. Following installation, the connector should be positioned below the lower edge and be easy to rotate. The obturator can then be placed in the tracheostomy tube.

We recommend rubbing the tube with a drop of stoma oil for care purposes.

No liability is accepted in cases of failure to comply with the cleaning instructions.

## **Sterilisation**

PRIMED® tracheostomy tubes are delivered sterile. They do not need to be resterilised by the user.

## **Period of use, storage**

The period of use of PRIMED® tracheostomy tubes is influenced by a number of factors. For example, the composition of the secretion, the thoroughness of the cleaning and other aspects can be of decisive importance.

We recommend replacing the tube after 29 days of use at the latest.

Replaced tracheostomy tubes must be disposed of properly immediately.

### **Regular checks avoid risks!**

To ensure continuous treatment, it is urgently recommended to have at least one replacement tube available at all times.

Tubes which have not yet been used should be stored in a dry environment, protected from direct sunlight and/or heat.

## **Scope of delivery**

PRIMED® tracheostomy tubes for neonates and children are packed individually and sterile and are available in different sizes and models as per the table.

Package content:

- PRIMED® tracheostomy tube
- Obturator
- Disconnecting wedge
- Neck strap
- Instructions for use

All the important data concerning the tracheostomy tube are provided on the label on the external, printed box.

## **Guarantee**

We guarantee that the delivered object is free from defects in accordance with the scope of our general terms and conditions (GTCs). A copy of our GTCs can be found among our commercial papers and on our website [www.primed-halberstadt.de](http://www.primed-halberstadt.de), for example.

Repairs and other work on PRIMED® tracheostomy tubes may only be performed by Primed Halberstadt Medizintechnik GmbH and demonstrably instructed specialists. Failure to comply with these criteria will void the guarantee!


## **Catalogue / accessories**

Our company specialises in supplying tracheotomy and laryngectomy patients with medical supplies. An up-to-date product overview can be found in our product catalogue, which you can request from us free of charge or view at [www.primed-halberstadt.de](http://www.primed-halberstadt.de).


### Sizes / lengths

The sizes of the PRIMED® tracheostomy tubes are based on the DIN EN ISO 5366-3 European standard. The specified sizes describe the internal diameter of the tracheostomy tube at the tip of the tube in mm. The requisite information concerning the individual product can be found on the packaging and the neck flange.

#### PRIMED® cuffless tracheostomy tubes for neonates and children (NEO / PED / PEDL)

#### PRIMED® cuffed tracheostomy tubes for children (PEDC / PEDCL)

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PRIMED® tracheostomy tubes are a Class IIb medical device in accordance with Annex IX of the European Medical Device directive 93/42/EEC.

The manufacturer reserves the right to make changes to the products at any time.

PRIMED® is a registered trademark of **Primed Halberstadt Medizintechnik GmbH**.

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