
PEAK FLOW METER

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Liebe Patientin, lieber Patient!

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Vertrieb

PARI GmbH
Postfach 1551
D-82317 Starnberg

Vitalograph Ireland Ltd
Gort Road Business Park
Ennis
Co. Clare
Ireland

Please read carefully!

Dear patient,

Your doctor has prescribed the PARI Peak Flow Meter for you for measuring the strength of your breathing. This brochure contains detailed information on how you can learn to be more accurate in estimating the current state of your respiratory tracts and any changes in them. Your doctor is drawing up an action plan for you that will help you to react quickly and optimally to changes in your illness.

What is a peak flow?

The peak flow measures the force or speed of your exhalation. Your peak flow result is higher when your health is better, and lower when your respiratory tracts have narrowed. Peak flow values are also the basis on which asthma is classified into different degrees of severity.

The aim of this is to achieve the most stable peak flow possible. This means that the differences between morning and evening measurement results (that also occur in healthy people) are as low as possible.

Why are peak flow measurements taken?


Both you and your doctor can already detect small changes in the state of your lungs through regular peak flow measurements. By taking medication in good time, you can get your asthma under control or prevent a serious asthma attack. Your peak flow log lets you assess whether you are responding to the medication prescribed for you by your doctor.

Measurement


How is the peak flow measurement carried out?

Measure your peak flow values at least twice a day, and if possible at the same time each day, e.g. in the morning on getting up and in the evening before going to bed. If you are breathless or have a bronchial infection, you should take additional readings.


- Push the pointer to the lower end of the measurement scale.


- Hold the device horizontally in front of your mouth. Please make sure that your fingers are not covering the pointer slit or the openings at the end of the device.


- Take a deep breath. Firmly surround the mouthpiece using your lips and teeth.
- Blow as **hard** as possible into the peak flow meter using just **one** breath.

- Read off the value from the pointer position. Take a note of the value and push the pointer back to the lower end of the scale.

- Take two more readings and then work out the colour zone of the highest of the three readings.

- Enter the highest of the three readings (personal best) into the enclosed log table.

Action Plan

How can I get my personal asthma Action Plan?

Plan for you, stating which measures will have to be taken if the state of your respiratory tract worsens, i.e. if the peak flow value falls below certain limits. This Action Plan (for the self-management of your illness) is based on the peak flow values measured over a certain period and taken from your table.

Green	Everything is O.K. Your peak flow value is more than L/Min (> 80 % of your personal best). Take the following continuous daily medication:
Amber	Caution! Measures must be taken. Your peak flow value is below L/Min (< 80 % of your personal best) Increase your medication dose as follows:
Red	Emergency program. Your peak flow value is below L/Min (< 50 % of your personal best). Take the following emergency medication immediately: If your peak flow value does not improve by more than at least 50 units, contact a doctor on emergency call immediately.

Your peak flow results normally fall into one of three coloured zones – green, amber or red. Follow the Action Plan instructions for the colour zone in which the highest of your measurements falls.

Stick the enclosed „Action Plan“ filled out by your doctor onto the rear side of your peak flow meter.

How often should the Action Plan be updated?

You should see your doctor at least twice a year to have your asthma checked and for possible changes to be made to your Action Plan.

Peak flow log for four weeks

Please take sufficient copies of the Action Plan before filling it out for the first time.

You can also obtain new copies from the PARI Internet page on www.pari.de.

Name of the patient:										Prescribed medication:										Month:											
																				1 2 3 4 5 6 7 8 9 10 11 12											
Date:	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
L/Min.	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀	☀
800																															
750																															
700																															
650																															
600																															
550																															
500																															
450																															
400																															
350																															
300																															
250																															
200																															
150																															
100																															
50																															
Cough																															
Breathless																															
Sputum																															

☀ = morning reading

☾ = evening reading

Important

Care and hygiene in cleaning your Peak Flow Meter

In order to prevent a contaminated Peak Flow Meter becoming a health hazard, e.g. from an infection, the following hygiene regulations must always be followed.

Your PARI Peak Flow Meter will usually give reliable results for up to 3 years. At the end of this time, you must ask your doctor for a new device.

The device will only function correctly if all moving parts move freely, otherwise a replacement will be required. Do not apply oil or other lubricants!

Please make sure that the device is not damaged and keep it clean and free from dust. If you suspect that the device is damaged or is giving incorrect results, contact your doctor immediately.

The mouthpiece of the Peak Flow Meter should be cleaned once a week with a soft sponge that has been immersed in a solution of washing up liquid diluted with water. The mouthpiece should then be rinsed thoroughly with fresh water.

The device should only ever be used by one patient. This is particularly important in the case of contagious illnesses.

Please do not ever try to dismantle the device. This can cause your Peak Flow Meter to give false measurements!

As your doctor is the only person who can draw up or amend your Action Plan, take your PARI Peak Flow Meter and the peak flow log with you when you visit the doctor.

The coloured bars on the left next to the measurement scale on your PARI Peak Flow Meter (for setting your individual limit values) can be moved.

Never change the position of the coloured bars yourself!

Should you inadvertently move the coloured bars, please ask your doctor to reset them. Fixing them using transparent adhesive tape should prevent them being moved inadvertently.

Guarantee

Your PARI Peak Flow Meter is supplied with a one-year guarantee.

If a defect occurs, contact your specialist dealer.

Technical Data

Material:	Recyclable ABS plastic
Accuracy:	Better than ± 10 l/min or $\pm 10\%$ of the measured value
Reproducibility:	Better than ± 10 l/min or $\pm 5\%$ of the measured value
Changes due to altitude:	Reduces measured values by about 5% per 1000 m (lower air density increases PEF by about 5% per 1000 m)
Max. flow resistance:	0.00384 kPa/l/min at 720 l/min
Standard range:	50 – 800 l/min BTPS
Storage:	Temperature: 0 – 50°C
Relative humidity:	10 % – 95 %
Frequency response:	Profile difference between A and B less than 15 l/min/15% (EN ISO 23747:2007)

Date of issue of information: August 2010

Sales

PARI GmbH
Postfach 1551
D-82317 Starnberg
Germany

Vitalograph Ireland Ltd
Gort Road Business Park
Ennis
Co. Clare
Ireland

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