# RespiCare Sp/ST

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Dräger Medical b.v.

## For Your Safety and that of Your Patients

### Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observation of these instructions. The apparatus is only to be used for purposes specified here.

### Maintenance

The apparatus must be inspected and serviced regularly by trained service personnel at twelve monthly intervals (and a record kept). Repair and general overhaul of the apparatus may only be carried out by trained service personnel. We recommend that a service contract be obtained with a Dräger Medical dealer and that all repairs also be carried out by them. Only authentic Dräger Medical spare parts may be used for maintenance. Observe chapter "Maintenance Intervals".

### Accessories

Do not use accessory parts other than those in the order list.

### Not for use in areas of explosion hazard

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Do not use this apparatus in areas with inflammable anesthetic agents.

### Safe connection with other electrical equipment

Electrical connections to equipment which is not listed in these Instructions for Use should only be made following consultations with the respective manufacturers or an expert.

### Liability for proper function or damage

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus is serviced or repaired by personnel not employed or authorized by Dräger Medical or if the apparatus is used in a manner not conforming to its intended use.

Dräger Medical cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Dräger Medical are likewise not modified by the recommendations given above.

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## Intended Use

RespiCare® – respiratory support unit for patients who need pressure-controlled ventilation support for therapy reasons. For use only with patients who can maintain spontaneous breathing.

The device is not intended for unconscious or non-spontaneously breathing patients.

RespiCare should only be used with a non-invasive ventilation connector such as nasal mask, mouth mask or mouth-and-nose mask.

For sleep apnea therapy,  
for weaning long-term ventilation patients and  
for disorders of the respiratory muscles.

For use  
in hospitals,  
in sleep laboratories,  
at home,  
for stationary or mobile use in vehicles and wheelchairs.

RespiCare is available in two models with specific ventilation modes:

### RespiCare Sp

CPAP – Spontaneous breathing with continuous positive pressure in the airways.

CPAP/ASB – Spontaneous breathing with synchronized pressure support.

With disconnection/apnea alarm.

### RespiCare ST

CPAP – Breathing with continuous positive pressure in the airways.

CPAP/ASB – Spontaneous breathing with synchronized pressure support and apnea ventilation.

With disconnection/apnea alarm.

**The unit may only be used as prescribed by a doctor and with the settings specified by the doctor!**

**The unit must not be used for life-saving functions.**

**Mobile telephones must not be used within a distance of 10 metres from the unit!**

Mobile telephones can impair the correct functioning of electro-mechanical equipment and can therefore put patients at risk.

**The unit must not be used in explosion-hazard areas!**

Risk of fire and explosion.

In case the device is used in an ambient temperature above 30 °C:

- make sure that the air going to the patient does not exceed 40 °C.

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
## Preparing for Use

### Positioning

- The RespiCare must be set up on a level, horizontal surface, such as a bedside table, wall shelf or trolley. Ventilation slits in the base of the housing must remain clear! Ensure that the unit cannot fall down!
  - The unit must not be set up near a heater and must not be exposed to direct sunlight, otherwise it may overheat! Ensure that air can circulate freely around the unit.
  - Filter and silencer on the rear of the unit must not be blocked by curtains, blankets, etc. otherwise air cannot be drawn into the unit!
  - Ensure that the unit cannot fall into full bathtubs, basins, etc. as this can cause electric shocks!
- When used out doors, e.g. in a wheelchair:
- Ensure that the unit is protected from wetness!

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### Insert coarse dust filter

- 1 Plug new coarse dust filter into insert on rear of unit and turn clockwise as far as possible = lock into place.
- 2 Plug silencer into filter housing until it engages.  
The coarse dust filter traps coarse particles, such as fluff.  
Allergic patients can use a fine dust filter to trap household dust, pollen, etc. See order list, page 32.

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
## Connect ventilation hose

Only use original ventilation hoses!

- 1 Firmly plug the ventilation hose into the socket on the unit.
  - 2 Plug the other end of the ventilation hose into the SilentFlow measurement hoses pointing towards ventilation hose.
  - 3 Plug measurement hoses into the sockets on the unit – note different diameters.
  - Connect the SilentFlow with the port on the nasal mask.
- A nasal mask without outflow port must be used!**
- CO<sub>2</sub> is purged from the system via the SilentFlow.
- Hoses and anti-static or electrically conductive tubes must not be used.

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
The RespiCare has an input for 24 V DC.

To operate it on the mains:

- 1 Connect the external AC/DC adapter to the 24 V DC input.
- 2 Plug the mains cable into the AC/DC adapter and into the mains socket.

The AC/DC adapter can be connected to any voltage between 100 and 240 V AC.

To connect it to 12 V DC source:

- Use the optional DC/DC converter, see order list, page 32.

To connect it to 24 V DC source:

- Use the 24 V DC cable, see order list, page 32.
- Do not connect to any power supply that might supply a too high voltage in case of any fault. The power at the input of the device may never exceed 24 V.

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

Every basic RespiCare unit can be set as required with the aid of the control panel.

**The control panel may only be used by a doctor or by someone who has been duly authorized by the doctor!**

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- 1 Insert the control panel in the magnetic holder on the front of the unit. Lay the unit on its side, press the lead into the duct in the base of the unit and

- 2 plug the control panel lead into the socket on the rear.
- Return the unit to an upright position.  
Remote operation of the control panel:  
● Connect the extension lead to the control panel plug and plug the extension lead into the socket on the rear of the unit.  
Only use the original extension lead.  
Other components must not be connected.

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
### Use of analog output

To connect the RespiCare to an external device with analog signal inputs, the optional control unit with analog outputs for pressure and flow can be used. (see order list, page 32)


- Disconnect the RespiControl from the RespiCare.
- Connect the "RespiControl Analog" to the RespiCare.
- Connect the analog data cable to the "RespiControl Analog" and to the external device.
- Observe Technical Data, page 29.

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

## Switch on unit

- Press on/off switch on rear of unit – until it engages = ON.
- The RespiCare unit starts a self-test for approx. 1 second to test the internal program memories. All LEDs briefly light up during this time.  
The audio alarm is sounded for a moment.
- Make sure the audio alarm works before using the device.

- The green LED in the button remains alight.

On the control panel:

- All pixels appear on the display and the display is illuminated for approx. 1 second. The model name then appears in the first line, together with the last ventilation mode selected.  
This line remains on display as the status line.  
Example:  
**CPAP/ASB**

## Set ventilation modes / parameters

This is done with the aid of the control panel.

**The control panel may only be used by a doctor or by someone who has been duly authorized by the doctor!**

- Display - example:  
↓ CPAP ↓  
Vertical arrows identify the selected line and the direction in which additional lines can be displayed.
  - Press button or to select a line.
  - Press button .
- The two arrows flash and point towards the text. Example:  
→ Paw = 20 mbar ←

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The mode or parameter can now be set:

- 2 Use the buttons  $\downarrow$  or  $\uparrow$  to set the required values. **The settings are immediately adopted by the unit.**

- 3 Press button  $\bullet$  to confirm.

The display is automatically illuminated for two minutes after pressing one of the three buttons.

The settings are saved even when the unit is switched off.

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**CPAP – Continuous Positive Airway Pressure**

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## Set CPAP

**CPAP – Continuous Positive Airway Pressure**

Spontaneous breathing with continuous positive pressure in the airways.

- Press  $\uparrow$  to select the first line and activate by pressing  $\bullet$
- Use  $\downarrow$  to set **CPAP** mode and press  $\bullet$  to confirm.
- Press  $\downarrow$  to select **Paw = \_ \_ mbar** and activate by pressing  $\bullet$
- Press  $\downarrow$  or  $\uparrow$  to set **Paw = \_ \_ mbar** and press  $\bullet$  to confirm. Adjustment range: 5 to 20 mbar in increments of 0.5 mbar

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
**CPAP – Continuous Positive Airway Pressure**

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## Set CPAP/ASB

**CPAP – Continuous Positive Airway Pressure**

Spontaneous breathing with continuous positive pressure in the airways.

**ASB – Assisted Spontaneous Breathing**  
Spontaneous breathing with synchronised assistance. With every spontaneous inspiration, the unit provides assistance until the pressure level **Pinsp** is reached and then switches over to pressure level **Pexp** when expiration commences or after not more than 4 seconds.

- Press  $\uparrow$  to select the first line and activate by pressing  $\bullet$
- Use  $\downarrow$  or  $\uparrow$  to set **CPAP/ASB** mode and press  $\bullet$  to confirm.
- Press  $\downarrow$  to select **Pinsp = \_ \_ mbar** and activate by pressing  $\bullet$   
**Pinsp:** inspiratory pressure.

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- Use  $\downarrow$  or  $\uparrow$  to set **Pinsp = \_ \_ mbar** and press  $\bullet$  to confirm.
- Press  $\downarrow$  to select **Pexp = \_ \_ mbar** and activate by pressing  $\bullet$   
Pexp: expiratory pressure.
- Use  $\downarrow$  or  $\uparrow$  to set **Pexp = \_ \_ mbar** and press  $\bullet$  to confirm.

After switching on or disconnecting the unit, wait for a few breaths until the system has stabilized before setting the trigger.

- Press  $\downarrow$  to select **Trigg.I = \_ of 5** and activate by pressing  $\bullet$   
Trigg.I: Sensitivity for synchronising spontaneous breathing and inspiration.  
1 = Maximum sensitivity
- Use  $\downarrow$  or  $\uparrow$  to set **Trigg.I = \_ of 5** and press  $\bullet$  to confirm.
- Press  $\downarrow$  to select **Trigg.E = \_ of 5** and activate by pressing  $\bullet$   
Trigg.E: Sensitivity for synchronising spontaneous breathing and expiration.  
1 = Maximum sensitivity
- Use  $\downarrow$  or  $\uparrow$  to set **Trigg.E = \_ of 5** and press  $\bullet$  to confirm.

Each inspiration triggered by spontaneous breathing will be indicated by a star in the first line of the control panel display.

- Press  $\downarrow$  to select **Ramp = 0.\_ s** and activate by pressing  $\bullet$   
Ramp: Time for pressure increase from Pexp to Pinsp.  
0.0 s = Steepest pressure increase
- Use  $\downarrow$  or  $\uparrow$  to set **Ramp = 0.\_ s** and press  $\bullet$  to confirm.

#### Setting ranges for CPAP/ASB

Parameter	Range	Increments
Pinsp	5 to 20 mbar	0.5 mbar
Pexp	2* to 18 mbar	0.5 mbar
Trigg.I	1 to 5**	1 step
Trigg.E	1 to 5**	1 step
Ramp	0 to 1.0 s	0.1 s

\* At extreme settings (Pinsp = 20 mbar, Pexp = 2 mbar)  
Pexp. = 2 mbar may not be attained in some cases.

\*\* Trigger sensitivity: see Technical Data, page 28

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## Apnea ventilation

### Only for RespiCare ST

CPAP/ASB mode can optionally be combined with apnea ventilation. In this case, the unit switches over to time controlled ventilation if it does not detect any spontaneous breathing within a definable time. When apnea ventilation is selected, ApneaVent appears in the first line. The unit stops time controlled ventilation when spontaneous breathing is restored, and at the latest after one minute.

- Press **↓** to select **ApneaVent = \_ \_**, and activate by pressing **●**
- Press **↓** to set **ApneaVent = On**, and press **●** to confirm.
  
- Press **↓** to select **AV-Delay = \_ \_ s**, and activate by pressing **●**  
AV = Apnea Ventilation  
AV-Delay = Time after end of inspiration phase.
- Press **↓** or **↑** to set **AV-Delay = \_ \_ s**, and press **●** to confirm.
  
- Press **↓** to select **AV-Freq = \_ \_ /min**, and activate by pressing **●**  
AV-Freq = Ventilation frequency during apnea ventilation
- Press **↓** or **↑** to set **AV-Freq = \_ \_ /min**, and press **●** to confirm.
  
- Press **↓** to select **AV-Insp = \_ %**, and activate by pressing **●**  
AV-Insp = Inspiration as a proportion of the overall breathing phase in %
- Press **↓** or **↑** to set **AV-Insp = \_ %**, and press **●** to confirm.

### Adjustment ranges for apnea ventilation

Parameter	Range	Increments
AV-Delay	2 to 30 s	1 s
AV Freq	6 to 40 /min	1 /min
AV Insp	20 to 50 %	2 %

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## Activate / deactivate alarm

### disconnection / apnea alarm

- Press **↑** or **↓** to select **Alarm** and activate by pressing **●**
- Select required status **Alarm = active** by pressing **↑** or **Alarm = inactive** by pressing **↓** and press **●** to confirm.

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
- Press or to select Patient = \_\_\_h.

Reset application time to 0:

- Press or to select Patient = \_\_\_h, press to confirm and
- press to reset display to 0 and press to confirm.

## When oxygen is applied

Oxygen (O<sub>2</sub>) actively supports combustion.

**No smoking, no naked flames!**

Use the oxygen adapter, see order list, page 32.

Other adapters must not be used!

The oxygen adapter prevents oxygen leaking into the unit when the fan is defective (no airway pressure).

**Oxygen must not be directed into the unit or intake port, due to the risk of fire!**

- 1 Plug the oxygen adapter onto the socket of the RespiCare. **Ensure that the valve housing points upwards.**
- 2 Plug the ventilation hose onto the socket of the oxygen adapter.
- 3 Plug the oxygen metering hose onto the socket of the adapter.
- Max. permissible oxygen flow: 10 L/min.

Switch on the RespiCare:

- Press on/off switch
- press button , the green LED goes out.
- 4 The ball in the oxygen adapter must move upwards when the fan is switched on
- Turn on the oxygen supply.

After use:

- **Turn off the oxygen supply before switching off the RespiCare.**

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## **When a humidifier is used**


In case of irritation or dryness of the upper airway a humidifier can be used in order to prevent this. For installation of a humidifier and its breathing tube system, carefully follow the instructions in the user manual of the breathing systems (5730161/5730176).

**When using a humidifier it is not allowed to use the E-Vent (2M86522).**

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

## Operation

### Start operation

- 1 Press on/off switch on rear of unit – until it engages = ON.  
The RespiCare performs a self-test.

- 2 The green LED in the button lights up and the unit is now on standby.
- 3 The min. (minutes) LED lights up while the unit is on standby.
- 4 The bar graph indicates the set time for the "soft start" function:  
Each gradation corresponds to 2 minutes.  
To set the soft start function, see page 14.
- 2 Press the button. The green LED goes out.
- The soft start time display goes out, and the blower starts delivering air.  
**After starting up the unit suppresses any alarms for 1 minute.**


### Using the soft start function

#### Can only be set in standby mode.

For slowly and "softly" increasing the set airway pressure over a variable time of 0 to 30 minutes.

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- 1 Press  to switch to standby. The green LED in the button lights up.
- 2 The **min.** LED lights up.
- 3 The time (minutes) elapsed before reaching the airway pressure set on the control panel is indicated on the bar graph. Each gradation corresponds to 2 minutes.
- 4 Press  or  to set the required time. The time changes by 2 minutes whenever one of the two buttons is pressed.
- 1 Press  to switch to operating mode. The green LED goes out and the set time becomes effective.
- 3 The current airway pressure is now displayed on the bar graph. The soft start function is effective whenever changing over from standby to operation and is also retained after switching off the unit.

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## Nasal mask

- For RespiCare, use only nasal masks without outflow port. CO<sub>2</sub> is purged via the SilentFlow.
- Put on nasal mask in accordance with separate Instructions for Use.

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
## In the event of a power failure

a continuous tone sounds for at least 2 minutes. This acoustic alarm cannot be switched off.

- Remove the nasal mask to avoid any accumulation of CO<sub>2</sub> and a shortage of oxygen!

Further information for trouble-shooting can be found in the table "Fault – cause – remedy" on page 19.

## In the event of an alarm

Information for trouble-shooting can be found in the table "Fault – cause – remedy" on page 19.

### 1 Audio alarm\* and flashing red alarm LED ⚠:

Disconnection or apnea.

The audio alarm can be suppressed for 2 minutes:

### 2 Press ⚠

- Eliminate the cause of the disconnection.

The flashing red alarm LED ⚠ goes out and the audio alarm stops when the cause of the fault has been remedied.

The following message flashes in the first line when the control panel is connected:

**Disconnect**

or

**ApnoeVent**

- Remedy the cause of the alarm. The display will clear automatically.

Parameters can be set even while an alarm is displayed.

### 1 Continuous tone and red alarm LED ⚠ lights up steadily:

Equipment fault.

Individual segments of the bar graph light up constantly.

The alarm tone cannot be switched off. The fan switches off.

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\* Activate / deactivate alarm, see page 11.

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- Remove the nasal mask to avoid any accumulation of CO<sub>2</sub> and a shortage of oxygen!
- Read off the scale value for the 2 illuminated gradations on the bar graph (example: 4 and 12) and report these to your service partner. The numbers refer to the equipment fault.

The following message appears on the display of the control panel when the latter is connected:

**Equipment fault, technical defect**  
with the same numbers as on the bar graph.

#### Continuous tone without visual indication


draws attention to a power failure. This alarm tone cannot be switched off.

- Remove the nasal mask to avoid any accumulation of CO<sub>2</sub> and a shortage of oxygen!

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#### Service indicator

The required inspection and maintenance intervals are monitored and signalled automatically by the RespiCare unit:

- 1 Yellow LED lights up :  
The unit must be inspected and serviced by experts within the next two months. The LED does not go out until the unit has been inspected and serviced.

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

## End of operation

When oxygen has been applied:

- **Turn off the oxygen supply first** and then switch the unit to standby or off.

In the event of brief interruptions up to approx. 15 minutes:

- 1 Press  and the green LED lights up. The unit remains on standby.

For longer breaks:

- 1 Press  and the green LED lights up. The unit is on standby.
- 2 Press switch on rear of unit completely = OFF.  
The unit is now switched off.

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## Fault – Cause – Remedy

If a fault occurs, the following table can help to identify the underlying cause and remedy it.

Fault	Cause	Remedy
Intermittent tone, red alarm LED  flashes.	Disconnection	Ensure that ventilation hose, test hoses, nasal mask are fitted without leaks or kinks.
Airway pressure is not reached during inspiration.	Coarse dust filter/ fine dust filter clogged.	Clean coarse dust filter, page 24, or replace it. Replace fine dust filter.
Airway pressure has increased significantly.	Leak in pressure measuring hose or hose has slipped off socket.	Ensure that there are no leaks or kinks in the test hoses.
Unit does not switch over to pressure assistance in CPAP/ASB mode.	Leak in test hoses or hoses buckled.	Ensure that there are no leaks or kinks in the test hoses.
Trigger function or pressure control defective.	Water residue in test hoses.	Fully dry out the test hoses.
Set pressure parameters too low after switching on.	Soft start function effective.	Check setting of soft start function, page 14.
Continuous tone, red alarm LED  lights up continuously.	Equipment fault.	Call Dräger Medical or an authorized expert.
Continuous tone, all lamps off.	No mains power.	Check connection of power cords.
Airway pressure not displayed on bar graph.	Bar graph illumination deactivated.	Press  to activate bar graph illumination.

If the faults cannot be remedied:

- Use a spare unit and call an authorized expert.

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

Otherwise, risk of electric shock.

## Stripping down

- 1 Remove SilentFlow from ventilation hose.
- 2 Pull ventilation hose off its socket and unplug the connector to the oxygen adapter if applicable.
- 3 Remove measurement hoses from sockets.

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Otherwise, risk of electric shock.

### Control panel, basic unit, cable

- Remove any dirt with a disposable cloth. Liquid must not be allowed to leak into the connections

### Breathing tube with integrated measuring tubes, SilentFlow and oxygen adapter

Specific cleaning and disinfecting procedure information is found in the user manuals of the accessories.

### General cleaning instructions

- Remove all visible deposits.
- Clean the accessories with warm soapy water.
- Rinse with boiled lukewarm water.
- Dry parts thoroughly.  
Water in the test hoses or in the oxygen adapter may impair the operation of the unit.  
Water in the ventilation hose or in the SilentFlow and measuring adapter can reduce the air flow through the test outlet and allow the CO<sub>2</sub> content of the inspiratory air to rise.

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## Cleaning / disinfection in the hospital

**The disinfectants mentioned here are not suitable to be used at home.**

Use a surface disinfectant, such as a product based on:

- Aldehydes
  - Alcohol
  - Quaternary ammonium compounds.
- The following are not suitable:
- Compounds containing phenol
  - Compounds which give off halogens
  - Strong organic acids
  - Compounds which give off oxygen.

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Users in the Federal Republic of Germany are advised to use the disinfectants specified in the current DGHM\* list. This list also specifies the active agent in each disinfectant. Products based on the active agents specified above are recommended for users in countries in which the DGHM list is not available.

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In case the device is used for another patient:

- make sure the device is properly disinfected, according the cleaning instructions, see Order List at page 32.

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#### **Control panel, basic unit, cable**

- Wipe-disinfecting  
The following disinfectants can e.g. be recommended:  
Incidur®, Frekanol®, Gigasept FF®  
Note the manufacturer's directions. Liquid must not be allowed to leak into the connections.

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#### **Breathing tube with integrated measuring tubes, SilentFlow, oxygen adapter**

- Disinfect in a bath, e.g. with Gigasept FF®  
The parts must subsequently be rinsed in clear boiled water. Thoroughly shake off all remaining water. Prevent liquid entering the measurement tubes. Allow the parts to dry completely. Water in the test tubes may impair the unit's correct operation.

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#### **Nasal mask**

- Must be cleaned in accordance with the separate Instructions for Use.

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


02237286

## Check Readiness for Operation

The unit's readiness must be checked whenever it has been cleaned and disinfected.

After assembling the unit ready for operation:

- 1 Switch on unit:  
Press button on rear of unit – until it engages = ON.
- 2 The green LED in the button lights up.
- 3 LED **min.** lights up.
- 4 The time for the soft start is shown on the bar graph in minutes.
- 5 Press and set the time for the soft start function to zero.
- Put on a suitable nasal mask and
- 2 press to switch on the fan.

02237286

- 6 LED **mbar** lights up.
- 7 Press :
- 8 The airway pressure must appear on the bar graph.
- 9 Press , the green LED lights up and the fan is switched off.
- Remove nasal mask.
- 10 Use the buttons and to set the time for the soft start.

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## Maintenance Intervals

### Coarse dust filter

Must be cleaned every three weeks.  
The filter must be replaced if damaged,  
at the latest every six months.  
It can be disposed of as household  
refuse.

### Fine dust filter


Must be replaced after approx.  
500 hours of application.  
It must be replaced at the latest after  
six months.  
The filter can be disposed of as house-  
hold refuse.

### Inspection and maintenance

As indicated by the service indicator  
 : see page 17.  
At least once per year.

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### Clean coarse dust filter

- Every three weeks.
- Draw foam insert out of filter housing.
- Wash foam insert in soapy water, rinse thoroughly in clear water and dry completely.
- Press foam insert back into filter housing.
- Damaged foam inserts must be replaced.

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### Disposal of the unit

- when it reaches the end of its useful life.
- RespiCare must be appropriately disposed of in agreement with the competent specialized waste disposal company.  
All applicable legal regulations must be obeyed.

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## Customer Service Mode

For selecting the language in which messages are displayed on the control panel.

For testing the alarm tone and lamps.

For indicating the hours of operation.

For indicating the software version.

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## Select customer service mode

- Disconnect RespiCare from the patient and connect the control panel.
  - Press on/off switch completely = OFF.
  - Press on/off switch – until it engages = ON simultaneously
  - press  for at least 3 seconds.
  - The RespiCare unit runs a self-test before the following message appears:  
**CSM - Customer Service Mode**
  - Release the  button.
- To leave customer service mode:
- Press on/off switch completely = OFF.

## Set language

- Select customer service mode  
Display:  
**Adjust language**  
    

To activate:

- Briefly press  and the set language appears on the second line.  
Display (example):  
**Adjust language**  
**deutsch**
- The following languages can be set as required using the  and  buttons:  
**English**  
German (deutsch)  
French (français)  
Spanish (español)  
Italian (italiano)  
Portuguese (português)  
Dutch (nederlands)
- Confirm language = press 

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## Display operating hours

To call up operation hours:

- Select customer service mode, see page 25.
- Press to select.

Display:

### Operation hours

- 003 +

Press to select:

Tot: Total hours of operation

or

Pat: Period of use by patient

or

Service: Hours of operation since unit was last serviced

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## Display software version

To display software version:

- Select customer service mode, see page 25.
- Press to select.

Display:

### Firmware version

- 004 +

Press to select:

Main device: Software version of the basic unit

Remote device: Software version of the control panel

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


## What's what

### Front view

- 1 Socket for measuring tube
- 2 Socket for measuring tube
- 3 Socket for ventilation hose
- 4 Yellow LEDs for bar graph
  - min. LED lights up: time for soft start (minutes)
  - mbar LED lights up: airway pressure
- 5 Bar graph for airway pressure (mbar) or time for soft start (min.)
- 6 Two buttons for setting the time for the soft start  
 $\downarrow\circlearrowleft$  = decrease,  $\uparrow\circlearrowright$  = increase

- 7 Red alarm LED  $\Delta$
- 8 Button  $\triangleleft$  to suppress the alarm tone for approx. 2 minutes
- 9 Yellow LED for service required  $\wrench$
- 10 Button  $\square$  to switch the fan on and off (standby)
- 11 Button  $\odot$  to switch the bar graph illumination on and off
- 12 Control panel for settings and data
- 13 Button  $\bullet$  to activate and confirm settings
- 14 Two buttons  $\downarrow$  and  $\uparrow$  for making selections and settings
- 15 Display for settings, data and messages

### Rear view

- 16 Handle for carrying the unit
- 17 On/off switch
- 18 24 V DC input
- 19 Connection for control panel
- 20 Coarse dust filter / fine dust filter
- 21 Silencer
- Rating plate on underside (not illustrated)

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### Ambient conditions

In operation

Temperature	5 to 35 °C
Air pressure	700 to 1100 hPa
Rel. humidity	10 to 95 % no condensation

In storage

Temperature	-20 to 70 °C
Air pressure	600 to 1200 hPa
Rel. humidity	10 to 95 % no condensation

### Settings

#### Ventilation mode CPAP

Airway pressure	5 to 20 mbar
<b>Ventilation mode CPAP/ASB</b>	
Inspiration pressure Pinsp	5 to 20 mbar *
Expiration pressure Pexp	2 to 18 mbar
Trigger sensitivity Trigg.I	5 settings (approx. 10 to 50 L/min)
Trigger sensitivity Trigg.E	5 settings (approx. 15 to 75 L/min below peak flow)
Ramp	11 settings (0 to 1.0 s)
apnea Ventilation	ON / OFF
Latency time	2 to 30 s
AV-Frequency	6 to 40 /min
AV-Insp.	20 to 50 %
Time for soft start function	0 to 30 min

### Device parameters

Device compliance (with patient hose system)	≤0.6 mL/mbar
Maximum pressure limit	40 mbar (on patient)
Inpiration and expiration resistance	max. 6 mbar at 60 L/min
Airway pressure display	
Range	0 to 30 mbar
Resolution	2 mbar
Accuracy	±2 mbar

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**Monitors**

## Disconnection/apnea

Triggered when alarm limit is exceeded for at least 5 seconds.  
Alarm limit for airway pressures >10 mbar: actual value at least 5 mbar below set value.  
Alarm limit for airway pressures ≤10 mbar: actual value less than 50 % of set value.  
The alarm is suppressed for approx. 1 minute after switching on the fan and when changing the ventilation mode.

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24 V DC ± 5 %  
max. 50 VA  
Schutzklasse II

**Operating parameters**

Rated voltage  
Power consumption  
Classification as per  
EN 60601-1

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AC/DC converter  
Rated voltage  
Power consumption  
Mean sound pressure at  
10 mbar  
Free-field measurement over a  
reflecting surface, distance 1 m  
  
Dimensions (W x H x D)

24 V DC ±5 %  
max. 50 VA  
Protection class II  
  
100 to 240 VAC ±10 %  
max. 60 VA  
max. 30 dB(A)

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255 x 176 x 225 mm  
150 x 65 x 75 mm  
3,5 kg  
0,25 kg  
ø22, nicht konisch

Basic unit  
(incl. hose sockets)  
Control panel  
  
Weight  
Basic unit  
Control panel  
  
Connection for ventilation hose ø22, not tapered

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3,5 mm Stereo-Klinken-  
stecker  
  
Outer connector  
Inner connector  
Shell

**Analog Output**

Connector type  
Wiring  
Outer connector  
Inner connector  
Shell

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Electromagnetic compatibility  
EMC

Tested to EN 60 601-1-2

Classification as per EC  
Directive 93/42/EEC Annex IX

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11-634

UMDNS-Code  
Universal Medical Device  
Nomenclature System –  
Nomenclature for medical  
products

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## Abbreviations and Symbols

ApneaVent	Pressure controlled, time controlled ventilation in the event of apnea
AV-Delay	Selectable delay time between the end of the inspiration phase and the start of apnea ventilation
AV-Freq	Ventilation frequency during apnea ventilation
AV-Insp.	Inpiration phase as a proportion of the overall respiration phase during apnea ventilation
ASB	<b>A</b> sisted <b>S</b> pontaneous <b>B</b> reathing Spontaneous breathing with synchronized assistance
Op. time	Operating time counter for the total operating time of the unit in hours
CPAP	<b>C</b> ontinuous <b>P</b> ositive <b>A</b> irway <b>P</b> ressure <b>B</b> reathing Spontaneous breathing with continuous positive pressure in the airways
Freq	Ventilation frequency
Insp %	Percentage inspiration in the respiratory phase
mbar	Bar graph indication of the airway pressure (mbar)
min.	Bar graph indication of the time (minutes) for the soft start function
Patient	Resettable display of the duration of application in hours
Paw	Airway pressure
Pexp	Expiration pressure
Pinsp	Inspiration pressure
Trigg.E	Expiratory flow trigger
Trigg.I	Inspiratory flow trigger
	Standby button for switching the fan on and off
	Setting the time for the soft start
	Switches the bar graph illumination on and off
	Making selections and settings
	Activating and confirming settings
	Suppresses the alarm tone for approx. 2 minutes
	Red alarm LED
	Yellow LED for Service
	Note the Instructions for Use!
	On/off switch
	Type BF = body floating

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<b>Designation/Description</b>	<b>Part No.</b>
RespiCare, consisting of: Basic unit, coarse dust filter, power cord, AC/DC converter, ventilation hose, SilentFlow, measurement hose system, Instructions for Use and carrying bag.	
RespiControl	57 30 310
Tube system:	
Tube system, disp.	57 02 176
Tube system, reus.	57 30 161
SilentFlow	57 30 003
Tube humid disp.	57 02 178
Tube humid reus.	57 02 177
Measuring adapter	57 02 200
Mask connection, set of 5, standard	57 30 315
Mask connection, set of 5, ISO 15	57 30 316
Coarse dust filter	57 30 317
Fine dust filter	57 30 461
Carrying bag	84 13 369
Silencer	57 30 320
Mains cable euro	57 30 459
Mains cable GB	57 30 458
<b>Accessories for use with the basic unit</b>	
Control panel with analog outputs	57 30 462
RespiLink software	57 01 902
Extension cable for control panel, 30 m	84 13 357
AC/DC converter	57 30 457
DC/DC converter (12 V DC to 24 V DC)	57 30 460
Cable for analog output	57 30 465
Oxygen adapter	84 13 397
Trolley	2M 86 029
Table for trolley	2M 86 030
Instructions for Use de/en/fr/nl/es/pt/tr/it	57 30 440
Cleaning instruction	57 02 068

## Patientenpass

Name \_\_\_\_\_

Vorname \_\_\_\_\_

Geburtstag \_\_\_\_\_

Anschrift \_\_\_\_\_

Telefon \_\_\_\_\_

Maske (Typ/Größe) \_\_\_\_\_

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Zubehör \_\_\_\_\_

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Alter \_\_\_\_\_ Jahre

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Datum \_\_\_\_\_

Gewicht \_\_\_\_\_

kg

Blutdruck \_\_\_\_\_

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CPAP/ASB		
Pinsp:		mbar
Pexsp:		mbar
Trigg.l:		
Trigg.E:		
Rampe:	0,	s
Alarm		aktiv/inaktiv*
ApnoeVent		AUS/EIN*
AV-Latenz:		s
AV-Freq:		min
AV-Insp:		%

\* streichen

## Patient Passport

Surname	First name	Date of birth
Address		Tel.
Mask (model/size)	Hood/harness size* Gr.	Accessories

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© Dräger Medical b.v.