# RespiCare VV

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### Home ventilator

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## Contents

<b>For Your Safety and that of Your Patients</b> .....	49	<b>Maintenance intervals</b> .....	75
Precautions .....	49	Battery maintenance .....	76
<b>Intended Use</b> .....	50	Disposal of internal battery .....	77
<b>Preparing for Use</b> .....	51	Disposal of the device .....	77
Connect power supply .....	51	<b>Fault – Cause – Remedy</b> .....	78
Connect hose system .....	53	<b>Preparation for renewed use</b> .....	80
Using a disposable bacteria filter (optional) .....	53	Dismantling .....	80
Using a heater/humidifier (optional) .....	54	Cleaning .....	80
O <sub>2</sub> enrichment .....	54	Disinfecting .....	81
Connecting an O <sub>2</sub> monitor (optional) .....	55	<b>What's what</b> .....	82
Connecting a PEEP valve (optional) .....	56	<b>Technical Data</b> .....	84
Connecting a remote alarm (optional) .....	56	<b>Symbols and abbreviations used</b> .....	88
<b>Before use</b> .....	57	<b>Operating principle</b> .....	90
Quick check .....	57	IPPV/IPPVASS mode .....	90
Preparing to set up (SET UP) .....	57	IPPVPRS mode .....	91
<b>Operation</b> .....	61	SIMV .....	91
IPPV mode .....	61	Maximum safe pressure .....	92
Selecting IPPV ASS mode .....	65	<b>Order List</b> .....	93
Selecting a different ventilation mode .....	65		
IPPVPRS mode .....	65		
SIMV .....	68		
Locking settings .....	71		
Operation with internal battery .....	71		
Switching off the device .....	72		
In the event of an alarm .....	73		
Maintenance menu .....	73		
Display operating hours .....	74		
Set contrast .....	74		

## For Your Safety and that of Your Patients

### Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observation of these instructions. The apparatus is only to be used for purposes specified here.

### Maintenance

The apparatus must be inspected and serviced yearly by trained service personnel (and a record kept).

Repair and general overhaul of the apparatus may only be carried out by trained service personnel.

We recommend that a service contract be taken out with Dräger Medical Service or a Dräger Medical authorised dealer and that all repairs also be carried out by them.

Only original Dräger Medical parts should be used for maintenance.

Observe chapter "Maintenance Intervals".

### Accessories

Do not use accessory parts other than those in the order list.

### Not for use in areas of explosion hazard

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

### Safe connection with other electrical equipment

Electrical connections to equipment which is not listed in these Instructions for Use should only be made following consultations with the respective manufacturers or an expert.

### Liability for proper function or damage

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus is serviced or repaired by personnel not employed or authorized by Dräger Medical or if the apparatus is used in a manner not conforming to its intended use.

Dräger Medical cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of Dräger Medical's terms of sale and delivery are likewise not modified by the recommendations given above.

Dräger Medical  
a division of Dräger Nederland B. V.

## Precautions

**This device should only be used as prescribed by a doctor and with the settings and monitoring parameters specified by the doctor.**

**During assisting ventilation, watch for possible risks of hyperventilation!**

**This device must not be used with flammable gases or anaesthetic agents. Risk of fire!**

**Avoid pollutants in the ambient air. RespiCare VV uses ambient air for ventilation. Pollutants would be inhaled by the patient and could be harmful to the patient.**

**Before using the device, always check the internal battery charge level, see "Quick check", page 57.**

## Intended Use

Time-cycled, constant-volume ventilator for home ventilation.  
Suitable for patients with a body weight of approximately 10 kg upwards for

- controlled ventilation IPPV MODE
- controlled assisting ventilation IPPVASS MODE
- controlled assisting ventilation with PEEP IPPVASS MODE/ PEEP\*
- controlled assisting ventilation with pressure control IPPVPRS MODE
- controlled assisting ventilation with pressure control and PEEP IPPVPRS MODE/PEEP\*
- synchronized intermittent mandatory ventilation SIMV
- synchronized intermittent mandatory ventilation with PEEP SIMV/PEEP\*

With monitoring of:

- Airway pressure

RespiCare VV is suitable for non-invasive and invasive ventilation.

Applications:

- Long-term ventilation at home
- Hospital use, e.g. in preparation for home ventilation
- For mobile use, e.g. in a vehicle or wheelchair

---

\* Optional PEEP valve required

## Preparing for Use

- Place RespiCare VV upright on a level surface, such as a bedside table or wall shelf. Ensure that the device cannot fall down.  
or place on trolley ("Accessories" on the Order List on page 93). Secure the trolley so that it cannot roll away.  
In specific situations, RespiCare VV can also be used in other positions, e.g. on its side.
  - Do not use the device straight away after it has been stored or transported outside the range of prescribed ambient conditions. Wait for the device to reach room temperature.
- **Never place objects on the device.**
  - **Protect the device from damp and liquids.**


## Connect power supply

RespiCare VV can operate with various power sources:

- mains power 110 to 230 V AC
- DC 12 to 36 V DC max. 4 A
- external battery (optional)

### Mains power

- For stationary use, connect the device to the mains.  
Voltage range: 110 to 230 V AC.
- 1 Connect power lead to socket at rear of device and plug into a properly installed shockproof mains socket.


00137178

### DC supply

Instead of an AC supply, the device can also be connected to a DC supply, for instance for use when in a wheelchair or car.

- 2 Insert DC supply cable ("Accessories" on the Order List on page 93) at the rear of the device and connect up to the DC supply, 12 V to 36 V DC, maximum 4 A.


The internal battery is not intended to provide a power supply during transport or when in mobile use. It is designed to be an emergency power supply in case of mains failure or when changing over from one type of supply to the other.

### External battery (optional)

Instead of mains power, RespiCare VV can also be powered by an optional battery pack.

- Follow the separate instructions for use of the battery pack and charger.
- 1 Connect battery lead to DC input at the rear of the RespiCare VV and
  - 2 insert into charger port.
  - 3 Connect cable between charger and battery pack.
  - 4 Whenever possible, always connect the charger to the mains.


## Connect hose system

The hose system supplied is shown here.

- Use a new disposable hose system.
- Other hose systems should only be used if they are of the same diameter and length as this hose system and a manufacturer's certificate of compatibility for use with RespiCare VV is supplied.

Hose systems of different diameters and lengths can impair the performance of the device.

- 1 Attach elbow connector to inspiration nozzle of hose system and to cone on device.
- 2 Connect control line of breathing valve to nozzle on device.


## Using a disposable bacteria filter (optional)

See "Consumables" on the Order List on page 93

- Connect filter to cone on device and attach elbow connector of breathing valve to filter cone.
- Replace disposable bacterial filter if the airway pressure rises sharply during inspiration. This indicates that the filter is clogged.
- Replace the disposable bacterial filter at least after every 800 hours' operation or every 60 days.

A dirty or clogged filter could impair the operation of the device and put the patient at risk!

- Used filters can be disposed of as domestic waste.


## Using a heater/humidifier (optional)

The humidifier moistens and warms up the air breathed.

**System compliance increases when a heater/humidifier is used.**

- **Check the tidal volume.**
  - Use the heater/humidifier indicated under "Accessories" on the Order List on page 93.
  - Prepare the heater/humidifier according to the separate instructions for use and place it near the RespiCare VV.
- 1 Attach elbow connector to inspiration nozzle of hose system and to cone on device.
  - 2 Attach breathing valve control line to nozzle on device.
  - 3 Detach inspiration hose from inspiration nozzle and Y-piece.
  - 3 Connect the heater/humidifier connecting hoses to the inspiration nozzle and the Y-piece.
- **When using a heater/humidifier, do not use a humidity/heat exchanger. There is a risk of increased condensation, impairment of device operation and putting the patient at risk.**

01437356


## O<sub>2</sub> enrichment

When concentrations between 25 and 50 % by volume are necessary, an O<sub>2</sub> enrichment set should be used:

- 4 Insert filter disc (disposable microbial filter) into the housing of the O<sub>2</sub> enrichment set, with the rough side of the filter disc towards the RespiCare VV housing.
- Replace the filter:
    - every 500 hours in outdoor use,
    - every 1200 hours in indoor use.


**A dirty filter could impair the operation of the device and put the patient at risk.**

- Used filters can be disposed of as domestic waste.
- 5 Press O<sub>2</sub> enrichment set housing firmly on to the air inlet.
  - 6 Attach O<sub>2</sub> supply hose to angled nozzle on housing and connect to an O<sub>2</sub> supply, such as an O<sub>2</sub> cylinder with a flow meter and flow valve, or to an O<sub>2</sub> concentrator, such as the Dräger Permax SC.
- **Optional O<sub>2</sub> monitoring with an alarm device should be used when using an O<sub>2</sub> enrichment set.**

00537356

- Using an optional O<sub>2</sub> monitor, see page 55.
- If O<sub>2</sub> enrichment is not necessary:
- 1 Fit silicone cap supplied to the angled nozzle.


### Setting the O<sub>2</sub> concentration

The O<sub>2</sub> flow necessary for a desired O<sub>2</sub> concentration depends on the minute volume applied.

- To set, refer to the mixing diagram on the O<sub>2</sub> enrichment set housing.
- Example:  
 Desired O<sub>2</sub> concentration: 35 % by vol.  
 Minute volume applied: 11 L/min  
 Necessary O<sub>2</sub> flow: 2 L/min

This table is valid for an O<sub>2</sub> supply with 100 % by vol.  
 When using the Permax O<sub>2</sub> concentrator, lower values apply because Permax delivers max. 95 % by vol. O<sub>2</sub>.

**Oxygen (O<sub>2</sub>) is a fire accelerant.**

- No smoking, no naked lights.


- Keep O<sub>2</sub> fittings free from oil and grease. Fire hazard.

**If the mask is leaking, the O<sub>2</sub> concentration set may not be reached.**

### Connecting an O<sub>2</sub> monitor (optional)

When using increased O<sub>2</sub> concentrations


- Use the O<sub>2</sub> monitor indicated under "Accessories" on the Order List on page 93.
- Detach the breathing valve and elbow connector.
- Fit O<sub>2</sub> monitor to the cone.
- Attach elbow connector to the O<sub>2</sub> monitor cone.
- Use the O<sub>2</sub> monitor according to the manufacturer's instructions for use.

### Connecting a PEEP valve (optional)

To use a positive end-expiratory pressure PEEP

- 1 Attach elbow connector to outlet nozzle of breathing valve and
  - 2 push PEEP valve into the elbow connector.
- Follow instructions for use of the PEEP valve.

### Connecting a remote alarm (optional)


This provides an alarm elsewhere if a fault occurs.

Use the remote alarm indicated under "Accessories" on the Order List on page 93.

- Switch on RespiCare VV.
- 3 Connect the connecting cable to the remote alarm and to the rear of RespiCare VV.
  - 4 Switch remote alarm to I. It is now on:  
If an alarm occurs on RespiCare VV, the audible alarm sounds and the red alarm lamp flashes. Disconnection of the connecting cable also triggers a remote alarm.

If an alarm occurs on the remote alarm:

- Immediately investigate and remedy the cause at the RespiCare VV.

## Before use

- Do not connect patient yet.
- Leave patient connection of breathing system open.

### Quick check


- Set ON/OFF switch on the rear to  position = ON
- Start up alarm is activated.
- 1 Unkink pressure measuring line:
- High pressure alarm is activated.
- Detach pressure measuring line: after about 15 seconds the low pressure alarm is activated.
- Reconnect pressure measuring line to the apparatus.

#### Disconnect power supply:

- Remove mains plug or disconnect DC supply:
- The "no power" alarm is activated and the internal battery charge illustration is displayed. See page 71, "Operation with internal battery".
- Reconnect power supply.
- Press and hold »●« key for 3 seconds, then
- turn ON/OFF switch to  = OFF.

#### If the battery is discharged:


- Before using the device, charge the internal battery, see "Battery maintenance", page 76.

01137356

## Preparing to set up (SET UP)

- Press and hold down the »↓« and »↑« keys at the same time and turn the ON/OFF switch to  = ON,
- hold down the »↓« and »↑« keys until the **SET UP** menu appears on the display.
- Use the »↓« or »↑« keys to select the parameters. The little box to the left of the parameter is highlighted.
- Use the »+« or »-« keys to set the values. Ventilation mode IPPV and the SIGH VOLUME and SIGH FREQ parameters additionally need to be confirmed by pressing the »●« key.

02437178

### Set language (LANGUAGE)

The following are available:


- French (FRANCAIS)
- English (ENGLISH)
- Italian (ITALIANO)
- German (DEUTSCH)
- Dutch (NEDERLANDS)
- Spanish (ESPAÑOL)
- Portuguese (PORTUGUESE)

After selecting, the messages will appear in the selected language.

### Select display mode (SCREEN MODE)

- AUTO** The screen display disappears automatically after one minute, irrespective of the ventilation mode (in the event of an alarm it appears automatically, or by pressing any key).
- ON** The screen display remains visible all the time

### Store alarm (MEMORY ALARM)

- NO** The device does not store the alarm. Once a fault is remedied, the audible alarm is silent, the alarm message disappears and the alarm indicator lamp goes out.
- YES** Once a fault is remedied, the audible alarm is silent and the alarm indicator lamp goes out. Alarm messages have to be reset manually by pressing and holding down the »  « key for 3 seconds.


### Setting the frequency alarm (FREQUENCY ALARM)

This is used to monitor the frequency rise during controlled, assisting ventilation IPPVASS or controlled, assisting ventilation with pressure control IPPVPRS in order to detect such things as hyperventilation.

The following settings are available:

- NO** No monitoring
- +5** An alarm is activated if the current frequency is 5 %, 10 %, 15 % or 20 % higher than
- +10** the set frequency.
- +15**
- +20**

The frequency alarm +10 is automatically activated if a positive trigger is selected (TRIGGER LEVEL + YES).


02497356

### Selecting the pressure curve (PRESSURE CURVE)

The pressure curve can be displayed with the pressure distribution over a complete breathing cycle in the lower field of the screen provided for alarm messages. The display is refreshed at each new breathing cycle.

The horizontal dashed line represents the zero line.

The maximum amplitude is 60 mbar. The curve is displayed in the form of separate bars.

02137178

**YES** The pressure curve is displayed in the alarm message field. An alarm occurring overwrites the pressure curve.

**NO** The pressure curve is not displayed.

### Setting the sigh (SIGH)

The sigh function is operative in the ventilation modes IPPV/IPPVASS. The flow pattern is rectangular. To preserve the set ratio I/E, the device extends the complete breathing cycle time.

**YES** The sigh is activated in the menu for ventilation mode IPPV. If the word "SIGH" appears light on a dark background, the function is operative at the time.

**NO** The sigh is switched off.

The following sigh parameters can be set:

The sigh volume can be set within a range from the single tidal volume  $V_T \times 1.0$  to twice the tidal volume  $V_T \times 2.0$  in increments of 0.1.

The sigh frequency can be set within a range of 50 cycles to 950 cycles in increments of 50 between two sighs.

### Selecting the ventilation mode (IPPV)

IPPV/IPPVASS, IPPVPRS and SIMV are available. When a selection is made and confirmed, the appropriate ventilation menu appears. The operating hours to date (TOTAL HOUR) are displayed briefly. The device ventilates using the last ventilation parameters set, see pages 61 onwards.

The SET UP setting remains active until a new setting is carried out.

Changing to ventilation mode:

- Select »IPPV« line,
- Set ventilation mode  
**VOL CTRL**  
or  
**PRS CTRL**  
or  
**SIMV**

**SET UP**

LANGUAGE      **ENGLISH**

SCREEN MODE      **AUTO**

MEMORY ALARM      **NO**

FREQUENCY ALARM      **NO**

PRESSURE CURVE      **YES**

SIGH      **YES**

VOLUME       $V_T \times 1.0$

SIGH FREQ.      50

IPPV      **VOL CTRL**

TRIG LEVEL +      **YES**

MAINTENANCE

03037566

**Selecting positive trigger (TRIG LEVEL +)**

To set a positive trigger pressure in the ventilation menus.

NO Positive trigger pressure cannot be set

YES Positive trigger pressure can be set  
+1 to +5 mbar in IPPV MODE  
+1 to +8 mbar in IPPVPRS MODE

**Maintenance**

This function allows access to the Maintenance menu, see page 73.

The maintenance menu is for service use only.

## Operation

- Switch ON/OFF switch on rear of device to  position = ON

The menu for the selected ventilation mode appears on the display with the last ventilation parameters set. The operating hours to date (TOTAL HOUR) are displayed briefly.

The little box in the **OFF** line is black (cursor position).

If a positive trigger level was previously selected, the little black box appears in the **ON** line and the message **TRIGGER LEVEL +** is displayed as a reminder.

- Press the »●« key to start and confirm operation with a positive trigger level.
- The device ventilates using the ventilation parameters displayed.

If no settings are made for one minute, the display disappears in "SCREEN MODE – AUTO" (for the convenience of the patient, for instance at night).

To switch on the display:


- Press any key.

In the event of an alarm, the display is activated automatically.

## IPPV mode

The following parameters are displayed on the screen for the IPPV ventilation mode:

- the ventilation parameters necessary for IPPV mode
- the alarm limits for high pressure alarm and low pressure alarm
- the pressure curve (if selected in the SET UP menu)
- alarm messages in the event of an alarm
- bar display of airway pressure with superimposed alarm limits for high pressure alarm and low pressure alarm (black triangles) and superimposed trigger threshold (white triangle)
- measured PEEP value below the bar display of airway pressure
- measured max. airway pressure above the bar display
- graphic display of internal battery charge level (when device is being powered by the battery).

To set the ventilation parameters:

- Press »↓« or »↑« key to select the parameter.  
The little box to the left of the parameter turns black.
- Press »+« or »-« keys to set the value.
- Press »●« to confirm.

During setting, the device calculates the speed necessary to inflate and deflate the built-in bellows. While this is taking place, for a duration of 3 to 15 seconds, the display reads:

**CALCULATING**

The device will not accept a setting which is not compatible with the parameters available.

In no case is ventilation interrupted.

**Frequency**

Setting range: 8 to 50 1/min

In the ACV ventilation mode, the frequency set operates as a minimal frequency. When the patient's efforts to breathe trigger the device, the frequency rises. The actual frequency is displayed. If a frequency alarm was set in the SET UP menu, an alarm is activated if the set alarm frequency is exceeded.

**Minute volume**


Setting range: 3 to 25 L/min

If a minute volume has been set, it is not possible to set the tidal volume. While ACV mode is operative, the minute volume set operates as a minimal value. This can be increased by means of the frequency rise triggered by the patient. The device displays the minute volume currently delivered.

**Tidal volume**

Setting range: 0.1 to 1.6 L

If the tidal volume has been set, it is not possible to set the minute volume.

02637/256

### Flow ramp

The following flow shapes can be selected:

- rectangular (R)
- sinusoidal (S)
- accelerating (A)
- decelerating (D)


The flow shape does not affect the frequency, tidal volume or ratio I/E.

The abbreviation and symbol for the selected flow shape are displayed on the screen.

Example:

R

for a rectangular flow shape.

### RATIO I/E

Setting range: 1:1 to 1:3

In the ventilation mode ACV, the efforts to breathe can trigger the device and alter the set RATIO I/E towards a ratio of 1:1.

### High pressure alarm (HIGH PRESSURE)

Setting range: 12 to 60 mbar

The set value is additionally displayed as a black triangle next to the bar display scale for airway pressure. If the airway pressure reaches this value, an alarm is activated and the bellows cease inflating and deflating in order to avoid unacceptably high pressure in the lungs.

### Low pressure alarm (LOW PRESSURE)

Setting range: 2 to 50 mbar

The set value is additionally displayed as a black triangle next to the bar display scale for airway pressure. If the airway pressure does not pass through this value during a breathing cycle, an alarm is activated.

### Trigger level

Setting range:

-1, -2, -3 mbar, OFF

Additionally, if TRIGGER LEVEL + was selected in the SETUP menu:

+1, +2, +3, +4, +5 mbar

-1 or +1 mbar = max. sensitivity

The set value is additionally displayed as a white triangle next to the bar display scale for airway pressure. Settings greater than -3 mbar cause the display "OFF" to appear and the ventilation mode changes from ACV MODE to CV MODE.

The device displays each ventilation stroke triggered by the patient's efforts to breathe by showing the trigger value light against a dark background.

- **Set a positive trigger level if a positive end-expiratory pressure PEEP is activated. If a positive trigger level is set without PEEP being active, the device triggers of its own accord and can cause hyperventilation.**

Instead of the selection in the SETUP menu, a positive trigger level can also be selected in the ventilation menu:

- Press »↓« key to select the »TRIGGER LEVEL« line,
- press »-« key to select the value -1,
- hold down the »-« key for 6 seconds to set a value between + 1 and + 8 mbar,
- press »●« to confirm.

### When PEEP is in use, please note:

- **Set the trigger level at least 3 mbar below the measured PEEP value to avoid self-triggering of the device and to enable the PEEP valve to be correctly set.**

### Selecting IPPV ASS mode

The following trigger levels may be selected:

-1, -2, -3 mbar,

or

+1, +2, +3, +4, +5 mbar if TRIGGER LEVEL + was selected.

### Switching off IPPV ASS mode

- Select »TRIGGER LEVEL« line,
- press »+« key to select **OFF**.

The device operates in ventilation mode IPPV without trigger.


### Selecting a different ventilation mode

- Press »↓« key to select »OFF« line,
- hold down »●« key for 3 seconds
- ON/OFF switch to  position = OFF.
- Press and hold down »↓« and »↑« keys at the same time,
- ON/OFF switch to  position = ON,
- hold down »↓« and »↑« keys until SET UP menu appears.
- Press »↓« key to select »MODE« line,
- press »+« or »-« key to select the desired ventilation mode,
- press »●« to confirm.

### IPPVPRS mode

The following parameters are displayed on the screen for the IPPVPRS ventilation mode:

- the ventilation parameters necessary for IPPVPRS mode
- the alarm limit for the low pressure alarm
- alarm messages in the event of an alarm
- pressure curve (if selected in the SET UP menu)
- graphic display of internal battery charge level (when device is being powered by the battery)
- bar display of airway pressure with superimposed alarm limit for low pressure alarm (black triangle), inspiration pressure value (black triangle) and superimposed trigger threshold (white triangle)
- measured PEEP value below the bar display of the airway pressure
- measured max. airway pressure above the bar display
- at the top of the pressure curve field the device displays the current tidal volume  $V_T$  and the current I/E ratio.

031372956

Setting the parameters

- Press »↓« or »↑« key to select the parameter.  
The little box to the left of the parameter turns black.
- Press »+« or »-« keys to set the value,
- press »●« to confirm.

**Minimal frequency**

Setting range: 4 to 20 1/min

The frequency set is a minimal frequency. The patient's efforts to breathe can trigger the device, thus raising the frequency. The current frequency is displayed.

**Inspiratory flow (FLOW INSP)**

Setting range: 0.4 to 1.5 L/s

**Inspiratory pressure (PRESSURE)**

Setting range: 7 to 40 mbar

The set value is additionally displayed as a black triangle next to the bar display of airway pressure.


When the ventilation pressure reaches the set inspiration pressure, the device stops inspiration and changes over to expiration (at plateau setting = 0 s).

**Plateau**

Setting range: 0 to 2 s

When plateau is operative, after reaching the set inspiratory pressure the device extends inspiration for the time set.

The plateau is inoperative if the inspiratory pressure is not reached.

02657336

**Low pressure alarm (PAW LOW)**

Setting range: 2 to 35 mbar

The low pressure alarm is additionally displayed by a black triangle next to the bar display of airway pressure.

If the airway pressure does not pass through the lower alarm limit in each breathing cycle, an alarm is activated after 15 seconds.

**Trigger level**

Setting range:

-1, -2, -3 mbar, OFF

Additionally, if TRIGGER LEVEL + was selected in the SETUP menu:

+1, +2, +3, +4, +5, +6, +7, +8 mbar

-1 or +1 mbar = max. sensitivity

The set value is additionally displayed as a white triangle next to the bar display scale for airway pressure. Settings greater than -3 mbar cause the display "OFF" to appear.

- **Set a positive trigger level if a positive end-expiratory pressure PEEP is activated. If a positive trigger level is set without PEEP being active, the device triggers of its own accord.**

Instead of the selection in the SETUP menu, a positive trigger level can also be selected in the ventilation menu:

- Press »↓« key to select the »TRIGGER LEVEL« line,
- press »-« key to select the value -1,
- hold down the »-« key for 6 seconds to set a value between +1 and +8 mbar,
- press »●« to confirm.

**When PEEP is in use, please note:**

- **Set the trigger level at least 3 mbar below the measured PEEP value to avoid self-triggering of the device and to enable the PEEP valve to be correctly set.**

**Vt MINI (MINIMAL VOLUME)**

Setting range: 0.05 to 1.2 L

In this setting the device monitors whether a specified minimal tidal volume is delivered. If this tidal volume is not achieved in 3 successive ventilation strokes, an alarm is activated and the display for the measured tidal volume flashes.


## SIMV

### Synchronized Intermittent Mandatory Ventilation

A "weaning" method for patients who can breathe spontaneously

The following parameters are displayed on screen in SIMV ventilation mode:

- the ventilation parameters necessary for SIMV
- alarm limit for high pressure alarm
- alarm messages in the event of an alarm
- pressure curve (if selected in the SETUP menu)
- graphic display of internal battery charge level (if powered by the internal battery)
- backup tidal volume  $V_{TS}$  MINI which the patient should breathe during spontaneous breathing
- bar display of airway pressure with superimposed alarm limit for high pressure alarm (black triangle)
- measured PEEP value below bar display of airway pressure
- measured max. airway pressure above bar display
- at the top of the pressure curve field the device displays the current tidal volume  $V_T$  and the current I/E ratio.

To set the parameters

- Press »↓« key or »↑« key to select the parameter. The little box to the left of the parameter turns black.
- Press »+« key or »-« key to set the value,
- press »●« to confirm.

### Minimum frequency (FREQUENCY)

Setting range: 4 to 20 1/min

The set frequency is a minimum frequency. The patient's efforts to breathe may trigger the mandatory ventilation stroke ahead of time and thus increase the frequency. The display shows the current frequency.

### Tidal volume $V_T$

Setting range: 0.1 to 1.6 L

### Inspiratory flow (FLOW INSP.)

Setting range: 0.4 to 1.5 L/s

**Backup rate (BACKUP RATE)**

Setting range: 6 to 22 1/min

To set a minimum frequency for ventilation in case of apnea.

- Set the backup rate at least 2 1/min higher than the minimum frequency for SIMV.

If this ventilation mode is automatically activated in the event of apnea, the message "VOL CTRL" appears and the VT display is shown in white against a black background during each stroke.

The device stops this "apnea ventilation" mode when spontaneous breathing is detected again.

**High pressure alarm (PAW HIGH)**

Setting range: 14 to 60 mbar

The set value is additionally displayed as a black triangle next to the bar display for airway pressure. If the ventilation pressure reaches the set value, the device activates a visual and an audible alarm. The maximum airway pressure is displayed as a figure above the bar display for airway pressure. The device updates this figure at each ventilation stroke.

**Inspiratory pressure support (PRESSURE)**

Setting range: OFF, 5 to 40 mbar

For inspiratory pressure support during spontaneous breathing between two mandatory strokes. Pressure support is activated by the trigger.

- Set a value of not less than 10 mbar for pressure support. Pressure support is relative to the trigger. A positive trigger value does not reduce the relative pressure support value.

**Trigger level**

Setting range:

-1, -2 mbar, OFF

Additionally, if TRIGGER LEVEL + was selected in the SETUP menu:

+1, +2, +3, +4, +5 mbar

-1 or +1 mbar = max. sensitivity

The set value is additionally displayed as a white triangle next to the bar display scale for airway pressure. Settings greater than -3 mbar cause the display "OFF" to appear.

- **Set a positive trigger level if a positive end-expiratory pressure PEEP is activated. If a positive trigger level is set without PEEP being active, the device triggers of its own accord.**

Instead of the selection in the SETUP menu, a positive trigger level can also be selected in the ventilation menu:

- Press »↓« key to select the »TRIGGER LEVEL« line,
- press »-« key to select the value -1,
- hold down the »-« key for 6 seconds to set a value between +1 and +5 mbar,
- press »●« to confirm.

The white triangle next to the bar display scale for airway pressure disappears.

**When PEEP is in use, please note:**

- **Set the trigger level at least 3 mbar below the measured PEEP value to avoid self-triggering of the device and to enable the PEEP valve to be correctly set.**

In combination with the positive trigger level set, the frequency alarm is automatically active during SIMV. This alarm indicates that the current frequency is 10 % higher than the set frequency. The default value 10 % can be set within a range between 5 % and 20 %.

During the spontaneous breathing phases, the device automatically sets the trigger level to -1 or +1 mbar.

**Apnea time delay (APNEA)**

Setting range: 5 to 40 s

To set a time delay when apnea occurs.

After this apnea time delay, the device ventilates with the set backup rate (BACKUP RATE) and the set tidal volume VT. The message "VOL CTRL" appears and the VT display appears white on a black background during every stroke.

A pressure drop to 3 mbar triggers an audible and a visual alarm 3 seconds after the apnea time delay.

### Backup tidal volume (VtS MINI)

Setting range: 0.05 to 1.5 L

The backup tidal volume indicates the minimum volume which the patient should deliver during spontaneous breathing. If the measured tidal volume is less than the backup tidal volume after 3 successive breaths, an intermittent alarm sounds and the display for the measured tidal volume flashes.


### Locking settings

The settings can be locked to prevent alteration:

- Press »↓« key to select »OFF« line.
- Press and hold down »+« and »-« keys at the same time for approximately 6 seconds until the key symbol  appears in the top left of the display.
- In this setting, only the OFF function in the bottom line can now be selected.

#### To release lock

- Press »↓« key to select »OFF« line.
- Press and hold down »+« and »-« keys at the same time for approximately 6 seconds until the key symbol  disappears.


02/23/356

### Operation with internal battery


If the external power supply fails, the internal battery automatically takes over the device power supply.

A battery symbol next to the alarm field shows the battery charge level. In the example: battery 50 % charged. The message "ALIM. EXT. FAIL" appears in the pressure curve field.

02/23/356

Possible battery charge displays:

The device activates an alarm on changing over to internal battery power, for instance in case of mains power failure.

- The red alarm lamp lights up, accompanied by an intermittent audible alarm.

This alarm can be suppressed:

- Press »**A**« key. The intermittent audible alarm is silent and the red alarm lamp remains lit.

If the "battery completely discharged" symbol appears, the device displays the alarm message »**BATTERY EMPTY**«. The remaining operating time is then only 5 to 15 minutes, depending on the age of the battery. The intermittent audible alarm sounds and in this state it cannot be suppressed.

- Connect to an external power supply immediately.

## Switching off the device

- First, disconnect the patient.
- Press »**↓**« key to select the »**OFF**« line,
- press and hold down »**●**« key for 3 seconds (a beep will sound 3 times),
- the message »**USER STOP**« appears in the display.
- Then turn the ON/OFF switch at the rear to  = OFF.

Switching off at the ON/OFF switch without following the steps described above is interpreted by the device as an unintentional stop and an alarm is activated.

To suppress this alarm:

- Turn the ON/OFF switch to  = ON and back again to  = OFF.

### In the event of an alarm

- 1 If an alarm occurs, the red alarm lamp lights, accompanied by an intermittent or continuous audible alarm.
- 2 The alarm message is displayed in the alarm field on the screen.
- 3 The audible alarm can be suppressed for 2 minutes by pressing the »A« key.


Only the "SAFETY STOP" alarm (internal device failure) cannot be suppressed, because of its significance.

Depending on the SET UP setting for "Store alarm", page 58, after the fault has been remedied either

- the audible alarm is muted and the message disappears from the display, or
- the audible alarm is muted but the message is stored.

To acknowledge the message:

- Press and hold down the »A« key for 3 seconds.
- To remedy device errors, refer to the Fault – Cause – Remedy list on page 78.

### Maintenance menu


The maintenance menu is intended for service purposes only.

With the device switched off:

- Press and hold down the »↓« and »↑« keys,
- turn the ON/OFF switch to = ON.
- Hold down the »↓« and »↑« keys until the SET UP menu appears on screen.

In the SET UP menu:

- Press »↓« to select the »MAINTENANCE« line and
- press »●« to confirm.

The maintenance menu is displayed:


This shows the following values:

- date and time
- operating hours (TOTAL HOUR)
- software version

#### To set

- Use »↓« and »↑« keys to select
- use »+« or »-« keys to set the value
- press »●« to confirm.

The »EVEVT DOWNLOAD« line takes the user into the service menu (for service purposes only).

#### Closing the maintenance menu

- Use »↓« and »↑« keys to select the »SET UP« line and press »●« to confirm.

### Display operating hours

When the device is switched on, the total number of operating hours (TOTAL HOUR) are displayed in the pressure curve field for approximately 10 seconds.

The TOTAL HOUR also appear in the maintenance menu.

To view the display during operation:

- Press »↓« to select the »OFF« line and
- Press »+«.

The operating hours are displayed for approximately 10 seconds.

Operating hours cannot be displayed whilst an alarm is occurring.

### Set contrast

During operation:

- Press and hold down »↓« and »↑« keys at the same time
- press »+« for more contrast or
- press »-« for less contrast.

## Maintenance intervals

**The device and its components should be disinfected and cleaned before any maintenance work is carried out. This includes when sending back to manufacturer for repair.**

Disposable bacterial filter for hose system	Replace immediately if resistance increases and at least every 800 hours' operation or every 60 days. The used filter can be disposed of as domestic waste.
Disposable bacterial filter for pressure measurement	Replace every 800 hours' operation or every 60 days. The used filter can be disposed of as domestic waste.
Disposable hose system with breathing valve	Replace for each new patient or if damaged. The old hose system can be disposed of as domestic waste.
Disposable microbial filter on the O <sub>2</sub> enrichment set	Replace filter <ul style="list-style-type: none"> <li>— every 500 hours in outdoor use</li> <li>— every 1200 hours in indoor use</li> </ul> The old filter can be disposed of as domestic waste.
Screen at bellows inlet	Replace if dirty. This should only be replaced by Dräger Medical Service or a Dräger Medical authorised dealer.
Internal battery	Replace every 3 years. This should only be replaced by Dräger Medical Service or a Dräger Medical authorised dealer.
Inspection and service	Annually, by authorised Dräger Medical Service personnel.

## Battery maintenance

### Charge internal battery

- Before using for the first time and
- if the device has not been used for about 15 days:
- Connect device to mains supply. There is no need to switch on the device.
- 1 The indicator lamp for mains power lights up green. The internal battery is being charged.
- Leave the device connected to the mains for about 17 hours. Even a battery which was completely discharged will then be fully charged.

Running the device on mains power keeps a constant trickle charge to the battery.

**The internal battery is not charged when operating on an external DC supply (12 to 36 V DC).**

### To ensure long battery life:

- Approximately every 2 months, run the device on the internal battery for 2 hours.
- This can be done by disconnecting from the mains supply. The device automatically switches over to internal battery power.
- Connect device to mains supply again after 2 hours.

The device operating concept is based on partial battery discharge. So:


- **Do not discharge the internal battery completely, as this will shorten its service life.**

If the device is not used for more than 2 months:

- Remove the internal battery to avoid any leakage. This should only be done by Dräger Medical Service.
- Observe battery storage conditions:  
The permitted storage time at the corresponding storage temperatures is shown in the table below.

Temperature °C	Storage time (months)
0 to 20	12
21 to 30	9
31 to 40	5
41 to 50	2.5

- Follow the manufacturer's instructions for maintenance of the (optional) external battery.

0163/356

### **Disposal of internal battery**

The battery in RespiCare VV contains toxic substances.  
It should be disposed of in accordance with the regulations applicable in the country of use.

### **Disposal of the device**

– at the end of its life

RespiCare VV should be properly disposed of as advised by an authorised disposal contractor and in accordance with the legal provisions in force.

## Fault – Cause – Remedy

The table below is a guide to identifying faults and remedying them.

The messages appear in combination with an intermittent or continuous audible alarm. The messages are listed below in alphabetical order.

Message	Cause	Remedy
<b>ALARM PAW HIGH</b> Intermittent audible alarm	Alarm limit for high pressure alarm is set too low	Correct setting
	Hollow needle blocked	Clean hollow needle
	Patient coughing	Switch off alarm
<b>ALARM PAW LOW</b> Continuous audible alarm	Alarm limit for low pressure alarm set too high	Correct setting
	Hose system or control line disconnected or damaged	Reconnect or replace hose system or control line
<b>ALIM EXT. FAIL</b> Intermittent audible alarm	No mains power	Check mains supply, if necessary use a different socket
	Power lead faulty	Replace power lead
Device is powered by the internal battery	No DC voltage	Connect DC supply: change from 12 – 36 V or external supply to mains power
	DC supply cable faulty	Replace cable
<b>APNEA</b> Continuous audible alarm	Hose system or pressure measuring hose is disconnected	Connect or replace hose system or pressure measuring hose
	No spontaneous breathing	Change over to CV mode
<b>BATTERY EMPTY</b> Intermittent audible alarm	Battery is discharged	Connect device to an external power supply, mains or DC
<b>EPROM FAIL</b> Continuous audible alarm	Internal device failure	The first time it occurs, switch off the alarm. If the alarm occurs again: Replace device and call Dräger Medical Service or a Dräger Medical authorised dealer
<b>FILLING</b> Continuous audible alarm	Air inlet is blocked	Clean air inlet. Replace filter, if fitted
<b>FREQUENCY ALARM</b> Intermittent audible alarm	Trigger level set too high	Correct setting
	Frequency alarm set too low	Correct setting

Message	Cause	Remedy
<b>MINIMAL VOLUME</b> Intermittent audible alarm	Alarm level for minimal volume too high	Correct setting
	Hollow needle blocked	Clean hollow needle
<b>OVERPRESSURE</b> Continuous audible alarm	Disposable bacterial filter for hose system is blocked	Replace disposable bacterial filter
<b>SAFETY STOP</b> Continuous audible alarm	Internal device failure	Replace device: call Dräger Medical Service or a Dräger Medical authorised dealer
<b>VOL CTRL</b> Intermittent audible alarm	VTS MINI setting too high	Correct the setting
	Disposable microbial filter in O2 enrichment set is dirty	Replace disposable microbial filter
	Hollow needle displaced	Clean the hollow needle
<b>SLIDING</b> Continuous audible alarm	Bellows motor inoperative	Replace device: call Dräger Medical Service or a Dräger Medical authorised dealer

If the fault cannot be remedied:

- Replace the device. Call Dräger Medical Service or a Dräger Medical authorised dealer.

Technical faults are displayed by a flashing number at the bottom right of the screen:

- 1 Fault in audible alarm
  - 2 Malfunction between keyboard and PCB
  - 3 Default setting not possible
- Notify the number displayed to Dräger Medical Service or the Dräger Medical authorised dealer. Switching the device off and on again will cause the fault number to disappear.

## Preparation for renewed use

### Dismantling

- **Disconnect mains plug to avoid risk of electric shock!**

- Detach the disposable hose system from the device and from the heater/humidifier, if installed, and dispose of as domestic waste.
- Detach O<sub>2</sub> monitor from the device.
- Remove disposable bacterial filter (if fitted) from the device.

**Do not allow any fluid to enter the filter. This would increase filter resistance and the operation of the device could be impaired!**

- Prepare the O<sub>2</sub> monitor, heater/humidifier and PEEP valve according to the separate instructions for use.

### For each new patient

The following should be carried out by authorised Dräger Medical Service personnel:

- Disinfect the bellows
- Replace all filters
- Dispose of disposable hose system and use a new one
- Disinfect other parts contaminated by patient breathing, such as the PEEP valve, according to the separate instructions for use.

### Cleaning

**Do not allow any fluid to enter the device or its connections. This could impair the operation of the device!**

### Device, cable and control line

- Remove coarse dirt with a disposable cloth.
- Wipe down device, cable and control line with a damp disposable cloth.
- Allow parts to dry thoroughly.

## Disinfecting

**Do not allow any fluid to enter the device or its connections. This could impair the operation of the device.**

### Device, O<sub>2</sub> enrichment set and cable

- Wipe-disinfect external surfaces using an aldehyde-based disinfectant.

For users in the Federal Republic of Germany we basically recommend the use of disinfectants included in the current DGHM list (DGHM: Deutsche Gesellschaft für Hygiene und Mikrobiologie, German Hygiene and Microbiology Association). The DGHM list (published by mhp-Verlag GmbH, Wiesbaden) also lists the active ingredient in each disinfectant.

In countries where the DGHM list is not available, we recommended the type of active ingredient stated above.

- Following manufacturer's directions for use.
- Allow parts to dry thoroughly.

### For each new patient


The following should be carried out by authorised Dräger Medical Service personnel:

- Sterilise the bellows
- Replace all filters
- Dispose of disposable hose system and use a new one. Disinfect other parts contaminated by patient breathing, such as the PEEP valve, according to the separate instructions for use.

## What's what

### Front view


- 1 Pressure measuring line
- 2 O<sub>2</sub> connection cap
- 3 O<sub>2</sub> connection angled nozzle
- 4 O<sub>2</sub> enrichment set with disposable bacterial filter
- 5 Hose system cone
- 6 Elbow connector
- 7 PEEP valve (optional)
- 8 Control line
- 9 Breathing valve, disposable item
- 10 Breathing valve control line nozzle
- 11 Disposable hose system

00837356

### Setting panel

- 12 Red alarm lamp
- 13 »A« key to suppress the audible alarm for 2 minutes or to acknowledge a message
- 14 »+« key to increase values
- 15 »-« key to decrease values
- 16 »●« key to confirm a setting
- 17 Mains power supply indicator lamp
- 18 DC supply indicator lamp
- 19 Display
- 20 »↓« key to select parameters (scroll up)
- 21 »↑« key to select parameters (scroll down)

00737356

**Rear view**

**22** Handle

**23** Trolley mounting point (x 4)


**24** RS 232 interface

**25** Remote alarm connection

**26** Device socket for DC  $\text{---}$

**27** Device socket for mains power  $\text{~}$  with fuse

**28** ON/OFF switch  $\text{●/○}$

## Technical Data

### Ambient conditions

In operation:

Temperature	1 to 30 °C
Air pressure	600 to 1060 hPa
Rel. humidity	30 to 70 %

In storage:

Temperature	0 to 60 °C
Air pressure	600 to 1060 hPa
Rel. humidity	10 to 85 % no dewing

### Performance characteristics:

Working principle Volume metering, time-cycled constant volume or pressure-controlled

Volumetric ventilation mode IPPV/IPPVASS

Frequency	8 to 50 1/min Adjustable in 1/min increments
Minute volume	3 to 25 L/min
Tidal volume	0.1 to 1.6 L Adjustable in 0.01 L increments
Flow ramp	Rectangular, sinusoidal, accelerating, decelerating
Ratio I/E	1:1 to 1:3
Trigger level	-3 to +5 mbar* Adjustable in 1 mbar increments
Sigh volume	VT x 1.0 to VT x 2.0 Adjustable in 0.1 increments
Sigh frequency	Between 50 and 950 strokes Adjustable in 50-stroke increments
Paw high alarm	12 to 60 mbar Adjustable in 1 mbar increments
Paw low alarm	2 to 50 mbar Adjustable in 1 mbar increments

Press-controlled ventilation mode IPPVPRS:

Minimum frequency	4 to 20 1/min Adjustable in 1/min increments
Flow insp.	0.4 to 1.5 L/s Adjustable in 0.1 L/s increments
Pressure insp.	7 to 40 mbar Adjustable in 1 mbar increments
Plateau	0 to 2 s Adjustable in 0.1 s increments

\* 1 mbar = 0,1 kPa

Trigger level	-3 to +8 mbar Adjustable in 1 mbar increments
Min. volume	0.1 to 1.6 L
Paw low alarm	12 to 50 mbar Adjustable in 1 mbar increments

Synchronized intermittent mandatory ventilation mode SIMV:

Frequency	4 to 20 1/min Adjustable in 1/min increments
Tidal volume	0.1 to 1.6 L Adjustable in 0.01 L increments
Backup rate	6 to 22 1/min Adjustable in 1/min increments
Flow	0.4 to 1.5 L/s Adjustable in 0.1 L/s increments
High pressure	14 to 60 mbar Adjustable in 1 mbar increments
Inspiratory pressure support	OFF, 5 to 40 mbar Adjustable in 1 mbar increments
Trigger level	-2 to +5 mbar Adjustable in 1 mbar increments
Apnea time delay	5 to 40 s Adjustable in 1s increments
Backup tidal volume	0.05 to 1.5 L Adjustable in 0.05 L increments

**Alarms**

Low pressure alarm	2 to 35 mbar Adjustable in 1 mbar increments
High pressure alarm	12 to 60 mbar Adjustable in 1 mbar increments
Mains power failure alarm	Alarm in the event of mains power failure or DC supply failure
Battery empty alarm	Alarm if the "battery empty" symbol appears in the display
Frequency alarm	NO, 5 %, 10 %, 15 %, 20 % Alarm if the current frequency is 5 % (10 %, 15 %, 20 %) higher than the set frequency
High pressure alarm	Alarm if the maximum safe pressure of 75 mbar is reached
Inflation volume alarm	Alarm if there is no air intake
Slip alarm	Alarm if the bellows motor fails
Safety stop alarm	Alarm in the event of device failure
EPROM inop.	Alarm in the event of an internal device error

Minimal volume alarm (in IPPVPRS mode and SIMV)	0.1 L to 1.6 L Alarm if the set tidal volume is not reached in 3 successive strokes
Apnea alarm (in SIMV)	Alarm if the device does not detect any spontaneous breathing
Lock settings	to prevent inadvertent altering of the ventilation parameters


**Device characteristics**

Compliance of bellows	0.5 L/mbar
Compliance of hose system supplied with device	0.9 mL/mbar
Device resistance with hose system:	
inspiratory at 60 L/min	5.5 mbar
expiratory at 60 L/min	3 mbar
Pressure display	-3 to 80 mbar
AC power supply	110 V to 230 V ± 10 %, 50/60 Hz
Fuse	F1.6H 250 V IEC 127-2/1
Power consumption	70 VA
DC supply	12 to 36 V DC
Fuse	5 A slow blow
Power consumption	36 VA
Internal battery	Comprises 2 lead-gel units, airtight, rechargeable from the AC power supply, designed for 2 hours' operation.
Protection class	Class I,  BF
Electromagnetic compatibility	Tested to EN 60601-1-2
Dimensions H x W x D	340 x 240 x 254 mm
Weight	12 kg
Sound level (corresponding to free field measurement over a reflective surface)	43 dB (A) in typical ventilation at a distance of 1 m
Classification according to EC Directive 93/42/EC Appendix IX	Class IIb
UMDNS code Universal Medical Device Nomenclature System	14-355

**Remote alarm**

Battery	9 V alkaline 6F22, quantity 1
Current consumption	
Standby	0.1 mA
Alarm	3.0 mA
Service life	approximately 1 year approximately 150 alarm hours
Dimensions	114 x 72 x 72 mm
Weight	150 g

## Symbols and abbreviations used

	Suppress audible alarm for 2 minutes
	Select parameters
	Set higher value
	Set lower value
	Confirm value
	ON/OFF switch ON
	ON/OFF switch OFF
	Setting locked
	Internal battery charge level (example: battery completely discharged)
	Caution, refer to Instructions for Use
	Protect interface connections at rear of device from electromagnetic discharge
	Remote alarm connection
	Mains power AC
	DC supply DC
AC 	Alternating current
DC 	Direct current
	Earthed conductor
	Safety class BF, bodyfloating
I/E	Inspiration : Expiration ratio
IPPV	Controlled ventilation
IPPVASS	Controlled assisting ventilation
IPPVPRS	Controlled pressure-controlled ventilation
PEEP	Positive end-expiratory pressure
R	Rectangular flow shape
S	Sinusoidal flow shape
A	Accelerating flow shape
D	Decelerating flow shape

SIMV	Synchronized intermittent mandatory ventilation
UMDNS	Universal Medical Device Nomenclature System
V <sub>T</sub>	Tidal volume
V <sub>TS</sub>	Backup tidal volume

## Operating principle


The ventilation principle is based on the inflation and deflation of a set of bellows which delivers a volume of air to the patient's lungs. The volume of air is controlled by directional valves during inflation and deflation.

The movement of the bellows is generated by a bellows drive actuated by a stepper motor.

The upper part of the bellows with the directional valves is fixed in position, whereas the lower part is free to move.

A magnet is attached to the mobile part. Magnetic position sensors are placed at the sides of the bellows for position determination.

- 1 Power supply
- 2 Setting panel
- 3 Electronic control
- 4 Pressure sensor
- 5 Bellows drive
- 6 Stepper motor
- 7 Bellows
- 8 Directional valve – inflation
- 9 Directional valve – deflation
- 10 Hose system
- 11 Upper magnetic position sensor
- 12 Lower magnetic position sensor
- 13 Breathing valve

### IPPV/IPPVASS mode

When the device is switched on, the bellows 7 are deflated via directional valve 9. Once the upper position sensor 11 is reached, the bellows 7 start to open and air is drawn in via directional valve 8 until the set tidal volume is reached.

#### Inspiration

During the inspiration phase (determined by the setting for frequency and I/E) the bellows are deflated.

The motor speed depends upon the selected tidal volume and flow shape. During the inspiration phase the tidal volume is delivered in the selected flow shape through directional valve 9, the hose system 10 and the breathing valve 13 to the patient's lungs. Directional valve 8 is closed. The expiration aperture of the breathing valve 13 is closed.

### Expiration

During the expiration phase the bellows **7** start to open again as far as the set tidal volume, and ambient air is taken in via directional valve **8**. The patient exhales into the ambient air through the open expiration aperture of the patient valve **13**. In ACV mode, the patient can trigger the next ventilation stroke by his efforts to breathe.

### IPPVPRS mode

When the device is switched on, the bellows **7** are inflated and ambient air is taken in via directional valve **8**. Directional valve **9** is closed. When the lower position is reached (max. volume 1.6 L), it remains on stand-by.


### Inspiration

The bellows are deflated either under the control of the minimal frequency set or triggered by the patient's efforts to breathe. The bellows are deflated and deliver a volume of air via directional valve **9** to the patient until the set inspiratory pressure is reached. Directional valve **8** is closed.

Inspiration can be extended by setting the plateau time.

### Expiration

In the expiration phase, the bellows **7** are inflated again as far as the lower position (max. volume 1.6 L). The pressure in the bellows **7** decreases, breathing valve **13** opens and the patient exhales through the open expiration aperture of the patient valve **13**.

### SIMV

When the device is switched on, the bellows **7** are inflated and ambient air is taken in via directional valve **8**. When the lower position is reached (max. volume 1.6 L), it remains on stand-by.

### Inspiration

At the start of a mandatory breathing stroke, either started by the minimal frequency set or triggered by the patient's efforts to breathe, the bellows are deflated until the set tidal volume has been delivered through directional valve **9** with the set inspiratory flow.

#### Spontaneous breathing phase

During the remaining cycle time, the bellows are deflated in rhythm with the spontaneous rate, stopped and travel back for gentle synchronisation with the patient's spontaneous breathing. An apnoea monitor activates an alarm if the spontaneous breathing necessary for the trigger is not present.

### **Maximum safe pressure**

#### **IPPV/IPPVASS mode**

When the set high pressure alarm is reached, bellows compression stops for the inspiration time. The bellows then return to the starting position.

#### **IPPVPrs mode**

The maximum inspiratory pressure is reached when inspiration stops and the system immediately changes over to the expiration phase (plateau time = 0 s).

#### **SIMV**

During the mandatory ventilation stroke, the backup pressure acts in the same way as in IPPV mode/IPPVASS mode.

During the spontaneous breathing phase, a very low pressure defines the stopping and return travel of the bellows to allow the patient to breathe out.

## Order List

Designation	Part No.	Designation	Part No.
<b>RespiCare VV basic device</b>		Fisher & Paykel type heater/humidifier MR 410 AGU (German)	57 30 494
RespiCare VV standard set	57 30 437	Fisher & Paykel type heater/humidifier HC 500 AEK (UK)	57 30 498
Euro power cable	57 30 042	Fisher & Paykel type heater/humidifier HC 500 AEU (Europe)	57 30 499
CH power cable	57 30 043	Fisher & Paykel type heater/humidifier HC 500 AGU (German)	57 30 500
Italy power cable	57 30 255	Water trap	57 30 503
UK power cable	57 30 256	Fisher & Paykel MR 370	57 30 504
Hose connection for pressure measurement	57 30 087	Accessories holder	57 30 081
Hose system with instructions for use, expiration valve	57 01 982		
O <sub>2</sub> enrichment set comprising:	57 30 078	<b>Consumables</b>	
– 5 filter discs		Hose system with breathing valve	57 01 982
– instructions for use		Fuse F1.6H 250 V IEC 127-2/I, set of 4	57 30 086
– Silicone cap		Disposable bacterial filter for hose system	57 30 084
– Disposable microbial filter, set of 5		Disposable microbial filter for O <sub>2</sub> enrichment set, set of 5	57 30 085
Silicone cap for O <sub>2</sub> enrichment set	57 30 088		
Fuse F1.6H 250 V IEC 127-2/I Set of 4	57 30 086	Instructions for Use de/en	57 01 981
<b>Accessories</b>			
Carrying case	57 30 074		
Cable for DC supply 12 to 36 V DC	57 30 075		
Trolley, base and pole (boxed)	57 30 076		
Swivel arm to hold the hose system (for fitting to a rail)	57 30 077		
O <sub>2</sub> enrichment set comprising:	57 30 078		
– 5 filter discs			
– instructions for use			
– Silicone cap			
– Disposable microbial filter, set of 5			
Remote alarm with connecting cable and cable holder	57 30 079		
O <sub>2</sub> monitor	57 30 080		
PEEP valve	84 07 475		
External battery	57 01 983		
Humidifier HME booster UK	57 01 990		
Humidifier HME booster EU	57 01 991		
T connector HME	57 01 992		
Fisher & Paykel type heater/humidifier MR 410 AEK (UK)	57 30 492		
Fisher & Paykel type heater/humidifier MR 410 AEU (Europe)	57 30 493		

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Directive 93/42/EEC  
concerning Medical Devices

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**57 01 981**  
90 37 356 - GA 5665.151 de/en  
© Dräger Medical  
3rd edition - October 2001  
Subject to alteration