TRACOE medical GmbH

Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany

Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200

info@tracoe.com, www.tracoe.com




TRACOE® *twist* plus Tracheostomy Tubes






TRACOE® *twist* plus Tracheostomiekannülen

**REF 311, REF 312, REF 313,
REF 314, REF 316, REF 888-316**

EN	Instructions for Use
DE	Gebrauchsanweisung
FR	Instructions d'utilisation
IT	Istruzioni d'uso
ES	Instrucciones de uso
PT	Instruções de uso
DA	Brugsanvisning
FI	Käyttöohje
NB	Bruksanvisning
HU	Használati utasítás
SV	Bruksanvisning
NL	Gebruiksaanwijzing
HR	Upute za uporabu
SL	Navodila za uporabo
CS	Návod k použití
RO	Instrucțiuni
TR	Kullanma talimatı
RU	инструкция по применению
EL	Οδηγίες χρήσης
KO	사용 설명서
ZH	用于新生儿
AR	دليل استعمال

	REF 311	REF 312	REF 313	REF 314	REF 316	REF 888- 316
	—	1	—	—	—	—
	1	—	—	—	—	—
	—	—	—	1	—	—
	—	—	1	—	—	—
	—	—	—	—	—	1
	—	—	—	—	1	—
	—	1	—	1	—	1
	2	1	2	1	2	1
	1	1	1	1	1	1
 ⑩	—	1	—	1	—	1
 ⑪	—	—	—	—	1	1
 ⑫	1	1	1	1	1	1

EN / Warning. Please read Instructions for Use · **DE** / Achtung. Gebrauchsanweisung beachten · **FR** / Attention. Respecter le mode d'emploi · **IT** / Attenzione. Consultare le istruzioni per l'uso · **ES** / Atención. Siga las instrucciones de uso · **PT** / Atenção. Ter em consideração as instruções de utilização. · **DA** / Vigtigt: Følg brugsanvisningen · **FI** / Huomio. Noudata käyttöohjetta · **NB** / OBS. Følg bruksanvisningen · **HU** / Figyelem. Olvassa el a használati útmutatót · **SV** / OBS! Se bruksanvisningen · **NL** / Let op. Rekening houden met de gebruiksaanwijzing · **HR** / Upozorenje. Pridržavajte se uputa za uporabu · **SL** / Pozor. Upoštevajte navodila za uporabo · **CS** / Pozor. Respektujte návod k použití · **RO** / Atenție. Respectați instrucțiunile de utilizare · **TR** / Dikkat. Kullanma talimatına dikkat ediniz · **RU** / Внимание! Соблюдать инструкцию по применению · **EL** / Προσοχή. Ακολουθήστε τις οδηγίες χρήσης · **KO** / 주의. 사용 설명서를 준수하십시오. · **ZH** / 注意. 注意使用说明书 · **AR** / تنبيه. يجب مراعاة دليل الاستعمال

EN / For single patient use only · **DE** / Nur für einen Patienten · **FR** / Pour un seul patient · **IT** / Prodotto destinato ad un solo paziente · **ES** / Sólo para un paciente · **PT** / Só para um paciente · **DA** / Må kun bruges til én patient · **FI** / Vain yhden potilaan käyttöön · **NB** / Kun for én pasient · **HU** / Csak egy betegnél alkalmazható · **SV** / Endast för en patient · **NL** / Slechts voor één patiënt bestemd · **HR** / Samo za jednog pacijenta · **SL** / Samo za enega bolnika · **CS** / Pouze pro jednoho pacienta · **RO** / Pentru un singur pacient · **TR** / Sadece tek bir hasta için · **RU** / Только для одного пациента · **EL** / Μόνο για έναν ασθενή · **KO** / 오직 일인 환자 용입니다 · **ZH** / 只供一位病人使用 · **AR** / تستعمل لمريض واحد فقط

EN / Item number · **DE** / Artikelnummer · **FR** / Référence · **IT** / Codice prodotto · **ES** / Número de artículo · **PT** / Número do artigo · **DA** / Artikelnummer · **FI** / Tuotenumero · **NB** / Artikkelnummer · **HU** / Cikkszám · **SV** / Katalognummer · **NL** / Artikelnummer · **HR** / Broj artikla · **SL** / Številka artikla · **CS** / Číslo výrobku · **RO** / Numărul de articol · **TR** / Ürün numarası · **RU** / Номер артикула · **EL** / Αριθμός είδους · **KO** / 품목 번호 · **ZH** / 产品编号 · **AR** / رقم السلعة

EN / Batch code · **DE** / Chargencode · **FR** / Code de lot · **IT** / Numero di lotto · **ES** / Código de lote · **PT** / Código do lote · **DA** / Batchkode · **FI** / Eräköodi · **NB** / Chargekode · **HU** / Tételkód · **SV** / Satsnummer · **NL** / Lotcode · **HR** / Oznaka serije · **SL** / Koda šarže · **CS** / Kód šarže · **RO** / Cod serie de producție · **TR** / Parti kodu · **RU** / Код партии · **EL** / Κωδικός παρτίδας · **KO** / 생산 단위 · **ZH** / 批次編碼 · **AR** / رمز الشحن

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Instructions for Use for TRACOE® *twist* plus Tracheostomy Tubes

NOTE:

- Please read the instructions for use carefully. They are part of the described product and must be available at all times. For your own safety and that of your patients, please observe the following safety notes.

- The illustrations to which the text refers can be found on the (foldout) image pages at the beginning of these instructions. The numbers and letters in brackets refer to the respective illustrations and product components of the tracheostomy tube. The employed symbols or icons are explained on pages 1 to 5.

1. General Information

CAUTION: The product may only be applied by physicians and other instructed personnel in charge of the patients' care! The physician in charge must determine the correct size and type.

Intended Purpose: This medical device is a tracheostomy tube.

Functional Description: The tracheostomy tube maintains a passage for the breathing air through the tracheostoma.

- The blocked cuff separates the upper airways from the lower airways so that the air cannot flow from the lungs into the mouth nor the other way around. This way, the patient can only breathe through the tube. In order to block the cuff, it is inflated with air via the inflation line. In its blocked state, the cuff seals the space between the trachea and the outer cannula wall. Furthermore it allows efficient ventilation with ventilation equipment and, at the same time, prevents subglottic secretion from getting into the lower airways.

- Fenestrated cannulas allow a part of the air, which is needed for speaking, to reach the upper airways through the fenestration. In addition, the airway resistance in the upper airways is reduced.

- For a more detailed functional description, please see the following text.

Maximum Useful Life of the Tracheostomy Tube:

- With cuff 29 - days,
- without cuff - 5 weeks,

from first use (see chapter 5 "General Precautions"). This maximum period of use also includes all times, during which the tracheostomy tube or inner cannula are not in use.

Product for Single-Patient Use: Multiple use by one and the same patient during the period of use is permitted.

Product Info Card: This product is supplied with a product info card with two detachable labels. The labels show product-specific data. The info card should be kept separately as it makes reordering easier. The detachable label can be attached, e.g. to the patient's record.

2. General Description

The product consists of an outer cannula (1) with a cuff (2) or without a cuff and an inner cannula with a 15 mm connector (3), which is already inserted in the outer cannula and is locked with a locking ring (4).

The cuff is inflated or deflated via the inflation line (2a).

The outer cannula is made of a radiopaque polyurethane and is attached to a neck flange (5) that can swivel on 2 axes (cardanic suspension).

The perforated obturator (6) is already inserted in the inner cannula and allows the use of a Seldinger wire ($\varnothing 1.27 \pm 0.04$ mm, e.g. REF 518) for cannulation.

3. Indications

The products are suitable for all patients where access to the airways via a tracheostoma is required.

- Tubes without a cuff are particularly suited for patients who breathe spontaneously and need to use a tube.
- Tubes with a cuff are suitable for all patients where access to the airways via a tracheostoma with sealing of the trachea is required.
- Fenestrated tubes make it easier for patients with an intact larynx to speak.
- Tubes with a suction line are used for patients where the suction of subglottic secretion is indicated.

4. Contraindications

4.1 Absolute contraindications

- The occlusion cap (10) / speaking valve must not be used for laryngectomised patients (patients without a larynx) - danger of suffocation!

4.2 Relative contraindications

Relative contraindications are cases where the risks and benefits of the procedure must be carefully evaluated:

- Applications in paediatrics.

5. General Precautions

- It is strongly recommend that a ready-to-use spare tube and several spare inner cannulas are kept at the patient's bedside. These must be stored under clean and dry conditions.

- The tube must be checked for integrity and correct functioning before each application or insertion. This check includes the free lumen, cuff leakage test, correct and stable fit of the inner cannula inside the outer cannula, absence of kinks, a stable connection between the tube and neck flange, etc. The cuff material must not be brittle. If this product is damaged, it must be replaced by a new product.

- No force may be used on the tracheostomy tube as, otherwise, there is danger of damaging it (e.g. breaking it). If the connections at the 15 mm connector are overly tight or jammed, always use a disconnecting wedge that is approved for tracheostomy tubes.

- The 15 mm connector (3) should be kept clean and dry.

- When changing the inner cannula, always make sure that the inflation line (2a) of the cuff is not positioned between the inner and outer cannulas because it can get jammed and damaged.

- When repositioning the patient in bed, make sure that the patient does not lie on the pilot balloon (2a), as this would lead to a strong increase of the cuff pressure and can possibly damage the trachea.

- During mechanical ventilation and frequent change of the patient's position or manipulation of the tube, the inner cannula may become separated from the outer cannula.

- The use of tracheostomy tubes can lead to pressure points, necrosis of the skin on the neck and skin irritation (e.g. moisture). In order to avoid this, we recommend that a dressing is placed underneath the neck flange (e.g. TRACOE purofoam REF 958 or REF 959 or TRACOE softpad REF 969).

- In order to avoid damage to the material, the cuff should not be in contact with lidocaine-containing aerosols or any ointments.

- Thin cuffs have a certain permeability for water vapour. Therefore, condensate may accumulate in the cuff. For small quantities, this is irrelevant. However, if bigger quantities of condensate get suctioned into the inflation line by accident, the cuff pressure can no longer be measured and adjusted correctly and the tube needs to be replaced. Before the tube is removed, as much as possible of the air or water that is contained in the cuff must be removed using a syringe.
- For the joint application of further medical devices, the respective instructions for use must be observed. If in doubt, please contact the manufacturer.

6. Warnings

- Do not use if the sterile packaging is damaged.
- Use the locking ring (4) for loosening the inner cannula (C). The 15 mm connector must not be turned because the inner cannula inside the outer cannula may rotate and the air supply may be interrupted unnoticed.
- Never use fenestrated inner cannulas for ventilation.
- Only use the occlusion cap with fenestrated tubes (outer + inner cannula) and an unblocked cuff.
- When using laser or electrosurgical treatment equipment, make sure that there is sufficient distance to the tracheostomy tube. There is a risk of fire, also toxic gases may form and the tube may get damaged.
- Select the tube so that the fenestration (if applicable) is positioned at sufficient distance to the stoma canal. If the above warning is not observed, there is a risk of emphysema for ventilated patients. Also, granulation tissue may form or the airway resistance may increase if speaking valves or occlusion caps are used.
- If a fenestrated tube is used, there is a risk of increased formation of granulation tissue.
- The cuff pressure can change, for example due to altitude changes (e.g. in an aircraft), or if laughing gas is used as an anaesthetic or during the connection and disconnection of a hand-held manometer.
- Excessive cuff pressure poses a risk of permanent damage to the trachea.
- If the cuff pressure is insufficient, there is a risk of aspiration.
- Make sure that the correct Luer connectors are used (e.g. 2a + 9).
- During the insertion and removal of the tube, irritation, a need to cough or bleeding may occur.

7. Adverse Side Effects

Pressure points, necroses, skin irritation, granulation tissue, need to cough and bleeding.

8. Inserting the Tube

8.1 Preparing the Tube

1. Check the pack content for completeness (D).
2. For tubes with a cuff, check the cuff (2) for leakage by inflating it. For this purpose, inflate the cuff with a hand-held manometer, for example the TRACOE cuff pressure monitor REF 720, to a pressure of 50 cm H₂O ($\approx$ 36.78 mm Hg) and watch for 1 minute to see whether the cuff becomes deflated. If the cuff is leak-tight, remove all of the air from the cuff using a syringe. Push the cuff upwards in the direction of the neck flange (5). This makes the insertion of the tube easier.

8.2 Preparing the Patient

- The patient should be well pre-oxygenated before the cannulation or recannulation.
- Slightly overstretch the patient's neck if possible in order to make the insertion of the tube easier.
- In the event that the cannulation or recannulation leads to complications, safety measures should be taken (e.g. tracheal spreader, tube with a smaller diameter), which allow the physician to carry out short-term ventilation via translaryngeal intubation or a larynx mask.

8.3 Inserting the Tube

Carry out the following steps:

1. When inserting the tube, make sure that the obturator is not pushed back into the tube by accident. This can easily be avoided by holding both the tube and obturator with one hand. To aid insertion, a thin layer of water-soluble lubricant may be applied to the protruding part of the obturator at the patient end and to this area of the tube including the cuff.
2. After the tube has been inserted into the patient's trachea, the obturator has to be removed immediately.
3. For the ventilation, a connection between the standardised 15 mm connector (3) of the inner cannula and the ventilation machine needs to be established.
4. Inflate the cuff via the pilot balloon (2a). Adjust the cuff pressure to the individual ventilation therapy and check it at regular intervals. Typically, the pressure should be between 20 cm H₂O ($\approx$ 15 mm Hg) and 30 cm H₂O ($\approx$ 22 mm Hg).
5. The neck strap (12) is attached to the neck flange (5) to

keep the tube in place.

6. Furthermore, the position of the tube has to be checked as well as its correct function (see chapter 6 "Warnings"). If required, trachea and tube should be *carefully* cleared by means of suctioning in order to ensure that the airways are free.

7. Now check the cuff pressure in order to make sure that the cuff has not been damaged during the insertion.

8.4 Removing the Tube

The following preparations need to be made to remove the tube:

- Bend the patient's head slightly backwards and
- completely unblock the cuff of cuffed tubes (see chapter 8.5).

8.5 Unblocking the Cuff

Before the cuff is unblocked, certain preparations need to be made to ensure that as little secretion as possible gets into the bronchia.

- This can be done by subglottic suctioning.
- For awake patients, we recommend that the secretion is removed during the unblocking process with the help of a suction catheter that is inserted through the tube.

• If tubes with subglottic suction are used, see chapter 9.5. For unblocking, connect a syringe to the pilot balloon (**2a**). Completely remove all air / condensate if applicable. (See chapter 5. "General Precautions").

9. Handling

9.1 Changing the Inner Cannulas

If viscous secretion collects in the inner cannula, which cannot be removed by means of suction and obstructs the airflow, the inner cannula needs to be replaced with a new or cleaned inner cannula.

By turning (**C**) the locking ring (**4**) the inner cannula is loosened and can be removed. After a new inner cannula has been inserted into the outer cannula, the inner cannula is locked by turning (**C**) the locking ring (**4**).

9.2 Fenestrated Tracheostomy Tubes

Unfenestrated inner cannulas have white locking rings and are used, amongst other things, during the ventilation of a patient.

To allow the patient to speak, a **fenestrated inner cannula** (blue locking ring) is inserted in **the fenestrated outer cannula**. After the cuff has been unblocked, a speaking valve can be placed on the 15 mm connector of this inner cannula. The instructions for

use of the respective speaking valve must be observed.

9.3 Occlusion Cap of the Fenestrated Tube for Weaning

The fenestrated inner cannula can be closed by attaching the enclosed occlusion cap (10) to the 15 mm connector.

CAUTION: When preparing the tube, it has to be ensured that the upper airways of the patient are clear. The upper airways are cleared by coughing up the secretion or removing it by means of suction. The cuff needs to be unblocked. When closing the tube, the patient's breathing and vital signs must be monitored by a nursing professional. If any signs of breathing difficulties are detected, the occlusion cap must be removed immediately.

9.4 Keeping the Fenestration Open

If the tube is used over a long period of time, it has to be ensured at regular intervals that the fenestration is not blocked by secretion, incrustation or in-growing tissue. If required, the tube may need to be changed.

9.5 Application of REF 316, REF 888-316 with Subglottic Suction Line

These tubes are fitted with a flat suction line (9) that is fixed to the outer bend of the outer cannula and ends with two openings immediately above the cuff (2). The suction line has a female Luer connector, through which any secretion can be removed using a syringe. As an alternative, a special suction device with a vacuum controller in connection with the enclosed connectors (11) can be used. After the suction process has been completed, the Luer connector needs to be closed.

CAUTION:

- During suction it has to be ensured that no excessive vacuum is used for a prolonged period of time (maximum of 200 mbar).
- In order to prevent the subglottic area from drying out, we recommend the use of intermittent suction.
- The suction might get blocked due to accumulated secretion or the adhesion of the tracheal mucosa. If rinsing the suction line is planned (e.g. with a small amount of air or physiological saline solution), it has to be ensured that the cuff is sufficiently blocked (risk of aspiration).
- The suction line may cause pressure points in the area of the tracheostoma or may cause the tracheostoma to become irregular. In such cases, the physician needs to decide whether this type of tube can still be used.

10. Cleaning, Disinfection and Storage

10.1 Cleaning

The tube (1 + 4) needs to be cleaned before each repeated use. Lukewarm water of drinking water quality can be used for cleaning. TRACOE medical recommends the use of the cleaning products offered by TRACOE for cleaning. After the tubes have been cleaned, they must be rinsed with water of drinking water quality and air-dried. When cleaning the outer cannulas, make sure not to damage the cuff.

CAUTION:

- Aggressive domestic cleaning agents, cleaning agents for dentures or solvents such as agents with a high alcohol content must never be used for cleaning as these may have a lasting adverse effect on the function.
- The tubes should not be heated to temperatures above 65 °C, since then the safety of the products cannot be guaranteed anymore.

10.2 Disinfection

Should disinfection appear necessary TRACOE medical can provide, on request, a list of tested disinfectants. This list is also available at www.tracoe.com. The use of unsuitable disinfectants can result in damage to the product.

After the tubes have been disinfected, they must be rinsed with sterile water and air-dried.

10.3 Storage

The clean tube has to be stored in a dry place.

11. Product Modifications

Modifications to the TRACOE products may only be carried out by employees or agents of TRACOE medical GmbH.

12. Disposal

These products may only be disposed of in accordance with the valid national regulations governing waste materials.

13. Returns

Returned used products will only be accepted subject to previous request and if a completed decontamination certificate is enclosed. This form is available either directly from TRACOE medical or via the website www.tracoe.com

14. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be effected exclusively on the basis of the valid General Terms and Conditions, which are available either from TRACOE medical GmbH or on the website www.tracoe.com.

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