







HEADstrap

HEADstrap for JOYCE/JOYCE Full Face WM 26360
HEADstrap with EasyClip WM 26365




Instructions for Use



JOYCE / JOYCE Full Face



SOMNOmask / SOMNOplus



Wöchentlich, nur Handwäsche!

Weekly, hand wash only!

Toutes les semaines,
lavage à la main uniquement!

Wekelijks, alleen handwas!

Settimanalmente, solo lavaggio
a mano!

Ugentligt, skal håndvaskes!

Ukentlig, kun håndvask

Viikoittain, vain käsinpesu!

Semanalmente, sólo lavado a mano!

En gång i veckan, bara handtvätt!

Gebrauchsdauer: 12 Monate

Working life: 12 months

Durée
d'utilisation: 12 mois

Gebruiksduur: 12 maanden

Durata utile: 12 mesi

Kan anvendes: 12 måneder

Brukstid: 12 måneder

Käyttöikä: 12 kuukautta

Vida útil: 12 meses

Livsläng: 12 månader

Manufactured by:

Weinmann

Geräte für Medizin GmbH+Co. KG
P.O. Box 54 02 68 · 22502 Hamburg
Germany

Phone +49/40/5 47 02-0

Fax +49/40/54 70 24 61

E-mail int.sales@weinmann.de

Internet www.weinmann.de

WM 16984a - 07.08